

# A NEW LC-MS APPROACH FOR SYNTHETIC PEPTIDE CHARACTERIZATION AND IMPURITY PROFILING

Waters

THE SCIENCE OF WHAT'S POSSIBLE.™


Nilini Ranbaduge, Ying Qing Yu, and Asish Chakraborty  
Waters Corporation, Milford, MA 01757

To learn more about synthetic peptide analysis, visit [www.waters.com/syntheticpeptides](http://www.waters.com/syntheticpeptides)

## INTRODUCTION

- Peptide therapeutics are an emerging group of pharmaceuticals applicable to a wide range of medical challenges.
- These peptide drugs exhibit relatively low toxicity and high biological activity compared to most conventional drugs.<sup>1,2</sup>
- Chemical synthesis is one of the most common method of production for synthetic peptides
- While the low toxicity of these drug products makes them more appealing as therapeutics, the impurities introduced during chemical synthesis of these peptides demand thorough characterization procedures to maintain efficacy and safety of the drug.
- The conventional synthetic peptide impurity profiling methods are mostly LC optical based assays relying on chromatographic separation of impurities.
- Even at optimal chromatographic performance, obtaining baseline resolved peaks for optical detection of low abundance impurities is a challenge. Adding HRMS to the analysis resolves this by improving the detection and identification of low level impurities based on the accurate mass and MS/MS information.
- In this study, we have developed a single LC-HRMS based analytical workflow for characterization and impurity profile monitoring of synthetic peptide drugs.

## METHODS



## RESULTS





Figure 3. A) TUV, TIC and zoomed-in TIC spectra are shown for the Eledoisin API peak. B) The API peak has co-eluting impurity peptides detected by high-resolution mass spectrometry. The most abundant peaks are shown on the spectrum.



## Screening Workflow: Synthetic Peptide/Impurity profiling

Impurity profiling using UNIFI screening workflow

Component name	Identification status	Neutral mass (Da)	Observed m/z	Mass error (ppm)	Observed RT (min)	RT Check	UV-Purity	% MS-Purity	% Impurity/API
1 API Pyroglutamic Acid E N-TERM, Amidation C-TERM	Identified	1187.6068	594.8055	1.4	15.41	Pass	94.705	64.155	100.000
2 EPKDAFIGLM +Fmoc -Ala -Gly -Met	Identified	1170.5730	586.2935	3.2	15.41	Pass	26.437	41.209	20.000
3 EPKDAFIGLM +Fmoc +Glu -Lys +t-Butyl -Phe -Met	Identified	1209.8052	605.7996	3.4	15.41	Pass	8.024	12.508	7.737
4 Pyroglutamic Acid E N-TERM, Amidation C-TERM, +Ser	Identified	1274.63291	638.3243	0.9	14.15	Pass	0.604	0.473	0.517
5 EPKDAFIGLM Pyroglutamic Acid E N-TERM, Amidation C-TERM, -Ala	Identified	1116.56777	559.2896	0.8	14.92	Pass	0.577	0.401	0.425
6 EPKDAFIGLM Pyroglutamic Acid E N-TERM, Amidation C-TERM, Oxidation M	Identified	1203.59579	602.8959	0.1	5.11	Pass	3.885	0.215	0.335
7 EPKDAFIGLM Pyroglutamic Acid E N-TERM, Amidation C-TERM, -Gly	Identified	1130.57942	566.2977	1.3	14.92	Pass	0.260	0.209	0.325
8 EPKDAFIGLM Pyroglutamic Acid E N-TERM, Amidation C-TERM, -Pro	Identified	1090.54812	546.2814	0.1	12.28	Pass	0.100	0.049	0.077
9 EPKDAFIGLM Amidation C-TERM	Identified	1205.61144	603.8130	0.0	10.75	Pass	0.032	0.051	0.051
10 EPKDAFIGLM Pyroglutamic Acid E N-TERM, Amidation C-TERM, +Lys	Identified	1315.69584	658.8540	-1.8	3.87	Pass	0.069	0.006	0.009



## CONCLUSION

The analytical workflow here demonstrates how a UPLC-QToF platform controlled by UNIFI Scientific Information System can be used to identify the API and its impurities of a synthetic therapeutic peptide sample. The key benefits of the workflow are:

- The instrument-informatics platform is a compliant-ready system that can be validated and implemented in regulated laboratory environments.
- The automated peptide mapping workflow identifies impurities based on accurate mass and validates the assignments using MS/MS data.
- The workflow can be used to characterize API and impurity sequences that are linear or cyclic.
- The scientific library can be used to a built custom impurity library for any selected peptide sequence.
- The screening workflow is utilized in robust high throughput profiling therapeutic peptide impurities.

## References:

- Mason, J. M., Design and development of peptides and peptide mimetics as antagonists for therapeutic intervention. *Future Medicinal Chemistry* 2010, 2 (12), 1813-1822.
- Pernot, M.; Vandersesse, R.; Frochet, C.; Guillemin, F.; Barberi-Heyob, M., Stability of peptides and therapeutic success in cancer. *Expert Opinion on Drug Metabolism & Toxicology* 2011, 7 (7), 793-802.

## METHODS

**Sample preparation**  
The synthetic peptide *Eledoisin* (pEPKDAFIGLM-NH<sub>2</sub>) was purchased from New England Peptides Ins (Gardner, MS) and *Calcitonin* (*Salmon*) (CSNLTCLVGLKLSQELHKLQTYPRNTGSGTP-NH<sub>2</sub>) was purchased from Bachem (Torrance, CA). A 2 µg/µL stock solution in water was further diluted to a final concentration of 0.2 µg/µL for the analysis.

**LC condition**  
LC system: ACQUITY UPLC H-Class Bio System  
Column: ACQUITY UPLC Peptide CSH C<sub>18</sub> 130 Å, 1.7 µm, 2.1 mm x 100 mm  
Column temperature: 65 °C  
Mobile phase: A: 0.1% FA in H<sub>2</sub>O and B: 0.1% FA in Acetonitrile  
Gradient: Eledoisin: 16-24% B over 30 min, Calcitonin: 14-34% B over 20 min

**MS condition**  
MS instrumentation: *Vion IMS QTof mass spectrometer*  
Capillary voltage: 2.8 kV  
Cone voltage: 50 °C, Source offset: 50 °C, Source temperature: 80 °C  
Desolvation Temp 300 °C, Cone gas: 20 L/hr, Desolvation gas: 500 L/hr

**Informatics**  
UNIFI Scientific Information System 1.9.2 (Peptide mapping, screening workflows)




Figure 1. The schematic shows the analytical workflow used in synthetic peptide impurity characterization and profiling. The UNIFI software platform controls data acquisition, processing and reporting.

## Validation of API/impurity using MS/MS




Figure 4. The fragmentation spectra for A) Eledoisin API sequence containing two modifications; N-terminal pyroglutamic acid and C-terminal amidation. B) fragmentation spectrum of oxidized Eledoisin API sequence C) fragmentation spectrum of Calcitonin API peak containing a disulfide bond (in orange). The UNIFI informatics platform can identify, label and verify fragment ions for each peptide sequence.

## Automated data processing and limit checks with UNIFI peptide mapping workflow

Component name	Protein name	Fragment label	Peptide	Modifiers	Observed m/z	Charge	Observed mass (Da)	Mass error (ppm)	Observed RT (min)	Matched 1st Gen Primary Ions	MS Response	Pass/Fail
1	Eledoisin	1-F1-11R	EPKDAFIGLM	Pyroglutamic Acid E N-TERM [1], Amidation C-TERM [11]	594.8077	2	1188.6081	0.0	16.21	1273189920	Pass	
2	Eledoisin	1-F1-11R	EPKDAFIGLM	Pyroglutamic Acid E N-TERM [1]	595.2987	2	1189.5902	-1.6	18.74	3639293	Fail	
3	Eledoisin	1-F1-11R	EPKDAFIGLM	Pyroglutamic Acid E N-TERM [1], Amidation C-TERM [1], Oxidation M [11]	602.8049	2	1204.6025	-0.5	5.65	1644417	Fail	
4	Eledoisin	1-F1-11R	EPKDAFIGLM	Pyroglutamic Acid E N-TERM [1], Amidation C-TERM [1], -Ala [6]	559.2889	2	1117.5705	-0.5	15.86	2995809	Fail	
5	Eledoisin	1-F1-11R	EPKDAFIGLM	Pyroglutamic Acid E N-TERM [1], Amidation C-TERM [1], -Gly [9]	566.2969	2	1131.5865	-0.2	15.47	1760597	Fail	
6	Eledoisin	1-F6-11R	AFGLM	Amidation C-TERM [6], -Ala [1]	579.3327	1	579.3327	0.6	16.22	781641	Fail	
7	Eledoisin	1-F1-11R	EPKDAFIGLM	Pyroglutamic Acid E N-TERM [1], Amidation C-TERM [1], -Ser [3]	551.2911	2	1101.5749	-1.1	17.38	1476992	Fail	
8	Eledoisin	1-F1-11R	EPKDAFIGLM	Pyroglutamic Acid E N-TERM [1], Amidation C-TERM [1], -Ser [3], Oxidation...	559.2889	2	1117.5706	-0.4	17.36	43053	Pass	
9	Eledoisin	1-F1-11R	EPKDAFIGLM	Pyroglutamic Acid E N-TERM [1], Amidation C-TERM [1]	594.8068	2	1188.6063	-1.6	17.56	724457	Fail	

Figure 5. The UNIFI review panel displays API and impurities identified. The MS response is used to determine the pass/fail status based on the relative% response compared to the MS response of the API.

TO DOWNLOAD A COPY OF THIS POSTER, VISIT [WWW.WATERS.COM/POSTERS](http://WWW.WATERS.COM/POSTERS)



## INTRODUCTION

- Peptide therapeutics are an emerging group of pharmaceuticals applicable to a wide range of medical challenges.
- These peptide drugs exhibit relatively low toxicity and high biological activity compared to most conventional drugs.<sup>1,2</sup>
- Chemical synthesis is one of the most common method of production for synthetic peptides
- While the low toxicity of these drug products makes them more appealing as therapeutics, the impurities introduced during chemical synthesis of these peptides demand through characterization procedures to maintain efficacy and safety of the drug.
- The conventional synthetic peptide impurity profiling methods are mostly LC-optical-based assays relying on chromatographic separation of impurities.
- Even at optimal chromatographic performance, obtaining baseline resolved peaks for optical detection of low abundance impurities is a challenge. Adding HRMS to the analysis resolves this by improving the detection and identification of low level impurities based on the accurate mass and MS/MS information.
- LC-MS based monitoring of impurity profile provides both accurate mass and relative abundance information that benefits process development and quality control.
- In this study, we have developed LC-HRMS-based analytical workflows for characterization and impurity profile monitoring of synthetic peptide drugs.