

From the Leader in Spectral Data

Checklist

Checklist	Status	Fix?
Noise	●	<input type="checkbox"/>
Contaminants	●	<input type="checkbox"/>
Technique	●	<input type="checkbox"/> ATR-IR

Optimized Corrections

Baseline of Query, Vertical Clipping, Vertical Offset

New Search

Search Status

1-Component Results: Top Hit: 78.1%

2-Component Results: Top Hit: 95.5%

3-Component Results: Click to Continue Searching

Classifications: Top Hit: 99.9%

Peak Results: Top Hit: 59.3%

Mixture of Two Steroids - ATR

Score | Info | Weight | Name | Chemical Structure | Spectrum

1-Component Results	2-Component Results	Classifications	Peak Results	Functional Groups
1	95.48	N.A.	Composite Spectrum	
	<input type="checkbox"/> 0.62	Ethisterone		
	<input type="checkbox"/> 0.38	Epiandrosterone		
		N.A.	Residual Spectrum	

KnowItAll Analytical Edition

Solutions to Accelerate IR, MS, NMR, Raman, & UV-Vis Analysis

Powerful Software. Quality Data. Results You Can Rely On

WILEY

Whether you use one or more techniques, KnowItAll has the right solutions for your lab!

Wiley's KnowItAll Analytical Edition offers solutions to **identify, analyze, and manage analytical data.**

This vendor neutral environment **supports multiple instrument vendor file formats** and techniques to streamline your laboratory workflow and work the way you do.

KnowItAll eliminates the need for multiple software packages with powerful tools integrated into a **single, easy-to-use interface.** We continually add spectral intelligence to our software, which also includes patented tools not available in other packages.

Combined with the **world's largest spectral reference databases***—including the renowned Sadtler™ libraries and spectra from trusted partners—KnowItAll gives chemists the most advanced technology available for fast, accurate analysis!



Accessibility features like keyboard access to menus, audio narration for icons, and tool tips.

Available in English, Japanese, Chinese, French, and German

KnowItAll and ChemWindow are trademarks of Wiley in certain jurisdictions.

**Subscription required to KnowItAll Spectral Libraries.*

Key Features

Basic Spectrum Analysis

Advanced Spectrum Search & Mixture Analysis

Database Building / Management

Structure Drawing & Reporting (ChemWindow)

Spectrum Processing

Comprehensive KnowItAll IR, MS, NMR, Raman, UV-Vis spectral library subscriptions available*

Access additional technique-specific tools such as: IR/Raman functional group analysis, MS reverse & adaptive search (patent pending), automated MS deconvolution/analysis, NMR prediction & much more! Plus 12K bonus IR spectra.

Data Types

IR, MS, NIR, NMR, Raman, UV-Vis

Chromatograms

Structures

By combining all the tools and data you need into one system, you have greater ability to extract meaningful knowledge from your data.

The KnowItAll interface is designed so the user can transfer information from one tool to another, and move from one task to the next, without having to leave the main interface or open another program. Multiple tasks are performed using logically grouped "toolboxes."

Because all the tools are located in a single, integrated environment, using this system will invariably **save time, improve workflow, and increase your ability to reach conclusions** from your data.

Powerful tools integrated into a single, easy-to-use interface. Customizable toolboxes.

Move seamlessly between tasks: search, process, manage data, draw structures, and more!

The screenshot displays the KnowItAll Informatics System interface. The main window shows a 1H NMR spectrum for 'DEMOX #10; Phenol' with peaks at 6.31, 6.83, 6.90, and 7.18 ppm. The interface includes a menu bar, a toolbar with various analysis tools, and a sidebar with icons for ID Expert, SearchIt, Minelt/Create Database, QC Expert, and AssignIt NMR. A table below the spectrum lists search results for Phenol, Poly(styrene), Cyclohexane, Benzene, and Carbon Tetrachloride. The Phenol entry is highlighted, and its chemical structure is shown. To the right, a 'Structure/Properties' panel displays the chemical structure of Phenol and a table of its properties.

Substructs		Sel. Substructs		Original Data Files	
All Properties		Attachments		Preferred Properties	
Name	Value				
Name	Phenol				
Boiling Point	181.8 °C				
CAS Registry Number	108-95-2				
Comments	Used in manufacturing many industrial compounds such as phenol-formaldehyde resins, bisphenol A, alkylphenols and certain dyes. Somewhat soluble in water; very soluble in alcohol, chloroform, ether, carbon disulfide. Highly toxic and caustic. A general disinfectant. NIOSH= SJ33250				
Density	20C=1.0767; 25C=1.132 G/ML				
Dielectric Constant	9.78 (60C)				
Flash Point	175F (CC)				
Formula	C6 H6 O				

What's Included?

Snapshot of the tools and applications included.

Data Toolbox

ID Expert™	One-click "first-pass" spectral identification tool
SearchIt™	Advanced database searching
Minelt™	Spectral database building and data mining
QC Expert™	QC spectrum comparison
AssignIt™ NMR	Add assignments to databases

Spectral Analysis Toolbox

AnalyzeIt™	Functional group analysis - IR, Raman, IR Polymer, Vapor Phase IR Synthetic Cannabinoids
PredictIt™ NMR	NMR chemical shift prediction

Spectral Processing Toolbox

ProcessIt™	Spectrum processing for IR, MS, NMR, Raman, MS
MS Expert™	Combines GC-MS deconvolution with automatic reference database search to analyze components

Basics Toolbox

ChemWindow	2D structure drawing
ReportIt™	Publish professional reports, with structures, spectra, chromatograms, etc.
BrowseIt™	Web portal with links to KnowItAll training resources and product news

Spectral Libraries*

KnowItAll Libraries	Subscribe to the world's largest collection of IR, MS, NMR, Raman, UV-Vis
---------------------	---

*Subscription required to KnowItAll Spectral Libraries.

Explore the KnowItAll Tools In-Depth

- Data Toolbox
- Spectral Analysis Toolbox
- Spectral Processing Toolbox
- Basics Toolbox

Data Toolbox



First Pass Spectral Identification Tool

When it comes to identifying unknown spectra, it's difficult to figure out where to begin. Wiley's KnowItAll ID Expert offers both novices and experts the perfect place to start. It provides fast, reliable answers to scientists identifying unknown spectra by matching against reference spectra.


The spectral intelligence built into KnowItAll ID Expert when combined with Wiley's high-quality KnowItAll Spectral Libraries*—the world's largest collection—makes this a **quick first pass analysis tool**.

How Does It Work?

- The user simply opens an unknown spectrum and KnowItAll ID Expert automatically performs a series of basic analyses—**single and multi-component search, peak search, functional group (if IR/Raman)**—and summarizes the results to give a complete overview of possibilities.
- If there are problems with the unknown spectrum or the reference spectra, ID Expert has the spectral intelligence to identify and fix some of these issues.
- Once the user has identified the unknown spectrum, a PDF report can be generated with a single click.

Includes Patented Optimized Corrections Technology to ensure best search results possible.

Also includes quick classification tool for IR spectra of designer drugs using Wiley's new IR Designer Drug Models.



*Subscription required to KnowItAll Spectral Libraries.




For Advanced Spectral Database Searching

Wiley offers the most powerful tools to search spectral data, built using the same technologies Wiley uses to analyze its vast datasets. With fast search speeds, powerful algorithms, and patented technologies, Wiley is able to deliver results you can rely on.

Key Features

- Import sample spectrum and search against user-generated databases or comprehensive Wiley KnowItAll spectral references libraries for IR, MS, NMR, Raman, UV-Vis*
- Searches are fully customizable and driven by powerful algorithms
- Optimized for speed and performance
- Search by name, structure, substructure, properties, spectra (full spectrum or selected range), and peak—in any combination
- Include or exclude regions from your search
- Perform mixture analysis for multiple components
- Include or exclude known components in mixtures to narrow results
- Manually select peaks or use automated peak picking capability
- Algorithms include Euclidean Distance, First and Second Derivative Euclidean Distance, Correlation, and Baseline Correction; For MS: Dot Product (Cosine), Wiley Dot Product (Cosine), Composite P1 and P3
- Perform simultaneous multi-technique search with spectra from various techniques to orthogonally validate for more confident analyses
- Includes patented Optimized Corrections technology to ensure best search results
- Includes MS Adaptive Search (patent pending) and Reverse Search
- Easily compare spectra using various views: overlay, offset, stack, butterfly, subtraction, etc.
- Supports multiple instrument types and vendor formats (see full list at www.knowitall.com/formats)



The Analytical Edition also includes these bonus databases:

- IR – Sadtler Polymers, Hummel – Wiley (1,907 spectra)
- IR – Sadtler Standards (Organic & Polymeric Compounds Subset) – Wiley (9,996 spectra)
- ATR-IR - Sadtler Solvents (629 spectra)
- **Microplastic Classifications Database** - This database is now included with the Analytical Edition as a fast and cost-effective way to classify microplastics of IR / ATR-IR spectra. For identification of microplastic samples, we recommend users subscribe to the comprehensive KnowItAll IR Spectral Library.

Powerful Spectral Search Tools that Set KnowItAll Apart

Wiley is committed to taking spectral analysis to the next level! We are continually adding spectrally intelligent tools to our portfolio that accelerate analysis. Here's a closer look at some of the unique and powerful solutions in **KnowItAll's SearchIt** that make Wiley a leader in spectral informatics!

Mixture Analysis: This industry leading capability for IR and Raman, is now available for MS! One of the most powerful features is KnowItAll's ability to analyze mixtures. When searching an unknown against a reference database, you can choose to search for multiple components. The result is a series of composite spectra, each accompanied by the individual component spectra that comprise the composite spectrum, as well as the residual spectrum (difference between the query and the composite spectrum). Composite spectra are then ranked by how closely they resemble the query spectrum.


Patented Optimized Corrections Technology: Searching is not always a straight-forward process. What if there are problems with your query spectrum, library reference spectra, or both? If there is, you may never find the right match—even if it is in the library. Wiley offers a unique patented solution to solve this complex problem and lead you to the best results.

Wiley's Optimized Corrections is a spectrally intelligent solution built into KnowItAll's ID Expert and SearchIt applications. It performs a computationally complex set of corrections on all query and reference spectra in a search to find the optimal match between the query and each reference spectrum. Multiple corrections are applied automatically to compensate for differences between spectra caused by the variability of different instruments and accessories as well as other factors, including human error. Corrections include: baseline correction, clipping, horizontal shift, vertical shift, intensity distortion, and ATR correction.

Multi-Technique "Simultaneous" Spectral Searching: KnowItAll is the world's first search system capable of searching spectra in multiple analytical techniques simultaneously from one or many databases. For example, query an NMR spectrum in one database and a mass spectrum in another database at the same time to find the most relevant hits from each database linked to one another by chemical structure.

Patent-Pending Adaptive Search for MS: When matching an unknown mass spectrum against a reference, this technology finds spectral matches that are similar to the unknown but have additional or missing selective fragment(s). It then suggests what might be causing the differences, where possible. This feature provides tremendous insight into structural possibilities to explore when there is no exact match. Ultimately, this may lead to more intelligent and confident identification and confirmation.

MS Reverse Search Feature: This search finds matches for peaks in the reference spectrum and ignores peaks which are only in the unknown spectrum.



Spectral Database Building and Management

Chemists and spectroscopists produce valuable data every day within their organizations. Because Wiley Science Solution's primary business is spectral databases, KnowItAll is built through years of experience in doing just that—building databases.


Build Single or Multi-Technique Databases from Various Vendors

Researchers can build searchable databases that include one or more analytical techniques (IR, MS, Near IR, NMR, Raman, UV-Vis), chemical structures, and other metadata. So even if a laboratory's instruments come from multiple manufacturers, KnowItAll can archive the data.

Key Features

Build Databases with Spectral & Chromatographic Data

- Build databases with one or more analytical techniques
- Build databases with multiple spectra in the same record
- Import analytical data even if laboratory has instruments from multiple vendors
- One-click import of common instrument file formats or *.csv format (spreadsheet)
- Enhance each record with peak information, structures, and properties, such as source of sample, boiling point, etc.
- Import multiple structure formats (with stereo-chemical bonds and identifiers)
- Use "Batch Import and Export" for efficient handling of spectra, structures, and property files
- Supports unrestricted spectral range and resolution - Store spectra at the precise range and resolution at which each spectrum was measured rather than being forced to conform to a fixed range and resolution
- Make database more powerful by attaching spreadsheets, MSDS, and other documents or adding hyperlinks to web pages
- Create cross-reference from record to data from another technique; i.e., an NMR spectrum can be linked to an IR spectrum
- NEW! Categorize chemical structures of controlled substances according to *Drug Enforcement Agency* regulations (batch processing available)
- Property calculators for single or batch calculations for entire datasets - formula, molecular weight, C-13 NMR prediction, bad baseline indicator, baseline analysis: area difference, SPLASH ID, various masses (average, exact, nominal)



- MS tools to calculate elemental composition and isotope distribution
- Quickly add properties and structures from PubChem to your database

Customize Databases

- Databases can be customized to meet laboratory specifications
- Create custom fields to support associated meta data relevant to their work
- Choose from three types of property fields: text, numeric, hyperlink
- Generate "preferred property" forms so users enter properties consistently
- Set spectral parameters such as x- and y-resolution

Extract the Most Information from Your Data

- Fully integrated with other KnowItAll applications for processing, database searching/mining, structure drawing, processing, reporting, and more

Multi-Technique Viewing & Mining


With Minelt, users can view reference databases, user-created databases, or search results. Access databases containing many types of data, such as spectra chromatograms, structures, physical properties, and more. Since analytical databases can contain one or more techniques in the same record, this tool is ideal for accessing databases of reference spectra.

Advanced Datamining Capabilities

Compare any two variables from a database using a scatter plot diagram to separate data that follow a desired trend from that which does not. Select any point on the scatter plot to display the compounds associated with that record.


Patented Overlap Density Heatmaps Technology

Traditionally, the visualization of multiple spectra takes place in an overlay, offset, or stacked plot. These traditional plotting methods, however, obscure trends when viewing large amounts of data. With Overlap Density Heatmaps, users can visualize trends and assess similarities and dissimilarities in massive amounts of data. Specifically, this technology allows the user to see common features of overlapped objects (such as spectra) by color coding spectral areas from highest to lowest overlap.




For example: This plotting feature is useful in the analysis of spectral searches for samples run in multiple techniques by plotting the quality of database search results (Hit Quality Indices - HQIs) against each other (e.g., IR HQI versus Raman HQI).

Overlap Density Heatmaps: An Example




Traditional Stacked Display

Thirty-one IR spectra of alkanes are shown. While some trends appear, the extent of the trends is obscured.



OD Heatmap
OD Level = 100

An Overlap Density Heatmap showing only those areas of overlap common to all spectra.



OD Heatmap
OD Level = 0

An Overlap Density Heatmap of the thirty-one alkanes shown revealing all overlap levels. High levels of overlap are displayed in red; low levels are displayed in violet.



IR, Raman, Chromatogram Quality Control Comparison

Wiley's KnowItAll QC Expert software performs a rapid quality check of a sample IR or Raman spectrum or chromatogram against a "gold standard" user spectrum or chromatogram to verify that a material meets control specifications.

Key Features

- Perform QC comparison of a sample to a selected reference file
- Validate results by also comparing the sample to a reference database to ensure the sample not only matches the selected reference, but that it also does not match anything else in the reference databases*
- Define user privileges, reference data, and other settings to ensure technicians follow set protocols and focus on output
- Identify problems with the sample spectrum - QC Expert's built-in spectral intelligence identifies issues and suggests ways to fix them

The screenshot displays the QC Expert software interface. The main window shows a comparison of a sample spectrum (red) against a selected reference spectrum (green). The x-axis represents wavenumber in cm⁻¹, ranging from 4000 to 500. The y-axis represents intensity, ranging from 0 to 1.0. A 'Sample Spectrum Status' table is visible, showing a 'Pass at 95% Threshold' and a 'Match: 100.0%'. Below this, a 'Match % of All Reference Spectra' bar chart shows a match of 83.1% for the next closest reference, D-GLUCOSE. The 'QC Comparison Status' section provides a detailed list of reference spectrum properties for D-GLUCOSE, including database code, name, formula, and source.

Checklist	Status	Fix?
Baseline	✓	ⓘ
Noise	✓	ⓘ
Contaminants	✓	ⓘ
Technique	✓	ⓘ

QC Comparison Status
Sample File Name: ATR_Glucose.lrf
Sample Title: DWX #238; D-(+)-Glucose
Pass at 95% Threshold
Match: 100.0%
Match % of All Reference Spectra: 83.1%
Next Closest: D-GLUCOSE: 83.1%

Selected Reference Spectrum Properties	
Name	Value
Database Code	DWX
Database Title	ATR-IR - Controlled & Prescription Drugs 1 - Bio-Rad
Sadtler	
Record ID	238
Name	D-(+)-Glucose
Correlation	
Algorithm	95 %
Threshold	
Proximity Warning	2 %
CAS Registry Number	50-99-7
Catalog Number	G7528
Formula	C6H12O6
Instrument Name	Bio-Rad P15
Lot Number	12700801
Mol.Weight	180.158 g/mol
Polar Surface Area Predicted	110.4
Source of Sample	Sigma-Aldrich Inc.
Source of Spectrum	Forensic Spectral Research

Structure: OC[C@H]1O[C@@H](O)[C@H](O)[C@@H](O)[C@H]1O

*Subscription required to KnowItAll Spectral Libraries.




AssignIt™ NMR

Create Fully Assigned NMR Databases

AssignIt NMR allows users to add NMR assignments to the structures in ^1H , ^{13}C , ^{19}F , ^{31}P , ^{15}N , ^{17}O , ^{11}B , and ^{29}Si NMR databases. AssignIt's easy-to-use interface allows quick database information input, such as peak shift assignments, intensities, coupling constants, and multiplicities—all linked to chemical structure.

Key Features

- Import of a wide variety of NMR formats
- Assign atoms to peaks in the experimental spectrum
- Interactive coupling calculation tool
- Automated calculation of J value within a multiplet signal
- "Find signals with same J" feature to find similar splitting within a spectrum
- Intuitive interface with summary view and data-entry forms to add/edit assignments
- Automated and manual peak picking tools
- Full integration with ChemWindow (structure drawing) and MineIt™



Spectral Analysis Toolbox



Advanced Functional Group Analysis

for IR, Raman, IR (Polymers), & Vapor Phase IR (Synthetic Cannabinoids)

Interpret a Spectrum: Simply load a spectrum and click a peak of interest; Analyzelt then lists all functional groups possible at that peak position. Compare peak regions for each group by overlaying with the spectrum and narrow results by tagging the “most likely” candidates.

Correlate a Structure with a Spectrum: This powerful feature helps determine if a proposed structure matches an observed spectrum. Just draw or import a structure to view its component functional groups. Then compare peak regions for each group by overlaying with the spectrum.

Classification Knowledgebases



IR & Raman	Over 200 functional groups and hundreds of interpretation frequencies
IR Polymer	Over 100 functional groups and hundreds of interpretation frequencies
NEW! Vapor Phase IR - Synthetic Cannabinoids (Patented)	Over 60 bridge carbonyl groups
Or build your own knowledgebases to improve interpretations	

Benefits

- Useful in the classification of spectra of unknown compounds
- Supplemental to other methods of spectral interpretation

Key Features

- Import and peak analysis of spectra
- Intelligent “Suggest a Peak” feature
- Determine if a structure matches a spectrum
- Browse knowledgebase by chemical class
- Tag and summarize negative or positive interpretations
- Peak overlay display for easy comparison
- Display/highlight structural bonds involved in vibrational frequency
- Link to additional data in Sadtler Handbook (AnalyzeIt IR only)
- View notes for functional groups when available





NMR Spectrum Prediction

With PredictIt NMR, perform database-based NMR spectrum predictions for ^{13}C , ^1H , and other nuclei.

Predictions are performed automatically when users open a structure in PredictIt NMR. To make predictions, this tool examines databases of substructures that have ^1H , ^{13}C or other shifts assigned to them. The substructures are defined by the number of shells that represent atoms within n bonds of the central atom.

For example, a shell of four would include the central carbon atom and all atoms within four bonds of this atom. After looking for exact matches, PredictIt NMR looks for matching shells for each atom in the structure, starting with shell four and proceeding to smaller shells until matches are found.


The tool searches the database(s) for specific chemical environments, which are described by a modified HOSE (Hierarchically Ordered Spherical description of Environment) code, a topology code used to describe the chemical surroundings of an atom in a molecular structure. The original structure and results are displayed in PredictIt NMR's main window. Each atom's average shift (and standard deviation) is displayed at the top level of the tree control.

Solvent-Specific Prediction for Improved Accuracy

KnowItAll offers the first solvent-specific NMR chemical shift prediction on the market. Users can choose from a list of common solvents such as chloroform, acetone, and dimethyl sulfoxide and KnowItAll will automatically recalculate all chemical shifts for that solvent.

More Than Just the Spectral Data

Predicted peak shifts are not the only piece of information that NMR spectroscopists need. PredictIt NMR not only allows easy retrieval of the real spectral data used to build the prediction, but also access to available information related to the reference spectrum, such as sample source, solvent, conditions of production, equipment, and properties for the molecule.



Spectral Processing Toolbox



IR Spectrum Processing


ProcessIt provides a variety of tools to process IR spectra and improve the quality of archived data and search results. It can also be used in conjunction with other KnowItAll tools. For example, a spectrum can be transferred from SearchIt to ProcessIt to correct potential searching problems and transferred back.

Processing Capabilities Include:

- Flatline
- Truncation / Padding
- Normalization
- Smoothing (Quad-Cubic Savitsky Golay, Fourier methods)
- Baseline Correction (spline, linear, and polynomial methods)
- ATR Correction
- Reverse ATR Correction
- Kubelka-Munk Transform
- Spectral Subtraction and Spectral Addition
- Average Spectra
- Peak Picking

Analysis Capabilities Include:

- Area Under the Curve (AUC)



Raman Spectrum Processing

ProcessIt provides a number of tools to process Raman spectra and improve the quality of archived data and search results. It can also be used in conjunction with other KnowItAll tools. For example, a spectrum can be transferred from SearchIt to ProcessIt to correct potential searching problems and transferred back.

Processing Capabilities Include:

- Flatline
- Truncation / Padding
- Normalization
- Smoothing (Quad-Cubic Savitsky Golay, Fourier methods)
- Baseline Correction (spline, linear, and polynomial methods)
- ATR Correction
- Reverse ATR Correction
- Kubelka-Munk Transform
- Spectral Subtraction and Spectral Addition
- Average Spectra
- Peak Picking

Analysis Capabilities Include:



- Area Under the Curve (AUC)

GC-MS Data Processing

With ProcessIt, users can view and select a mass spectrum or the average of mass spectra and define the background. The selected mass spectrum or extracted spectra can be transferred for database searching or imported into a user database.

Selected Ion Chromatograms (SICs) - ProcessIt allows the display of a selected ion chromatogram. Multiple ion chromatograms can be displayed in the top pane.

Spectral Subtraction - This feature allows calculation of the average mass spectrum and also allows the elimination of background noise via manual background subtraction. Single or multiple ranges for either process can be specified. As shown in the screenshot to the right (*bottom*), the lower pane holds either selected MS, selected MS - background (*red bar*), average MS (*green bar*), or average MS (*green bar*) - background (*red bar*).



NMR Spectrum Processing


With ProcessIt, import and process NMR spectra from various sources to improve the quality of archived data and search results. This tool is easy to use, yet offers a comprehensive set of processing features to correct experimental artifacts and improve the appearance of your spectra.


Chemists and spectroscopists can use ProcessIt at their own desktops to process and re-process experimental data. In addition to being more convenient for the user, ProcessIt also saves valuable processor time at the instrument, thereby improving sample throughput.

Because ProcessIt is fully integrated in the KnowItAll informatics environment, processed spectra can be transferred to other KnowItAll tools with a single click.

Key Features:

- Import 1D processed or FID spectra from multiple formats
- Processing features: zero filling, interactive window functions, and Fourier transform
- Automatic and manual phase correction
- Automatic and manual baseline correction, includes polynomial, spline and linear algorithms
- Automatic and manual peak picking
- Automatic and manual integration
- Addition and subtraction of spectra
- Overlay multiple spectra for easy comparison
- Macro capability for quick and efficient processing
- Export in JCAMP format
- Spectrum handling tools, such as horizontal zoom, box zoom, hand cursor, and scaling
- Integrated with Minelt for archiving of processed spectra, ReportIt to create reports containing spectra, peak, and integral tables, and SearchIt for spectral searches





MS Expert™

Automated Non-Targeted GC-MS Analysis


GC-MS data analysis can be time consuming, especially when examining complex analytes. This application automatically processes, deconvolutes, and analyzes GC-MS data. Combined with KnowItAll's fast database search, it suggests matches to knowns and allows unknowns to be further examined.

When it is necessary to send coeluted components to the sophisticated mixture analysis application in KnowItAll, MS Expert facilitates this seamlessly. Novel compounds' structural characteristics can be further deduced by applying the patent-pending MS Adaptive Search (in SearchIt) that uses fragmentation and structural information to propose probable structural details.

Key Features

- Automatic GC-MS data deconvolution
- Automatic search of component's extracted MS spectra to match against comprehensive KnowItAll MS spectral references libraries*
- Also allows user databases to be included in the search
- Simultaneous regular spectral search and reverse search
- Adjustable match score calculations based on regular spectral search and reverse search HQI values
- Can analyze unit mass data as well as accurate mass data
- Display presents TIC, component profile, extracted spectrum vs raw spectrum, extracted spectrum vs. matched reference spectrum graphs, as well as the matched reference data structure
- Allows user to manually pick additional peaks from the TIC
- Select TIC regions for analysis
- Adjustable analysis sensitivity parameters
- Adjustable input data resolution set-up for accurate mass data
- Seamlessly transfer components which are not matched well to reference spectra to SearchIt for manual examinations using either Adaptive Search or mixture analysis
- Generate reports with results
- Supports multiple instrument types and vendor formats (www.knowitall.com/formats)
- Optimized for speed and performance

*Subscription required to KnowItAll Spectral Libraries.



How Does It Work?

Simply import a GC-MS data file and the software automatically deconvolutes the TIC into components. All extracted component MS spectra are then automatically searched against reference libraries* to find matches. Results are shown as a hit list for each component.

Unidentified components or components with low match scores can then be sent to KnowItAll's SearchIt tool for manual investigation. In SearchIt, one can use adaptive search to match similar components or use mixture analysis to separate coeluted components.

Basics Toolbox



A Full-Featured 2D Structure Drawing Program

ChemWindow is the software chemists worldwide choose for chemical structure drawing. It provides an advanced set of drawing tools that's easy to use— just click and drag to draw any chemical structure. Access the most comprehensive set of tools to draw rings, bonds, atoms, electrons, charges, chains, arrows, and more.


Key Features

- Customizable toolbars with tools to draw chemical structures, including bonds, rings, atom labels, charges, etc.
- Chemical recognition features such as hot keys, chemical syntax checker
- Advanced stereo-chemical recognition—using technology not available in other packages.
- OLE (Object Linking and Embedding) technology for in-place editing in word processing and presentation software
- Tools to calculate mass and formula, MS tools to calculate elemental composition and isotope distribution
- Predefined styles for captions and structures
- Links to OPSIN Name2Structure to convert a name to a structure
- Easily import existing structures from multiple file formats (ChemDraw - *.cdx, CML - *.xml, Hampden - *.hsf, InChI - *.txt, JCAMP - *.dx, *.jdx, BIOVIA/MDL - *.mol, *.rxn, Smiles - *.smi, XYX - *.xyz, etc.)
- Supports reaction files including RInChI, as well as CDX and CDXM files



Web Training Resources

Browselt is a web browser built into the KnowItAll software with links to KnowItAll tutorial videos and other resources for KnowItAll users.






ReportIt™

A Full-Featured Publishing Program

With ReportIt, create standard reports, design papers, presentations, and web publications that include annotations, tables of data, spectra, 2D and 3D structures, and more.

Key Features

- Custom templates to create uniform reports for enterprise-wide format standardization
- Customizable toolbars to draw chemical reactions and other reports, including arrows, text boxes, shapes, etc.
- Clip art libraries with hundreds of laboratory glassware drawings and engineering symbols
- OLE technology (Object Linking and Embedding) for in-place editing in word processing and presentation software
- MS fragmentation tool to display a mass for each fragment
- Advanced editing options to align, space, center graphics, and rotate captions
- Predefined styles for captions and structures
- 3D structure visualization for high-quality, realistic 3D drawings
- Table tool to enter and organize your data
- Spectrum / chromatogram import in common native file formats
- Multi-spectrum display in three display modes: overlay, stack, and offset
- Advanced spectrum display editing features to customize the appearance of spectra and chromatograms, including axes, colors, labels, etc.
- Custom annotation tool to link objects like spectral peaks to text graphics or chemical structure captions



Accelerate Your Analyses with KnowItAll Spectral Libraries

From the Leader in Spectral Data

Spectral analysis software is not complete without a high-quality database of reference spectra. And with KnowItAll, Wiley offers the best of both!

Wiley is the leading producer and publisher of spectral databases, with a collection that contains over 2 million spectra (IR, MS, NMR, Raman, and UV-Vis) covering pure compounds and a broad range of commercial products.

Analyze Samples Faster

The KnowItAll software when combined with Wiley's KnowItAll spectral library subscriptions offers users an unparalleled solution for spectral identification. With access to massive collections of high-quality reference spectra, the likelihood and speed of analysis increases. Your lab can ultimately analyze samples faster and save valuable research time. It's that simple.

Broad Compound Coverage

The KnowItAll collection is an essential tool for the identification, classification, and verification of unknown compounds in a wide range of applications such as polymer/materials, environmental, forensics/toxicology, pharmaceutical, biotech, automotive/aerospace, food/cosmetics, and many more.

Trusted Data from a Trusted Source

Wiley is an authoritative source for spectral data. Their renowned databases are processed according to rigorous protocols to ensure they are of the highest quality. These qualification procedures start at data acquisition and continue throughout the database development process. Any data acquired from trusted partners is thoroughly vetted before inclusion in our collections.



Our experts can help you determine the best data mix for your laboratory:

- KnowItAll IR Spectral Library
- KnowItAll Mass Spectral Library
- KnowItAll NMR Spectral Library
- KnowItAll Raman Spectral Library
- KnowItAll UV-Vis Spectral Library

With KnowItAll annual subscriptions, researchers can get access to comprehensive collections of spectra—plus any updates.

Powerful Software. Quality Data. Results You Can Rely On.

sciencesolutions.wiley.com