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

Waters™

Application Note

Analysis of Per- and Polyfluoroalkyl Substances in Groundwater by Direct Injection Using the Benchtop Multi-Reflecting Time-of-Flight Xevo™ MRT Mass Spectrometer

Hania Khoury-Hollins, Sara Stead, Jayne Kirk, Richard Lock

Waters Corporation, United States

Published on November 20, 2025

Abstract

The increasing regulatory scrutiny and scientific interest in per- and polyfluoroalkyl substances (PFAS) are driving a shift in analytical workflows from traditional tandem quadrupole (TQ) methods toward high-resolution mass spectrometry (HRMS) platforms. While TQ systems are the preferred choice for absolute quantification, they are restricted by the availability of commercialized standards and the limited number of compounds that can be identified and quantified within a single method. In contrast, HRMS enables both absolute quantification and non-targeted screening (NTS) workflows simultaneously, offering a more holistic understanding of PFAS profiles in complex environmental matrices.

To demonstrate the capabilities of HRMS, the Xevo MRT Mass Spectrometer was evaluated for sensitivity and dynamic range by direct injection of groundwater samples and standards. 30 PFAS compounds were detected



Analysis of Per- and Polyfluoroalkyl Substances in Groundwater by Direct Injection Using the Benchtop Multi-Reflecting Time-of-Flight Xevo™ MRT Mass Spectrometer

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




Lisa Reid, Matthew E. Daly, Lee A. Gettings, Jayne Krikorian
Waters Corporation

Abstract

We demonstrate the ability of data processing to identify and pathway profile a typical metabolomic study acquired using a Xevo™ MRT Mass Spectrometer coupled with an ACQUITY™ Premier LC System. The resulting data shows a clear unsupervised PCA separation between the three study groups with exogenous or environmental markers being the most significant differences from the different volunteer lifestyles. The Xevo MRT MS delivers rapid, robust data acquisition for large cohort metabolomic studies, capable of providing a high mass resolution of up to 100,000 full width half max (FWHM) and consistent sub-ppm mass accuracy, ensuring the analyst has increased confidence in putative feature identification.

Benefits

Simple and robust workflow delivering excellent data reproducibility and confident results

A powerful benchtop mass spectrometer routinely providing high mass resolution for both precursor and product ions, enabling separation of compounds with near identical molecular weights e.g. isomers

Delivering a consistent mass accuracy <1 ppm RMS, increasing confidence in putative feature identification

Ability to acquire data at extremely fast acquisition rates, maintaining mass resolution and mass accuracy benefits. Allowing compatibility with shorter chromatographic methods, increasing sample throughput

Metabonomics Workflow using a Xevo[®] MRT Mass Spectrometer

Xevo MRT– toxikologický screening

Waters™

LC-MS Conditions


LC/MS methods with ACQUITY UPLC I-Class Plus & Xevo MRT

Ionization Mode	ESI Positive	ESI Negative
Run time	15 min gradient elution	7.5 min gradient elution
Mobile Phase A	5mM ammonium formate pH 3.0	Water with 0.001% formic acid
Mobile Phase B	Acetonitrile with 0.1% formic acid	Acetonitrile with 0.001% formic acid
Acquisition range	m/z 50-1000	
MS ^E conditions	Collision energy ramp 10-40 eV	
ToF MRM conditions	23 MRM transitions for compounds of interest	
Column	ACQUITY UPLC® HSS C18 Column	
Flow rate	400 µL/min	
Injection vol.	5 µL	

Quantitation Accuracy: ToF MRM vs. MSe

Measured quantitation deviation (%)* of a wide range of analytes at different concentrations


The ToF MRM method demonstrates greater quantitation accuracy compared to MSe.

* Deviation(%) = $\frac{\text{Measured} - \text{Expected}}{\text{Expected}} \times 100$, ng/mL (mira odchylky mezi naměřenou a skutečnou koncentrací)


Quantitation Sensitivity: ToF MRM vs. MSe

Ratio of S/N measured between ToF MRM and MSe acquisition @ the same level

ToF MRM is approximately 9 times more sensitive than MSe for targeted quantitative analysis, although the actual improvement varies by analyte.


Ratio of S/N = Analyte S/N (ToF MRM) / Analyte S/N (MSe)



Amphetamine and oxazepam are excluded, as their S/N in MSe mode is <3

THC-COOH is not detectable within the positive mode

Tramadol (25.6-720 ng/mL)

The ToF MRM method demonstrates greater linearity than MSe in the quantitation of tramadol.

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
- Amphetamines
- Anaesthetic agents
- Analgesics
- Antiarrhythmics
- Anticonvulsants
- Antiepileptic
- Antimalarial
- Antipsychotics
- Anti-seizure
- Antidepressants
- Antihypertensive
- Benzodiazepines
- Cocaine + metabolites
- Opioids
- Fentanyl/Fentanyl analogues
- Novel Psychoactive Substances
- Cathinone's
- Hallucinogen's
- Nitazene's
- Synthetic Cannabinoid's
- Pesticides (136)
- Piperazines
- SSRIs /TCAs

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

Name of the library	NPS library ¹	CHU Life Developed library ²	Wiley Library Maurer, Meyer, Heffer, Weber: LC-HRMS/MS Library of Drugs, Poisons, and Their Metabolites
# of analytes	2373	1,483	Over 3,200

Xevo MRT – lipidy



Waters™

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Conventional – 12 min

High-throughput – 5 min

1. FFA/Lysophospholipids; 2. Glycerophospholipids; 3. Glycerolipids

Waters™

Application Note

High Sensitivity, High Throughput LC-MS
Analysis of Eicosanoids Using the Xevo™ MRT
QToF

Nyasha Munjoma, Gemma Molyneux, Lee A. Gethings, Jayne Kirk, Richard Lock

Waters Corporation

Waters™

Application Note

The Achievability of Single Cell Lipidomics
Utilizing the Novel Xevo™ MRT Mass
Spectrometer

Scarlet Ferrinho, Nyasha Munjoma, Preeti Mourya, Shazneel Briones, Lee A. Gethings, Robert S. Plumb, Oliver

Census, Paul A. Townsend

Waters Corporation, Faculty of Health and Medical Sciences, University of Surrey

[PRODUCT SOLUTION]

A High-Throughput Lipodomic Workflow Using the Waters Xevo MRT Mass Spectrometer

INTRODUCTION

The Xevo™ MRT Mass Spectrometer delivers an unrivaled combination of performance at speed, providing high resolution time-of-flight data independent of spectral acquisition rate. The MS specification of 100,000 (FWHM) resolution at scan rates of up to 100 Hz combined with <500 ppb mass accuracy enables confident identification of analytes at throughput levels that enable large cohort experiments.

The cutting-edge multi-reflecting ToF technology that Waters™ Corporation pioneered in the SELECT SERIES™ MRT, has been scaled into a compact benchtop platform to solve the most challenging problems in biomedical research and epidemiological studies.

In large scale metabolic phenotyping experiments performance is critical, from instrument stability over large studies to precision of measurement, generating high quality experimental outcomes. The novel acquisition system, with dual gain analogue-to-digital converter, delivers long term system stability over a high dynamic range, ensuring exceptional quality data with consistent robustness and reproducibility.

Here we describe, using an exemplar sample set, the use of the Xevo MRT Mass Spectrometer for the analysis of colorectal cancer human serum samples. This combined with market leading chemistry and separations, using premier technology, and informatics including LipoStar2 from Mass Analytica™ (MassAnalytica.com) describes a complete lipidomic workflow.

EXPERIMENTAL CONDITIONS

SAMPLES

Samples	
	6 Healthy control plasma
	4 Colon cancer plasma
Colorectal cancer (CRC) human serum samples	2 Rectum cancer plasma
	NIST SRM 1950 plasma
	Study Reference/pooled (QC)
	A 100x dilution of EquiSPLASH

LC conditions

LC System	Waters ACQUITY™ Premier
Column	ACQUITY Premier UPLC™ CSH™ C ₁₈ 2.1 x 50 mm, 1.7 µm
Column Temperature	55 °C
Sample temperature	8 °C
Injection volume	0.5 µL
Mobile phase A	50:25:25 H ₂ O:IPA:MeCN w/5 mM ammonium acetate and 0.05% acetic acid
Mobile phase B	50:50 IPA:MeCN w/5 mM ammonium acetate and 0.05% acetic acid
Flow rate	0.4 mL/min

Xevo MRT – DESI XS (2025)

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DESI XS with Xevo MRT Mass Spectrometer for MS Imaging

Accelerating imaging. Amplifying sensitivity. Advancing visualization.

INTRODUCTION

This product solution highlights the capabilities of the Xevo[™] MRT mass spectrometer equipped with the DESI[™] (Desorption Electrospray Ionization) XS source (Figure 1), demonstrating its ability to provide rapid acquisitions while maintaining high sensitivity for MS imaging. With combined high mass resolution and mass accuracy, your lab can achieve greater spectral quality, improved analyte detection, and higher confidence in compound identification. The data presented includes data acquired in both positive and negative ionization modes and compares various acquisition speeds.



THE POWER OF DESI XS

Analysis on the DESI XS source offers a powerful and versatile approach to mass spectrometry imaging, delivering a range of benefits that enhance a wide range of applications, including pharma, chemical materials, and pre-clinical research applications. One of the most significant advantages of DESI is its ability to perform ambient image acquisitions without the need for extensive sample preparation or matrix application. The direct sample analysis without the application of matrices makes DESI ionization particularly valuable for high spatial resolution acquisitions where application of a matrix could lead to analyte movement. Whilst simple, rapid sample preparation makes it ideal for real-time analysis applications.


The DESI XS source supports high spatial and high mass resolution imaging, currently offering analyses with pixel sizes as small as 5 μm , enabling single-cell image resolution. This level of detail allows researchers to visualize molecular distributions with exceptional fidelity, crucial for understanding complex biochemical processes and tissue heterogeneity. The integration of Microscope Mode within High-Definition Imaging (HDI[™]) Software further enhances efficiency by allowing rapid survey scans followed by targeted high-resolution re-acquisition of regions of interest [1].

Figure 1. Next-generation solution for MS imaging: DESI XS with Xevo MRT mass spectrometer.


Waters[™]

Xevo MRT – mAb, proteiny, peptidy, oligonukleotidy

Waters™

Waters™

Application Note

Enhanced Host Cell Protein Identification and Quantitation with Multi-Reflective LC-MS

Jonathan Fox, Malcolm Anderson, Ying Qing Yu, Scott J Berger, Nick Pittman, Laetitia D

Waters Corporation, United States

Published on August 07, 2025



Waters™

Application Note

Applying Peptide Mapping and Multi-Attribute Method (MAM) Workflow for Biosimilar mAb Drug Products Comparison on the Xevo™ G3 QToF Platform

Kellen DeLaney, Samantha Ippoliti, Lisa Reid, Owen Cornwell, Ying Qing Yu, Emma Harry, Mark Towers

Waters Corporation

Waters™

 waters_connect™