

Agilent EasiCal Polystyrene Standards for Individual Pore Size GPC/SEC **Columns**

Technical Overview

Introduction

The Agilent EasiCal system for organic solvents consists of two different combs, each with 10 detachable spatulas supporting a mixture of five polymer standards. The thin film of polymer (approximately 5 mg) on the tip of the PTFE spatulas rapidly dissolves when immersed in eluent to provide two calibration solutions for gel permeation chromatography/size exclusion chromatography. A single pack provides 10 spatulas of each type, with molecular weights selected to provide equidistant calibration points for greater accuracy.

Conditions

Calibrants Agilent EasiCal PS-1 Polystyrene Standards

Columns Agilent PLgel 10 μ m 10⁶Å, 300 × 7.5 mm (p/n PL1210-6160) Agilent PLgel 10 μ m 10⁴Å, 300 × 7.5 mm (p/n PL1210-6140)

Agilent PLgel 10 μ m 500Å, 300 × 7.5 mm (p/n PL1210-6120)

Eluent Flow rate 1.0 mL/min Conc 0.05% Detector UV. 254 nm System Agilent PL-GPC 50

GPC/SEC Columns and Calibrants from Agilent

Agilent offers a comprehensive portfolio of GPC/SEC columns and calibrants for highperformance separations based on molecular size in solution. Agilent delivers leading solutions for characterizing and separating polymers by GPC/SEC, and manufactures all components for accurate polymer analysis.

Look at the Agilent Literature Library on www.agilent.com/chem/gpc-sec for a comprehensive range of application notes and technical overviews to help you get the best from your Agilent GPC/SEC columns and instruments.








Figure 1. Separation of low molecular weight polystyrene standards on Agilen. PLgel 10 µm columns.

Table 1. Agilent EasiCal Standards

Comb A MW mixture	Comb B MW mixture
8,500,000	3,040,000
1,030,000	330,000
156,000	66,000
28,500	9,200
3,250	580

www.agilent.com/chem

Agilent shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material.

Information, descriptions, and specifications in this publication are subject to change without notice.

© Agilent Technologies, Inc., 2015 Printed in the USA April 30, 2015 5990-8224EN

