

REVERSED-PHASE LC PRIMARY STRUCTURE CHARACTERIZATION WORKFLOW

AGILENT LC ADVANCEBIO RP-mAb COLUMNS



In this document Agilent applications chemists share their recommendations for an optimum LC system and its configuration for characterizing biomolecules. They also offer guidance on a generic method to get you started, and how this method can be further optimized to meet your specific separation goals.

Additional application information is available at www.agilent.com/chem/advancebio

Agilent 1260 Infinity Bio-Inert LC System

Guidelines

Bonded phase: C4 and SB-C8 are routinely used. For alternative selectivity use the Diphenyl column.

Gradient: IPA:ACN:water with 0.1% TFA or 0.1% FA to elute all components of interest.

Sample solubility: Mix with starting mobile phase.

Temperature: Higher column temperature can dramatically improve resolution and recovery of proteins.

Resolution/selectivity: A blend of IPA:ACN provides better resolution. Other organic solvent substitutions can be used for different selectivity. Use of IPA results in increased pressure, which can be managed by column temperature and flow rate.

LC/MS

Desalt protein samples before injection.

Small id columns (e.g. 2.1 mm) are often the best choice with a 0.5 to 1.0 mL/min flow rate.

FA provides better MS signals. Use less TFA in the eluent to enhance MS signals, or use AcOH.

IPA with ACN in the mobile phase provides sharper TIC peaks.

Mobile phases

Eluent A: 0.1% TFA

Eluent B: IPA, ACN, and water with 0.09% TFA

Pump (G5611A)

Typical flow rate for 2.1 mm id columns is 1.0 mL/min

Sample injection (G5667A)

1 to 5 μ L injection for samples containing 1 to 5 mg/mL of mAb. Samples can be dissolved in water or eluent A.

Column compartment (G1316C)

60 to 90 °C is a typical temperature for good separation

Detection (G1315D)

UV (210, 280 nm) with a 10 mm bio-inert standard flow cell

**BIO
inert**



Column selection


Column	C4	SB-C8	Diphenyl
2.1 x 50 mm, 3.5 μ m	799775-904	789775-906	799775-944
2.1 x 75 mm, 3.5 μ m	797775-904	787775-906	797775-944
2.1 x 100 mm, 3.5 μ m	795775-904	785775-906	795775-944
2.1 x 150 mm, 3.5 μ m	793775-904	783775-906	793775-944
4.6 x 50 mm, 3.5 μ m	799975-904	789975-906	799975-944
4.6 x 100 mm, 3.5 μ m	795975-904	785975-906	795975-944
4.6 x 150 mm, 3.5 μ m	793975-904	783975-906	793975-944



Agilent Technologies

Intact mAb Analysis


Fast and high-resolution separation



Columns: 2.1 x 100 mm, 3.5 μ m
 Eluent A: 0.1% TFA in water:IPA (98:2)
 Eluent B: IPA:ACN:Eluent A (70:20:10)
 Flow rate: 1.0 mL/min
 Gradient: 10-58% B in 4 min, 1 min wash at 95% B,
 1 min re-equilibration at 10% B
 Injection: 5 μ L (1 mg/mL)
 Sample: Herceptin IgG1 variant
 Temperature: 80 $^{\circ}$ C
 Detection: UV, 254 nm


mAb Fragment Analysis

Chemical digestion – heavy chain/light chain



Column: AdvanceBio RP-mAb Diphenyl,
 2.1 x 50 mm, 3.5 μ m
 Eluent A: 0.1% TFA in water
 Eluent B: IPA:ACN:water (70:20:10)
 + 0.09 % TFA
 Flow rate: 1 mL/min
 Gradient: 0 min, 15% B; 0.5 min, 25% B;
 1.5 min, 35% B; 1.51 min, 35% B;
 3 min, 60% B; 4 min, 60% B
 Injection: 1 μ L (1 mg/mL) (TECP reduction)
 Sample: Rituximab innovator
 Temperature: 80 $^{\circ}$ C
 Detection: UV, 220 nm

Enzymatic digestion – Fab/Fc regions




Column: AdvanceBio RP-mAb C4,
 2.1 x 100 mm, 3.5 μ m
 Eluent A: 0.1% TFA in water
 Eluent B: n-Propanol:ACN:eluent A (80:10:10)
 Flow rate: 0.8 mL/min
 Gradient: 5-40% B in 5 min, 1 min wash at 95% B,
 1 min re-equilibration at 10% B
 Injection: 1 μ L (2 mg/mL)
 Sample: Herceptin IgG1 variant – papain digested
 Temperature: 60 $^{\circ}$ C
 Detection: UV, 220 nm


LC/MS Analysis of Intact mAbs

Fast chromatography with excellent peak shape and MS data using formic acid mobile phase





AdvanceBio RP-mAb C4,
 2.1 x 50 mm, 3.5 μ m



AdvanceBio RP-mAb Diphenyl,
 2.1 x 50 mm, 3.5 μ m



Eluent A: 0.1% FA in water
 Eluent B: IPA:ACN:water (80:10:9.9) + 0.1% FA
 Gradient: 0 min, 20% B; 4 min, 20% B;
 5 min, 40% B; 10 min, 70% B;
 11 min, 90% B; 11.1 min, 20% B
 Flow rate: 0.6 mL/min
 Temperature: 80 $^{\circ}$ C
 Injection: 1 μ L (1 μ g/ μ L)
 Sample: Innovator Herceptin
 Detection: Agilent 6530 Accurate-Mass
 Q-TOF LC/MS



This information is subject
 to change without notice.

© Agilent Technologies, Inc. 2015
 Printed in USA, October 6, 2015
 Publication number 5991-6321EN



Agilent Technologies