

# Multiclass Residue Analysis of Veterinary Drugs in Animal Tissues and Eggs Using Liquid Chromatography Coupled With Tandem Mass Spectrometry (LC/MS/MS)

## Application Note

### Authors

Tony Zhang, Dan-Hui Dorothy Yang,  
Tarun Anumol, Jianzhong Li,  
Andy Zhai, and Joan Stevens  
Agilent Technologies Inc.  
412 Yinglun Road  
Shanghai 200131  
China

### Abstract

This application note describes a workflow solution that was developed for the screening and quantification of multiclass veterinary drugs in different animal-derived food samples, from sample preparation to quantification of results. These veterinary drugs, belonging to 32 different chemical classes, included macrolides, sulfonamides, tetracyclines,  $\beta$ -lactams,  $\beta$ -agonists, chloramphenicols, nitroimidazoles, cepheids, avermectins, benzimidazoles, glucocorticosteroids, nonsteroid anti-inflammatory drugs (NSAIDs), peptides, anticoccidiosis drugs, quinolones, imidazoles, androgens, polyethers, triphenylmethanes, phenothiazines, quinoxalines, trematocides, antiviral, pesticides, dapsones, organic acids, nitros, tranquilizers, and so forth. The animal-derived food samples included muscle and liver of swine, cow, and chicken, as well as fish and hen eggs. The sample preparation involved a rapid and efficient protein precipitation extraction by acidified acetonitrile, followed by Agilent EMR—Lipid dSPE (p/n 5982-1010) and a polish kit (p/n 5982-0102) for further cleanup. The workflow solution supports the mid-to-high-end of Agilent tandem mass spectrometers (Agilent 6460 Triple Quadrupole LC/MS, Agilent 6470 Triple Quadrupole LC/MS, and Agilent 6495 Triple Quadrupole LC/MS) coupled with the Agilent Jet Stream Ionization Source.



**Agilent Technologies**

## Introduction

Veterinary drugs (VDs) have been used widely in veterinary practice to treat and prevent diseases, or to enhance growth and feed efficiency. If veterinary drugs are not administered correctly, it can lead to the presence of VD residues in foods of animal origin for human consumption, which can pose a health risk. With prolonged exposure, some antimicrobials, even at very low levels in foodstuff, may lead to an increase in antimicrobial resistance in human populations.

To better control the use of VDs and protect public health, the residue levels of VDs in meat and other foods are regulated by federal regulatory agencies. In China, Regulation No. 235, issued by the ministry of agriculture (MOA), separates veterinary drugs into four categories: Drugs in category I are permitted without maximum residue levels (MRLs) defined. Drugs in Category II are permitted for disease treatment with MRLs required. Drugs in Category III are permitted for disease treatment, however, no detectable levels are allowed in samples for human consumption. Drugs in Category IV are prohibited with no detection permitted in samples. The Chinese government annually sets up a series of monitoring plans to control the use of Category II, III, and IV veterinary drugs. Currently, the published veterinary drug residue analysis standard methods in China focus on only one chemical class or several chemical classes of drugs (for example, sulfonamides, quinolones, tetracyclines, and so forth). It can take days to finish testing and get the reports. Such a workflow is not only time-consuming, labor-intensive, and costly, but is also inefficient when dealing with a large number of samples. The aim of our workflow solution is to enable laboratories to analyze >180 commonly monitored VDs across multiple chemical classes simultaneously. The list of VDs in our method were based on general notice No. 235 (MOA, China).

Animal-derived food samples contain high amounts of proteins and lipids, and pose significant challenges with regards to sample preparation. In this study, acidified acetonitrile (5 % formic acid) was selected as an extraction solvent capable of efficient protein precipitation to extract analytes into the organic phase. Agilent Bond Elut EMR—Lipid dispersive SPE (EMR—Lipid) was used for cleanup with a polishing step. EMR—Lipid is a novel sorbent material that selectively removes major lipid classes from the sample without affecting analyte recoveries.

## Experimental

### Standards and reagents

Veterinary drug standards were purchased from Dr. Ehrenstorfer GmbH, WITEGA laboratorien Berlin-Adlershof GmbH, Toronto Research Chemicals (TRC) or AccuStandard, Inc. Ultrapure water (>18.2 M $\Omega$ , ELGA VEOLIA PureLab Chrous system), acetonitrile (LC/MS grade, Fluka), and formic acid (~98 %, for mass spectrometry, Fluka) were used for mobile phase preparation. Acetonitrile (ACN, HPLC grade, Sigma-Aldrich), dimethyl sulfoxide (DMSO, >99.9 %, Aldrich), ammonium acetate ( $\geq$ 98 %, Sigma-Aldrich), and formic acid (~98 %, HPLC grade, Fluka) were used during sample preparation.

### Sample preparation

A quick and simple sample preparation method was developed. The homogenized sample was extracted using acidified acetonitrile, which in essence served as a protein precipitation step. The extract was further cleaned using Agilent EMR—Lipid to selectively remove lipids while not trapping contaminants of interest, followed by a final polishing step for further cleanup.

The sample preparation method suits many types of animal-derived food samples, including muscle and liver of swine, cow, and chicken, as well as fish and hen eggs.

### Instrumentation

Analysis was performed on an Agilent 1290 Infinity II LC system consisting of:

- Agilent 1290 Infinity II Binary Pump with a 35  $\mu$ L Jet Weaver (G7120A)
- Agilent1290 Multisampler with a 20- $\mu$ L loop (G7167B)
- Agilent1290 Infinity II Multicolumn Thermostat (G7116B)

Three models of tandem mass spectrometers were tested:

- Agilent 6460 Triple Quadrupole LC/MS with Jet Stream technology (G6460CA)
- Agilent 6470 Triple Quadrupole LC/MS with Jet Stream technology (G6470AA)
- Agilent 6495 Triple Quadrupole LC/MS with the iFunnel and Jet Stream technology (G6495A)

## Method development and validation

The following performance characteristics were considered during method development and validation: dynamic range, limits of detection (LODs), limits of quantification (LOQs), recovery, and repeatability.

The full validations were carried out on an Agilent 6495 Triple Quadrupole LC/MS with five matrices (pork, swine liver, eel, chicken, and egg) at three spiking levels (1 ng/g, 5 ng/g, and 20 ng/g). The full validation was carried out on

an Agilent 6470 Triple Quadrupole LC/MS with pork matrix at three spiking levels (1 ng/g, 5 ng/g, and 20 ng/g). The full validations were carried out on an Agilent 6460 Triple Quadrupole LC/MS with pork and swine liver matrices at three spiking levels (1 ng/g, 5 ng/g, and 20 ng/g).

## Results and Discussion

Table 1 lists the compounds covered in this solution.

Table 1. Veterinary Drugs Analyzed

Classification	English name	CAS	Classification	English name	CAS
β-Agonists	Cimaterol	54239-37-1		Thiabendazole-5-hydroxy	948-71-0
	Clenbuterol	21898-19-1		Triclabendazole	68786-66-3
	Clorprenaline	3811-25-4	Cephems	Cefaclor	53994-73-3
	Penbutolol	38363-40-5		Cefamandole	34444-01-4/ 58648-57-0
	Propranolol	318-98-9		Cefapirin	21593-23-7
	Ractopamine	90274-24-1		Cefazolin	25953-19-9
	Salbutamol	18559-94-9		Cefetamet pivoxyl	65243-33-6
	Terbutaline	23031-32-5		Cefoperazone	62893-19-0/ 62893-20-3
	Tulobuterol	56776-01-3		Cefotaxime	63527-52-6
Androgens	Mengestrol Acetate	595-33-5		Ceftiofur	80370-57-6
	Methyltestosterone	58-18-4		Cephalexin	15686-71-2
	Nadrolone/19-Nortestosterone	434-22-0	Cephalonium	5575-21-3	
	Testosterone	58-22-0	Cephradine	38821-53-3	
	β-trenbolone	10161-33-8	Chloramphenicols	Chloramphenicol	56-75-7
Avermectins	Avermectin B1a	65195-55-3		Florfenicol	73231-34-2
	Doramectin	117704-25-3		Thiamphenicol	15318-45-3
	Eprinomectin	123997-26-2	Glucocorticosteroids	Beclomethasone	4419-39-0
	Ivermectin	70288-86-7		Betamethasone	378-44-9
Benzimidazoles	2-Aminoflubendazole	82050-13-3		Cortisone	53-06-5
	5-Hydroxymebendazole	60254-95-7		Dexamethasone	50-02-2
	Albendazole	54965-21-8		Flumethasone	2135-17-3
	Albendazole sulfone	75184-71-3		Hydrocortisone	50-23-7
	Albendazole sulfoxide	54029-12-8		Methylprednisolone	83-43-2
	Albendazole-2-aminosulfone	80983-34-2		Prednisolone	50-24-8
	Cambendazole	26097-80-3		Prednisone	53-03-2
	Fenbantel	58306-30-2		Triamcinolone	124-94-7
	Fenbendazole	43210-67-9	Triamcinolone acetonide	76-25-5	
	Flubendazole	31430-15-6	Macrolides	Acetylisovaleryltalosin/Tylvalosin	63409-12-1
	Mebendazole	31431-39-7		Erythromycin	59319-72-1
	Mebendazole-amine	52329-60-9		Kitasamycin/Leucomycin	1392-21-8
	Oxfendazole	53716-50-0		Oleandomycin	7060-74-4
	Oxfendazole sulfone/ Fenbendazole sulfone	54029-20-8		Tilmicosin	108050-54-0
	Oxibendazole	20559-55-1		Tylosin	74610-55-2
	Thiabendazole	148-79-8			

Table 1. Veterinary Drugs Analyzed (continued)

Classification	English name	CAS	Classification	English name	CAS
Nitroimidazoles	Dimetridazole	551-92-8		Sulfamoxole	729-99-7
	Dimetridazole- OH (HMMNI)	936-05-0		Sulfantran	122-16-7
	Metronidazole	443-48-1		Sulfaphenazole	526-08-9
	Metronidazole-OH	4812-40-2		Sulfapyridine	144-83-2
	Ronidazole	7681-76-7		Sulfaquinoxaline	59-40-5
Peptides	Bacitracin A	1405-87-4	Sulfathiazole	72-14-0	
	Virginiamycin M1	211411-53-0	Sulfisomidine	515-64-0	
Polyethers	Lasalocid A	25999-20-6	Sulfisoxazole/Sulfafurazole	127-69-5	
	Maduramicin	61991-54-6	Trimethoprim	738-70-5	
	Monensin	22373-78-0	Tetracyclines	Chlortetracycline	57-62-5
Quinolones	Ciprofloxacin	85721-33-1		Doxycycline	564-25-0
	Danofloxacin	112398-08-0		Oxytetracycline	6153-64-6
	Difloxacin	98106-17-3		Tetracycline	60-54-8
	Enoxacin	74011-58-8	β-lactams	Amoxicillin	26787-78-0
	Enrofloxacin	93106-60-6		Ampicillin	69-53-4
	Fleroxacin	79660-72-3		Benzylpenicillin/Procaine	61-33-6
	Flumequine	42835-25-6		Benzylpenicillin/Penicillin G	
	Lomefloxacin	98079-51-7	Oxacillin	66-79-5	
	Nalidixic acid	389-08-2	Sulbactam	68373-14-8	
	Norfloxacin	70458-96-7	Triphenylmethanes	Crystal Violet/Basic violet 3	548-62-9
	Ofloxacin	82419-36-1		Leucomalachite green	129-73-7
	Orbifloxacin	113617-63-3		Malachite green	569-64-2
	Oxolinic acid	14698-29-4	Diterpene	Valnemulin	133868-46-9
	Pefloxacin	70458-92-3	Phenothiazines	Chlorpromazine	50-53-3
	Sarafloxacin	98105-99-8		Xylazine	7361-61-7
	Sparfloxacin	110871-86-8	Furans	Nifurstyrenate	54992-23-3
Sulfonamides	Sulfabenzamide	127-71-9	Quinoxalines	Olaquinox	23696-28-8
	Sulfacetamide	144-80-9	Trematocides	Closantel	57808-65-8
	Sulfachloropyridazine	80-32-0		Nitroxynil	1689-89-0
	Sulfaclozine	102-65-8		Rafoxanide	22662-39-1
	Sulfadiazine	68-35-9	Anticoccidiosis	Clopidol	2971-90-6
	Sulfadimethoxine	122-11-2		Decoquinat	18507-89-6
	Sulfadoxine	2447-57-6		Diclazuril	101831-37-2
	Sulfaguanidine	57-67-0		Ethopabate	59-06-3
	Sulfamerazine	127-79-7		Halofuginone	55837-20-2
	Sulfameter/Sulfamethoxydiazine	651-06-9		Nequinat	13997-19-8
	Sulfamethazine/Sulfadimidine	57-68-1		Nicarbazine	330-95-0
	Sulfamethizole	144-82-1		Robenidine	25875-50-7
	Sulfamethoxazole	723-46-6		Toltrazuril	69004-03-1
	Sulfamethoxy-pyridazine	80-35-3		Toltrazuril sulfone	69004-04-2
	Sulfamonomethoxine	1220-83-3	Toltrazuril sulfoxide	69004-15-5	
			Zoalene	148-01-6	

Table 1. Veterinary Drugs Analyzed (continued)

<b>Classification</b>	<b>English name</b>	<b>CAS</b>
Antivirus	Amantadine	768-94-5
Pesticides	Carbofuran	1563-66-2
	Chlordimeform	6164-98-3
	Coumaphos	56-72-4
	Deltamethrin	52918-63-5
	Diazinon	333-41-5
	Dichlorvos	62-73-7
	Fenthion sulfoxide	3761-41-9
	Malathion	121-75-5
	Pentachlorophenol (PCP)	131-52-2
	Propetamphos	31218-83-4
	Trichlorfon	52-68-6
Lincosamide	Lincomycin	7179-49-9
Dapsones	Dapson	80-08-0
	N-Acetyl dapsone	565-20-8
Organic acids	Heptadecafluorooctanesulfonic acid (PFOS)	1763-23-1
	Perfluorooctanoic acid (PFOA)	335-67-1
Nitros	3-Amino-5-nitro- <i>o</i> -toluamide (ANOT)	3572-44-9
	4,6-Dinitro- <i>o</i> -cresol (DNOC)	534-52-1
	Nitrovin	2315-20-0
	Sodium nitrophenolate/ 4-nitrophenol	100-02-7/ 63317-67-9
Tranquilizer	Azaperol	2804-05-9
	Azaperone	1649-18-9
Non-steroid anti-inflammatory drugs (NSAIDs)	4-Acetylamino antipyrine	83-15-8
	4-Formylaminoantipyrine	1672-58-8
	Carprofen	53716-49-7
	Diclofenac/Diclofenac acid	15307-86-5
	Flufenamic acid	530-78-9
	Flunixin	42461-84-7
	Indomethacin	53-86-1
	Indoprofen	31842-01-0
	Ketoprofen	22071-15-4
	Mefenamic acid	61-68-7
	Meloxicam	71125-38-7
	Piroxicam	36322-90-4
	Sulindac	38194-50-2
	Tenoxicam	59804-37-4
Tolfenamic acid	13710-19-5	
Tolmetin	26171-23-3	
Parasiticide	Levamisole	14769-73-4

A set of nine matrix spiked calibration standards (0.1 ng/g, 0.2 ng/g, 0.5 ng/g, 1.0 ng/g, 2.0 ng/g, 5.0 ng/g, 10 ng/g, 20 ng/g, and 40 ng/g) were analyzed consecutively, and linear fittings were generated with coefficient of correlation values ( $R^2$ ). Figures 1A and 1B show calibration curves obtained in pork on a 6495 Triple Quadrupole LC/MS for 12 compounds that belong to different chemical classes: ractopamine, testosterone, albendazole, cefapirin, chloramphenicol, prednisolone, oleandomycin, metronidazole, ciprofloxacin, sulfabenzamide, amoxicillin, and carprofen.

The LODs and LOQs were calculated by Agilent MassHunter Quantitative Analysis software (Ver. B.06.00 or above) under *Replicate Injection MDL-LOQ-LOD Calculation* using data files from seven replicates of each level (1 ng/g, 5 ng/g, and 20 ng/g). The level chosen for LOD and LOQ calculation should meet the following criteria:

- The concentration level chosen was within the calibration range.
- The %RSD of replicates was less than or equal to 20.
- The concentration level should be as low as possible.









Figure 1A. Calibration curves of ractopamine (A), testosterone (B), albendazole (C), cefapirin (D), prednisolone (E), and oleandomycin (F) from 0.1 ng/g to 40 ng/g in pork on an Agilent 6495 Triple Quadrupole LC/MS.









Figure 1B. Calibration curves of metronidazole (G), ciprofloxacin (H), sulfabenzamide (I), amoxicillin (J), carprofen (K), and chloramphenicol (L) from 0.1 ng/g to 40 ng/g in pork on an Agilent 6495 Triple Quadrupole LC/MS.

The validation results in pork on a 6495 Triple Quadrupole LC/MS are summarized to demonstrate the performance of our workflow solution:

- Good linearity results were achieved, with 98 % of compounds having  $R^2 \geq 0.99$ .
- The recovery data at three different levels were grouped into four categories: % recovery between 50–79, % recovery between 80–120, % recovery between 121–150, and % recovery greater than 150. The results shown in Figure 2 demonstrate good recoveries achieved through matrix spiked calibration. For example, at 5 ng/g spiking level, the recoveries of 92 % compounds were within 80–120.
- The repeatability at three different levels were also grouped in four categories:
  - % RSD between 0–10
  - %RSD between 11–15
  - %RSD between 16–20
  - %RSD greater than 20

Figure 3 shows the results for all veterinary drugs tested in this method. Only a few compounds had %RSD greater than 20 %. For example, at 5 ng/g spiking level, 11 compounds (sulfanitran, cephradine, diclazuril, rafoxanide, closantel, tylvalosin, oxacillin, toltrazuril-sulfoxide, sodium nifurstyrenate, amoxicillin, and doramectin) had RSDs greater than 20 %, that could be caused either by weak signal response or greater variance in signal response for unstable analytes.

In this application note, an additional nitrogen concentration step was applied for the 6460 Triple Quadrupole LC/MS and the 6470 Triple Quadrupole LC/MS to match the detection sensitivity of the 6495 Triple Quadrupole LC/MS. The validation results showed that similar performance results could be achieved through this additional step. Table 2 lists validation data for some representative compounds on different tandem mass spectrometers. Note that some unstable compounds, such as penicillins, undergo decomposition during the nitrogen concentration step. The following nine compounds failed to be detected on the 6460 Triple Quadrupole LC/MS and the 6470 Triple

Quadrupole LC/MS due to decomposition under the nitrogen concentration process: carprofen, closantel, deltamethrin, diazinon, sodium nifurstyrenate, oxacillin, rafoxanide, sulbactam, and  $\beta$ -trenbolone. In the final reconstitution step, using membrane filtration (not centrifugation) also caused great losses of the following compounds: avermectin B1a, cephradine, decoquinatate, dichlorvos, diclazuril, doramectin, eprinomectin, flufenamic acid, heptadecafluorooctanesulfonic acid (PFOS), ivermectin, lasalocid, maduramicin, mefenamic acid, monensin, nequinatate, pentachlorophenol (PCP), perfluorooctanoic acid (PFOA), robenidine, and tolfenamic acid.




Figure 2. Analyte recoveries at different spiking levels in pork on an Agilent 6495 Triple Quadrupole LC/MS.




Figure 3. Analysis repeatability at different spiking levels in pork on an Agilent 6495 Triple Quadrupole LC/MS.



Table 2. Representative Validation Data in Pork on Different Tandem Mass Spectrometers (Agilent 6460 Triple Quadrupole LC/MS, Agilent 6470 Triple Quadrupole LC/MS, and Agilent 6495 Triple Quadrupole LC/MS)

Name	Model	R <sup>2</sup>	Dynamic range (ng/g)	Spiking level						LOQ (ng/g)	LOD (ng/g)
				1 ng/g		5 ng/g		20 ng/g			
				% Recovery*	%RSD	% Recovery*	%RSD	% Recovery*	%RSD		
Albendazole	G6460CA	0.998	0.1–40	97	3	99	3	109	4	0.3	0.1
	G6470AA	0.999	0.2–40	101	7	102	3	102	5	0.7	0.2
	G6495A	0.999	0.1–40	98	6	105	5	108	4	0.6	0.2
Amoxicillin	G6460CA	1.000	5–40	103	21	87	24	108	10	22.0	7.3
	G6470AA	0.999	2–40	78	33	103	11	107	5	5.9	1.8
	G6495A	1.000	1–40	154	73	80	21	117	8	8.5	2.6
Cefapirin	G6460CA	0.999	0.5–40	121	44	107	16	123	7	8.6	2.9
	G6470AA	0.995	2–40	148	39	130	10	101	7	6.7	2.0
	G6495A	0.999	1–40	101	43	150	10	110	19	7.3	2.2
Chloramphenicol	G6460CA	0.999	2–40	211	36	111	23	85	13	12.9	4.3
	G6470AA	0.998	1–40	92	16	90	13	106	9	1.4	0.4
	G6495A	0.997	0.2–40	84	15	96	11	95	5	1.3	0.4
Ciprofloxacin	G6460CA	0.999	1–40	191	27	93	19	101	9	8.9	3.0
	G6470AA	0.989	1–40	189	66	94	9	96	8	1.4	0.4
	G6495A	0.998	0.5–40	91	15	133	8	103	7	1.4	0.4
Metronidazole	G6460CA	0.999	0.1–40	92	9	104	6	102	3	0.8	0.3
	G6470AA	0.999	0.5–40	101	11	104	7	103	6	1.1	0.3
	G6495A	0.996	0.5–40	87	7	105	4	113	12	0.6	0.2
Oleandomycin	G6460CA	0.998	0.1–40	105	5	100	2	104	4	0.6	0.2
	G6470AA	0.995	0.2–40	119	7	111	8	97	6	0.9	0.3
	G6495A	1.000	0.1–40	120	8	100	8	95	4	0.9	0.3
Prednisolone	G6460CA	1.000	0.5–40	113	8	91	3	101	2	0.9	0.3
	G6470AA	1.000	0.2–40	88	12	112	3	104	7	1.0	0.3
	G6495A	1.000	0.2–40	101	3	100	5	91	6	0.3	0.1
Ractopamine	G6460CA	0.997	0.1–40	102	6	103	7	100	5	0.7	0.2
	G6470AA	0.999	0.2–40	78	16	116	8	99	5	1.3	0.4
	G6495A	0.999	0.1–40	94	11	93	4	95	11	1.0	0.3
Sulfabenzamide	G6460CA	0.998	0.5–40	100	6	87	3	105	4	0.6	0.2
	G6470AA	0.999	1–40	139	23	107	8	98	5	4.4	1.3
	G6495A	1.000	0.1–40	104	4	96	3	115	10	0.4	0.1
Testosterone	G6460CA	0.999	1–40	91	11	110	13	111	6	1.0	0.3
	G6470AA	0.999	1–40	62	13	105	13	101	5	0.8	0.2
	G6495A	0.999	0.5–40	112	28	98	5	102	10	2.5	0.8

\* Average recovery of seven replicates.

## Conclusions

The China multiclass residue veterinary drug workflow solution is developed to screen and quantify multiclass veterinary drugs in different animal-derived food samples. It supports all Agilent liquid chromatography tandem mass spectrometers coupled with an Agilent Jet Stream Ionization Source. We have demonstrated that our total workflow solution, from sample preparation to results, can be applied to the analyses of most types of animal-derived food samples with high sensitivity, and good accuracy and precision. With >180 veterinary drugs spanning 32 different classes in one method, our workflow solution would significantly help increase sample throughput without sacrificing data quality, robustness, or sensitivity.

## Acknowledgements

The authors sincerely thank Chunxiao Wang (Project leader), Vivian Xianyu Chen, Jianguo Ji, Zhixu Zhang, Chunye Sun, Jimmy Chan, Jerry Zweigenbaum, and Thomas Glauner (Agilent Technologies, Inc.) for their contributions.

## References

1. Announcement No. 235, issued by the Ministry of Agriculture. *Maximum residue limits in animal derived food samples*. China (2002).
2. L. Zhao, D. Lucas. *Multiresidue analysis of veterinary drugs in bovine liver by LC/MS/MS*, *Agilent Technologies Application Note*, publication number 5991-6096 (2015).
3. T. Anumol, J. M. Stevens, J. Zweigenbaum. *Analysis of Veterinary Drugs in Meat with the Agilent 6495 Triple Quadrupole LC/MS*. *Agilent Technologies Application Note*, publication number 5991-7895EN (2017).

## For More Information

These data represent typical results. For more information on our products and services, visit our Web site at [www.agilent.com/chem](http://www.agilent.com/chem).

[www.agilent.com/chem](http://www.agilent.com/chem)

Agilent shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material.

Information, descriptions, and specifications in this publication are subject to change without notice.

© Agilent Technologies, Inc., 2017  
Printed in the USA  
May 18, 2017  
5991-8084EN



**Agilent Technologies**