


## RapiGest SF Protein Digestion Surfactant

Waters patented\* *RapiGest*™ SF (surfactant) radically enhances protein enzymatic digestions in terms of speed and percent recovery.


*RapiGest* is a patented anionic surfactant that accelerates the in-solution production of peptides generated by proteases, such as trypsin, Asp-N, Glu-C, and Lys-C. Many hydrophobic proteins are resistant to proteolysis because their cleavage sites are inaccessible to endoproteases. *RapiGest*, a mild denaturant, helps solubilize and unfold proteins making them more amenable to cleavage without denaturing or inhibiting common proteolytic enzymes.

- Improves solubility of hydrophobic proteins for improved enzymatic digest
- Compatible with various enzymes and does not cause protein modifications
- Does not inhibit enzyme activity, unlike conventional denaturants
- Reduces the digestion time and requires less enzyme to achieve optimum digestion
- Improves the digestion of enzyme resistant proteins such as membrane proteins
- Acid labile and the degradation products from *RapiGest* do not interfere with LC-MS or MALDI MS analysis


### 1 Hour Proteolysis of Myoglobin Using Various Endoproteases



### How *RapiGest* Works



### UPLC/MS Analysis of Trypsin Digested mAb in the Presence of 0.1% *RapiGest*




A human mAb IgG sample was solubilized in 0.1% *RapiGest* (w/v) during the reduction (DTT) and alkylation (iodoacetamide) steps that preceded overnight trypsin digestion. *RapiGest* enhances the protein solubilization, therefore, improving the overall sequence coverage. Since *RapiGest* is rapidly degraded by mixing the digested sample with 1% formic acid (1:1 volume), direct sample injection onto a Waters Peptide Separation Technology  $C_{18}$  Column is possible without affecting the quality of the LC/MS analysis.

Hundreds of published scientific articles and presentations have documented the benefits of using *RapiGest* for improving the protein sequence coverage and reducing the sample preparation time. *RapiGest*'s application areas are diverse, ranging from proteomic research to therapeutic protein characterization, and it is effective for in-solution digestion protocols.

**RapiGest shows significant advantage over other denaturants since it is not disruptive to endoprotease activity**


Trypsin <sup>1</sup> Solution	Trypsin Activity (%) <sup>2</sup>	Trypsin <sup>1</sup> Solution	Trypsin Activity (%) <sup>2</sup>
No additive	<b>100</b>	0.1% SDS/0.1% RapiGest	67
0.1% RapiGest	<b>100</b>	50% Methanol	59
0.5% RapiGest	<b>100</b>	50% Acetonitrile	87
0.1% SDS	24	1M Urea	97
0.5% SDS	1	2M Urea	83

1 0.5 µg of trypsin in 50 mM ammonium bicarbonate, pH 7.9; 0.2 mM of BEAA.  
 2 Measured as delta BEAA absorbance @ 253 nm (slope within 5 min).



RapiGest hydrolyzes under acidic condition (half life time = 7.6 minutes at pH 2). Therefore, it is compatible with LC/MS and MALDI MS analysis.

**Use of RapiGest SF to Assist in Protein Deglycosylation**



**ORDERING INFORMATION**

Description	Part No.
RapiGest SF 1 mg vial	186001860
RapiGest SF 1 mg vial (5 pack)	186001861
RapiGest SF 10 mg vial	186002123
RapiGest SF 50 mg vial	186002122
RapiGest SF Custom	186002118

\* US 7,229,539 and US 8,580,533

**Waters**

THE SCIENCE OF WHAT'S POSSIBLE.®



Waters and The Science of What's Possible are registered trademarks of Waters Corporation. RapiGest is a trademark of Waters Corporation. All other trademarks are property of their respective owners.

©2014 Waters Corporation. Produced in the U.S.A.  
 February 2014 720000468EN SC-PDF

**Waters Corporation**  
 34 Maple Street  
 Milford, MA 01757 U.S.A.  
 T: 1 508 478 2000  
 F: 1 508 872 1990  
 www.waters.com