

[ DISQUE DISPERSIVE SAMPLE PREPARATION EXTRACTION PRODUCTS ]

QuEChERS Simplified

**DisQUE**<sup>TM</sup>  
Dispersive Sample Preparation

Waters  
THE SCIENCE OF WHAT'S POSSIBLE.®



The QuEChERS method (an acronym for **Q**uick, **E**asy, **C**heap, **E**ffective, **R**ugged, and **S**afe) followed by dispersive solid-phase extraction (d-SPE), is a simple and straightforward sample preparation technique ideal for multi-residue analysis for pesticides, veterinary drugs, and mycotoxins in a wide variety of food and agricultural products. DisQuE Dispersive Sample Preparation Extraction Products are conveniently packaged with pre-weighed sorbents and buffers in pouches and tubes as described in regulatory methods and protocols.



## The Science Behind DisQuE



To overcome the deficiencies of traditional methods, an improved extraction procedure was developed (Anastassiades et al.)\* to extract pesticides from fruits and vegetables. The procedure is based on an initial single-phase extraction using 10 or 15 g sample with acetonitrile at 1 mL acetonitrile per 1 g of sample. A liquid-liquid partition is created by adding excess salts and buffers to the extract. After centrifugation, the acetonitrile layer, containing the pesticide, is collected. The matrix can be further cleaned and the excess water removed with a single d-SPE step by mixing acetonitrile extract with anhydrous  $MgSO_4$  and primary secondary amine (PSA) sorbents. The biggest advantage, besides time and effort, is that the final extract can be analyzed directly by either GC-MS or by LC-MS/MS technology with a simple dilution.

Some QuEChERS methods add buffers during the extraction step. Today, the two most common methods are the European Committee for Standardization (CEN) Method 15662 which uses citrate buffer for extraction and the Association of Analytical Communities (AOAC) Official Method 2007.01 using acetic acid buffer.

\*Anastassiades, M., Lehotay, S.J., Tajnbaheer, D., & Schenck, F.J. (2003) *J. AOAC. Int.* **86**, 412-431

Method	Sample Size	Solvent	Tube Content
CEN Method 15662 Citrate Buffer	10 g	10 mL acetonitrile	4 g $MgSO_4$ , 1 g NaCl 1.5 g sodium citrate
AOAC Method 2007.01 Acetate Buffer	15 g	15 mL 1% acetic acid in acetonitrile	6 g $MgSO_4$ 1.5 g sodium acetate

*Kit choice based on preferred method*

### Low Water Content Commodity Extraction


QuEChERS extraction is designed for multi-residue pesticide analysis of fruits and vegetables with high water content. For QuEChERS analysis of commodities with low-water content, additional water is added to optimize the extraction. CEN method guidelines for adding water are shown below.

Sample Type	Sample Weight	Water Added	Note
Fruits and vegetables > 80% water content	10 g	–	
Fruits and vegetables 25-80% water content	10 g	X g	X = 10 g - water amount in 10 g sample
Cereals	5 g	10 g	
Dried fruits	5 g	7.5 g	Water can be added during comminution step
Honey	5 g	10 g	
Spices	2 g	10 g	

*CEN Method Guidelines for Adding Water into Commodities with Low Water Content*



## Simplified QuEChERS Method Sample Preparation Protocol



### 402 Pesticide Residues at 10 ppb ng/g In One 10 Minute Run

#### UPLC® Conditions


LC system: Waters ACQUITY UPLC® System  
 Column: ACQUITY UPLC BEH C<sub>18</sub>  
 Part number: 186002352  
 Column temp.: 40 °C  
 Sample temp.: 4 °C  
 Flow rate: 0.450 mL/min  
 Mobile phase A: 98:2 water: methanol + 0.1% formic acid  
 Mobile phase B: Methanol + 0.1% formic acid  
 Gradient:

Time	A%	B%
0.00	90	10
0.25	90	10
7.75	0	100
8.50	0	100
8.51	90	10


Total run time: 10 min  
 Injection volume: 20 µL, full loop injection

#### MS Conditions

MS system: Waters ACQUITY® TQ Detector  
 Ionization mode: ESI positive polarity  
 Capillary voltage: 1 kV  
 Desolvation gas: Nitrogen, 800 L/hr, 400 °C  
 Cone gas: Nitrogen, 5 L/hr  
 Source temp.: 120 °C  
 Acquisition: Multiple Reaction Monitoring (MRM)  
 Collision gas: Argon at 3.5 mBar



### Recovery Data for Three Types of Sample Matrices Fortified at 10 ng/g



## Selection of d-SPE Product for Sample Cleanup

Commodity Type	Cleanup Provided	Size	AOAC Method 2007.01	CEN Method 15662
General Fruits and Vegetables (Celery, Head Lettuce, Melon)	Removes polar organic acids, some sugars and lipids	2 mL Tubes	50 mg PSA, 150 mg MgSO <sub>4</sub> Part #186004572	25 mg PSA, 150 mg MgSO <sub>4</sub> Part #186004831
		15 mL Tubes	400 mg PSA, 1200 mg MgSO <sub>4</sub> Part #186008072	150 mg PSA, 900 mg MgSO <sub>4</sub> Part #186004833
Fruits and Vegetables with Fats and Waxes (Cereals, Nuts, Dairy, Avocado)	Removes polar organic acids, some sugars, more lipids and sterols	2 mL Tubes	50 mg PSA, 50mg C <sub>18</sub> , 150 mg MgSO <sub>4</sub> Part #186004830	25 mg PSA, 25 mg C <sub>18</sub> , 150 mg MgSO <sub>4</sub> Part #186004832
		15 mL Tubes	400 mg PSA, 400 mg C <sub>18</sub> , 1200 mg MgSO <sub>4</sub> Part # 186008073	150 mg PSA, 150 mg C <sub>18</sub> , 900mg MgSO <sub>4</sub> Part #186004834

### Cleaning Up Pigmented Fruits and Vegetables

Some fruits and vegetables, such as spinach, red sweet pepper, and carrots, have high content of non-polar pigments, such as carotenoids or chlorophyll. QuEChERS effectively removes some common matrix constituents, such as fatty acids and sugars; however, samples that contain high levels of pigment require additional treatment. Graphitized carbon black (GCB) effectively removes these pigments, and can be added to the clean-up tube

with PSA and MgSO<sub>4</sub> during the d-SPE cleanup step. The clean-up of pigments is especially critical for gas chromatographs because pigments often accumulate in the injection port, quickly fouling the injection liner.



The effect of GCB amounts on the de-coloring of spinach extracts.







Successful laboratories demand timely, high quality results. Waters scientists carefully evaluate product performance and can recommend methods for sample preparation, analytical separation, and instrument optimization for a wide array of food adulterants. This expert advice is readily available in our Food Testing Application Notebook.

For your copy, visit [www.waters.com/FT](http://www.waters.com/FT)






## Ordering Information


### Kits

Name	Description	Qty	Part No.	
DisQuE Dispersive SPE Kit—Pouch Format (AOAC Method 2007.01)	Pouch: 1.5 g sodium acetate and 6 g MgSO <sub>4</sub> 50 mL Tube: Empty 2 mL Tube: 150 mg MgSO <sub>4</sub> and 50 mg PSA	100	176002922	
DisQuE Dispersive SPE Kit—Pouch Format (CEN Method 15662)	Pouch: 1 g trisodium citrate dihydrate, 0.5 g disodium hydrogencitrate sesquihydrate, 1 NaCl and 4 g MgSO <sub>4</sub> 50 mL Tube: Empty 2 mL Tube: 150 mg MgSO <sub>4</sub> , 25 mg PSA, and 25 mg C <sub>18</sub>	100	176002923	
DisQuE Dispersive SPE Kit (AOAC Method 2007.01)	50 mL Tube: 1.5 g sodium acetate, 6 g MgSO <sub>4</sub> 2 mL Tube: 150 mg MgSO <sub>4</sub> , 50 mg PSA	100	176001676	
DisQuE Dispersive SPE Kit (CEN Method 15662)	50 mL Tube: 4 g MgSO <sub>4</sub> , 1 g NaCl, 1 g trisodium citrate dehydrate, 0.5 g disodium hydrogencitrate sesquihydrate 2 mL Tube: 150 mg MgSO <sub>4</sub> , 25 mg PSA, 25 mg C <sub>18</sub>	100	176001903	

### Extraction Tubes


Name	Description	Qty	Part No.	
50 mL Empty Tube for QuEChERS Extraction	—	50	186006814	
50 mL Tube—Blue Cap (AOAC Method)	1.5 g sodium acetate, 6 g MgSO <sub>4</sub>	100	186004571	
50 mL Tube—White Cap (CEN Method)	4 g MgSO <sub>4</sub> , 1 g NaCl, 1 g trisodium citrate dehydrate, 0.5 g disodium hydrogencitrate sesquihydrate	100	186004837	

### Sorbent Pouches

Name	Contents	Qty	Part No.	
DisQuE Pouch (AOAC Method)	1.5 g sodium acetate, 6 g MgSO <sub>4</sub>	50	186006812	
DisQuE Pouch (CEN Method)	4 g MgSO <sub>4</sub> , 1 g NaCl, 1 g trisodium citrate dehydrate, 0.5 g disodium hydrogencitrate sesquihydrate	50	186006813	

## Clean-Up Tubes

Name	Contents	Qty	Part No.
2 mL Tube–White Cap	150 mg MgSO <sub>4</sub> , 50 mg PSA	100	186004572
2 mL Tube–Blue Cap	150 mg MgSO <sub>4</sub> , 50 mg PSA, 25 mg C <sub>18</sub>	100	186004830
2 mL Tube–Yellow Cap	150 mg MgSO <sub>4</sub> , 25 mg PSA	100	186004831
2 mL Tube–Green Cap	150 mg MgSO <sub>4</sub> , 25 mg PSA, 25 mg C <sub>18</sub>	100	186004832
2 mL Tube–Blue Cap	150 mg MgSO <sub>4</sub> , 25 mg PSA, 25 mg C <sub>18</sub> , 7 mg GCB	100	186008071
2 mL Tube–Yellow Cap	150 mg MgSO <sub>4</sub> , 50 mg C <sub>18</sub>	100	186008075
2 mL Tube–Blue Cap	150 mg MgSO <sub>4</sub> , 25 mg PSA, 2.5 mg GCB	100	186008076
2 mL Tube–Blue Cap	150 mg MgSO <sub>4</sub> , 50 mg PSA, 30 mg C <sub>18</sub> , 30 mg Al-N	100	186008081
15 mL Tube–White Cap	900 mg MgSO <sub>4</sub> , 150 mg PSA	50	186004833
15 mL Tube–Blue Cap	900 mg MgSO <sub>4</sub> , 150 mg PSA, 150 mg C <sub>18</sub>	50	186004834
15 mL Tube–Clear Cap	1200 mg MgSO <sub>4</sub> , 400 mg PSA	50	186008072
15 mL Tube–Blue Cap	1200 mg MgSO <sub>4</sub> , 400 mg PSA, 400 mg C <sub>18</sub>	50	186008073
15 mL Tube–Blue Cap	1200 mg MgSO <sub>4</sub> , 400 mg PSA, 400 mg C <sub>18</sub> , 400 mg GCB	50	186008074
15 mL Tube–Clear Cap	900 mg MgSO <sub>4</sub> , 300 mg PSA	50	186008077
15 mL Tube–Blue Cap	900 mg MgSO <sub>4</sub> , 300 mg PSA, 300 mg C <sub>18</sub>	50	186008078
15 mL Tube–Blue Cap	900 mg MgSO <sub>4</sub> , 450 mg PSA, 300 mg C <sub>18</sub> , 50 mg GCB	50	186008079
15 mL Tube–Blue Cap	750 mg MgSO <sub>4</sub> , 250 mg PSA, 150 mg C <sub>18</sub> , 150 mg Al-N	50	186008080



Al-N = Alumina N

## Bulk Sorbent

Name	Qty	Part No.
Graphitized Carbon Black	25 g	186004835
Trifunctionally-bonded C <sub>18</sub> Silica	100 g	WAT035672



**SALES OFFICES:**

Austria 43 1 877 18 07  
Australia 61 2 9933 1777  
Belgium and Luxembourg 32 2 726 1000  
Brazil 55 11 4134 3788  
Canada 1 800 252 4752  
China 86 21 6156 2666  
Czech Republic 420 2 617 11384  
Denmark 45 46 59 8080  
Finland 358 9 5659 6288  
France 33 1 30 48 72 00  
Germany 49 6196 400 600  
Hong Kong 852 2964 1800  
Hungary 36 1 350 5086  
India 91 080 49292200 03  
Ireland 353 1 448 1500  
Israel 9723 3731391  
Italy 39 02 265 0983  
Japan 81 3 3471 7191  
Korea 82 2 6300 9200  
Malaysia 603 7841 1688  
Mexico 52 55 52 00 1860  
The Netherlands 31 76 508 7200  
Norway 47 6 384 6050  
Poland 48 22 101 5900  
Portugal 351 21 893 61 77  
Puerto Rico 1 787 747 8445  
Singapore 65 6593 7100  
Spain 34 93 600 9300  
Sweden 46 8 555 115 00  
Switzerland 41 56 676 7000  
Taiwan 886 2 2508 5500  
UAE 971 4 214 62 38  
UK 44 208 238 6100  
US 1 800 252 4752

[www.waters.com/disque](http://www.waters.com/disque)

**Waters**

**THE SCIENCE OF WHAT'S POSSIBLE.®**

**Waters Corporation**

34 Maple Street  
Milford, MA 01757 U.S.A.  
T: 508 478 2000  
F: 508 872 1990  
[www.waters.com](http://www.waters.com)

Waters, The Science of What's Possible, UPLC, ACQUITY, and ACQUITY UPLC are registered trademarks of Waters Corporation. DisQuE is a trademark of Waters Corporation. All other trademarks are the property of their respective owners.

©2015 Waters Corporation. Printed in the U.S.A.  
July 2016 720003048EN AW-SIG