6.1015.000 / 6.01015.030 Metrosep C Trap 1 - xx0/4.0

6.1015.000 Metrosep C Trap 1 - 100/4.0

6.01015.030 Metrosep C Trap 1 - 30/4.0

DE

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EN

Column material

Polystyrene-divinylbenzene with sulfonic acid groups

Dimensions

6.1015.000: 100 × 4.0 mm

6.01015.030: 30 × 4.0 mm

pH range

1 - 14

Maximum pressure

25 MPa (250 bar)

Application

Elimination of interfering cations from eluents and reagent solutions.

The Metrosep C Trap 1 - 30/4.0 can also be used to eliminate polyvalent cations in anion analysis. Insert the trap column between the injector and the analytical column for this.

Preparation

- Elimination of interfering cations from the eluent:

Rinse the trap column with eluent separately in advance. Insert the trap column directly in front of the injection valve into the eluent path.

- Elimination of interfering cations from a rinsing solution:
Rinse the trap column with the rinsing solution separately in advance. Insert the trap column directly in front of port 1 of the Dosino.
- Elimination of polyvalent cations in ion analysis:
Rinse the trap column with eluent separately in advance. Insert the trap column between the injector and the analytical column or the guard column.

Storage

Store the trap column in ultrapure water.

Regeneration

- Rinse the trap column against the flow direction with 20 mL 5% H₂SO₄.
- Rinse the trap column against the flow direction with 20 mL ultrapure water.
- Additionally rinse the trap column with 20 mL eluent if it is being used in anion analysis.

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