

# Hydrocarbons, $C_6 - C_7$ , aromatic hydrocarbons, $C_6 - C_8$

## Analysis of impurities in benzene

### Application Note

Materials Testing & Research

#### Authors

Agilent Technologies, Inc.

#### Introduction

The fast analysis of impurities in ethanol can be done very well using 0.15 mm id capillary columns in combination with split injection and FID detection. The 0.15 mm capillary offers a high plate number (ca. 150,000/25 m), can be operated at practical pressures (150 - 250 kPa) and can be used for a wide range of different applications. The problem of reduced loadability has been fully overcome by using a thick film Agilent CP-Sil 5 CB (1.2  $\mu$ m) column. This allows accurate trace analysis of impurities in many chemical products.

The reproducibility of the analysis is within 3% standard deviation, even for compounds which are present at 5 - 500 ppm. Despite split injection, impurities can be measured at 1 - 5 ppm. Typically, a 100% method is used for integration, with hydrogen as the carrier gas as analysis time. Helium can also be used. Typical analysis times are within 10 - 15 minutes.




**Agilent Technologies**

## Conditions

Technique : GC-capillary  
Column : Agilent CP-Sil 5 CB, 0.15 mm x 25 m, fused silica  
WCOT (df = 1.2 µm) (Part no. CP7693)  
Temperature : 70 °C (2 min) → 200 °C, 20 °C/min  
200 °C (5 min)  
Carrier Gas : H<sub>2</sub>, 150 kPa (1.5 bar, 21 psi)  
Injector : Split, T = 250 °C  
Detector : FID, T = 250 °C  
Sample Size : 2.0 µL  
Concentration Range : 1-200 ppm  
Solvent Sample : benzene balance

## Peak identification

1. impurity	5 ppm
2. methylcyclopentane	37 ppm
3. benzene	
4. cyclohexane	96 ppm
5. methylcyclohexane	427 ppm
6. toluene	49 ppm
7. ethylbenzene	11 ppm



[www.agilent.com/chem](http://www.agilent.com/chem)

This information is subject to change without notice.

© Agilent Technologies, Inc. 2011

Printed in the USA

31 October, 2011

First published prior to 11 May, 2010

A01394



**Agilent Technologies**