
Application Bulletin

Of interest to: Metals, electroplating

A 10

Potentiometric analysis of anodizing baths

Summary

This bulletin describes potentiometric titration methods for checking sulfuric acid and chromic acid anodizing baths. In addition to the main components aluminum, sulfuric acid and chromic acid, chloride, oxalic acid and sulfate are also determined.

Instruments

- 702 SET/MET Titrino, 716 DMS Titrino, 736 GP Titrino, 751 GPD Titrino or 785 DMP Titrino or 726 or 796 Titroprocessor with 700 Dosino or 685 Dosimat
- 2.728.0040 Magnetic Stirrer

The necessary accessories are listed under the different baths.

1. Sulfuric acid anodizing bath

Accessories

- 6.3014.223 Exchange Units (for sodium hydroxide solution possibly with 6.1608.040 reagent bottle made of polyethylene)
 - 6.0222.100 combined LL pH glass electrode with 6.2104.020 electrode cable
 - 6.0430.100 Ag Titrode with Ag₂S coating
 - 6.0431.100 Pt Titrode
-

Reagents

- Titrant: sodium hydroxide solution, c(NaOH) = 1 mol/L
- Titrant: silver nitrate solution, c(AgNO₃) = 0.1 mol/L
- Titrant: potassium permanganate solution c(KMnO₄) = 0.02 mol/L (0.1 N)
- Sulfuric acid, w(H₂SO₄) = 10%
- Manganese(II) sulfate monohydrate MnSO₄ * H₂O, p.a.

Analysis**a) Determination of sulfuric acid and aluminum**

Pour approx. 50 mL dist. water into the glass beaker and add 2.0 mL of the anodizing bath sample. Using the combined pH glass electrode titrate with $c(\text{NaOH}) = 1 \text{ mol/L}$ just past the second equivalence point.

Calculation

The first equivalence point corresponds to the H_2SO_4 content, the difference between the second and the first equivalence point to the Al^{3+} content.

1 mL $c(\text{NaOH}) = 1 \text{ mol/L}$ corresponds to 49.04 mg H_2SO_4 or 8.994 mg Al^{3+}

$$\% \text{H}_2\text{SO}_4 = \text{EP1} * \text{C01} * \text{C03} / \text{C00}$$

$$\% \text{Al}^{3+} = (\text{EP2} - \text{EP1}) * \text{C02} * \text{C03} / \text{C00}$$

EP1 = titrant consumption to reach the first EP in mL

EP2 = titrant consumption to reach the second EP in mL

C00 = 2.0 (sample volume in mL)

C01 = 49.04

C02 = 8.994

C03 = 0.1 (conversion factor for %)

b) Determination of chloride

Dilute 50.0 mL of the anodizing bath sample with approx. 50 mL dist. water and titrate with $c(\text{AgNO}_3) = 0.1 \text{ mol/L}$ using the Ag Titrode with Ag_2S coating.

Calculation

1 mL $c(\text{AgNO}_3) = 0.1 \text{ mol/L}$ corresponds to 3.545 mg Cl^-

$$\% \text{Cl}^- = \text{EP1} * \text{C01} * \text{C02} / \text{C00}$$

EP1 = titrant consumption in mL

C00 = 50.0 (sample volume in mL)

C01 = 3.545

C02 = 0.1 (conversion factor for %)

c) Determination of oxalic acid

Pipet 50.0 mL of the anodizing bath sample into the glass beaker, add 10 mL $w(\text{H}_2\text{SO}_4) = 10\%$ and approx. 0.5 g $\text{MnSO}_4 \cdot \text{H}_2\text{O}$ and titrate with $c(\text{KMnO}_4) = 0.02 \text{ mol/L}$ using the Pt Titrode.

Calculation

1 mL $c(\text{KMnO}_4) = 0.02 \text{ mol/L}$ corresponds to 4.502 mg oxalic acid

% oxalic acid = $\text{EP1} \cdot \text{C01} \cdot \text{C02} / \text{C00}$

EP1 = titrant consumption in mL

C00 = 50.0 (sample volume in mL)

C01 = 4.502

C02 = 0.1 (conversion factor for %)

2. Chromic acid anodizing bath**Accessories**

- 6.3014.223 Exchange Units
- 6.0431.100 Pt Titrode with 6.2104.020 electrode cable
- 6.0430.100 Ag Titrode with Ag_2S coating
- 6.1248.050 W electrode rod with 6.1241.030 electrode shaft and 6.2114.000 electrode cable as well as 6.1248.000 Pt electrode rod with 6.1241.030 electrode shaft and 6.2106.020 electrode cable

Reagents

- Titrant: sodium thiosulfate solution, $c(\text{Na}_2\text{S}_2\text{O}_3) = 0.1 \text{ mol/L}$
- Sulfuric acid, $w(\text{H}_2\text{SO}_4) = 96\%$
- Potassium iodide, p.a.
- Ethanol
- Titrant: silver nitrate solution, $c(\text{AgNO}_3) = 0.01 \text{ mol/L}$
- Barium chloride solution, $c(\text{BaCl}_2) = 0.05 \text{ mol/L}$:
Dissolve 12.214 g $\text{BaCl}_2 \cdot 2 \text{H}_2\text{O}$ in $c(\text{HCl}) = 0.1 \text{ mol/L}$ and make up to 1 L.
- Titrant: $c(\text{EGTA}) = 0.05 \text{ mol/L}$:
Dissolve 19.02 g ethylene glycol-bis-(2-aminoethyl)-tetraacetic acid in 250 mL $c(\text{NaOH}) = 1 \text{ mol/L}$ and, after cooling down, make up to 1 L with dist. water.
- Buffer solution pH = 10.5:
Dissolve 10 g NH_4Cl and 60 mL $w(\text{NH}_3) = 25\%$ in dist. water and make up to 1 L.
- Strong acidic cation exchanger (e.g. Dowex 50)

Analysis**a) Determination of chromic acid**

Pipet 10.0 mL of the anodizing bath sample into a 100 mL volumetric flask, fill to the mark with dist. water and mix thoroughly. Pipet 2.0 mL of this dilution (corresponding to 0.2 mL of the original sample) into a glass beaker, add 40 mL dist. water, 1 mL w(H₂SO₄) = 96% and approx. 1 g potassium iodide and titrate with c(Na₂S₂O₃) = 0.1 mol/L using the Pt Titrode.

Calculation

1 mL c(Na₂S₂O₃) = 0.1 mol/L corresponds to 3.333 mg CrO₃

$$\% \text{CrO}_3 = \text{EP1} * \text{C01} * \text{C02} / \text{C00}$$

EP1 = titrant consumption in mL

C00 = 0.2 (sample volume used in mL original sample)

C01 = 3.333

C02 = 0.1 (conversion factor for %)

b) Determination of chloride

Pipet 5.0 mL of the anodizing bath sample into a glass beaker. Add 20 mL dist. water, 20 mL ethanol and 1 mL w(H₂SO₄) = 96% and boil for 5 min in order to reduce all Cr(VI) present to Cr(III) (use the fume cupboard!). After cooling down, titrate with c(AgNO₃) = 0.01 mol/L using the Ag Titrode with Ag₂S coating.

Calculation

1 mL c(AgNO₃) = 0.01 mol/L corresponds to 0.355 mg Cl⁻ or 0.5844 mg NaCl

$$\text{mg/L Cl}^- = \text{EP1} * \text{C01} * \text{C03} / \text{C00}$$

$$\text{mg/L NaCl} = \text{EP1} * \text{C02} * \text{C03} / \text{C00}$$

EP1 = titrant consumption in mL

C00 = 5.0 (sample volume in mL)

C01 = 0.355

C02 = 0.5844

C03 = 1000 (conversion factor in mL/L)

c) Determination of sulfate

Pipet 1.0 mL of the anodizing bath sample of the bath into a glass beaker. Add 20 mL dist. water, 20 mL ethanol and 1 mL conc. HCl and boil for 5 min in order to reduce all Cr(VI) present to Cr(III) (use the fume cupboard!). After cooling down, percolate through a cation exchanger column into a 100 mL volumetric flask, rinse thoroughly with dist. water, fill to the mark and mix.

Pipet 20 ... 50 mL of the prepared bath solution (corresponding to 0.2 ... 0.5 mL of the original sample) into a glass beaker, add 5.00 mL $c(\text{BaCl}_2) = 0,05 \text{ mol/L}$ and allow to react for 3 min under stirring. Then add 10 mL buffer solution $\text{pH} = 10.5$ and titrate back the excess of Ba^{2+} with $c(\text{EGTA}) = 0.05 \text{ mol/L}$ using the W and Pt electrode (MET mode).

Calculation

1 mL $c(\text{EGTA}) = 0.05 \text{ mol/L}$ corresponds to 4.803 mg SO_4^{2-} or 4.904 mg H_2SO_4

$$\text{g/L SO}_4^{2-} = (\text{C30} - \text{EP1}) * \text{C01} / \text{C00}$$

$$\text{g/L H}_2\text{SO}_4 = (\text{C30} - \text{EP1}) * \text{C02} / \text{C00}$$

EP1 = titrant consumption in mL

C00 = 0.2 ... 0.5 (sample volume used in mL original sample)

C01 = 4.803

C02 = 4.904

C30 = 5.00 [added volume of $c(\text{BaCl}_2) = 0.05 \text{ mol/L}$ in mL]

Literature

- Metrohm Application Bulletin No. 130
Chloride titrations with potentiometric indication
Metrohm Ltd., Herisau.
- Metrohm Application Bulletin No. 140
Analytical determination of sulfate
Metrohm Ltd., Herisau.
- P. W. Wild
Moderne Analysen für die Galvanotechnik
Eugen Leuze Verlag, Saulgau, 1972.

Figures

```
'fr
785 DMP Titrimo      02287  785.0010
user                 th
date 2000-02-07     time 16:32      5
card label:         785
U(init)             322 mV MET U      AB89 a
smp1 size           2.00 ml
EP1                  4.203 ml          222 mV
EP2                  6.995 ml          19 mV
EP3                  8.684 ml         -218 mV
H2SO4                10.31 %
Al                   1.26 %
stop V reached
=====
```




Fig. 1: Result block and titration curve for the determination of sulfuric acid and aluminum in an anodizing bath.

```
'pa
785 DMP Titrino      02287  785.0010
user                th
date 2000-02-08     time 09:29      6
MET U              AB89 c
parameters
>titration parameters
V step              0.1 ml
dos.rate            max. ml/min
signal drift        40 mV/min
equilibr.time       28 s
start V:            OFF
pause               30 s
meas.input:         1
temperature         25.0 °C
>stop conditions
stop V:             abs.
stop V              3 ml
stop U              OFF mV
stop EP             9
filling rate        max. ml/min
>statistics
status:             OFF
>evaluation
EPC                 30 mV
EP recognition:     greatest
fix EP1 at U        OFF mV
pK/HNP:             OFF
>preselections
req.ident:          OFF
req.smpl size:      OFF
limit smpl size:    OFF
activate pulse:     OFF
=====
```

Fig. 2: Parameter settings for the determination of oxalic acid.

```
'fr
785 DMP Titrino      02287  785.0010
user                th
date 2000-02-08     time 09:29      6
card label:         785
U(init)             -134 mV MET U    AB89 c
smpl size           50.0 ml
EP1                  1.762 ml          -563 mV
Oxalic a             0.0159 %
stop V reached
=====
```




Fig. 3: Result block and titration curve for the determination of oxalic acid in an anodizing bath.

```
'pa
785 DMP Titrino      02287  785.0010
user                th
date 2000-02-08    time 11:46    13
MET U              AB89 A
parameters
>titration parameters
V step              0.1 ml
dos.rate            max. ml/min
signal drift        30 mV/min
equilibr.time       32 s
start V:            OFF
pause               0 s
meas.input:         1
temperature          25.0 °C
>stop conditions
stop V:             abs.
stop V              11 ml
stop U              OFF mV
stop EP             9
filling rate        max. ml/min
>statistics
status:             OFF
>evaluation
EPC                 30 mV
EP recognition:     greatest
fix EP1 at U        OFF mV
pK/HNP:             OFF
>preselections
req.ident:          OFF
req.smpl size:      OFF
limit smpl size:    OFF
activate pulse:     OFF
=====
```

Fig. 4: Parameter settings for the determination of chromic acid.

```
'fr
785 DMP Titrino      02287  785.0010
user                th
date 2000-02-08    time 11:46    13
card label:         785
U(init)             9 mV MET U    AB89 A
smpl size           0.20 ml
EP1                 8.883 ml      159 mV
CrO3                14.80 %
stop V reached
=====
```




Fig. 5: Result block and titration curve for the determination of chromic acid in an anodizing bath.


```
'pa
785 DMP Titrino      02287  785.0010
user                th
date 2000-02-09    time 14:29      7
MET U              AB89 C
parameters
>titration parameters
V step             0.05 ml
dos.rate           max. ml/min
signal drift       50 mV/min
equilibr.time      26 s
start V:           OFF
pause              30 s
meas.input:        1
temperature        25.0 °C
>stop conditions
stop V:            abs.
stop V             6 ml
stop U             OFF mV
stop EP           9
filling rate       max. ml/min
>statistics
status:            OFF
>evaluation
EPC                10 mV
EP recognition:    greatest
fix EP1 at U      OFF mV
pK/HNP:           OFF
>preselections
req.ident:         OFF
req.smpl size:     OFF
limit smpl size:  OFF
activate pulse:    OFF
=====
```

Fig. 6: Parameter settings for the determination of sulfate.

```
'fr
785 DMP Titrino      02287  785.0010
user                th
date 2000-02-09    time 14:29      7
card label:        785
U(init)            380 mV MET U      AB89 C
smpl size          0.2 ml
EP1                4.846 ml          397 mV
SO4 2-             9.99 g/l
H2SO4              10.20 g/l
stop V reached
=====
```




Fig. 7: Result block and titration curve for the determination of sulfate in an anodizing bath.