

# HPLC Analysis of Twenty One Preservative Compounds Found in Cosmetics Using a Thermo Scientific Hypersil GOLD Phenyl Column

Joanna Freeke, Thermo Fisher Scientific, Runcorn, Cheshire, UK

## Abstract

This application note demonstrates the use of a Thermo Scientific Hypersil GOLD Phenyl HPLC column for the analysis of twenty one preservative compounds found in cosmetics. All analytes were retained, with thirteen showing full baseline resolution and four pairs being partially resolved. The analysis can be performed on standard HPLC instrumentation with solution conditions suitable for MS detection.

## Introduction

The Hypersil GOLD® range of HPLC columns were developed to give reproducible and reliable chromatography analysis with excellent peak shape. The Hypersil GOLD Phenyl HPLC column offers excellent peak shape and alternative selectivity to the standard Hypersil GOLD HPLC column. In particular the presence of the phenyl group enhances the separation of aromatic compounds and moderately polar compounds by providing opportunity for a second mode of interaction through the phenyl ring.

Preservatives are added to cosmetic formulations to ensure that they do not carry pathogenic microorganisms and to eliminate the growth of microbes. In the USA the FDA requires that all cosmetics have to be adequately preserved and proof of this has to be prepared for every shipment. In Europe there is a list of allowed and restricted preservatives in Cosmetic Directive 76/768/EEC. It is therefore necessary for cosmetic product manufacturers to develop methods to screen cosmetics for these compounds. This analysis has become more important in recent years as there has been a rising public awareness of the additives in found in cosmetics and their health effects. In this application twenty one preservative compounds found in cosmetics are analyzed by HPLC. Some of these are aromatic and the Hypersil GOLD Phenyl HPLC column was selected to offer the best performance in separating these, in many cases, closely related compounds.



## Key Words

- Hypersil GOLD Phenyl
- Cosmetics
- Preservatives
- Parabens

## Experimental Details

Chemicals and Reagents	Part Number
Fisher Scientific LC-MS grade water	W/0112/17
Thermo Scientific Pierce LC-MS grade acetonitrile	TS-511001
Fisher Scientific ammonium acetate	A/3440/50
Fisher Scientific acetic acid	A/0400/PB08

## Sample Handling Equipment

Liquid handling hardware: FinnPipette Kit 1;	4700870
Includes 3 Thermo Scientific Finn pipette F2 Pipettors from 1 to 1000 µL	
Vials and closures: NSC Mass Spec Certified 2 mL clear vial with blue bonded PTFE silicone cap	MSCERT4000-34W

Separation Conditions	Part Number	
Instrumentation:	Thermo Scientific HPLC system	
Column(s):	Hypersil GOLD Phenyl 3 µm, 100 x 2.1 mm	25903-102130
Mobile phase:	(A) 5:95 acetonitrile:20 mM ammonium acetate, pH 4.45 (B) acetonitrile	
Gradient:	100 – 15% A over 20 minutes	
Flow rate:	0.5 mL/min	
Run time:	30 minutes (including equilibration time)	
Column temperature	30 °C	
Injection details:	2.0 µL	
UV detector wavelength:	254 nm and 214 nm	

## Solutions

Standard preparation: A 1 mg/mL solution of each of the preservative compounds in methanol was diluted with water to give a final concentration of 12 µg/mL.

## Data Processing

Software:	Thermo Scientific Chromatography Software
-----------	---

## Results

The analysis of twenty one preservative compounds commonly found in cosmetics was carried out on a Hypersil GOLD Phenyl HPLC column and the chromatograms are given in Figures 1 and 2 for detection at 254 and 214 nm respectively. Detection was carried out at these two wavelengths because some of the compounds do not absorb light at 254 nm. The chromatogram for detection at 214 nm (Figure 2) has had the baseline subtracted to correct for the large change in absorbance over the course of the gradient program. The list of compounds analyzed is given in Table 1 and the retention times and method precision for three selected compounds in six injections are given in Table 2.




Figure 1: Chromatogram for twenty one preservative compounds found in cosmetics separated on an Hypersil GOLD Phenyl 3 µm, 100 x 2.1 mm column at 254 nm (compounds 6, 8 and 13 do not absorb at 254 nm)

Compound	Compound
1 4-hydroxybenzoic acid	12 Ethylparaben
2 Salicylic acid	13 Dichlorobenzyl alcohol
3 Benzoic acid	14 Isopropylparaben
4 Benzyl alcohol	15 Propylparaben
5 Sorbic acid	16 Chlorhexidin
6 Phenoxyethanol	17 Isobutylparaben
7 p-Anisic acid	18 Butylparaben
8 Dehydroacetic acid	19 o-Phenylphenol
9 Methylparaben	20 Triclocarban
10 Chlorphenisic acid	21 Triclosan
11 Hexamidine 2-hydroxyethansulfonate	

Table 1: List of preservative compounds analyzed

[www.thermoscientific.com/chromatography](http://www.thermoscientific.com/chromatography)

© 2011 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

## Conclusions

- Twenty one preservative compounds commonly found in cosmetics are identified and separated by a Hypersil GOLD Phenyl HPLC column.
- For thirteen of the compounds baseline resolution was achieved and four pairs of compounds were partially resolved (6 and 7, 8 and 9, 11 and 12, 13 and 14).
- The analysis is reproducible with low % RSD and the mobile phase components are suitable for both MS and UV detection.




Figure 2: A chromatogram for twenty one preservative compounds found in cosmetics separated on a Hypersil GOLD Phenyl 3 µm, 100 x 2.1 mm column at 214 nm (all the compounds can be detected at this wavelength). A blank chromatogram is subtracted from the data and the raised baseline is due to absorption from mobile phase components

Peak	t <sub>R</sub> /min	% RSD
5	3.73	0.00
19	11.03	0.13
20	14.51	0.08

Table 2: Average retention time and method precision (% RSD) for six replicate injections for three of the twenty one compounds

In addition to these offices, Thermo Fisher Scientific maintains a network of representative organizations throughout the world.

**North America  
USA and Canada**  
+1 800 332 3331

**Europe  
France**  
+33 (0)1 60 92 48 34

**Germany**  
+49 (0) 2423 9431 -20  
or -21

**United Kingdom**  
+44 1928 534110

**Asia  
Japan**  
+81 3 5826 1615

**China**  
+86-21-68654588  
or +86-10-84193588  
800-810-5118

**India**  
+91-22-6742 9494

**Thermo Fisher  
Scientific Australia  
Pty Ltd**  
1300 735 292 (free call  
domestic)

**Thermo Fisher  
Scientific New  
Zealand Ltd**  
0800 933 966 (free call  
domestic)

**All Other Enquiries**  
+44 (0) 1928 534 050

## Technical Support

**North America**  
800 332 3331

**Outside North  
America**  
+44 (0) 1928 534 440

ANCCSHYPGOLDCOS 1211