

SOLUTIONS BY



# Automated Determination of Steroids in Urine of Food-Producing Animals

Dr. Susan Ehlers, Timo Horny  
Chemisches und Veterinäruntersuchungsamt  
Münsterlande-Emscher-Lippe (MEL)

## Contents

1	Introduction .....	3
1.1	Chemicals.....	3
2	Method Development.....	4
2.1	Material .....	4
2.2	Sample Material.....	4
2.3	Sample Treatment - Extraction .....	4
2.4	Solid Phase Extraction (manual).....	4
2.5	Solid Phase Extraction (automated by LCTech).....	5
2.6	Measurement using HPLC-MS/MS.....	7
2.7	Configuration FREESTYLE.....	9
3	Results .....	10
3.1	Recovery .....	10
3.2	Reproducibility.....	10
4	Conclusions .....	11
5	Literature .....	12

## 1 Introduction

In the European Union the use of steroids as growth-promotor for food-producing animals is strictly forbidden. In National Residue Control Plans (Directive 96/23/EC) the analysis of steroids in urine of food-producing animals is prescribed to uncover illegal treatment. Some official laboratories in the European Union have to analyse 1000 samples and more per year. The matrix urine leads to some interference in final measurement using LC-MS/MS or HPLC-HRMS, respectively. Therefore an efficient cleaning of the sample solution using two different SPE columns is needed. Due to the number of samples to be monitored time saving is crucial, and thus an automated sample clean-up is recommendable.

### 1.1 Chemicals

The chemicals of interest are listed in Table 1. These substances are often named synthetic steroids, but in the literature it is described, that substances like boldenone and nandrolone can be generated endogenously in boars for example [1].

Chemicals	Internal Standard
Alpha-Boldenone	Beta-Boldenone-D3
Beta-Boldenone	Beta-Boldenone-D3
Alpha-Nortestosterone	Beta-Nortestosterone-D3
Beta-Nortestosterone	Beta-Nortestosterone-D3
Alpha-Trenbolone	Beta-Trenbolone-D3
Beta-Trenbolone	Beta-Trenbolone-D3
Methyl-Testosterone	Methyltestosterone-D3
Stanozolol	Stanozolol-D3
Hydroxystanozolol	Hydroxystanozolol-D3
Ethinylestradiol	Ethinylestradiol-D4

Table 1: Steroids of interest

## 2 Method Development

### 2.1 Material

- Formic acid (0.1 %) in acetonitrile (ULC/MS), 0.1 % formic acid in water (ULC/MS) both from Biosolve (the Netherlands)
- Sodium acetate, anhydrous, per analysis from Merck (Germany)
- Ethyl acetate (picograde), methanol (picograde) and acetonitrile (picograde) from Promochem (Germany)
- Glucuronidase/arylsulfatase solution from Sigma (USA)

### 2.2 Sample Material

The urine is centrifuged to separate any particulate matter.

### 2.3 Sample Treatment - Extraction

4 mL of centrifuged urine is diluted with 4 mL of sodium acetate buffer (50 mmol/l, pH 5.2). 20 µL of a glucuronidase/arylsulfatase solution and internal standard solution is added to the sample. The sample is incubated at 37 °C over night. After a centrifugation step, the sample solution is cleaned-up using two different SPE columns.

### 2.4 Solid Phase Extraction (manual)

Column 1	CHROMABOND® HR-X, 6 mL, 200 mg
Column conditioning	4 mL methanol 4 mL double distilled water
Loading step	The sample is completely applied onto the column.
Washing step	4 mL methanol/double distilled water (v/v, 5/95) after a drying step (5 min.) 1 mL methanol
Elution Step	3 mL methanol

Table 2: Conditions of the SPE method for the first column (manual)

Column 2	CHROMABOND® NH <sub>2</sub> , 3 mL, 500 mg
Column conditioning	5 mL methanol 5 mL methanol/ethyl acetate (v/v, 40/60)
Loading step	The eluate from column 1 is completely applied onto the column.
Elution Step	The solution running through the column is collected. After waiting 10 min, vacuum was put to the column for 20 sec.

Table 3: Conditions of the SPE method for the second column (manual)

## 2.5 Solid Phase Extraction (automated by LCTech)

The manual solid phase extraction can be easily transferred to the robotic system FREESTYLE SPE for automated sample preparation. The software of the system contains pre-defined wizards for every step making the method development fast and easy.

 <small>AUTOMATISIERTE PROBEVORBEREITUNG</small>				
LCTech FreeStyle - Report on Methods: DSM      Date: 22.03.2016    Time: 16:49:22				
Name: Steroids.dsm				
Column:	LCTech_6ml.col	Extension cannula:	no	
Conditioning 1:	ON			
Volume:	4 ml	Dispensing Speed:	6 ml / min	
Suction Speed:	40 ml / min	Waiting time:	0.1 min      Port : 1 MeOH	
Conditioning 2:	ON			
Volume:	4 ml	Dispensing Speed:	6 ml / min	
Suction Speed:	50 ml / min	Waiting time:	0.1 min      Port : 7 H2O	
Conditioning 3:	OFF			
Load 1:	ON			
Volume:	8 ml	Dispensing Speed:	2 ml / min	
Suction Speed:	10 ml / min	Input Vial Type:	Type1@15	
Load 1 into Waste				
Quantitativ Transfer	2 x 2ml			Port : 7 H2O
Drying 0:	OFF			
Washing 1:	ON	stay on actual position		
Volume:	1 ml	Dispensing Speed:	6 ml / min	
Suction Speed:	40 ml / min	Waiting time:	0.1 min      Port : 8 MeOH/H2O 5/95	
Washing 2:	OFF			
Drying 1:	ON	stay on actual position		
Time:	1 min			
Washing 3:	OFF			
Elution 1:	OFF			
Drying 2:	OFF			
Washing 4:	ON	into Waste		
Volume:	1.4 ml	Dispensing Speed:	1.5 ml / min	
Suction Speed:	50 ml / min	Waiting time:	00.1 min      Port : 1 MeOH	
Drying 3:	OFF			
Elution 2:	ON			
Volume:	2.5 ml	Dispensing Speed:	2 ml / min	
Suction Speed:	50 ml / min	Waiting time:	0 min	
		Vial Type:	Type3@14	
Drying 4:	ON	stay on actual position		
Drying by defined air volume		Volume: 10 ml	Dispensing Speed: 6 ml / min	
Cleaning cycle 1 :	ON			
Volume:	1 ml	Dispensing Speed:	1 ml / min	
Suction Speed:	1 ml / min			Port : 1 MeOH

Figure 1: Method report (part 1) for the dual SPE method (automated)

# APPLICATION NOTE | AN0013

Column SPE_2: Extension cannula:	LCTech_3ml.col no	Same solvent as SPE_1		
Conditioning 2.1: Volume: Suction Speed:	ON 5 ml 6 ml / min	Dispensing Speed: Waiting time:	40 ml / min 0.1 min	Port : 1 MeOH
Conditioning 2.2: Volume: Suction Speed:	ON 5 ml 50 ml / min	Dispensing Speed: Waiting time:	6 ml / min 0.1 min	Port : 9 40/60 MeOH/Ethylacetate
Conditioning 2.3:	OFF			
Load_2: Volume: Suction Speed: Quantitativ Transfer	ON 3.2 ml 10 ml / min 1 x 1 ml	From vial with Eluate 2 Dispensing Speed: Input Vial Type:	2 ml / min Type3@14	Port : 1 MeOH
Drying 2.0: Time:	ON 0.5 min	into Waste		
Washing 2.1:	OFF			
Washing 2.2:	OFF			
Drying 2.1:	OFF			
Washing 2.3:	OFF			
Elution 2.1:	OFF			
Drying 2.2:	OFF			
Washing 2.4:	OFF			
Drying 2.3:	OFF			
Elution 2.2:	OFF			
Drying 2.4:	OFF			
Washing 4:	OFF			
SETUP :				
Use pressure limitation function during loading and washing	ON			
Pressure limit for syringe pump	160 digits			
Maximum count of triggered samples in series	3 Sample /s			
System - rinsing and conditioning with solvent from port:	Port :	1 MeOH		

Figure 2: Software Wizard (part 2) for the dual SPE method (automated)



Figure 3: Reusable plungers from LCTech GmbH allow higher operating pressure, higher flow-rates and reduced dead volume.




Figure 4: Rack arrangement for the dual SPE method (automated)

The eluate (from manual and automated clean-up) is dried in a gentle stream of nitrogen until dryness using a TurboVap from Biotage. After this step the sample is resolved in 500  $\mu$ L of 0.1% formic acid in acetonitrile / 0.1 % formic acid in water (v/v, 10/90). The solution is measured using LC-MS/MS.

## 2.6 Measurement using HPLC-MS/MS

The HPLC-MS measurement can be performed on an Agilent 6495 combined with a 1290 Infinity Quat. Pump. The MRM-Transitions are shown in Table 4 and possible HPLC-MS/MS conditions are shown in Table 5. The sample extracts can also be measured using HPLC-HRMS (conditions not shown).

# APPLICATION NOTE | AN0013

Chemicals	[M+H] <sup>+</sup>	Q1 (Quantifier)	Q2 (Quantifier)
Alpha-Boldenone	287	134.9	120.8
Beta-Boldenone	287	134.9	120.8
Alpha-Nortestosterone	275	109	90.8
Beta-Nortestosterone	275	109	90.8
Alpha-Trenbolone	271	199	164.9
Beta-Trenbolone	271	199	164.9
Methyl-Testosterone	303	108.9	96.9
Stanozolol	329	95	81
Hydroxystanozolol	345.1	96.8	66.8
Ethinylestradiol	279	158.9	132.8

Table 4: MRM-Transitions

Column	Phenomenex Luna 3u C18(2) 100A, 150 x 2.00 mm
Eluent A	0.1 % formic acid in acetonitrile
Eluent B	0.1 % formic acid in water
Gradient	10 % B to 93 % B (13 min.); 97 % B (hold 3 min.), 97 % B to 10 % B (0.5 min.), equilibration 3.5 min Total run time 20 min.
Flow rate	0.2 mL/min
Injection Volume	10 µL
Column Temperature	30 °C
Detection	Agilent 6495 with 1290 Infinity Quat. Pump Ion source: ESI Ionization mode: positive Gas temperature: 200 °C Gas flow: 11 L/min. Nebulizer: 20 psi Sheath gas heater: 380 °C Sheath gas flow: 11 Capillary: 3000 V VCharging: 1500 Ion funnel pos high pressure RF 150 Ion funnel pos low pressure RF 60

Table 5: Possible HPLC-MS/MS conditions



## 2.7 Configuration FREESTYLE

1.	FREESTYLE BASIC	P/N 12663-12
2.	Liquid level switch for organic solvents	P/N 12709
3.	FREESTYLE SPE module	P/N 12668
4.	Rack for solvent delivery	P/N 13156
5.	Rack for SPE-columns (2 needed) for up to 18 columns	P/N 13946
6.	Tray (Elution for SPE) for 54 pcs / test tubes 100x16 mm	P/N 13948
7.	Sample rack for 52x15 mL tubes (Falcon/Saarstedt)	P/N 14445
8.	Plunger for 6 mL-cartridges (4 needed) (M&N, 10 pcs. / pck)	P/N 14596
9.	Plunger for 3 mL-cartridges (4 needed) (M&N, 10 pcs. / pck)	P/N 14597
10.	Column adapter for 6 mL SPE columns (4 needed)	P/N 14613
11.	Column adapter for 3 mL SPE columns (4 needed)	P/N 14612
12.	Plunger removal tool, 6 mL	P/N 14643
13.	Software upgrade "Extended functions"	P/N 14044
14.	FREESTYLE Software-Dongle	P/N 14111

## 3 Results

### 3.1 Recovery

Official laboratories in the European Union are obliged to detect 1 µg/kg of synthetic steroids in urine samples. This is easily achieved by using the manual as well as the automated application.

The recoveries achieved were better than the requirements (-50 % to + 20 % at concentrations < 1 µg/kg) set by the European law (decision 2002/657/EG) that means.

### 3.2 Reproducibility

Internal standards were added to ten different urine samples twice. The samples were purified manually using two SPE columns on the one hand. In the second approach the same samples were purified automatically using the robotic system FREESTYLE.

All 20 samples were measured using HPLC-HRMS (Thermo, Q-Exactive)

The absolute peak areas of the deuterated internal standards (concentration 3 µg/L urine) in the ten samples were compared.

Sample Results					
Sample Name	Excluded	Calculated Amt	Area	Expected RT	Actual RT
Δa		Δa	Δa	Δa	Δa
KAL 0 (0 ug/kg)	<input type="checkbox"/>	3.000	161260700	15.44	15.39
KAL 1 (0.05 ug/kg)	<input checked="" type="checkbox"/>	3.000	158490589	15.44	15.39
KAL 2 (0.1 ug/kg)	<input checked="" type="checkbox"/>	3.000	160560920	15.44	15.37
KAL 3 (0.5 ug/kg)	<input checked="" type="checkbox"/>	3.000	159467625	15.44	15.37
KAL 4 (1 ug/kg)	<input checked="" type="checkbox"/>	3.000	158529360	15.44	15.37
KAL 5 (2 ug/kg)	<input checked="" type="checkbox"/>	3.000	155093991	15.44	15.36
KAL 6 (4 ug/kg)	<input checked="" type="checkbox"/>	3.000	150667545	15.44	15.36
KAL 7 (10 ug/kg)	<input checked="" type="checkbox"/>	3.000	138354112	15.44	15.34
Blank	<input type="checkbox"/>	3.000	117887	15.44	15.34
2015MEL000703 b	<input type="checkbox"/>	3.000	42653979	15.44	15.36
2015MEL000704 b	<input type="checkbox"/>	3.000	43599181	15.44	15.34
2015MEL000705 b	<input type="checkbox"/>	3.000	53412056	15.44	15.34
2015MEL000706 b	<input checked="" type="checkbox"/>	3.000	36386635	15.44	15.33
2015MEL000710 b	<input type="checkbox"/>	3.000	39525642	15.44	15.32
2015MEL000877 b	<input type="checkbox"/>	3.000	61450110	15.44	15.31
2015MEL000878 b	<input type="checkbox"/>	3.000	58637548	15.44	15.60
2015MEL000879 b	<input type="checkbox"/>	3.000	54795571	15.44	15.31
2015MEL001057 b	<input type="checkbox"/>	3.000	29468347	15.44	15.30

Peak area mean: 46 658 785  
 Relative standard deviation: 23 %  
 Relative standard deviation based on the measurement system\*: 5 %

Figure 5: Table with trenbolone reproducibility measurement results (automated via FREESTYLE)

# APPLICATION NOTE | AN0013

Sample Results					
Sample Name	Excluded	Calculated Amt	Area	Expected RT	Actual RT
Δa		Δa	Δa	Δa	Δa
KAL 0 (0 ug/kg)	<input type="checkbox"/>	3.000	134932998	15.44	15.45
KAL 1 (0.05 ug/kg)	<input type="checkbox"/>	3.000	132422820	15.44	15.44
KAL 2 (0.1 ug/kg)	<input type="checkbox"/>	3.000	134256586	15.44	15.44
KAL 3 (0.5 ug/kg)	<input type="checkbox"/>	3.000	134872252	15.44	15.43
KAL 4 (1 ug/kg)	<input type="checkbox"/>	3.000	135563552	15.44	15.44
KAL 5 (2 ug/kg)	<input type="checkbox"/>	3.000	134129367	15.44	15.45
KAL 6 (4 ug/kg)	<input type="checkbox"/>	3.000	132651993	15.44	15.43
KAL 7 (10 ug/kg)	<input type="checkbox"/>	3.000	122441471	15.44	15.44
Blank	<input type="checkbox"/>	3.000	108184	15.44	15.45
Blank	<input type="checkbox"/>	N/F	N/F	15.44	N/F
2015MEL000703 h	<input type="checkbox"/>	3.000	31680924	15.44	15.45
2015MEL000704 h	<input type="checkbox"/>	3.000	38051752	15.44	15.44
2015MEL000705 h	<input type="checkbox"/>	3.000	59861816	15.44	15.45
2015MEL000706 h	<input type="checkbox"/>	3.000	43476228	15.44	15.45
2015MEL000710 h	<input type="checkbox"/>	3.000	40881054	15.44	15.42
2015MEL000877 h	<input type="checkbox"/>	3.000	58796317	15.44	15.45
2015MEL000878 h	<input type="checkbox"/>	3.000	98294831	15.44	15.32
2015MEL000879 h	<input type="checkbox"/>	3.000	84798351	15.44	15.43
2015MEL001057 h	<input type="checkbox"/>	3.000	28588292	15.44	15.46

Peak area mean: 53 825 507  
 Relative standard deviation: 45 %  
 Relative standard deviation based on the measurement system\*: 5 %

Figure 6: Table with trenbolone reproducibility measurement manually

\*Relative standard deviation based on the measurement system was tested using the measurement of 8 standard solutions

	FREESTYLE		Manually	
	peak area mean	Relative standard deviation	peak area mean	Relative standard deviation
Boldenone	17.718.489	21	18.656.051	32
Nortestosterone	87.055.332	18	95.172.536	29
Stanozolol	323.036.143	16	339.948.453	22
Methyl testosterone	161.171.697	19	184.326.218	21
Ethinylestradiol	488.678	18	647.134	39

Table 6: Reproducibility measurement results

## 4 Conclusions

A comparison of the results achieved manually and on the automated system, respectively, clearly shows that the sample clean-up for the analysis of synthetic steroids employing two SPE columns can easily be processed on the FREESTYLE SPE system.

The measured peak area mean values (Table 6) are very similar, and the relative standard deviations indicate a slightly better reproducibility for the automated processing.

## 5 Literature

Confirmatory Method for the Determination of Steroids in Bovine Urine by LC-MS/MS;  
Bundesamt für Verbraucherschutz und Lebensmittelsicherheit 2011

[1] S. Polemans, K. De Wasch, H. Noppe, N. Van Hoof, S. Van Cruchten, B. Le Bizec, Y. Deceunick, S. Sterk, H.J. Van Rossum, M.K. Hoffmann, H.F. Brabander; Endogenous occurrence of some anabolic steroids in swine matrices; Food Additives and Contaminants, September 2005; 22(9):808-815



## Contact

LCTech GmbH  
Daimlerstraße 4  
84419 Obertaufkirchen  
Germany

Tel.: +49 8082 2717-0  
Fax: +49 8082 2717-100  
E-Mail: [info@LCTech.de](mailto:info@LCTech.de)

[www.LCTech.de](http://www.LCTech.de)  
[www.LCTech-online.com](http://www.LCTech-online.com)

SOLUTIONS BY

