

Extraction and Clean-up of Aminoglycoside Antibiotics from Chicken Meat Using EVOLUTE® EXPRESS WCX Solid Phase Extraction Columns with Analysis by LC-MS/MS

Jessica Lance¹, Gloria Gaucin¹, Carrie Snyder¹, B. J. Bench¹, Frank Kero², Victor Vandell², Martin Cherrier², Elena Gairloch²

1. Tyson Foods, Inc. Springdale Corporate Laboratory, Springdale, AR 72762, USA
2. Biotage, 10430 Harris Oaks Blvd., Suite C, Charlotte, NC USA 28269, USA




Figure 1. Analyte structures.

Introduction

This application note describes a polymer based-based weak cation exchange mixed-mode SPE protocol for the extraction of the aminoglycoside antibiotics gentamicin, streptomycin, spectinomycin and neomycin from chicken meat prior to LC-MS/MS analysis.

The utility of aminoglycoside antimicrobial additives and foods has recently been challenged due to the effects of bioaccumulation in an animal host. The persistence of these selected residues may result in population resistance to antibiotic treatment and as a consequence, has been considered an issue of public health.

This application note describes an approach to extraction of these analytes from chicken meat, with optimized analyte recovery, minimal ion suppression, and acceptable method precision.

Analytes

Gentamicin, Streptomycin, Spectinomycin, Neomycin

Sample Preparation Procedure

Format

EVOLUTE® EXPRESS WCX 500 mg/6 mL (Tablets) SPE columns, part number 612-0050-CXG.

[Note: This application was originally developed using EVOLUTE® WCX 500 mg/6 mL SPE columns, part number 612-0050-C. This format has now been replaced by the EXPRESS format, which provides improved flow characteristics].

Sample Pre-treatment

Weigh 5 g of chicken meat into a 50 mL conical tube.

Spike with internal standard as required.

Add 20 mL of extraction solution (10 mM NH₄OAc, 0.4 mM EDTA, 0.5% NaCl and 2% trichloroacetic acid in water), vortex for 1 min and shake for 10 mins. After shaking, centrifuge for 10 mins at 4000rpm, and decant the supernatant into a clean tube taking care not to transfer any of the tissue.

Repeat the extraction procedure and combine the two extracts from each sample. Filter using a glass fibre filter.

Adjust the sample pH to 6.5+/-0.25 using 30% NaOH, 1 N NaOH and 1 N HCl.

Solid Phase Extraction

Conditioning

Condition each column with methanol (10 mL)

Equilibration

Equilibrate each column with water (10 mL)

Sample Loading

Load pre-treated sample (40 mL) using vacuum or positive pressure

Interference Wash

Elute interferences with water adjusted to pH 6.5 (10 mL)

Elution

Elute analytes with 25% formic acid in water (6 mL)

Post Extraction

Add HPLC grade water to the extract to give a total volume of 10 mL, and filter through a 0.2 µm disc.

UPLC Conditions

Instrument

Waters Acquity UPLC

Column

Waters Acquity UPLC BEH Amide 1.7 μm 2.1 x 50mm

Mobile Phase

A: 20 mM heptafluorobutyric acid (HFBA) in H_2O /acetonitrile (95/5, v/v)

B: 20 mM HFBA in acetonitrile

Flow Rate

0.2 mL/min

Table 1. Gradient Conditions.

Time (min)	Flow rate (mL/min)	% A	% B
Initial	0.2	100	0
0.5	0.2	80	20
1	0.2	80	20
2	0.2	60	40
2.05	0.2	10	90
2.5	0.2	10	90
2.55	0.2	100	0
3	0.2	100	0
3.05	0	100	0

Injection Volume

7.5 μL

Sample Temperature

10 $^{\circ}\text{C}$

Column Temperature

40 $^{\circ}\text{C}$

Mass Spectrometry Conditions

Instrument

Xevo TQD Triple Quad Mass Spec equipped with an electrospray ionization source operated in positive ion mode. The compound selective MRM transitions are detailed in **Table 3**.

Desolvation Temperature

350 $^{\circ}\text{C}$

Ion Source Temperature

150 $^{\circ}\text{C}$

Collision Cell Pressure

3.4×10^{-3} mbar

Positive ion acquired in the multiple reaction monitoring (MRM) mode.

Table 2. MRM transitions for selected analytes.

Analyte	MRM Transition	Cone Voltage (V)	Collision Energy (eV)
C1 Gentamicin	479>157,160,322	40	25
C1a Gentamicin	450>112,160,322	35	25
C2+C2a Gentamicin	464>160,163,322	35	20
Streptomycin	582>176,246,263	70	32
Spectinomycin	351>98,140,333	40	20
Neomycin	615>160,163,293	52	30

Results

A summary of the performance for this method is given in Table 3. The observed linearity for each analyte over the concentration range of interest was $r^2 > 0.990$. The reference range defined by USDA-FSIS guidelines was 100–500 ppb.

Table 3. Method performance criteria.

Analyte	Linearity	LOQ (ppb)	LOQ spec
Gentamicin	0.9915	33	100
Streptomycin	0.9986	382	500
Spectinomycin	0.9959	75	100
Neomycin	0.9993	361	500




Figure 2. Representative chromatogram of a fortified meat specimen processed using EVOLUTE® WCX.

Recovery

Relative recoveries of the selected analyte from fortified specimens were determined at 3 concentration levels.

The results are given in **Figure 3**.

Streptomycin and neomycin:

Level 1 = 400 ppb; Level 2 = 1000 ppb; Level 3 = 10000 ppb

Gentamicin and spectinomycin:

Level 1 = 50 ppb; Level 2 = 500 ppb; Level 3 = 1000 ppb




Figure 3. Relative recovery from fortified samples.

Conclusions

This method was demonstrated as a viable option for residue measurements over a relevant concentration range in food safety laboratory applications.

Additional Notes

- Buffer Preparation.** Aminoglycoside extraction solvent mixture (10mM NH_4OAc , 0.4 mM EDTA, 0.5% NaCl and 2% TCA in water): Add 1.54 g of NH_4OAc to 1.95 L of water. Adjust the pH of the solution to 4.0 with 1N HCl and/or 1 N NaOH. Add 0.3 g $\text{Na}_2\text{EDTA}\cdot 2\text{H}_2\text{O}$, 10 g of NaCl and 40 g TCA. Mix to ensure the salts dissolve and adjust final volume to 2 L with pure water.
- Processing Conditions.** Samples were loaded on the SPE cartridge at a flow rate of 2 mL/min.

Ordering Information

Part Number	Description	Quantity
612-0050-CXG	EVOLUTE® EXPRESS CX 500 mg/6 mL (Tablets) SPE Columns	30
PPM-48	Biotage® PRESSURE+ 48 Positive Pressure Manifold 48 Position	1

EUROPE

Main Office: +46 18 565900
Toll Free: +800 18 565710
Fax: +46 18 591922
Order Tel: +46 18 565710
Order Fax: +46 18 565705
order@biotage.com
Support Tel: +46 18 56 59 11
Support Fax: +46 18 56 57 11
eu-1-pointsupport@biotage.com

NORTH & LATIN AMERICA

Main Office: +1 704 654 4900
Toll Free: +1 800 446 4752
Fax: +1 704 654 4917
Order Tel: +1 704 654 4900
Order Fax: +1 434 296 8217
ordermailbox@biotage.com
Support Tel: +1 800 446 4752
Outside US: +1 704 654 4900
us-1-pointsupport@biotage.com

JAPAN

Tel: +81 3 5627 3123
Fax: +81 3 5627 3121
jp_order@biotage.com
jp-1-pointsupport@biotage.com

CHINA

Tel: +86 21 68162810
Fax: +86 21 68162829
cn_order@biotage.com
cn-1-pointsupport@biotage.com

KOREA

Tel: +82 31 706 8500
Fax: +82 31 706 8510
korea_info@biotage.com
kr-1-pointsupport@biotage.com

INDIA

Tel: +91 11 45653772
india@biotage.com

Distributors in other regions
are listed on www.biotage.com

Literature Number: AN833.V.1

© 2020 Biotage. All rights reserved. No material may be reproduced or published without the written permission of Biotage. Information in this document is subject to change without notice and does not represent any commitment from Biotage. E&OE. A list of all trademarks owned by Biotage AB is available at www.biotage.com/legal. Other product and company names mentioned herein may be trademarks or registered trademarks and/or service marks of their respective owners, and are used only for explanation and to the owners' benefit, without intent to infringe.