


Halogens and sulfur in residual solvent applying Combustion IC



Determination of fluoride, chloride, bromide, and sulfate in a residual solvent applying combustion as sample preparation with subsequent anion chromatography with conductivity detection after sequential suppression. Analysis is crucial for waste classification in halogen-free and halogen-containing solvents.

Results

	Mean [mg/kg]
Fluoride	2.5
Chloride	7.1
Nitrite	n.q.
Bromide	147.0
Nitrate	n.q.
Phosphate (internal standard, [mg/L])	1.0
Sulfate	631.2

Sample

Used methanol

Sample preparation

Combustion followed by Inline Matrix Elimination

Columns

Metrosep A Supp 5 - 150/4.0	6.1006.520
Metrosep A Supp 4/5 Guard/4.0	6.1006.500
Metrosep A PCC 1 HC/4.0	6.1006.310

Solutions

Eluent	3.2 mmol/L sodium carbonate 1.0 mmol/L sodium hydrogen carbonate
Suppressor regenerant	100 mmol/L sulfuric acid
Rinsing solution	Ultrapure water
Absorption solution	300 mg/L hydrogen peroxide 1.0 mg/L phosphate

Analysis

Conductivity after sequential suppression

Parameters

Flow rate	0.7 mL/min
Injection volume	100 µL
P _{max}	15 MPa
Recording time	17 min
Column temperature	30 °C

Combustion parameters

Oven temperature	
Inlet/outlet	900/1000 °C
Absorption solution	5.0 mL

Instrumentation

881 Compact IC pro – Anion – MCS	2.881.0030
IC Conductivity Detector	2.850.9010
800 Dosino	2.800.0010
Remote box	6.2148.010
Mitsubishi AQF-100	*
Mitsubishi ABC-100	*
Mitsubishi WS-100	*
Mitsubishi GA-100	*

* from local Mitsubishi distributor / not shown in system graphic below

