


# Fluorine in coal sample applying Metrohm Combustion IC according to ASTM D8247



Coal contains a certain amount of fluorine, chlorine, and sulfur compounds. During combustion of the coal, these components release corrosive acids (e.g., fluorine compounds form hydrofluoric acid). Thermal power plants therefore request low-fluorine coal to avoid massive hydrofluoric acid production. In this application note, fluorine content in coal is determined by ion chromatography after pyrohydrolysis.

## Results

|      | 1 Fluorine<br>[mg/kg] | 2 Chlorine<br>[mg/kg] | 3 Unknown | 4 Sulfur |
|------|-----------------------|-----------------------|-----------|----------|
| Coal | 87.9                  | 55.8                  | -         | n.q.     |

## Sample

Coal

## Sample preparation

The sample is analyzed by Metrohm Combustion IC with flame sensor technology and intelligent Partial Loop Injection Technique with Inline Matrix Elimination.

## Columns

|                               |            |
|-------------------------------|------------|
| Metrosep A Supp 5 - 150/4.0   | 6.1006.520 |
| Metrosep A Supp 4/5 Guard/4.0 | 6.1006.500 |
| Metrosep A PCC 1 HC/4.0       | 6.1006.310 |
| Metrosep A Trap 1 - 100/4.0   | 6.1014.000 |
| Metrosep I Trap 1 - 100/4.0   | 6.1014.200 |

## Solutions

|                       |   |
|-----------------------|---|
| Eluent                | 3.2 mmol/L sodium carbonate<br>1.0 mmol/L sodium hydrogen carbonate |
| Suppressor regenerant | 100 mmol/L sulfuric acid  |
| Rinsing solution      | STREAM  |
| Absorber solution     | 200 mg/L hydrogen peroxide  |

## Parameters

|                       |              |
|-----------------------|--------------|
| Flow rate             | 0.7 mL/min   |
| Injection volume (IC) | 50 µL (MiPT) |
| P <sub>max</sub>      | 15 MPa       |
| Recording time        | 12 min       |
| Column temperature    | 30 °C        |

## Combustion parameters

|                                       |            |
|---------------------------------------|------------|
| Argon                                 | 100 mL/min |
| Oxygen                                | 300 mL/min |
| Oven temperature                      | 1100 °C    |
| Post-combustion time                  | 400 s      |
| Initial volume of absorption solution | 2.0 mL     |
| Water inlet                           | 0.2 mL/min |



## Analysis

Conductivity after sequential suppression

## Instrumentation

|                                     |             |
|-------------------------------------|-------------|
| 930 Compact IC Flex Oven/SeS/PP/Deg | 2.930.2560* |
| IC Conductivity Detector            | 2.850.9010* |
| MSM Rotor A                         | 6.2832.000* |
| Adapter sleeve for Suppressor Vario | 6.2842.020* |
| 920 Absorber Module                 | 2.920.0010* |
| Combustion Module (oven and ABD)    | 2.136.0700* |
| Autosampler MMS 5000                | 2.136.0800  |
| Kit for solid sampling              | 6.7302.000  |

\* available as 930 Metrohm Combustion IC (2.930.9010)



[www.metrohm.com](http://www.metrohm.com)

 **Metrohm**