

Thermo. Titr. Application Note No. H-008

Title: Determination of phosphate by magnesium titration

Scope: Determination of soluble orthophosphate ion, for example soluble phosphate in fertilizers such as DAP.

Principle: Orthophosphate will react exothermically with magnesium ions in the presence of ammonia to form an insoluble magnesium ammonium phosphate, $\text{MgNH}_4\text{PO}_4 \cdot 6\text{H}_2\text{O}$

Reagents: 1 mol/L $\text{Mg}(\text{NO}_3)_2$ – standardized
 $\text{NH}_3/\text{NH}_4\text{Cl}$ buffer solution. Add 142mL concentrated NH_3 solution to 17.5g NH_4Cl and dilute to 250mL with D.I. water

Method:

Basic Experimental Parameters:

Data rate (per second)	10
Titrant delivery rate (mL/min.)	1
No. of endothermic endpoints	1
Data smoothing factor	35

Procedure:

Weigh accurately approximately 10g of finely ground DAP fertilizer, and dissolve with deionized water in a 200mL volumetric flask. Make to volume, mix well, and filter through a dry Whatman #1 filter. Pipette a 20mL aliquot into a titration vessel and add 5mL $\text{NH}_3/\text{NH}_4\text{Cl}$ buffer solution. Titrate to an exothermic endpoint.

Determine the method blank by titrating aliquots of 20, 15, 10 and 5mL of a similar sample solution. Perform a linear regression on the titration data (Y-axis = titres (mL), X-axis = aliquot volumes (mL)). The Y-axis intercept is the method blank in mL.

The titrant is standardized by thermometric titration against standard 0.2M EDTA solution.

Results: Analysis of di-ammonium phosphate fertilizer:
 Mean (n=5) = 17.96 ± 0.007% w/w di-ammonium phosphate as P

Calculation:

$$\%P = \frac{((\text{Titre, mL} - \text{blank, mL}) \times M \text{ Mg}^{2+} \times \text{FW P} \times 100)}{(\text{Sample mass, g} \times 1000)}$$

$$\%P = \frac{((5.545 - 0.0555) \times 1.05933 \times 30.9738 \times 100)}{(1.00340 \times 1000)}$$

$$= 17.95\%$$
