

# Thermo. Titr. Application Note No. H-021

|               |  |
|---------------|--|
| <b>Title:</b> | <b>Determination of Free Acid in Copper Refining Solutions</b> |
|---------------|--|

|               |   |
|---------------|---|
| <b>Scope:</b> | Determination of free acid in copper refining solutions |
|---------------|---|

|                   |  |
|-------------------|--|
| <b>Principle:</b> | Thermometric titration of “free acid” values in copper refining solutions with standard NaOH after complexation of Fe(III) with potassium oxalate solution. “Free acid” is expressed as g/L sulfuric acid. |
|-------------------|--|

|                  |   |
|------------------|---|
| <b>Reagents:</b> | Standard 1 mol/L NaOH solution<br>Saturated K <sub>2</sub> C <sub>2</sub> O <sub>4</sub> solution, approximately 300g/L |
|------------------|---|

|                                 |   |                        |    |                                 |   |                             |   |                       |    |
|---------------------------------|---|------------------------|----|---------------------------------|---|-----------------------------|---|-----------------------|----|
| <b>Method:</b>                  | <p>Basic Experimental Parameters:</p> <table> <tr> <td>Data rate (per second)</td> <td>10</td> </tr> <tr> <td>Titrant delivery rate (mL/min.)</td> <td>2</td> </tr> <tr> <td>No. of exothermic endpoints</td> <td>1</td> </tr> <tr> <td>Data smoothing factor</td> <td>50</td> </tr> </table> <p>Basic method. Depending on likely free acid strength, pipette between 1 – 10mL of test solution into a titration beaker. Add 15mL saturated potassium oxalate solution, and dilute to approximately 25mL with D.I. water. Titrate to the first exothermic endpoint with standard 1 mol/L NaOH.</p> | Data rate (per second) | 10 | Titrant delivery rate (mL/min.) | 2 | No. of exothermic endpoints | 1 | Data smoothing factor | 50 |
| Data rate (per second)          | 10  |                        |    |                                 |   |                             |   |                       |    |
| Titrant delivery rate (mL/min.) | 2   |                        |    |                                 |   |                             |   |                       |    |
| No. of exothermic endpoints     | 1   |                        |    |                                 |   |                             |   |                       |    |
| Data smoothing factor           | 50  |                        |    |                                 |   |                             |   |                       |    |

|                 |  |
|-----------------|--|
| <b>Results:</b> | Analysis of synthetic and actual copper leach and refining solutions |
|-----------------|--|

|  | Sample I.D.                        | Free Acid g/L H <sub>2</sub> SO <sub>4</sub> |
|--|------------------------------------|--|
|  | Synth. Soln. ~9 g/L Fe (III)       | 9.73, 9.80                                   |
|  | Synth. Soln. ~4.5 g/L Fe(III)      | 9.92, 9.94                                   |
|  | Synth. Soln. ~2.25 g/L Fe(III)     | 9.99, 9.99                                   |
|  | Synth. Soln. no Fe(III)            | 10.31, 10.33                                 |
|  | Synth. Electrolyte, 0.2g/L Fe(III) | 202.6, 203.0                                 |
|  | Leach Soln. #1                     | 1.86, 1.88                                   |
|  | Leach Soln. #2                     | 2.12, 2.20                                   |
|  | Leach Soln #3                      | 1.69, 1.65                                   |
|  | “Final” Soln.                      | 6.39, 6.34                                   |

**Calculation:**

$$g/L \text{ Free Acid} = \frac{((\text{Titre} - \text{blank}) \times \text{FW } H_2SO_4 \times M \text{ NaOH})}{(\text{sample vol, mL} \times 2)}$$

$$g/L \text{ Free Acid} = \frac{((4.067 - 0.031) \times 98.074 \times 1.0256)}{(1.00 \times 2)} = 203.0 \text{ g/L}$$

**Thermometric Titration Plots:**




Fig.1. Free acid in electrolyte solutions




Fig. 2. Free acid in "final" solutions




Fig. 3. Free acid in leach solutions

Legend:

Red = solution temperature curve

Black = second derivative curve