

Biomarker Identification in Nonalcoholic Fatty Liver Disease via GC/MS

Todd Richards¹, David E. Alonso¹, Xiang Zhang² and Joseph E. Binkley¹ | ¹LECO Corporation, St. Joseph, MI; ²Department of Chemistry, University of Louisville, Louisville, KY

Introduction

- Non-alcoholic fatty liver disease (NAFLD) is a major medical problem in developed countries
- NAFLD is associated with:
 - Increase obesity and insulin resistance
 - Metabolic syndrome
- NAFLD is influenced by lifestyle:
 - Exercise
 - Diet
- Research Samples & Study
 - 168 Rodent samples (original sample set, Pegasus 4D Study)
 - Investigate two critical factors: 1) Increased fructose consumption and 2) Inadequate Cu intake
 - Complementary investigation: Utilize a more sensitive instrument (20 fg) with extended linear dynamic range (10⁶) to further characterize the samples
- Analysis Objective: Identify potential NAFLD biomarkers in pooled rodent samples using high performance, benchtop GC-TOFMS (Figure 1)



Figure 1. Pegasus[®] BT GC-TOFMS.

Workflow & Methods

- Sample preparation
 - Dry: 1) Speed Vac & 2) Lyophilizer
 - Derivatize: 1) MEOX (30 µL, 20mg/mL) & 2) MTBSTFA (20 µL)
 - (For each step, heat at 60 °C for 1 hr)
- Analyze derivatized samples using high performance GC-TOFMS, ChromaTOF[®] brand software, and statistical processing software
- Instrument Parameters




Figure 2. Analysis Workflow.

Table 1. GC-TOFMS Parameters

Gas Chromatograph	Agilent 7890 with LECO L-Pal3 Autosampler
Injection	1 µL, Split 20:1; 270 °C
Carrier Gas	He @ 0.8 mL/min, Constant Flow
Column	Rxi-5 Sil MS, 20 m x 0.18 mm i.d. x 0.18 µm (Restek, Bellefonte, PA, USA)
Temperature Program	60 °C (0.50 min), ramped 36 °C/min to 320 °C (3 min)
Mass Spectrometer	LECO Pegasus BT
Ion Source Temperature	250 °C
Ionization Mode	El
Mass Range (m/z)	45-600
Acquisition Rate	10 spectra/s

GC-TOFMS Results




Figure 3. AICs for 6 pooled liver samples.




Figure 4. A) eXtracted Ion Chromatogram (XIC) & B) TIC of liver sample 6*.

Table 2. Representative analytes in liver sample 6

Name	Formula	R.T. (s)	Area	Similarity	Name	Formula	R.T. (s)	Area	Similarity
Lactic acid, 2TBDMs derivative	C ₁₀ H ₁₆ O ₄ Si ₂	273	3889358872	866	L-Serine, 3TBDMs derivative	C ₁₃ H ₂₀ N ₂ O ₄ Si ₃	366	1804043501	930
Glycolic acid, 2TBDMs derivative	C ₁₀ H ₁₆ O ₄ Si ₂	276	392513491	865	Thr, (3TBDMs)-	C ₁₃ H ₂₀ N ₂ O ₄ Si ₃	372	1907273568	870
Alanine diTBDMs	C ₁₃ H ₂₀ N ₂ O ₄ Si ₃	283	2627670348	894	Myristic acid, 2TBDMs derivative	C ₁₉ H ₃₄ O ₄ Si ₂	379	2626194161	906
Glycine, 2TBDMs derivative	C ₁₀ H ₁₆ O ₄ Si ₂	287	2093903802	909	Alanine, 2TBDMs derivative	C ₁₀ H ₁₆ O ₄ Si ₂	382	246632471	802
3-Hydroxybutyric acid, (R)-, 2TBDMs derivative	C ₁₃ H ₂₀ N ₂ O ₄ Si ₃	294	1338064236	911	L-Phenylalanine, 2TBDMs derivative	C ₁₇ H ₂₂ N ₂ O ₄ Si ₃	383	861540012	913
L-2-Aminobutyric acid, 2TBDMs derivative	C ₁₀ H ₁₆ O ₄ Si ₂	297	980488165	889	Malic acid, 3TBDMs derivative	C ₁₂ H ₁₈ O ₄ Si ₂	384	133024256	852
β-Alanine, 2TBDMs derivative	C ₁₀ H ₁₆ O ₄ Si ₂	303	811778683	856	Aspartic acid, triTBDMs	C ₁₄ H ₂₀ N ₂ O ₄ Si ₃	391	1444900099	897
Dimethylpentobarbital	C ₁₃ H ₂₀ N ₂ O ₄ Si ₃	306	445555160	900	L-Cysteine, 3TBDMs derivative	C ₁₃ H ₂₀ N ₂ O ₄ Si ₃	399	533411147	943
Valine diTBDMs	C ₁₃ H ₂₀ N ₂ O ₄ Si ₃	307	208564362	864	L-Glutamic acid, 3TBDMs derivative	C ₁₃ H ₂₀ N ₂ O ₄ Si ₃	408	1539429832	939
L-Leucine, 2TBDMs derivative	C ₁₃ H ₂₀ N ₂ O ₄ Si ₃	314	3168854906	937	Palmitic acid, TBDMs derivative	C ₂₁ H ₄₀ O ₂ Si	410	1491680293	833
Niacinamide, TBDMs derivative	C ₁₃ H ₂₀ N ₂ O ₄ Si ₃	318	2512779336	923	L-Asparagine, 3TBDMs derivative	C ₁₇ H ₂₄ N ₄ O ₄ Si ₃	414	1044773406	873
Isoleucine, 2TBDMs derivative	C ₁₃ H ₂₀ N ₂ O ₄ Si ₃	320	2056141263	943	L-Lysine, 3TBDMs derivative	C ₁₃ H ₂₀ N ₂ O ₄ Si ₃	424	1125874542	917
4-Aminobutanoic acid, 2TBDMs derivative	C ₁₃ H ₂₀ N ₂ O ₄ Si ₃	325	46004677	831	Vaccenic acid, (Z)-, TBDMs derivative	C ₂₂ H ₄₀ O ₂ Si	435	4203899325	837
Uracil, 2TBDMs derivative	C ₁₃ H ₂₀ N ₂ O ₄ Si ₃	327	3063305152	929	Histidine triTBDMs	C ₁₇ H ₂₂ N ₄ O ₄ Si ₃	453	1428415311	895
L-Proline, 2TBDMs derivative	C ₁₃ H ₂₀ N ₂ O ₄ Si ₃	328	2171908625	918	Arachidonic acid, TBDMs derivative	C ₂₀ H ₃₈ O ₂ Si	455	1921897209	899
Itaconic acid, 2TBDMs derivative	C ₁₀ H ₁₆ O ₄ Si ₂	331	211731605	864	L-Tyrosine, 3TBDMs derivative	C ₁₇ H ₂₂ N ₂ O ₄ Si ₃	459	1990728852	922
D-Pyroglutamic acid, 2TBDMs derivative	C ₁₃ H ₂₀ N ₂ O ₄ Si ₃	360	2010302971	854	Uric acid, 4TBDMs derivative	C ₁₂ H ₁₂ N ₄ O ₄ Si ₂	488	239470552	898
Met, 2TBDMs	C ₁₃ H ₂₀ N ₂ O ₄ Si ₃	363	795167581	868	L-Tryptophan, 3TBDMs derivative	C ₁₇ H ₂₂ N ₂ O ₄ Si ₃	495	1189857045	866
					Cholesterol, TBDMs derivative	C ₂₇ H ₄₆ O ₂ Si	630	179334215	823

Ave. = 886/1000

GC-TOFMS Results




Figure 5. A) ChromaTOF deconvolution results - liver sample 6. Peak True (deconvoluted) and library spectra for aconitic acid (B,C) and lysine (D,E).




Figure 6. Pegasus BT high quality spectral data: Peak True and library spectra for pyroglutamic acid (A,B) and threonine (C,D).




Figure 7. Pegasus BT high performance spectral data with mass delta values <0.10 Da: Peak True and library spectra for glutamic acid (A,B) and succinic acid (C,D).

Pegasus BT Data Export and Statistical Analysis

- Raw data was exported (netCDF) using ChromaTOF (Figure 8)
- Data was uploaded and processed using XCMS online (xcmsonline.scripps.edu)
- PCA and Cloud Plot data (Figure 9) were used to find features that differentiated samples
- Potential biomarker metabolites for NAFLD were confidently identified using ChromaTOF (Figure 10)
- Potential biomarkers characterized included: Amino acids (e.g., histidine, alanine, valine, etc.), fatty acids (e.g., arachidonic acid) and misc. compounds (e.g., uridine, 7-methylxanthine)




Figure 8. ChromaTOF raw data export feature.




Figure 9. XCMS online Cloud Plot (m/z vs retention time) showing a metabolite feature with a p-value of 4.84e-4, retention time of 7.55 min, and EIC m/z = 440.396.




Figure 10. ChromaTOF - A) XIC 440.40, B) Generated formulas for fragment and molecular ions, C, D) Peak True and library match data for histidine.

Summary

- The Pegasus BT instrument produced high quality data
- This data was exported into statistical processing software for sample differentiation and feature selection
- ChromaTOF brand software was used to fully characterize metabolites using spectral similarity searches against large, well-established databases and formula determinations using data with mass deltas <0.10 Da

