

Better Data Faster: MassHunter Quantitative Analysis



Jade C. Byrd

MassHunter Product Manager

July 29, 2015

Agenda

Brief Introduction

How to fully leverage the information in the raw data during data analysis

Tips and Tricks for fast data review of even thousands of complex samples

Ways to automate peak detection to eliminate the need for manual integration

New features to quickly and automatically create pdf reports

Integrating MassHunter Quantitative Analysis with other laboratory systems


Labs Have More Samples than Ever



Agilent MassHunter Suite


Acquisition

- Agilent LCs, GCs, ICP and MS
- dMRM, tMRM, All Ions, Ion Mobility




Quantitative Analysis

- High-throughput quantitation
- Easily visualize results (Review-by-Exception)
- Integration with LIMS/LIS
- Custom reports



Qualitative Analysis

- ID and characterize compounds in all markets
- Novel algorithms for feature finding



MassHunter for all your needs



MassHunter Software

One software for all Agilent Mass Spectrometers

Minimize the training costs and delays – A Single User Experience for all techniques



LC/MS



GC/MS




ICP-MS

MassHunter Quantative Analysis Workflow



MassHunter Quantative Analysis Workflow



Including Additional Sample Information about Samples

We have an electronic sample list with sample information (pH, company identifiers, etc.) and we wish this information to appear in the final reports. How can we do this?

- Worklist/Sequence has two columns, “Comments” and “Info.” that will get read by MassHunter Quantitative Analysis.
- 10 additional columns can be created that are labeled “UserDefined, UserDefined1, UserDefined2, etc.”
- Custom columns of other names will get stored in the .D folder but not read by MassHunter Quantitative Analysis to be included in the results for reporting.

The screenshot shows a dialog box titled "Add Columns". On the left, under "Column Type", a list includes "MS Parameter", "User Defined" (highlighted), "MFC", "Protein", "Sequences", "Modification Profiles", and "Masses". On the right, under "Column Information", the "User-defined information" section contains: "Column name:" with the text "UserDefined2", "Value:" with the text "12345", and "User defined data type" with radio buttons for "Text" and "Number" (the "Number" option is selected). At the bottom right are "OK" and "Cancel" buttons.

Remembering Worklist/Sequence Column Layout



We hide, add, or reorder columns but next time MassHunter Acquisition brings it up the original way. How do we get MassHunter to remember these changes?

For LC/MS


- After getting Worklist the way we want, save the Worklist to the MassHunter/Worklist Templates folder
- Rename the file from myworklist.wkl to myworklist.wkt
- To use this template, select New Worklist
- If we don't need more than one, we name it default.wkt

For GC/MS

- Use the Sequence Table Tools



MassHunter Quantative Analysis Workflow




Unified Method


One location, one name


- Defines acquisition method
- Defines analysis method
- Defines report method

Allows granular control
over analysis and reporting
by sample





MassHunter Quantative Analysis Workflow





Sample						
	Name	Data File	Type	Level	Acq. Date-Time	
	Blank-1	CMAMBlk_01.d	Blank		5/12/2006 1:48 PM	
	Calib-L1	CMAMCal_L1.d	Cal	L1	5/12/2006 1:51 PM	
	Calib-L2	CMAMCal_L2.d	Cal	L2	5/12/2006 1:54 PM	
	Calib-L3	CMAMCal_L3.d	Cal	L3	5/12/2006 1:57 PM	
	Calib-L4	CMAMCal_L4.d	Cal	L4	5/12/2006 2:00 PM	
	Calib-L5	CMAMCal_L5.d	Cal	L5	5/12/2006 2:03 PM	
	QC-L2	CMAMQC_L2.d	QC	L2	5/12/2006 2:06 PM	
	QC-L4	CMAMQC_L4.d	QC	L4	5/12/2006 2:09 PM	
	Sample-1	CMAMSam_01.d	Sample		5/12/2006 2:12 PM	
	Sample-2	CMAMSam_02.d	Sample		5/12/2006 2:15 PM	
	Sample-3	CMAMSam_03.d	Sample		5/12/2006 2:18 PM	



Agilent MassHunter Quantitative Analysis (for QQQ) - DrugsOfAbuse - NewBatch.batch.bin

File Edit View Analyze Method Update Library Report Tools Help

New Method from Acquired MRM Data...
 New Method from Acquired Scan Data...
 New Method from Acquired Scan Data with Library Search...
 New Method from Acquired Chromatographic Data...
 New Method using Manual Setup
 New Method from File...

L2	5/12/2006 1:54 PM
L3	5/12/2006 1:57 PM
L4	5/12/2006 2:00 PM
L5	5/12/2006 2:03 PM
L2	5/12/2006 2:06 PM
L4	5/12/2006 2:09 PM
	5/12/2006 2:12 PM
	5/12/2006 2:15 PM
	5/12/2006 2:18 PM

Compound Information

Calibration Curve

Type: Origin: Weight: ISTD QC CC

Chromatogram

Responses

Modified 11 Samples (11 total)

Agilent MassHunter Quantitative Analysis (for QQQ) - DrugsOfAbuse - NewBatch.batch.bin

File Edit View Analyze Method Update Library Report Tools Help

New Method from Acquired MRM Data...
 New Method from Acquired Scan Data...
 New Method from Acquired Scan Data with Library Search...
 New Method from Acquired Chromatographic Data...
 New Method using Manual Setup

ISTD:

New Method from Acquired Data

Look in: DrugsOfAbuse

Name	Date modified	Type	Size
CMAMBIK_01.d	1/6/2015 8:39 AM	File folder	
CMAMCaL_L1.d	1/6/2015 8:39 AM	File folder	
CMAMCaL_L2.d	1/6/2015 8:39 AM	File folder	
CMAMCaL_L3.d	1/6/2015 8:39 AM	File folder	
CMAMCaL_L4.d	1/6/2015 8:39 AM	File folder	
CMAMCaL_L5.d	1/6/2015 8:39 AM	File folder	
CMAMQC_L2.d	1/6/2015 8:39 AM	File folder	
CMAMQC_L4.d	1/6/2015 8:39 AM	File folder	
CMAMSam_01.d	1/6/2015 8:39 AM	File folder	
CMAMSam_02.d	1/6/2015 8:39 AM	File folder	
CMAMSam_03.d	1/6/2015 8:39 AM	File folder	
DOA_ReportMethod.m	11/19/2014 11:34 ...	File folder	
QuantReports	6/18/2015 10:11 AM	File folder	
QuantResults	7/28/2015 9:54 AM	File folder	
Background - DOA in OF poster ASMS.ppt	5/23/2006 8:01 PM	Microsoft PowerP...	1,855 KB

Object name:

Objects of type:

Open Cancel

Chromatogram

Compound Information

Modified 11 Samples (11 total)

Method Tasks

- New / Open Method
- Workflow
- Method Setup Tasks
 - MRM Compound Setup
 - Retention Time Setup
 - ISTD Setup
 - Concentration Setup
 - Qualifier Setup
 - Calibration Curve Setup
- Globals Setup
- Save / Exit
 - Validate
 - Save
 - Save As...
- Exit
- Manual Setup Tasks
- Outlier Setup Tasks
- Advanced Tasks

Method Table


Level Name Prefix: # of Levels: 10 Create Levels Time Segment: <All> Compound: Reset Table View

Sample	Name	Data File	Type	Level	Acq. Method File	Acq. Date-Time
▶	CMAMCal_L3.d	CMAMCal_L3.d				

Quantifier										
Name	TS	Transition	Scan	Type	ISTD Compound Name	CAS#	ISTD Flag	ISTD Conc.	Time Reference Flag	
Amphetamine	1	136.2 -> 91.4	MRM	Target	<None>		<input type="checkbox"/>		<input type="checkbox"/>	
Qualifier										
Precursor Ion	Product Ion	Transition	Rel. Resp.	Uncertainty						
136.2	119.4	136.2 -> 119.4	26.6	20.0						
Quantifier										
Name	TS	Transition	Scan	Type	ISTD Compound Name	CAS#	ISTD Flag	ISTD Conc.	Time Reference Flag	
Amphetamine-d5	1	141.1 -> 93.4	MRM	ISTD	<None>		<input checked="" type="checkbox"/>		<input type="checkbox"/>	
Qualifier										
Precursor Ion	Product Ion	Transition	Rel. Resp.	Uncertainty						
141.1	124.4	141.1 -> 124.4	26.2	20.0						
Quantifier										
Name	TS	Transition	Scan	Type	ISTD Compound Name	CAS#	ISTD Flag	ISTD Conc.	Time Reference Flag	
Cocaine	1	304.1 -> 182.0	MRM	Target	<None>		<input type="checkbox"/>		<input type="checkbox"/>	
Qualifier										
Precursor Ion	Product Ion	Transition	Rel. Resp.	Uncertainty						
304.1	82.0	304.1 -> 82.0	3.9	20.0						
Quantifier										
Name	TS	Transition	Scan	Type	ISTD Compound Name	CAS#	ISTD Flag	ISTD Conc.	Time Reference Flag	
Cocaine-d3	1	307.1 -> 185.0	MRM	ISTD	<None>		<input checked="" type="checkbox"/>		<input type="checkbox"/>	
Qualifier										
Precursor Ion	Product Ion	Transition	Rel. Resp.	Uncertainty						
307.1	85.0	307.1 -> 85.0	4.4	20.0						

Sample Information

Compound Information



4 Compounds (4 total) 4 ISTD (4 total)

Target, Internal Standard and Qualifier Ion Ratios learned from the Data

Quantifier					
Name	TS	Transition	Scan	Type	
Amp	1	136.2 -> 91.4	MRM	Target	
Qualifier					
Precursor Ion	Product Ion	Transition	Rel. Resp.	Uncertainty	
136.2	119.4	136.2 -> 119.4	26.6	20.0	

Quantifier					
Name	TS	Transition	Scan	Type	
Amp-d5	1	141.1 -> 93.4	MRM	ISTD	
Qualifier					
Precursor Ion	Product Ion	Transition	Rel. Resp.	Uncertainty	
141.1	124.4	141.1 -> 124.4	26.2	20.0	

Name	TS	RT	Transition	ISTD Flag	
13C-Ochratoxin A	1	6.425	424.2 -> 250.1	<input checked="" type="checkbox"/>	
13C-T2-Toxin	1	5.959	508.3 -> 322.1	<input checked="" type="checkbox"/>	
13C-Zearalenone	1	6.437	335.2 -> 290.0	<input checked="" type="checkbox"/>	
Aflatoxin B1	1	4.426	313.1 -> 241.0	<input type="checkbox"/>	
Aflatoxin B2	1	4.176	315.1 -> 258.9	<input type="checkbox"/>	
Aflatoxin G1	1	3.868	329.1 -> 243.0	<input type="checkbox"/>	
Aflatoxin G2	1	3.585	331.1 -> 313.0	<input type="checkbox"/>	
Deoxynivalenol	1	1.493	297.1 -> 249.0	<input type="checkbox"/>	

New in B.07.00 – SIM from GC/MSD works the same as MRM

Available for MSD, QQQ, TOF, and Q-TOF

Peak Detection

Each peak can be configured to use 1 of 5 peak detection algorithms

The screenshot displays the Agilent MassHunter Quantitative Analysis (for QQQ) - [New Method] interface. The main window shows the Method Table with the following data:

Sample	Name	Data File	Type	Level	Acq. Method File	Acq. Date-Time
...	CMAMCal_L3.d	CMAMCal_L3.d				

Quantifier							
Name	TS	Transition	Scan	Type	RT	Int.	Int. Params.
Amp	1	136.2 -> 91.4	MRM	Target	2.134	Agile2	...






Qualifier						
Precursor Ion	Product Ion	Transition	Rel. Resp.	Uncertainty	Int. Params.	
136.2	119.4	136.2 -> 119.4	26.6	20.0		

Quantifier							
Name	TS	Transition	Scan	Type	RT	Int.	Int. Params.
Amp-d5	1	141.1 -> 93.4	MRM	ISTD	2.122	Agile2	...

The Integration Parameters Setup dialog box is open, showing the Integrator dropdown menu with the following options:

- Agile2
- Agile
- MS-MS
- MS-MS (GC)
- General
- Universal
- Spectrum Summation

Agile and Agile2 – parameter-less integrators



Compound Group Sample & Compound based Review/Report

The screenshot shows the Agilent MassHunter Quantitative Analysis (for GCMS) interface. The main window displays a 'Method Table' for sample 'CAL_L03'. The table lists various compounds and their associated data files, scan types, and retention times. A red box highlights the 'Cmpd. Group' column, which contains values like 'A,B' and 'B' for different compounds.


Sample									
Name	Data File	Type	Level	Acq. Method File	Acq. Date-Time				
CAL_L03	CAL_L03.D	Cal	3	624A.M	6/20/2008 8:53...				
Quantifier									
Name	TS	Scan	Type	MZ	RT	Ion Polarity	Criteria	Cmpd. Group	
Tetrahydrofuran	1	Scan	Target	42.0	9.311	Positive	Close RT	A,B	
Acetone	1	Scan	Target	43.0	6.194	Positive	Close RT	B	
Vinyl Acetate	1	Scan	Target	43.0	8.064	Positive	Close RT	A,B	
2-Butanone (ME...	1	Scan	Target	43.0	8.437	Positive	Close RT	A,B	
4-Methyl-2-Pent...	1	Scan	Target	43.0	12.382	Positive	Close RT	A,B	
Carbon Disulfide	1	Scan	Target	44.0	7.177	Positive	Close RT	B	
Methylene Chlори...	1	Scan	Target	49.0	6.877	Positive	Close RT	B	
Chloromethane	1	Scan	Target	50.0	4.493	Positive	Close RT	B	
2-Hexanone	1	Scan	Target	58.0	13.562	Positive	Close RT	A,B	

Useful when there is one analysis method but wish to review and report different sets of the compound list based on individual samples

Also used for StreamSelect (2 or 4 streams)

Compound Math: group peaks for calibration and quantitation

Five peaks representing Aroclor PCBs mixture



Method Table

Compound: Total 1016-1 Reset Table View

Sample				
Name	Data File	Type	Level	Acq. M
1016/1260 @ 50...	GC7-1842.D	Cal	500	
Quantifier				
Name	Type	Cmpd. Group	Compound Math	
TCMX	Target	S		
Aroclor1016-1	Target	L1-1		
Aroclor1016-2	Target	L1-1		
Aroclor1016-3	Target	L1-1		
Aroclor1016-4	Target	L1-1		
Aroclor1016-5	Target	L1-1		
Aroclor1260-1	Target	L2-1		
Aroclor1260-2	Target	L2-1		
Aroclor1260-3	Target	L2-1		
Aroclor1260-4	Target	L2-1		
Aroclor1260-5	Target	L2-1		
DCB	Target	S		
▶ Total 1016-1	Target	L1-1	Concentration Sum	
Total 1260-1	Target	L2-1	Concentration Sum	

New Find Feature in Method Editor

Method Table

Sample	Name	Data File	Type	Level	Acq. Method File	Acq. Date-Time
Blank-1		CMAMBik_01.d	Blank		APClautotune.m	5/12/2006 1:48...

Quantifier	Name	TS	Transition	Scan	Type	Precursor Ion	Product Ion	RT	Ion Polarity
▶	Amp	1	136.2 -> 91.4	MRM	Target	136.2	91.4	2.102	Positive
	Amp-d5	1	141.1 -> 93.4	MRM	ISTD	141.1	93.4	2.078	Positive
	Meth	1	150.1 -> 119.3	MRM	Target	150.1	119.3	2.239	Positive
	Meth-d5	1	155.2 -> 92.3	MRM	ISTD	155.2	92.3	2.233	Positive
	MDMA	1	194.2 -> 163.3	MRM	Target	194.2	163.3	2.269	Positive
	MDMA-d5	1	199.2 -> 164.3	MRM	ISTD	199.2	164.3	2.269	Positive
	Cocaine	1	304.1 -> 182.0	MRM	Target	304.1	182.0	2.449	Positive
	Cocaine-d3	1	307.1 -> 185.0	MRM	ISTD	307.1	185.0	2.450	Positive

Find

Find: Quantifier

Column: Product Ion

Operator: Equals

Value: 135


Find Next Close

Operator: Equals

Is null
Is not null
Equals
Does not equal
Contains
Does not contain
Greater than
Greater than or equal to
Less than
Less than or equal to

Search quickly in methods with a large number of compounds

After applying the Analysis Method



2 Configurable Views: both take advantage of Outliers

Batch-at-a-Glance


Batch Table

Sample: Sample Type:

Sample							Amp Results		
?	▼	Name	Data File	Type	Level	Acq. Date-Time	RT	Final Conc.	Accuracy
▶	▼	Calib-L1	CMAMCal_L1.d	Cal	L1	5/12/2006 1:51 PM	2.140	3.2209	128.8
	▼	Calib-L2	CMAMCal_L2.d	Outlier(s)					15.0
		Calib-L3	CMAMCal_L3.d	Amp: Accuracy value = 128.8 is outside the allowed range [80.0, 120.0]					09.4
		Calib-L4	CMAMCal_L4.d	Cal	L4	5/12/2006 2:00 PM	2.022	26.7565	107.0
		Calib-L5	CMAMCal_L5.d	Cal	L5	5/12/2006 2:03 PM	2.101	124.4862	99.6
		QC-L2	CMAMQC_L2.d	QC	L2	5/12/2006 2:06 PM	2.142	5.2293	104.6
		QC-L4	CMAMQC_L4.d	QC	L4	5/12/2006 2:09 PM	2.135	27.8044	111.2
		Sample-2	CMAMSam_02.d	Sample		5/12/2006 2:15 PM	2.143	4.8978	
		Sample-3	CMAMSam_03.d	Sample		5/12/2006 2:18 PM	2.105	14.2185	

http://www.chem.agilent.com/en-US/products-services/Software-Informatics/MassHunter-Workstation-Software/Pages/Batch_at_a_glance_demo.aspx

Compounds-at-a-Glance



http://www.chem.agilent.com/en-US/products-services/Software-Informatics/MassHunter-Workstation-Software/Pages/quant_flat_table.aspx

Outlier Options: Configured in the Method

40+ Outliers can be configured (B.07.01):

1. Retention Time
2. Relative Retention Time
3. Peak Resolution
4. Peak Symmetry
5. Peak Purity
6. Signal-to-noise Ratio
7. Limit of Detection
8. Limit of Quantitation
9. Method Detection Limit
10. Qualifier Ratio
11. Internal Standard Response
12. Internal Standard % Deviation
13. Sample Amount
14. Sample % RSD
15. Blank Concentration
16. Blank Response
17. Accuracy
18. Average Response Factor
19. Average Response Factor %RSD
20. Curve Fit R²
21. Relative Response Factor
22. Response Factor
23. QC Max Deviation
24. QC Relative Standard Deviation
25. QC LCS Recovery
26. CC Average Response Factor
27. CC Internal Std Response Ratio
28. CC Relative Response Factor
29. CC Response Ratio
30. CC Retention Time
31. Matrix Spike
32. Matrix Spike % Difference
33. Matrix Spike Percent Recovery
34. Matrix Spike Group Recovery
35. Surrogate
36. Surrogate % Recovery
37. Response Check
38. Mass accuracy
39. Mass Match Score
40. Library Match Score
- 41. Custom Calculation**



Method Tasks

- New / Open Method
- Workflow
- Method Setup Tasks
- Save / Exit
- Manual Setup Tasks
- Outlier Setup Tasks**

- Retention Time
- Relative Retention Time
- Peak Resolution
- Peak Symmetry
- Peak Full Width Half Maximum
- Peak Purity
- Signal-to-Noise Ratio
- Limit Of Detection
- Limit Of Quantitation
- Method Detection Limit
- Qualifier Ratio
- ISTD Response
- ISTD Response Percent Deviation
- Sample Amount
- Sample RSD
- Blank Concentration
- Blank Response
- Accuracy
- Average Response Factor
- Average Response Factor RSD
- Curve Fit R²
- Relative Response Factor
- Response Factor
- QC

Method Table

Level Name Prefix: # of Levels: 10 Create Levels Time Segment: <All>



Sample	Name	Data File	Type	Level	Acq. Method File	Acq. Date-Time
	Calib-L1	CMAMCal_L1.d	Cal	L1	APCIautotune.m	5/12/2006 1:51...
Quantifier						
	Name	TS	Transition	Scan	Type	QC Max. Dev. /
	7-amino clonaz...	1	286.0 -> 121.0	MRM	Target	0.2
	7-amino flunitraz...	1	284.0 -> 226.0	MRM	Target	0.2
	D4-7-amino clon...	1	290.0 -> 226.0	MRM	ISTD	0.2
	zopiclone	2	389.1 -> 245.1	MRM	Target	0.2
	alpha-OH alpraz...	3	325.0 -> 297.0	MRM	Target	0.2
	clonazepam	3	316.0 -> 270.0	MRM	Target	0.2
	D4-clonazepam	3	320.0 -> 274.0	MRM	ISTD	0.2
	flunitrazepam	3	314.0 -> 268.0	MRM	Target	0.2
	alprazolam	4	309.0 -> 281.0	MRM	Target	0.2
	D5-alprazolam	4	314.0 -> 286.0	MRM	ISTD	0.2
	D5-oxazepam	4	292.0 -> 246.0	MRM	ISTD	0.2
	flurazepam	4	388.0 -> 315.0	MRM	Target	0.2
	lorazepam	4	321.0 -> 275.0	MRM	Target	0.2
	oxazepam	4	287.0 -> 241.0	MRM	Target	0.2
	triazolam	4	343.0 -> 308.0	MRM	Target	0.2
	D5-temazepam	5	306.0 -> 260.0	MRM	ISTD	0.2
	desalkylflurazep...	5	289.0 -> 226.0	MRM	Target	0.2
	temazepam	5	301.0 -> 255.0	MRM	Target	0.2
	chlordiazepoxide	6	300.0 -> 283.0	MRM	Target	0.2
	D5-nordiazepam	6	276.0 -> 140.0	MRM	ISTD	0.2
	nordiazepam	6	271.0 -> 140.0	MRM	Target	0.2
	diazepam	7	285.0 -> 222.0	MRM	Target	0.2
	midazolam	7	326.0 -> 291.0	MRM	Target	0.2


- Batch Table
- Compound Information
- Calibration Curve
- Sample Information
- Metrics Plot
- Method Table
- Method Development Tasks
- Method Error List
- Status Bar
- Compounds-at-a-Glance...
- Toolbars
- Window Layout
- Batch Table Layout
- Add/Remove Columns...
- Restore Default Columns
- Load Column Settings...
- Save Column Settings...
- Reset Sort
- Lock Sample/Compound Columns
- Expand All
- Collapse All
- Auto Review Samples Ctrl+1
- Auto Review Compounds Ctrl+2

Layout: Restore Default Layout

Sample Type: <All> Compound: Deoxyynivalenol ISTD: 13C-Deoxyynivalenol

Type	Level	Acq. Date-Time	Deoxyynivalenol Res.		Aflatoxin G2 Results		Aflatoxin G1 Results		Aflatoxin B2 Results		Aflatoxin B1 Results		HT2-Toxin Results		Fumonisin B1 Results		T2-Toxin Results		Ochratoxin A Results		Zearalenone Results		Fumonisin B2 Results	
			Final Conc.	Accuracy	Final Conc.	Accuracy	Final Conc.	Accuracy	Final Conc.	Accuracy	Final Conc.	Accuracy	Final Conc.	Accuracy	Final Conc.	Accuracy	Final Conc.	Accuracy	Final Conc.	Accuracy	Final Conc.	Accuracy	Final Conc.	Accuracy
Sample		8/3/2011 8:48 PM	95.3381		28.1892		5.1051		6.8315		8.2204		43.2515		41.4520		7.0834		13.5732		30.5107		51.4217	
Sample		8/3/2011 9:13 PM	0.0000		0.0014		0.0042		0.0035		0.0014		0.0257		0.0853		0.0005		0.0000		0.0000		0.0709	
Sample	1	8/3/2011 9:26 PM	0.0963	128.2	0.0092	122.2	0.0096	127.5	0.0096	128.9	0.0090	120.4	0.0851	126.6	0.0356	143.1	0.0102	134.8	0.0111	145.4	0.0000	0.0	0.0394	157.1
Sample	2	8/3/2011 9:38 PM	0.2454	108.9	0.0248	109.7	0.0229	101.4	0.0225	100.6	0.0214	95.7	0.1969	97.7	0.0767	102.6	0.0284	125.8	0.0228	99.3	0.0825	109.3	0.0516	68.7
Sample	3	8/3/2011 9:51 PM	0.6825	90.9	0.0706	93.7	0.0715	94.8	0.0691	92.6	0.0713	95.6	0.5869	87.3	0.2214	88.9	0.0657	87.2	0.0706	92.2	0.3043	120.9	0.3238	129.2
Sample	4	8/3/2011 10:03 PM	2.0200	89.7	0.2038	90.1	0.2098	92.8	0.2074	92.6	0.2110	94.2	1.9161	95.0	0.7461	99.9	0.1770	78.3	0.2008	87.5	0.6621	87.7	0.6111	81.3
Sample	5	8/3/2011 10:16 PM	6.9948	93.2	0.6943	92.1	0.7180	95.3	0.6972	93.4	0.7293	97.7	6.1391	91.4	2.1444	86.1	0.6900	91.5	0.6851	89.6	2.1540	85.6	2.0497	81.8
Sample	6	8/3/2011 10:29 PM	21.1500	93.9	2.1893	96.8	2.1354	94.4	1.8739	83.7	2.2344	99.8	20.0025	99.2	6.6171	88.6	2.0369	90.1	2.1378	93.1	6.9392	91.9	6.6411	88.4
Sample	7	8/3/2011 10:41 PM	72.9867	97.2	7.3218	97.1	7.2630	96.4	8.3728	112.2	7.4287	99.5	68.6467	102.2	22.4563	90.2	7.1423	94.8	7.1220	93.1	26.4649	105.1	23.1225	92.3
Sample	8	8/3/2011 10:54 PM	216.7979	96.3	21.8531	96.6	21.5742	95.4	21.3032	95.2	21.4249	95.7	203.5765	101.0	73.7492	98.7	21.5221	95.2	22.5482	98.2	75.2638	99.7	74.8952	99.7
Sample	9	8/3/2011 11:06 PM	763.4102	101.7	76.5047	101.5	76.8672	102.0	75.2327	100.8	75.6584	101.4	669.4873	99.6	253.6088	101.9	77.1990	102.4	77.6985	101.6	251.7291	100.0	254.0880	101.4
Sample		8/3/2011 11:19 PM	3286.0170		32.6179		36.3763		40.8611		55.0843		958.9777		113.7979		20.2239		38.1642		98.8875		120.9559	
Sample		8/3/2011 11:57 PM	0.1786		0.0024		0.0055		0.0052		0.0046		0.0137		0.0920		0.0015		0.0052		0.0000		0.0375	
Sample		8/4/2011 12:09 AM	5.0085		0.4949		0.4384		0.4498		0.4654		3.6190		1.5708		0.4315		0.4360		1.8909		1.3844	
Sample		8/4/2011 12:22 AM	15.7770		1.4740		1.4470		1.4434		1.4543		12.1565		5.0211		1.4902		1.6385		4.7109		3.5953	
Sample		8/4/2011 12:35 AM	43.9983		4.4282		4.5907		4.3826		4.7169		38.4992		14.0406		4.3725		3.8254		17.2780		10.9553	
Sample		8/4/2011 12:47 AM	149.3063		15.3387		15.8213		15.5714		16.1309		135.1231		46.6525		15.0457		13.9448		57.0075		42.6722	
Sample		8/4/2011 1:00 AM	451.1134		47.0283		47.6888		46.5171		47.5886		429.2127		141.6326		47.7655		42.0208		168.6467		131.0514	
Sample		8/4/2011 1:12 AM	1528.7403		148.9406		150.4875		148.1366		158.5599		1390.5888		490.9122		157.1570		143.3065		514.9689		478.2135	
Sample		8/4/2011 1:25 AM	38.5318		84.1637		5.7807		4.9510		13.8139		155.1843		0.0000		9.4868		15.8711		65.9222		44.1909	





X 3.33 Y 70826.171 | 20 Samples (20 total) | 22 Compounds (22 total)



The image displays four overlapping 'Setup Graphics' dialog boxes, each showing a different tab: 'Samples', 'Compounds', 'Organize', and 'Outlier'.

Setup Graphics - Samples Tab: Shows a list of samples with columns 'Name' and 'Data File'. The list includes Solvent Blank, ISTD-Blank, ISTD-10000, ISTD-3000, ISTD-1000, ISTD-300, and ISTD-100.


Setup Graphics - Compounds Tab: Shows a list of compounds with columns 'Name' and 'Type'. The list includes 13C-Deoxynivalenol (ISTD), Deoxynivalenol (Target), 13C-Aflatoxin G2 (ISTD), Aflatoxin G2 (Target), 13C-Aflatoxin G1 (ISTD), Aflatoxin G1 (Target), 13C-Aflatoxin B2 (ISTD), and Aflatoxin B2 (Target).

Setup Graphics - Organize Tab: Shows options for 'Organize Rows by' (Compounds, Samples) and 'Overlay' (None - target only, None - target and qualifiers, Qualifiers, ISTD, Matrix Spike, Compound Groups, Sample Groups, Compounds, Samples).

Setup Graphics - Outlier Tab: Shows a tree view of metrics and qualifiers. The 'Peak Result' section includes Peak Not Found, Peak Resolution Front, Peak Resolution Rear, Retention Time, Relative Retention Time, Integration Metric, Symmetry, Full Width Half Maximum, Purity, Signal To Noise Ratio, Limit Of Detection, Limit Of Quantitation, and Method Detection Limit. The 'Qualifier' section includes Qualifier Peak Not Found, Qualifier Peak Resolution Front, Qualifier Peak Resolution Rear, Qualifier Signal To Noise Ratio, Qualifier Ratio, Qualifier Integration Metric, Qualifier Peak Symmetry, and Qualifier Peak Full Width Half Maximum. Other sections include ISTD, Sample, Blank, Calibration, and Average Response Factor RSD. A 'Filter' section at the bottom allows selecting 'Show all panes', 'Show panes without outliers', or 'Show panes with outliers'.




NEW in B.07.01





MassHunter Quantative Analysis Workflow



Excel Reports: A Case Study




Reduced Cost for Pesticides Residues Analysis by GC/MS/MS Using Mini-QuEChERS and a High Efficiency Source

Melissa Churley¹ and Joan Stevens²

¹Agilent Technologies Inc., 5301 Stevens Creek Blvd., Santa Clara, CA 95051 USA

²Agilent Technologies Inc., 2850 Centerville Road, Wilmington, DE 19808 USA



Distribution of Average Recoveries (n=6) for 126 Pesticides Spiked at 1, 5, 10 and 50 ng/g in Carrot, Tomato and Celery



These graphs took just **3 minutes** to prepare!


<http://www.chem.agilent.com/Library/applications/5991-6069EN.pdf>

How reports work in MassHunter Quantitative Analysis




Significant Gains in Performance

Quant Report Generation Time in Excel



Performance Comparison




**Also benefit of modernization of code base

Automatically Generating Reports

Report Method

- Consolidates report settings
 - Templates, Graphic settings, post process
- Supports both Excel and PDF



Include Report Method in Unified Method


Dozens of PDF Report Templates Available at Installation



AuditTrail.report.xml	5/19/2015 5:42 PM	XML Document	1 KB
DrugAnalysis.report.xml	5/19/2015 5:42 PM	XML Document	1 KB
DrugAnalysis_DopingScreening.report.xml	5/19/2015 5:42 PM	XML Document	1 KB
Env_CC_Avg.report.xml	5/19/2015 5:42 PM	XML Document	1 KB
Env_CC_MidPoint.report.xml	5/19/2015 5:42 PM	XML Document	1 KB
Env_CC_Previous.report.xml	5/19/2015 5:42 PM	XML Document	1 KB
Env_DualGCResults.report.xml	5/19/2015 5:42 PM	XML Document	1 KB
Env_DualGCResults_TIC.report.xml	5/19/2015 5:42 PM	XML Document	1 KB
Env_InitialCal.report.xml	5/19/2015 5:42 PM	XML Document	1 KB
Env_LCSpSpike.report.xml	5/19/2015 5:42 PM	XML Document	1 KB
Env_MSD.report.xml	5/19/2015 5:42 PM	XML Document	1 KB
Env_QA_Check.report.xml	5/19/2015 5:42 PM	XML Document	1 KB
Env_Results.report.xml	5/19/2015 5:42 PM	XML Document	1 KB
Env_Results_withGraphics.report.xml	5/19/2015 5:42 PM	XML Document	1 KB
Env_TPH_Validation.report.xml	5/19/2015 5:42 PM	XML Document	1 KB
Gen_ByCompound.report.xml	5/19/2015 5:42 PM	XML Document	1 KB
Gen_BySample.report.xml	5/19/2015 5:42 PM	XML Document	1 KB
Gen_BySample_withSN.report.xml	5/19/2015 5:42 PM	XML Document	1 KB
Gen_Calibration.report.xml	5/19/2015 5:42 PM	XML Document	1 KB
Gen_Complete.report.xml	5/19/2015 5:42 PM	XML Document	1 KB
Gen_LIMsIntegration.report.xml	5/19/2015 5:42 PM	XML Document	1 KB
Gen_ResultsQualifierRatio.report.xml	5/19/2015 5:42 PM	XML Document	1 KB
Gen_ResultsSummary.report.xml	5/19/2015 5:42 PM	XML Document	1 KB
Gen_Samples.report.xml	5/19/2015 5:42 PM	XML Document	1 KB
Gen_Samples.template.xml	6/1/2015 3:45 PM	XML Document	84 KB
gen_samplesJade.template.xml	6/1/2015 4:21 PM	XML Document	84 KB
LIMsIntegration_ExportResults.report.xml	5/19/2015 5:42 PM	XML Document	1 KB
logo.tif	5/19/2015 5:42 PM	TIFF image	743 KB
OPCW.report.xml	5/19/2015 5:42 PM	XML Document	1 KB
Outlier-BlankConc.template.xml	5/19/2015 6:14 PM	XML Document	72 KB
PAH_BioMarkers.report.xml	5/19/2015 5:42 PM	XML Document	1 KB
Pesticide_Residues.report.xml	5/19/2015 5:42 PM	XML Document	1 KB
QuantCSV.report.xml	5/19/2015 5:42 PM	XML Document	1 KB
SIMScan.report.xml	5/19/2015 5:42 PM	XML Document	1 KB
TargetedDeconvolution.report.xml	5/19/2015 5:42 PM	XML Document	1 KB

PDF Logo File

Customize the PDF reports by replacing the shipping Agilent logo file with a desired graphic (jpeg, etc.)



2 mechanisms to create reports: via Excel or Directly (PDF reporting)

Excel

Fully supported

Uses Industry Standard Excel Plug-In

Available since initial version of MassHunter Quantitative Analysis (~10 years), 400+ templates

Excel file format can be useful for repurposing

- Word documents, Power Point Presentations

Rich content, feature-rich

PDF Reporting

Available since B.06.00 SP1

Significantly faster, especially for large data sets

Graphics are PDF direct content, virtually zero Graphics Device Interface resources consumed


Minimum Memory Usage

Logos, labels and graphics can be customized

Based on Python programming language

Foundation for LIMS/LIS fileless integration

MassHunter Quantative Analysis Workflow




Integration with LIMS/LIS/Intranet, etc.


Done via reporting method

Tutorial videos include details


- SQL Direct Insert is fileless transfer




Uploading Results




*.csv file



LIMS, LIS,
Intranet



“SQL Direct Insert”
SQL Insert Statements



LIMS, LIS,
Intranet


User Community Resources Highlights

2 min videos and 30-60 min webcasts on www.agilent.com

In online Help, Videos, What's new

Migration from ChemStation to MassHunter

30 min. to 1 hr. Technical Webinars



Library>Videos> Search "MassHunter"

Agilent MassHunter Quantitative Analysis Video Demos


What's New in Quant B.07.01

- Centralized Network Solutions
 - o [OpenLAB ECM MassHunter Compliance \(Starter\)](#) **#**
 - o [Network Attached Storage \(NAS\) \(Advanced\)](#) **#**
- GC-MSD Library Search
 - o [Library Method \(Starter\)](#) **#**
 - o [Library Search \(Starter\)](#) **#**
 - o [Library Search Results \(Advanced\)](#) **#**
 - o [Go To NIST MS Program \(Advanced\)](#) **#**
- GC-MSD Library Search Report
 - o [TIC Library Search report \(Starter\)](#) **#**
 - o [TIC Library Search report Technical walkthrough \(Advanced\)](#) **#**
- GC Workflows
 - o [GC Instrument flavor – Simplified presentation \(Starter\)](#) **#**
 - o [Chromatogram Information – GC/TIC signal overlay \(Starter\)](#) **#**
 - o [Compound Math – Aggregate peaks \(Advanced\)](#) **#**
 - o [GC PCB Aroclor Workflow \(Advanced\)](#) **#**
- GC-MSD Quant Workflows
 - o [GC-MSD Target Deconvolution Wizard \(Starter\)](#) **#**
- GC-MSD EnviroQuant Workflows
 - o [Total Petroleum Hydrocarbons \(TPH\) \(Starter\)](#) **#**
 - o [GC-MSD EnviroQuant Tune Evaluation enhancements \(Starter\)](#) **#**
 - o [Peak Annotation Addin – Manual Integration Reasons \(Advanced\)](#) **#**
- QQQ Workflows
 - o [GC-QQQ Dioxin PCBs \(Advanced\)](#) **#**
 - o [LC-QQQ StreamSelect \(Starter\)](#) **#**
- LIMs Integration
 - o [SQL Direct Insert using Python \(Advanced\)](#) **#**
 - o [SQL Direct Insert using Python – Technical walkthrough \(Expert\)](#) **#**



Publication Part Number: G6845-90031

http://www.chem.agilent.com/Library/user_manuals/Public/EQ_ChemStation.pdf



www.chem.agilent.com/en-US/Training-Events/eSeminars/236E/Pages/default.aspx

Feedback about MassHunter Quantitative Analysis

“We work with QC analyses that involves tracking 40 parameters, looking at a sequence of 10 samples at a time... before, we had to go sample by sample, data file by data file to evaluate the sequence... Now, with MassHunter, we have created one customized overview screen that enables us to scroll parameter by parameter all the QCs for the full sequence very quickly.”

DR. ARMAND VERBUEKEN, GC/MS USER, BAYER, ANTWERP, BELGIUM ENVIRONMENTAL TESTING LAB

“At the American games in Mexico, having an option like Compounds-at-a-Glance made it very easy to detect which compounds are positive. Also, having the capability to customize our reporting is very important because we are continuously adding new compounds to the methods and we need to incorporate them into the reporting template.”

DR. JOSEP MARCOS, LC/MS USER, THE HOSPITAL DEL MAR IN BARCELONA

