


Agilent J&W GC Columns for Energy and Chemical Applications

# POWER YOUR PRODUCTIVITY AND ENSURE TOP PRODUCT QUALITY

The Measure of Confidence



**Agilent Technologies**

## AGILENT J&W PLOT PT GC AND DB-SULFUR SCD COLUMNS

### *Ruggedness and reliability* to meet your lab's demands

From product quality and safety... to price-driven efficiency improvements... to environmental stewardship... the pressure to stay productive, profitable, and compliant is more intense than ever. You simply do not have room in your day for downtime caused by problems with your analytical column.

Porous layer open tubular (PLOT) columns are the GC columns of choice for the analysis of light gases and volatile compounds because of their high retentive character, and therefore, have been widely adopted in the hydrocarbon processing industry. However, the stationary phase layer in many PLOT columns is not mechanically stable – resulting in particle shedding – which can adversely affect your analysis.

Likewise, if you use Sulfur Chemiluminescence Detection (SCD) to analyze sulfur compounds in natural gas, petroleum products, and process streams, you may have experienced detector sensitivity issues. Column bleed from thick film polydimethylsiloxane (PDMS) phases, which are commonly used, can contribute to “fouling” phenomena on SCD ceramic burner tubes. This leads to frequent and costly tube replacement and associated downtime. In addition, since light sulfur gases are reactive, flow path inertness is critical to avoiding peak tailing.

But the *good news* is... Agilent has developed two new column technologies specifically designed to reduce GC system downtime and maintenance, while expanding your lab's analytical capabilities.






## Reduce maintenance costs and extend your capabilities with Agilent J&W PLOT PT GC columns

New integrated particle-trapping technology on *both ends* of Agilent PLOT PT GC columns minimizes particle shedding in all applications. That means you can worry less about signal spiking – and unnecessary downtime. Our particle-trapping design also takes the worry out of performing GC or GC/MS analysis with PLOT columns for more confident identification of unknowns.

## Greatly improve SCD stability with Agilent J&W DB-Sulfur SCD columns

These new columns are optimized for low bleed, which reduces fouling of SCD ceramic tubes – minimizing instrument downtime and operational costs. In addition, DB-Sulfur SCD columns provide excellent peak shape and extended detector stability for all GC SCD methods that utilize PDMS stationary phases – such as ASTM D 5623 and 5504.



*Integrated particle trap at column ends*

## Inside: our latest column technologies for overcoming productivity barriers

**Page 4**

**Agilent J&W PLOT PT GC columns** for analyzing light gases and volatile organic compounds in hydrocarbon process streams

**Page 6**

**Agilent J&W DB-Sulfur SCD columns** for detecting sulfur compounds in petroleum and petrochemical processing

**Page 7**

**Agilent Gas Clean Filters and Agilent Inert Flow Path solutions** for maintaining gas purity and flow path inertness

To learn more about Agilent J&W PLOT PT GC and Agilent J&W DB-Sulfur SCD columns, visit [agilent.com/chem/energyGC](https://www.agilent.com/chem/energyGC)

# AGILENT J&W PLOT PT GC COLUMNS

## Analyze light gases and VOCs without worrying about particle shedding


Most Porous Layer Open Tubular (PLOT) columns are prone to shedding stationary phase particles, which can flow into the detector and cause the signal to spike. This prevents your system software from accurately identifying and quantifying compounds.

You *could* simply live with the problem... or try to work around it, either by installing a particle trap on the end of the column, or by connecting an inline filter at valve and detector interfaces. However, these approaches can cause problems such as leaks, clogs, method adjustments, decreased system performance, increased downtime, and higher labor costs.


### Now, there's a *better* option: Agilent J&W PLOT PT GC columns

These unique columns virtually eliminate particle shedding through integrated particle-trapping technology on *both ends of the column* – giving you the advantages of:

- **No need for cumbersome (and leaky) unions or downstream filters:** The integrated particle trapping technology is *built directly into the column* as one continuous length of fused silica.
- **Worry-free GC/MS and switch valve applications:** PLOT PT columns protect your switch valves and detectors from PLOT phase particles, so you can perform reliable backflush/heart cut techniques and expand your analytical capability by using capillary flow technology (CFT) and mass spec detection.
- **Decreased system maintenance/downtime:** You won't waste time replacing filters and column switching valves, which can be damaged by particle shedding.
- **Easy upgrade:** No selectivity changes minimize method adjustments and revalidation.



Coal-to-chemical process gas analysis on an Agilent J&W HP-PLOT Q PT GC column. Built-in particle trapping delivers greater stability than conventional PLOT columns, while enabling reliable operation with MS detection and capillary flow technology.



Fixed-gas mix comparison between an Agilent J&W CP-Molsieve 5Å PT column and a standard Molsieve capillary column using TCD. The dual-ended, integrated particle traps eliminated spiking (blue trace).



## ORDERING INFORMATION

ASTM Method	Title	PLOT PT Phase	Dimensions*	Part Number
D 1945	Standard Test Method for the Analysis of Natural Gas by GC	HP-PLOT Q PT	15 m x 0.53 mm x 40 µm	19095P-Q03PT
D 1946	Standard Test Method for the Analysis of Reformed Gas by GC	HP-PLOT Q PT	15 m x 0.53 mm x 40 µm	19095P-Q03PT
D 2163	Standard Test Method for the Analysis of Liquefied Petroleum (LP) Gases and Propene Concentrates by GC	HP-PLOT Al <sub>2</sub> O <sub>3</sub> KCl PT	30 m x 0.53 mm x 15 µm	19095P-K23PT
		HP-PLOT Al <sub>2</sub> O <sub>3</sub> S PT	30 m x 0.53 mm x 15 µm	19095P-S23PT
D 2427	Standard Test Method for Determination of C2 through C5 Hydrocarbons in Gasoline by GC	CP-Al <sub>2</sub> O <sub>3</sub> /KCl PT	50 m x 0.53 mm x 10 µm	CP7518PT
		GS-Alumina PT	30 m x 0.53 mm	115-3532PT
D 2593	Standard Test Method for Butadiene Purity and Hydrocarbon Impurities by GC	CP-Al <sub>2</sub> O <sub>3</sub> /KCl PT	50 m x 0.32 mm x 5 µm	CP7515PT
		CP-Al <sub>2</sub> O <sub>3</sub> /KCl PT	50 m x 0.53 mm x 10 µm	CP7518PT
D 2712	Standard Test Method for Hydrocarbon Traces in Propylene Concentrates by GC	GS-Alumina PT	30 m x 0.53 mm	115-3532PT
		GS-Alumina PT	50 m x 0.53 mm	115-3552PT
D 3271	Standard Practice for Direct Injection of Solvent-Reducible Paints into a Gas Chromatograph for Solvent Analysis	PoraPLOT Q PT	25 m x 0.53 mm x 20 µm	CP7554PT
D 3749	Standard Test Method for Residual Vinyl Chloride Monomer in Poly(Vinyl Chloride) Resins by Gas Chromatographic Headspace Technique	PoraBOND Q PT	10 m x 0.53 mm x 10 µm	CP7353PT
D 3792	Standard Test Method for Water Content of Coatings by Direct Injection Into a Gas Chromatograph	PoraBOND Q PT	25 m x 0.32 mm x 5 µm	CP7351PT
		PoraBOND Q PT	25 m x 0.53 mm x 10 µm	CP7354PT
D 4322	Standard Test Method for Residual Acrylonitrile Monomer Styrene-Acrylonitrile Co-polymers and Nitrile Rubber by Headspace Gas Chromatography	PoraBOND Q PT	25 m x 0.53 mm x 10 µm	CP7354PT
D 4424	Standard Test Method for Butylene Analysis by GC	HP-PLOT Al <sub>2</sub> O <sub>3</sub> S PT	50 m x 0.53 mm x 15 µm	19095P-S25PT
D 4509	Standard Test Methods for Determining the 24-Hour Gas (AIR) Space Acetaldehyde Content of Freshly Blown PET Bottles	PoraBOND Q PT	25 m x 0.32 mm x 5 µm	CP7351PT
		PoraBOND Q PT	25 m x 0.53 mm x 10 µm	CP7354PT
D 4961	Standard Test Method for GC Analysis of Major Organic Impurities in Phenol Produced by the Cumene Process	HP-PLOT Q PT	15 m x 0.53 mm x 40 µm	19095P-Q03PT
D 5303	Standard Test Method for Trace Carbonyl Sulfide in Propylene by GC	HP-PLOT Q PT	30 m x 0.53 mm x 40 µm	19095P-Q04PT
D 5507	Standard Test Method for Determination of Trace Organic Impurities in Monomer Grade Vinyl Chloride by Capillary Column/Multi-dimensional GC	HP-PLOT Q PT	15 m x 0.53 mm x 40 µm	19095P-Q03PT
		HP-PLOT U PT	30 m x 0.53 mm x 20 µm	19095P-U04PT
D 5508	Standard Test Method for Determination of Residual Acrylonitrile Monomer in Styrene-Acrylonitrile Co-polymer Resins and Nitrile-Butadiene Rubber by Headspace Capillary GC	HP-PLOT Q PT	30 m x 0.53 mm x 40 µm	19095P-Q04PT
D 6159	Standard Test Method for Determination of Hydrocarbon Impurities in Ethylene by GC	HP-PLOT Al <sub>2</sub> O <sub>3</sub> KCl PT	50 m x 0.53 mm x 15 µm	19095P-K25PT
		GS-Alumina PT	50 m x 0.53 mm	115-3552PT

### Additional PLOT PT columns not associated with specific ASTM Methods

PLOT PT Phase	Dimensions*	Part Number	PLOT PT Phase	Dimensions*	Part Number
PoraBOND Q PT	25 m x 0.25 mm x 3 µm	CP7348PT	HP-PLOT Al <sub>2</sub> O <sub>3</sub> KCl PT	50 m x 0.32 mm x 8 µm	19091P-K15PT
PoraBOND Q PT	50 m x 0.32 mm x 5 µm	CP7352PT	CP-Al <sub>2</sub> O <sub>3</sub> /Na <sub>2</sub> SO <sub>4</sub> PT	50 m x 0.32 mm x 5 µm	CP7565PT
PoraPLOT Q PT	10 m x 0.32 mm x 10 µm	CP7551PT	CP-Al <sub>2</sub> O <sub>3</sub> /Na <sub>2</sub> SO <sub>4</sub> PT	50 m x 0.53 mm x 10 µm	CP7568PT
PoraPLOT Q PT	25 m x 0.32 mm x 10 µm	CP7551PT	HP-PLOT Al <sub>2</sub> O <sub>3</sub> S PT	25 m x 0.32 mm x 8 µm	19091P-S12PT
PoraPLOT Q-HT PT	25 m x 0.32 mm x 10 µm	CP7557PT	HP-PLOT Al <sub>2</sub> O <sub>3</sub> S PT	50 m x 0.32 mm x 8 µm	19091P-S15PT
GS-Q PT	30 m x 0.53 mm	115-3432PT	HP-PLOT Al <sub>2</sub> O <sub>3</sub> M PT	50 m x 0.53 mm x 15 µm	19095P-M25PT
HP-PLOT Q PT	15 m x 0.32 mm x 20 µm	19091P-Q03PT	CP-Molsieve 5Å PT	30 m x 0.32 mm x 10 µm	CP7534PT
HP-PLOT Q PT	30 m x 0.32 mm x 20 µm	19091P-Q04PT	CP-Molsieve 5Å PT	25 m x 0.32 mm x 30 µm	CP7536PT
PoraPLOT U PT	25 m x 0.53 mm x 20 µm	CP7584PT	CP-Molsieve 5Å PT	25 m x 0.53 mm x 50 µm	CP7538PT
CP-Al <sub>2</sub> O <sub>3</sub> /KCl PT	25 m x 0.53 mm x 10 µm	CP7517PT	CP-Molsieve 5Å PT	50 m x 0.53 mm x 50 µm	CP7539PT

\*PLOT PT columns have 2.5 m of integrated particle traps on both ends that extend the stated length of the column by 5 meters.

Need a different configuration? The Agilent custom column shop can help you determine the right Agilent DB-Sulfur SCD column configuration for your application.

Minimize particle shedding and maximize productivity. Learn more at [agilent.com/chem/energyGC](http://agilent.com/chem/energyGC)

# AGILENT J&W DB-SULFUR SCD COLUMNS FOR PETROLEUM/PETROCHEMICAL PROCESSING




## Improve your SCD signal stability and system performance for sulfur compounds

Thick-film PDMS columns commonly used with SCD are prone to excessive bleed at high temperatures. When this happens, column bleed components accumulate and foul (coke) onto SCD burner ceramic tubes, destabilizing detector response over time. This leads to system downtime — as well as expensive ceramic tube and column replacement.

## Reduce bleed, improve data reliability, and lower maintenance costs with Agilent J&W DB-Sulfur SCD columns

DB-Sulfur SCD columns are engineered for optimal SCD stability and lower bleed — minimizing ceramic tube fouling. Less fouling increases the life of the tube, reducing downtime and operational costs.



Excellent long-term SCD stability. Here, an average of 10 injections of sulfur compound test mix was run over five months. Graph courtesy of Jim Luong, Ronda Gras, and Myron Hawryluk of Dow Chemical Canada

## Good resolution of H<sub>2</sub>S and COS at room temperature

### 7890A GC Conditions

Columns: Agilent J&W DB-Sulfur SCD, 60 m x 0.32 mm, 4.2 μm (p/n G3903-63001)

Carrier: Helium, constant flow mode, 2.8 mL/min

Inlet: 275 °C, Split ratio 10:1 (Sulfiner-treated capillary inlet system)

Oven: 35 °C for 3 min, 35 °C- 250 °C at 10 °C/min, 250 °C for 10 min

### SCD conditions


Burner temperature: 800 °C

Vacuum of burner: 364 torr

Vacuum of reaction cell: 5 torr

H<sub>2</sub>: 40 mL/min

Air: 60 mL/min



This separation produced excellent peak shape for ASTM D 5623 target sulfur compounds, including baseline separation of H<sub>2</sub>S and COS at room temperature and resolution of Thiophene and 2-Methyl-1-propanethiol.

## ORDERING INFORMATION

Part Number	Description	Temperature Limit
G3903-63001	DB-Sulfur SCD 60 m x 0.32 mm x 4.2 μm	-60 to 250/270 °C
G3903-63002	DB-Sulfur SCD 40 m x 0.32 mm x 0.75 μm	-60 to 270/290 °C
G3903-63003	DB-Sulfur SCD 70 m x 0.53 mm x 4.3 μm	-60 to 250/270 °C
G3903-63004	DB-Sulfur SCD 40 m x 0.32 mm x 3.0 μm	-60 to 250/270 °C

Need a different configuration? The Agilent custom column shop can help you determine the right Agilent DB-Sulfur SCD column configuration for your application.



## Protect the purity of your gas and the integrity of your results

### Agilent Gas Clean Filters

Gas contaminants – particularly moisture – can jeopardize sensitivity and cause shifts in retention time, elution time, and elution order. An Agilent Gas Clean Filter reduces oxygen, moisture, and hydrocarbons, giving you better PLOT analysis results. These filters provide:

- **Increased GC/MS sensitivity** for trace-level analysis
- **Instrument and column protection:** The filter changes color when it needs replacing, keeping moisture content <0.1 ppm
- **Easy replacement** *without tools or gas shut-off*



### Agilent Inert Flow Path solutions

Gaseous sulfur compounds pose sampling and analytical problems because they are polar, reactive, and need to be determined at trace levels. An Agilent inert flow path reduces active sites on the metal and glass surfaces that come into contact with the sample – ensuring less compound adsorption and more reliable trace-level results of sulfur components. Other advantages include:

- **Reproducible sample vaporization:** highly inert glass wool in column liners creates an exceptional surface for vaporization and mixing *without losing active analytes*
- **Longer column life,** because nonvolatiles remain inside the inlet
- **Greater sensitivity,** allowing you to run more samples

Agilent Inert Flow Path components include:


- Ultra Inert inlet liners and CrossLab Ultra Inert liners
- Inert Flow Path split/splitless inlet
- Ultra Inert Gold Plated inlet seals
- UltiMetal Plus Flexible Metal ferrules
- Agilent J&W Ultra Inert GC columns
- Gas Clean purifier

To learn more about building an inert flow path, visit [agilent.com/chem/inert](http://agilent.com/chem/inert)



To learn more about preventing untimely burner failure, visit [agilent.com/chem/energyGC](http://agilent.com/chem/energyGC)





## No time for downtime? Count on these rugged Agilent columns for energy and chemical applications

- **Agilent J&W PLOT PT GC columns** feature integrated trapping technology that prevents stationary phase particles from exiting the column. This virtually eliminates detector spikes and helps to prevent damage to switching valves for reduced instrument maintenance.
- **Low-bleed Agilent J&W DB-Sulfur SCD columns** maximize the life of SCD ceramics on your Chemiluminescence Detector – increasing instrument uptime and reducing operational costs.

### For more information

To learn more about  
Agilent J&W PLOT PT GC  
and Agilent J&W DB-Sulfur SCD  
columns, visit  
[agilent.com/chem/energyGC](http://agilent.com/chem/energyGC)

Find an Agilent customer center  
in your country:  
[agilent.com/chem/contactus](http://agilent.com/chem/contactus)

Save time finding the right GC column at  
[agilent.com/chem/selectGC](http://agilent.com/chem/selectGC)

U.S. and Canada  
**1-800-227-9770**  
[agilent\\_inquiries@agilent.com](mailto:agilent_inquiries@agilent.com)

Europe  
[info\\_agilent@agilent.com](mailto:info_agilent@agilent.com)

Asia Pacific  
[inquiry\\_lsca@agilent.com](mailto:inquiry_lsca@agilent.com)

India  
[India-lsca\\_marketing@agilent.com](mailto:India-lsca_marketing@agilent.com)



This information is subject to change without notice.

©Agilent Technologies, Inc. 2013  
Printed in the USA, August 26, 2013  
5991-2977EN



**Agilent Technologies**