

# Chromium Speciation in Drinking Water using LC(IC)-ICP-MS

Mina Tanoshima, Tetsushi Sakai, and Craig Jones, Agilent Technologies

Winter Plasma Conference 2012  
Poster FP38



## Introduction

Chromium is a transition metal which may be present in the environment in various forms depending on the sample type and origin. Hexavalent chromium, Cr(VI) is well known to be highly toxic, while the other stable oxidation state, trivalent chromium, Cr(III) is an essential element for humans. Various regulations around the world exist to minimize the risk of exposure to dangerous levels of Cr(VI). More recently, permissible levels of Cr(VI) in drinking water have been re-evaluated and significantly lower limits have been proposed.

In this study, an isocratic separation of Cr(III) and Cr(VI) using HPLC coupled to ICP-MS is used to quantify ultra-trace levels of both Cr(III) and Cr(VI) in highly mineralized waters in less than 4 minutes. Data is presented to highlight the improved analytical capability for these species in waters and the applicability of the method to the determination of Cr species in food and environmental samples.

## Experimental



Figure 1. Agilent 7700 Series ICP-MS and HPLC

The separation and detection of the two Cr species is important because the total chromium concentration does not provide adequate information on toxicity. The anionic, hexavalent form of Cr is toxic, while in its cationic trivalent oxidation state, chromium(III) is an essential element for human nutrition. Hence separating the two (or species) is necessary before quantifying using ICP-MS as a detector. LC(IC)-ICP-MS enables Cr species to be separated and measured with high accuracy and good sensitivity.

**HPLC conditions:** An Agilent 1200 liquid chromatograph equipped with a binary HPLC pump, autosampler and vacuum degasser were used in this study. The HPLC system was connected to the ICP-MS using the Agilent LC connection kit. An anion exchange column (4.6 mm i.d. x 30 mm polyhydroxymethacrylate base resin) was used for separation. The column temperature was maintained at ambient for all experiments. The Agilent's bio-compatibility kit (Part # 5085-9972) was installed for sample delivery line. The details of the operating conditions are reported in Table 1.




Figure 2. Peak Height and Signal to Noise ratio vs He gas flow

**ICP-MS:** An Agilent 7700x ICP-MS was used for Cr detection. Instrument operating conditions are shown in Table 1. Initially, the removal of interferences on the primary Cr isotope at m/z 52 was evaluated. The use of collision/reaction cell technology with ICP-MS allows Cr to be measured with good accuracy and sensitivity, with removal of the primary matrix-based interferences due to ArO and ClOH. The relationship between He gas flow and analyte signal is shown in Figure 2. From this graph, 4mL/min flow rate was applied in this method. H<sub>2</sub> gas was also evaluated, however it showed no significant improvement compared with He.

HPLC Parameters	ICP-MS Parameters
Column	Agilent anion exchange column, G3268-80001: 4.6mm x 30mm id
Mobile phase	5mM EDTA (2Na) / 5mM NaH <sub>2</sub> PO <sub>4</sub> /15mM Na <sub>2</sub> SO <sub>4</sub> , pH=7.0 adjusted by NaOH
Flow rate	1.2mL/min
Temperature	Ambient
Injection volume	100µL
RF power	1550 W
Sample depth	8 mm
Carrier gas	1.05 L/min
Dwell time	0.5 sec
Isotope monitored	<sup>52</sup> Cr, <sup>51</sup> Cr
Cell gas	He
Flow rate of cell gas	4mL/min

Under the conditions described above, ICP-MS detection using He gas made yielded detection limits (DLs) of < 200 ng/L for both <sup>52</sup>Cr(III) and <sup>52</sup>Cr(VI) with injection volume of 100µL. The detection limits were calculated as three times the peak to peak height to height of noise as measured on standard chromatograms. However, increasing the injection volume should provide better DLs. The DLs with various injections volumes from 5µL-100µL are shown in Table 2.

## Results and Discussion

Table 2. DLs for Cr species to injection volume

Injection volume	Peak Height / counts		Area / counts		DL	
	<sup>52</sup> Cr(III)	<sup>52</sup> Cr(VI)	<sup>52</sup> Cr(III)	<sup>52</sup> Cr(VI)	<sup>52</sup> Cr(III)	<sup>52</sup> Cr(VI)
5 µL	32620.0	24233.0	514586	503774	<b>1.88</b>	<b>2.83</b>
20 µL	130764	97833.0	2101007	2007572	<b>0.72</b>	<b>0.96</b>
50 µL	323592	241948	5154324	4927731.2	<b>0.29</b>	<b>0.37</b>
100 µL	63267.0	475244	1E+07	9786463	<b>0.13</b>	<b>0.17</b>




Figure 3. 8 hours stability test using 5µg/L standard solution (n=100)

The long-term stability was evaluated using a Cr standard solution. Figure 3 shows good reproducibility over 8 hours (n=100) using 5µg/L standard solution. The RSDs for both <sup>52</sup>Cr(III)-EDTA and <sup>52</sup>Cr(VI) were less than 5%.

### Drinking Water Analysis

In order to test the suitability of the method for these real-world sample types, the method was applied to the determination of both Cr species in both spiked and unspiked mineral water samples. The three samples evaluated were a Japanese mineral water referred as Water A, and two French mineral waters referred as Water B and Water C. The drinking waters selected covered a wide range of typical mineral water compositions, including Water C which is at the extreme end of highly mineralized drinking water (over 450ppm Ca and over 1000ppm sulfate). The major element composition of the water samples is shown in Table 4.

Table 4. The major composition of the water samples

	Water A (ppm)	Water B (ppm)	Water C (ppm)
Na	4.5	11.6	3.4
Ca	9.7	11.5	469
Mg	1.5	8	74.5
K	2.8	0.2	2.8
Sulfate			1121




Figure 4. Chromatogram for spiked (blue) and unspiked (black) mineral water A

Table 5. Results of 10 µL mineral water stability test (8 hours, n=30/each water)

Sample	<sup>52</sup> Cr(III)-EDTA			<sup>52</sup> Cr(VI)		
	AVG	Area	Conc. (ppb)	AVG	Area	Conc. (ppb)
Water A	12879	196130	0.148	18745	103	0.211
Water B	850	4.4	4.4	2.1	2.1	2.1
Water C	103560	16.7	16.7	900154	10.1	10.1
Spiked A	8968	8968	0.103	21331	0.240	0.240
Spiked B	9	9	1.0	2.3	2.3	2.3
Spiked C	90072	10.1	878234	9.9		
STD	6868	6868	0.078	12499	0.140	0.140
Sample C RSD (%)	9.8	9.8	1.4	1.4		

### Applying bio-inert LC pump

The possibility to lower the background using Agilent 1260 Infinity Bio-inert pump was also evaluated. The metal-free components in the sample flow-path and the absence of iron and steel in solvent delivery was expected to help achieve lower detection limit when quantifying Chromium. Mobile phase and 1% HNO<sub>3</sub> (EL grade) and 0.1% HCl were quantified after passing through the bio-inert pump and also conventional pump. Ultra Pure Water (Millipore), 5mM EDTA, 5mM NaHP0<sub>4</sub>, and Na<sub>2</sub>SO<sub>4</sub> were also quantified for the reference of the background. The results are shown in Table 6.

Table 6. Results of Cr concentration and background

Mobile phase	52 Cr		52 Cr		52 Cr	
	Conc. (ppb)	CPS	Conc. (ppb)	CPS	Conc. (ppb)	CPS
1% HNO <sub>3</sub>	0.40	0.781	0.31	4658.5	<0.001	142
0.1% HCl	0.61	1163.4	2.70	3983.5	0.04	766
Ultra Pure Water	0.51	753.8	0.43	6405.8	0.23	312.5
5mM Na <sub>2</sub> SO <sub>4</sub>					<0.001	10.0

When comparing between conventional and bio-inert pump, the solution which passed through bio-inert pump showed slightly lower background. When quantified each component of the solution, the background of 5mM NaHP0<sub>4</sub> showed relatively high. In order to achieve lower background and lower detection limit, applying higher purified NaHP0<sub>4</sub> might help.

## Results and Discussion

### Quantification of Cr(VI) at ultra-trace levels

While we have developed the method to measure both Cr(III) and Cr(VI) in drinking water simultaneously, the State of California in the US has recently proposed (2009) setting a new "public health goal" for Cr(VI) in drinking water of 60ng/L. To meet this goal, another method for higher sensitivity and selectivity for Cr(VI) quantification was developed. The method is outlined in Table 7. The same column was used for separation but the injection volume was increased to 1500µL. In this method, Cr(III) and Cr(VI) were not quantified because with large injection volume, the water dip will cover the peak of Cr(III) and also in order to quantify Cr(III), EDTA will be added to the sample solution and this EDTA, anion, might help excess of the capacity of the column, causing peak shape change or retention time delay. From these reasons, only Cr(VI) was focused to measure.

Table 7. Operating parameters of ICP-MS and HPLC

HPLC Parameters	ICP-MS Parameters
Column	Agilent anion exchange column, G3268-80001: 4.6mm x 30mm id
Mobile phase	1mM EDTA (2Na) / 5mM NaH <sub>2</sub> PO <sub>4</sub> /15mM Na <sub>2</sub> SO <sub>4</sub> , pH=7.0 adjusted by NaOH
Flow rate	1.2mL/min
Temperature	Ambient
Injection volume	1500µL
RF power	1550 W
Sample depth	8 mm
Carrier gas	1.05 L/min
Dwell time	0.5 sec
Isotope monitored	<sup>52</sup> Cr, <sup>51</sup> Cr
Cell gas	He
Flow rate of cell gas	4mL/min

By applying a larger injection volume (1500 µL), DLs in the region of ng/L for Cr(VI) became achievable. Although other high concentration anions exist in drinking water besides Cr(VI), no peak shape change nor retention time delay occurred. The calibration linearity showed good correction coefficient which was better than 0.9995 in the low level Cr(VI) (1000-50ng/L). Detection limit calculated by S/N was about 8ng/L.




Figure 5. Chromatograms of Cr(VI) for standard solution and various waters from the State of California, USA. The left chromatogram shows 50ng/L standard solution and the right shows non-spiked chromatographs of three water samples.

Table 8 shows the spike recovery test for tap water from two different places and river water sample from California. Chromatograms are shown in Figure 5. While the amount of Cr(VI) for River Water A was calculated below this published health limits, the results for both Tap Water A and B are above the proposed goal.

Table 8. Results of 50ng/L spike recovery test for three different water samples

Tap Water A	Recovery		Tap Water B	Recovery		River Water A	Recovery	
	Non spiked	Spiked (%)		Non spiked	Spiked (%)		Non spiked	Spiked (%)
1	0.3440	0.622	93.6	0.1202	0.618	91.2	0.2011	0.5291
2	0.1772	0.6470	93.28	0.1281	0.6222	99.62	0.5423	0.5282
Average	0.1806	0.6403	93.53	0.1242	0.6210	93.33	0.5417	0.5256

## Conclusions

Cr(III)-EDTA and Cr(VI) were successfully separated and quantified using LC(IC)-ICP-MS in natural, high matrix water samples:

- with good detection limits for both species
- with good long stability (8 hours)
- with good reproducibility within different columns

Cr(VI) was successfully quantified at much lower levels, using high volume injection, achieving 8ng/L of DL. The applying the bio-inert LC pump may expect the better DL with lower background. Also applying higher purified NaHP0<sub>4</sub> might help.