


## Quaternary Surfactants on a Thermo Scientific™ Acclaim™ Surfactant Column with Suppressed Conductivity Detection



Column: Thermo Scientific™ Acclaim™  
Surfactant 5  $\mu$ m, 4.6  $\times$  150 mm  
Pump: Thermo Scientific™ Dionex™  
ICS-3000 DP  
Mobile Phases: (A) Water  
(B) 100 mM formic acid  
(C) 70:30 acetonitrile:water (v/v)  
Flow: 1.0 mL/min  
Gradient Times:     -5    0    7   12  
                  %A   65   65   30   30  
                  %B   15   15   15   15  
                  %C   20   20   55   55  
Temperature: Dionex ICS-3000 DC at 30 °C  
Injection: Dionex ICS-3000 AS autosampler; 15  $\mu$ L  
Detector: Suppressed conductivity  
Suppressor: Thermo Scientific™ Dionex™ CSRS™  
ULTRA II 4 mm, external water mode,  
1.0 mL/min, 44 mA

Peaks: (100  $\mu$ g/mL each)  
1. Dodecylpyridinium chloride  
2. Dodecyl-benzyl-dimethylammonium chloride  
3. (Octylphenoxy)ethoxyethyl-benzyl-  
dimethylammonium chloride  
4. Hexadecyl-trimethylammonium bromide  
5. Hexadecylpyridinium chloride

23629

Quaternary surfactants are used as microbiocides and preservatives in a wide variety of applications, from industrial to pharmaceutical. Ordinary HPLC columns don't cleanly resolve these compounds, producing peaks with severe tailing. The Acclaim Surfactant column resolves this class of surfactants well, with good peak shapes. Suppressed conductivity can detect as little as 15 ng of surfactant, even for substances with weak UV absorbtion.