

HPLC Cloud Software Set-up Manual

Manufactured by Newcrom Co. division of SIELC Technologies, Inc.

E: mail@sielc.com Ph. +1 (847) 229-2629

Introduction

HPLC.Cloud is the very first and only HPLC software that is accessible via the cloud. SIELC Technologies' Alltesta™ is an innovative instrument with a cutting-edge approach to instrument control, data collection, data manipulation, and data storage.

Each unit of the Alltesta™ is directly connected to a secure cloud server via a local Ethernet line or wireless network.

No installation or support is needed to gain access to the software as HPLC Cloud is a web-based application where only a username and password are needed. Simply sign up online to get started.

SIELC's **Power Tower** module provides the connection to the cloud. It is a chromatography data acquisition device that connects all the units on your account to the software. This device allows you to access and control your laboratory from anywhere in the world, at any time, on any device with only an internet connection needed.

DEVICES

(This set-up can be done via any device that has an internet connection: phone, tablet, laptop, etc.)

1. In order to set up and connect your devices to the cloud, start by creating an account on the website HPLC.Cloud. Click “**Sign Up**” as shown below:




Cloud-based Software for Chromatography

A free web app to run laboratory tests, control injections, record chromatograms, integrate peaks and estimate compounds.




Chromatography - A technique of separating a mixture into individual components. Can be made highly selective with variation in both mobile and stationary phases.

Cromite and Stream LC supported.




2. By clicking the **Sign Up** button you will be directed to the Sign-Up page:



The Sign Up page features the HPLC Cloud logo at the top. Below the logo are two links: "Sign Up" (with a user icon) and "Log In". There are two input fields: one for email (with an @ icon) and one for password (with a key icon and the text "password, 8 symbols minimum"). Below the input fields is a disclaimer: "By clicking Sign Up, you agree to our Terms of service, Cookie Policy, Data Policy." A large blue "Sign Up" button is positioned below the disclaimer. At the bottom, there is a "Log In with Google" button featuring the Google logo.

3. Enter in the required information, click **Sign Up** and your account set-up is complete.
4. You will receive an email shortly after confirming that your email has been registered with the software.
5. After that is complete, login to your account and you will be taken to a screen that should show this:



6. To add your device(s) to the software you will need to click the **Add Device** button shown above.

7. You are now going to be shown the screen below, which will ask for the serial number and key of your device. These can be found on the back of each power tower.

Serial


Key

Add device

Scan

Cancel

- **Note:** All units have been tested by SIELC. A message might pop up saying that “these units belong to sales@sielc.com.” Just confirm that you want to delete them from that account and transfer them to your own account.
8. Repeat steps 6-7 to add each Power Tower that you have to your HPLC.Cloud account.
 9. Once you have added each device to your account you will be able to see them all by clicking the **Devices** tab at the top of the screen.



10. Now you can select the system that you want to use by clicking on the StreamLC/Tower (you can rename each device by clicking on the serial number once you have the specific one opened):

- If your device is not plugged in, it will appear **Offline**.
- If it is plugged in, it should appear **Online**.



11. When you click on your device, it may look like the image below.

12. First, check that your firmware is up to date by clicking on the text starting with **fv**. If there is a new firmware version, click **Update Now**. Alternatively, click **Remind me later** to set an update reminder.



13. Once it is updated, click **OK**.

14. To connect equipment (Autosampler, Pump, Detector, etc.), click on “Settings”.

Example Alltesta

Stream LC T70I-A3LF fv1.106 • Online



Settings



No connected equipment

15. Select the devices from the appropriate drop-down menus (red arrows). The options are:

- Auto: This will automatically select the device.
- H(P/D/A)000000000: The letter in the 2nd position depends on the device (Pump, Detector, Autosampler).

16. Input the Volume of your Mobile Phase (A) and Waste Containers (blue arrows). A standard Alltesta container is 250 mL.

Settings

Pumps

Pump

Auto

Pump B (gradient)

Off

Injection

☐ Combine injection and fraction collector

Autosampler

Auto

Tray

auto

Valve loop

Off

Detection

Detector

Auto

☐ LED saving mode

☐ Analog signal

Phase

Volume of A container, mL

←

Add B container

←

Waste

Volume of waste container, mL


←

Save

17. If you are using a gradient system (i.e. a system with 2 pumps), select the 2nd pump from the Pump B drop-down menu and click “Add B Container” to enter the volume of the Mobile Phase (B) container (green arrows). Using the gradient system will be discussed in the ***Inject*** section.

18. Once you add all the component devices, you will be able to see and control each device (pump, detector, valve, etc.) that is associated with that respective Alltesta.

- **Note:** if your Alltesta™ Mini/standalone devices are not connected to the Power Tower, no devices will be shown.



19. In order to control the Detector manually (without a test/sample) there are 2 buttons highlighted above in blue: **Turn on LEDs** and **Zero**.


- **Turn on LEDs** refers to the 4 different colored diodes that the detector has. Each diode emits light with the maximum wavelength indicated. The image above shows a real-time plot of a signal of all four available wavelengths.
- **Zero** allows you to manually “zero” all four signals at the same time. **Note:** The system will always zero the detector automatically before an injection.

20. In order to control the Pump manually, you can click **Start** or **Purge**, highlighted in black above.

- **Start** will run the pump with the current set flow rate. This is best used to quickly replace the mobile phase in the pump.
- **Purge** will run the pump with a 4.00 mL/min flow rate. **This button should be used only after the purge valve has been opened on the pump.** This is used to remove any air that may have entered the pump barrel. If the valve is not opened prior to purging the system, the pressure will build up and some of the tubing connections may pop out of their connections. This can be fixed by refastening the tubing and fittings.
- **Note:** The **Stop** button will replace the **Start** button once the pump starts.


21. The two bottles highlighted in orange correspond to the Mobile Phase A and Waste Container bottles and will show track the current volume in each bottle. If you added a bottle for Mobile

Phase B in the settings for a gradient system, then a 3rd bottle labeled “B” will show as well. You can update the volume for each bottle by clicking on the respective bottle, typing in the updated value, and then clicking **Set**.




SAMPLES

1. In order to set-up a sample to begin your testing, you'll need to click on the **Calibrations** tab at the top of the screen.



2. Click on the **New** button to create your first sample:



3. This screen below will now be prompted in front of you where you will be asked to detail your sample. Fill in everything related to your sample and click **Create**.
 - Select the device and method you will be using for this sample.
 - We will explain in the “Injections” section how to create a “Method.”
 - In this example, we’ll be testing for caffeine in a sample of instant coffee. We’ll use the device “Example Alltesta” and a method using a mobile phase of 10% ethanol, 90% water, and a buffer solution of H_3PO_4 at .3%.

Name

Coffee

Description, optional

Instant Coffee

Expected compounds, separate with semicolon, optional

Caffeine;

Device, optional

Example Alltesta


Method

EtOH/H₂O/H₃PO₄ - 10/90/.3%

Use Ctrl+Enter to create

Create


4. Once you have filled out everything related to your sample, it will now be shown under the “**Sample**” tab of the application.



The screenshot shows the HPLC Cloud application interface. At the top, there is a navigation bar with the HPLC Cloud logo and several tabs: Devices, Samples (which is highlighted with a blue underline), Calibrations, Sequences, and Account. Below the navigation bar, the main heading is 'Samples'. To the right of the heading are three buttons: 'View Compare', 'View Archive', and a '+ New' button. Below these are three filter boxes: 'Keywords', 'Any device', and 'Any time'. A table displays the sample data with columns: Name, Sequence, Device, and Time. The table shows one entry for 'Coffee' on 'Friday, Mar 3', marked as 'New', with the device 'Example Alltesta' and the time '12:05 PM'.

Name	Sequence	Device	Time
Coffee		Example Alltesta	12:05 PM

5. Repeat steps 1-4 to add as many samples as you'd like. Samples can be added at any time and need not be done all at once.




This screenshot shows the HPLC Cloud application interface after adding more samples. The layout is identical to the previous screenshot, but the table now contains three entries: 'Energy Drink', 'Tea', and 'Coffee'. All three entries are marked as 'New' and were recorded on 'Friday, Mar 3' at '12:05 PM' using the 'Example Alltesta' device.

Name	Sequence	Device	Time
Energy Drink		Example Alltesta	12:05 PM
Tea		Example Alltesta	12:05 PM
Coffee		Example Alltesta	12:05 PM

COLUMNS

1. Click on the **Account** tab:



2. You will be automatically directed to the method presets section, so click on **My Columns** in order to get started:
3. Once here, click on **New** to set up your column.
4. You will now be prompted with a screen that will ask for all details about your column. There are 3 ways to input the column information.
 - a. **Manual Entry** : Manually enter the information found on the column label into the appropriate boxes
 - b. **Autofill via Serial:** Type in the serial number in the “Serial” box and click the arrow or hit “Enter” on your keyboard. This should autofill all of the column information into the appropriate boxes.
 - c. **Autofill via QR Code:** If you have a camera attached to your device, click the camera button and show the QR code on the column packaging.

Once the information is all filled out, click **Create**. **Note:** Make sure to input a max PSI of 2000.

Create Column

Column name

Column

Type

Length, mm

Serial

nr1rce7a8



Internal diameter, mm

Lot

Particle size, μm


Max pressure, psi

Pore size, \AA


Create

CALIBRATIONS

1. Calibrations are used to calculate the concentration of a compound in an unknown sample.
2. In order to set-up a Calibration click on the **Calibrations** tab at the top of the screen.



3. Click **New** to create a new Calibration.



4. Input the compound name you are setting the Calibration for and select the Device, Method, and Column, then click **Create**. In this example, we'll be using a Caffeine standard with a concentration of 1 mg/mL.

★

Compound name

Caffeine

≡

Description, optional

1 mg/mL

📄

Device

Example Alltesta

🔧

Method

EtOH/H₂O/H₃PO₄ - 10/90/.3%

🔍

Column


Primesep SB_Example

✓


Create

Use Ctrl+Enter to create

- Once the Calibration is created, it will appear under the Calibrations tab.




- Click on the Calibration's name (ex. **Caffeine** to open its page, seen below).




- To start a new Calibration Injection, click **New Injection**.

8. Input the Vial location and Concentration (in mg/mL) information, then click **Start**.
 - **Note:** The Device, Method, and Column information should automatically populate since we already input this info when creating the Calibration.


Injections
+ New Injection


1


Setup new injection



Description




A1



1

Device and configuration


Device



Example Alltesta


• Online

Method




EtOH/H2O/H3PO4 - 10/90/.3%

Column




Primesep SB_Example


 Start




9. After the scan is complete, select the peak of the analyte that you want to use for the calibration. To do this, hover over the row of the desired peak and click **Select**.



10. To see the Calibration, click "Injections" on the top left of the graph.




11. You should see 2 graphs: One graph (left) showing the linear regression used to calculate the calibration (peak area vs. amount) and the other graph (right) overlaying peaks of the various standards used for that Calibration (see below where multiple standards were used).




12. We will explain later in the “Sequences” and “Injections” sections how to utilize the Calibration to calculate the unknown concentration of the compound of interest in a Sample.

SEQUENCES

1. Sequences allow users to automate and group sample testing. While optimized for use with an Autosampler, it can be utilized with manual injections.
2. To create a sequence, click the **Sequences** tab at the top of the screen.



3. Click **New**.



4. Input all relevant information and then click **Create**.

Name

Caffeine Test

Description

Testing for Caffeine in Coffee, Tea, and an Energy Drink

Expected compounds, separate with semicolon

Caffeine;

Device

Example Alltesta


Method

EtOH/H₂O/H₃PO₄ - 10/90/.3%

Use Ctrl+Enter to create


Create

- Created sequences will then appear under the Sequences tab as shown below.




The screenshot shows the HPLC Cloud interface. At the top is a navigation bar with icons for HPLC Cloud, Devices, Samples, Calibrations, Sequences (selected), and Account. Below the navigation bar is the 'Sequences' section. It features a header with 'View Compare', 'View Archive', and a '+ New' button. There are three filter boxes: 'Keywords', 'Any device', and 'Any time'. Below these is a table with columns 'Name', 'Device', and 'Time'. A date separator 'Monday, Mar 6' is present. The table contains one entry: 'Caffeine Test' (marked as 'New') on 'Example Alltesta' at '10:22 AM'.

- Click on the newly created Sequence (in this example, **Caffeine Test**). This will open up the Sequence window.




The screenshot shows the 'Caffeine Test' sequence window. At the top is the title 'Caffeine Test' with a menu icon. Below it is a subtitle 'Sequence 3/6/2023 Example Alltesta ! New'. The 'Description' section reads 'Testing for Caffeine in Coffee, Tea, and an Energy Drink'. The 'Compounds' section has a search bar and a '+ Add' button. The 'Methods' section shows 'EtOH/H2O/H3PO4 - 10/90/.3%' and a '+ New' button. On the right is a 'Start' button. At the bottom is a 'Sequence Table' section with a checked 'Only Binded' checkbox. Below this is a 'Samples' section with a '+ Add sample or calibration' button, which is highlighted by a red arrow.


7. To add Samples and Calibrations to your Sequence, Click the **Add Sample or Calibration** button as shown above. Type in the name of the Sample. To make this entry a regular sample, click **Sample**; to make it a calibration, click **Calibration**. Repeat for each sample or calibration you want to include in your sequence.




8. To edit a particular Sample or Calibration, click on the sample name in the list.
9. To add an injection to the sequence, click **Add Instruction** for each Injection you want to do.



10. To select a Sample/Calibration for an Instruction, click on the Name of the Sample for the respective instruction and a drop-down list will appear of all available Samples and Calibrations.



11. Choose the Sample or Calibration you want to set for that instruction.
12. The Method you input when creating the Sequence should automatically be selected for each Instruction.
13. To edit the current method, click on the method name. The standard menu or method editing will appear.



14. At the bottom are 4 options:
 - a. **Save as Preset:** Save the current settings as a preset method.
 - b. **Delete:** Delete the current method.
 - c. **Save as New:** Save as a new method within the sequence.
 - d. **Save:** Save and update the settings of the current method. **Note:** This option can only be used if this particular method has not yet been used in an injection. Once an injection uses this method, the only way to edit the method will be to **Save as New**.



15. If you want to add additional methods, click **+ New**.

Caffeine Test

Sequence 3/6/2023 [Example Alltesta](#) ! New

Description

Testing for Caffeine in Coffee, Tea, and an Energy Drink

Methods

EtOH/H2O/H3PO4 - 10/90/.3%;

+ New

16. Select the Method for each Sample as you prefer. If you have more than one method in this sequence, you will be able to select the desired method from a drop down menu.

Caffeine Test

Sequence 3/6/2023 [Example Alltesta](#) ! New

Description

Testing for Caffeine in Coffee, Tea, and an Energy Drink

Methods

EtOH/H2O/H3PO4 - 10/90/.3%;

+ New

Sequence Table

☒ Only Binded

Start

Sample name	Method	Vial	Repeat	Run
★ Caffeine	EtOH/H2O/H3PO4 - 10/...	A1	1	<input checked="" type="checkbox"/>
☉ Coffee	EtOH/H2O/H3PO4 - 10/...	A1	1	<input checked="" type="checkbox"/>
☉ Tea	EtOH/H2O/H3PO4 - 10/...	A1	1	<input checked="" type="checkbox"/>
☉ Energy Drink	EtOH/H2O/H3PO4 - 10/...	A1	1	<input checked="" type="checkbox"/>

+ Add instruction

17. For each sample, select the vial location for the Autosampler to test each sample by clicking on the vial location on the map that pops up. Alternatively, you can type in the vial location.

18. To set the number of runs for that injection, type in the desired number in the Repeat column.
Note: Entering 1 will correspond with 1 run for that instruction.
19. To prevent a particular instruction from running during the sequence, uncheck the blue box in the Run.
20. To rearrange the order of the instructions, hover the cursor over the far left side of the Instructions until it is over the dots. Click and drag to the desired location and then let go of the mouse to release the instruction at its new location.
21. Verify that all of the Instructions you want to run in this sequence are correct and in the correct order.
22. Click **Start** to run the Sequence.

Sample name	Method	Vial	Repeat	Run
★ Caffeine	EtOH/H2O/H3PO4 - 10/...	A1	1	<input checked="" type="checkbox"/>
☉ Coffee	EtOH/H2O/H3PO4 - 10/...	A1	1	<input checked="" type="checkbox"/>
☉ Tea	EtOH/H2O/H3PO4 - 10/...	A1	1	<input checked="" type="checkbox"/>
☉ Energy Drink	EtOH/H2O/H3PO4 - 10/...	A1	1	<input checked="" type="checkbox"/>

23. To view the Injection results for a particular sample, click on the name of the desired Sample in the Sample list. Alternatively, click on the Eye (which has now replaced the blue box) in the Run column to view the results of that particular Instruction. If multiple injections were performed during a particular instruction, you can cycle through them by clicking on the corresponding number on the Sample page.

Injections	1	2
------------	---	---

24. To pause a Sequence, click **Pause**. This will pause the sequence in its current Instruction. If the sequence is resumed, it will repeat the current instruction.
- For example, in the Sequence shown below, if the Sequence was paused during the Caffeine instruction, resuming a Paused Sequence will repeat Injection 1 of the Caffeine instruction. If 3 injections were called for this particular instruction and the Sequence was paused during Injection 2, the sequence will resume the Sequence by repeating Injection 2.
 - Pausing and then Resuming a Sequence will result in an extra injection for that instruction.
 - Note:** This will not stop/pause the current injection.
25. To Stop a Sequence, click **Stop**. If the Sequence is resumed, it will start from the next instruction.
- For example, in the Sequence shown below, if the Sequence was stopped during the Caffeine instruction, the Sequence will resume with the Coffee instruction. If 3 injections were called for the Caffeine instruction and the Sequence was stopped during Injection 1, the sequence will skip the remaining 2 injections and start with the Coffee instruction.
 -

Caffeine Test



Sequence 3/6/2023 Example Alltesta ▶ **Executing**

Description

Testing for Caffeine in Coffee, Tea, and an Energy Drink

Methods

[EtOH/H2O/H3PO4 - 10/90/.3%;](#) New

Sequence Table



Sample name	Method	Vial	Repeat	Run
★ Caffeine	EtOH/H2O/H3PO4 - 10/90/...	A1	▶ 1 of 1	
☕ Coffee	EtOH/H2O/H3PO4 - 10/...	A1	1	

26. To stop the current injection, open the currently running Sample (either by clicking on the Sample name in the Sample list or by clicking the Eye in the Run column) and click **Abort**. If the Sequence was not already Stopped, then aborting will Pause the Sequence.

27. To resume the Sequence, click “Resume.”

Caffeine Test

Sequence 8/30/2021 [Example Alltests](#) **Paused**

Description
Testing for Caffeine in Coffee, Tea, and an Energy Drink

Methods
[EtOH/H2O/H3PO4 - 10/90/.3%](#) [+ New](#)

Sequence Table

Resume **Stop**

Name	Amount	Method	Vial	Repeat	Run
Coffee		EtOH/H2O/H3PO4 - 10/9...	F1	1	👁
Tea		EtOH/H2O/H3PO4 - 10/9...	F8	1	👁
Energy Drink		EtOH/H2O/H3PO4 - 10/9...	A8	1	👁
★ Caffeine	100 %	EtOH/H2O/H3PO4 - 10/9...	A8	1 of 1	👁
+ Add instruction					

INJECTIONS

Autosampler Injection

1. For users with an Autosampler, continue below. For users manually injecting their samples, skip to “*Manual Injections.*”
2. When you’re ready to begin an injection using the Autosampler, first click on the **desired sample** you would like to use (ex. **Coffee**):

Samples [View Compare](#) [View Archive](#) [+ New](#)

Keywords Any device Any time

Name	Sequence	Device	Time
Monday, Mar 6			
Coffee		! New	12:30 PM

3. Once the sample is opened, click on **New Injection** to begin:


Coffee [...](#)

Sample 3/6/2023 ! New < Share

Injections [+ New Injection](#)

No injections for this sample. Click **New Injection** to create an injection.

4. You will now be prompted with this screen that will ask you for all information related to the injection. **Select the device, method, column, and vial.** Once you have filled in all the info, click **Start**.
 - **Note:** Make sure you have the correct device selected on the first line before proceeding. You will not be able to continue without selecting a device.
 - **Note:** You have the option to add in a brief description of the sample as well.



Coffee

⋮

Sample 3/6/2023
! New
↩ Share

Injections
[+ New Injection](#)

~ 1

Setup new injection

≡

Description

⊕

Vial

Device and configuration

🔧

Device

• Online

▼

📄

Method


📏

Column

▶ Start

- Before you begin the injection, you will need to create the Method. Click on the **Method name** (above). Below is all the information that needs to be filled out. **Note:** These entries below are for this particular method/column/sample and can be modified for each experiment's individual needs.

Method

 Search and use preset

Method name

EtOH/H2O/H3PO4 - 10/90/.3%

- General
- Pump
- Autosampler

🕒 General

Run time, min

6.5

Comment

- Processing
- ⊗ Report
- ⊗ Peaks

➔ Pump

A phase

MeCN/H2O - 50/50

B phase

Flow rate, mL/min

0.5

Refill flow rate, mL/min

8

Pressure min, psi

100

Pressure drop time, s

30

Detector

Channel mode

All

Default channel for integration

UV

Injection

Sample amount, μL

5

Injection time, s

360

Injection delay, s

Needle depth, mm

Intake rate, $\mu\text{L}/\text{min}$

100

Shaking mode before injection

Shaking duration, s

Wash after injection

OFF

Processing

☒ Auto baseline correction

☒ Auto peaks

Saving **method** doesn't change **method presets** that you use.

Save as preset


Revert

Save


- To access a pre-set method, click on the **search bar** (above) at the top of the screen and select the desired method.

7. When you select a method preset, it will fill in all the respective fields. Some fields do not require an entry, including Pressure min, Pressure drop time, Injection delay, Needle depth, Intake rate, Shaking mode, and Shaking duration. If left blank, some of these fields will run on default settings (like Pressure min, Pressure drop, Needle depth, and Intake rate), and others aren't required at all. Fill in these fields and click **Save** to finish and use this method for your injection. Alternatively, click **Save as Preset** to save this method to your list of Preset methods.
 - The **Flow Rate** correlates to the internal pressure (higher flow rate = higher pressure, etc.).
 - Enabling **Save Pressure** will allow you to track the pump pressure throughout the test.
 - A longer **Injection Time** can ensure all of the collected **Sample Amount** is injected as well as helping to minimize noise in the data. You can set the Injection time (in s) to equal the **Run Time** to delay the valve switching until the end of the run.
 - **Injection Delay** can be useful for ensuring the detector reaches a proper baseline before the sample is introduced.
 - **Intake Rate** sets how fast the Autosampler will collect the sample. The default setting is 10,000 $\mu\text{L}/\text{min}$. Generally speaking, slower speeds will produce more accurate results. See the Index section "Autosampler" for recommended rates based on sample size.
 - Toggling on **Wash After Injection** is generally recommended
 - For most applications, you can leave **Needle Depth**, **Intake Rate**, **Shaking Mode before Injection**, and **Shaking Duration** blank.

8. Click **Start** and the chromatogram will automatically present itself on the screen.
 - **Note:** When injecting your sample, volumes between 1-10 μL (0.001-0.01 mL) should be enough for a quality result. Larger volumes injected may cause the sample to be too concentrated, which, in turn, may result in overloading the detector.




Coffee




Sample 3/6/2023
! New
Share


Injections
+ New Injection

 1

Setup new injection




Description



A1

Device and configuration


Device



Example Alltesta


• Online

Method





EtOH/H2O/H3PO4 - 10/90/.3%

Column




Primesep SB

 Start



9. The image below shows the window that will appear when the scan starts. The shaded bars labeled **A** and **W** (and **B** in gradient systems) correspond to the Mobile Phase volumes depicted on the device's home page. After a few seconds the page will update with the Pressure graph and UV/Vis detector data.



Coffee

⋮

Sample 3/6/2023 [Example Alltesta](#) ⌂ Running [Share](#)

Injections 1 + New Injection

~ 1 Stop


Sample injection 3/6/2023 2:55 PM ⌂ Preparing

Device	Example Alltesta	Vial A1
Method	EtOH/H2O/H3PO4 - 10/90/.3% ⌂ 🔗	
Column	Primesep SB 🔗	

A W

⌂ Injecting

10. The image below shows a completed scan of a Coffee sample using an Autosampler.



11. Now we will calculate the amount of the compound in the sample using the Calibration created earlier (in this case the amount of Caffeine in the Coffee sample).

12. Move your mouse over the peak's row and click the **Blue Star** that appears.



13. A popup will appear allowing you to select a Calibration from a drop-down menu.

Use calibration

Peak at 2.699 min
Area 12.152 mAU×min

Calibration

Estimate

14. Click “Estimate” after you have selected the Calibration you want to use.

15. The graph will now show the calculated amount of the compound in question for this particular injection based on the Calibration.


#	Time	Area	Compound	Amount
1	2.699	12.152	→ Caffeine	0.280989 mg/mL

16. You can also run an injection using the “Quick Start” method.

Example Alltesta

Stream LC T70I-A3LF fv1.239 • Online

Settings




Method EtOH/H₂O/H₃PO₄ - 10/90/.3%

Run time	10 min	Injection	5.0 µl in 360 sec
A phase	MeCN/H ₂ O - 50/50	Sample intake	at 100 µl/min
B phase		Flow rate	0.50 mL/min
Pressure	min. 100 psi reach in 30 sec		


Column Primesep SB

Serial	4UB	Length × ID	100 × 2.1 mm
Type	SB-21.100.0510	Particle, pore size	5 µm, 100 Å

17. First, ensure the Method and Column selected for your device are correct.
18. Then click **Quick Start** and select the Vial you want to test. The system will then start the Injection.



19. To see the chromatograph, click on the **Sample title** that appears at the top of the page.



20. As shown above, one can stop the Injection by clicking the **Stop** button.

21. If you are viewing the chromatograph, you can also stop the Injection by clicking **Stop**.



Sample 12964

⋮

Sample 3/6/2023 [Example Alltesta](#) Running [Share](#)

Injections 1 + New Injection


1 


Sample injection 3/6/2023 3:39 PM Running

Device	Example Alltesta	Vial A1
Method	EtOH/H2O/H3PO4 - 10/90/.3%	 
Column	Primesep SB	

A W

Data


Pressure psi

...


Stop


Manual Injection

1. For users manually injecting their samples, continue below.
2. When you're ready to begin a manual injection, first click on the **desired sample** you would like to use (ex. **Coffee**):



Name	Device	Time
Monday, Aug 30 Coffee	! New Example Alltesta	4:26 PM

3. Once the sample is opened, click on **New Injection** in order to begin:



← Samples

Coffee

Sample 8/30/2021 Example Alltesta ! New < Share


Description
instant Coffee

Compounds
Caffeine

Injections + New Injection

No injections for this sample. Click **New Injection** to create an injection.

4. You will now be prompted with this screen that will ask you for all information related to the injection. ****You will need to choose a method and a column.****
 - **Note:** Make sure you have the correct device selected on the first line before proceeding. You will not be able to begin your injection without doing so.
 - **Note:** You have the option to add in a brief description of the sample as well.



Coffee

Sample 8/30/2021 [Example Alltesta](#) < Share

Description
instant Coffee

Compounds
Caffeine

Injections 2 + New Injection

~ 3

Setup new injection

Description

Device and configuration

Device
Example Alltesta Online

Method
EtOH/H2O/H3PO4 - 10/90/.3%

Column
Primesep SB_Example

Start

5. Before you begin the injection, you will need to create the Method.
Click on the **Method name** (above). Below is all the information that needs to be filled out:

Method

Search and use preset

Method name

Method 1

- General
- Pump
- ☒ Thermostat

General

Run time, min

Comment

→ Pump

Mobile phase

Buffer

Flow rate, mL/min

Pressure min, psi

Pressure reach time, s

Save pressure

ON

Injection

Sample amount, µl

Injection time, s

Injection delay, s

Saving method doesn't change method presets that you use.

Revert

Save

6. Each system also has 3 methods that have been downloaded onto it. You can see them by clicking on the **“search bar”** (above) at the top of the screen and selecting which method you would like.

- **Note:** when manually injecting your sample, do not rely on the measurement from the syringe - it is imprecise. Manual injections will be “full loop” injections, meaning the injection volume will be based on the (fixed) volume of the capillary tube connecting **positions 1 & 4**. You can either calculate the volume of a custom-cut capillary (length x internal diameter) or purchase and use pre-cut capillaries with known volumes.
7. When you select a method preset, it will fill in everything automatically for you, except the “**buffer**”, “**pressure min, psi**”, and “**pressure reach time, sec**”. Fill these out and click **Save** to finish and use this method for your injection. Alternatively, you can fill in the information manually to customize the method. Below is an example of a custom method.
- Enabling “**Save Pressure**” will allow you to track the pump pressure throughout the test.
 - A longer “**Injection Time**” will give more time for the mobile phase to pass through the loop.
 - “**Injection Delay**” can be useful for ensuring the detector reaches a proper baseline before the sample reaches the detector.
 - For most applications, you can leave “**Needle Depth**,” “**Intake Rate**,” “**Shaking Mode before Injection**,” and “**Shaking Duration**” blank.

Method

Method name
EtOH/H2O/H3PO4 - 10/90/.3%

- General
- Pump
- ⊗ Thermostat

General

Run time, min
4

Comment

→ Pump

Mobile phase
EtOH/H2O - 10/90%

Buffer
H3PO4 - 0.1%

Flow rate, mL/min
0.5

Pressure min, psi
200

Pressure reach time, s
30

Save pressure
ON


⊕ Injection

Sample amount, µl
10


Injection time, s
10

Injection delay, s
30


8. As mentioned in the Alltesta™ Mini Analyzer Set Up Guide, your valve set-up should look like this:



- This set-up will allow you to easily swap out sample and wash vials between tests.
 - The only 'moving' part will be the sample vial; the syringe will remain fastened in place.
9. Before you start the Injection, raise the plunger of the syringe slightly so that you pull your sample through the **green capillary** loop (connecting **positions 1 & 4**). Once the sample reaches the syringe, you know the sample has filled the loop.
 - With the valve in the default **Sample Loading** mode, when you raise the plunger, the sample will take the following route from the vial (**position 5**) to the syringe (**position 6**): **5 > 4 > 1 > 6**.
 - With the valve in **Sample Injection** mode, the pump (**position 2**) will push the mobile phase through the loop, 'injecting' the sample into the column (**position 3**): **2 > 1 > 4 > 3**.



- Click **Start**; the chromatograph should automatically present itself on the screen (refresh the page if it does not).



Coffee

Sample 8/30/2021 [Example Alltesta](#) < Share

Description
instant Coffee

Compounds
Caffeine

Injections 2 + New Injection

~ 3

Setup new injection

Description

Device and configuration

Device

Example Alltesta • Online

Method


EtOH/H₂O/H₃PO₄ - 10/90/.3%

Column

Primesep SB_Example

Start

11. The image below shows the window that will appear when the scan starts.



Coffee


Sample 8/30/2021 [Example Alltesta](#) Running [Share](#)

Description
Instant Coffee

Compounds
Caffeine

Injections

2 3 Preparing + New Injection



3


Sample injection 9/13/2021 3:33 PM Preparing Stop

Description

Device [Example Alltesta](#) [Vial A1](#)

Method [EtOH/H2O/H3PO4 - 10/90/.3%](#)

Column [Primesep SB_Example](#)



No recorded data

12. Steps 13-23 are the same as when using an Autosampler.
13. The image below shows a completed scan of a Coffee sample.



Injections 2 + New Injection

~ 2

Sample injection 8/30/2021 5:06 PM

Device Example Altista Visi A1
Method EtOH/H2O/H3PO4 - 10/90/.3%
Column Primesep SB_Example

Data
UV

Peaks editing Baseline editing Tools

#	Time	Area	Compound	Amount
1	2.699	12.152		

15. A popup will appear allowing you to select a Calibration from a drop-down menu.

Use calibration


Peak at 2.699 min
Area 12.152 mAU×min

Calibration

Estimate

16. Click **Estimate** after you have selected the Calibration you want to use.

17. The graph will now show the calculated amount of the compound in question for this particular injection based on the Calibration.



#	Time	Area	Compound	Amount
1	2.699	12.152	→ Caffeine	0.280989 mg/mL

18. You can also run an injection using the “Quick Start” method.

Stream LC T701-A3LF fv1.106
Online

Settings

→ Pump
203938554752500300330042
Ready
0.00 mL/min
Flow rate
-3 psi 28.50 °C
Pressure Temperature
5.80 mL
Start Stop

Valve loop
HV00000000
Ready

Sample loading
Quick start

Detector
HD00000000
Ready

Red 630nm
Green 520nm
Blue 460nm
UV 275nm

☐ Turn on leds
Calibrate

Method EtOH/H2O/H3PO4 - 10/90/.3%

Run time 4 min
Mobile phase EtOH/H2O - 10/90%
Buffer H3PO4 - 0.1%
Flow rate 0.50 mL/min
Pressure min. 200 psi reach in 30 sec

Injection 10.0 µl in 10 sec
Temperature 27 °C

Column Primesep SB_Example

Serial SBQLPLGL
Type Primesep SB


Length × ID 100 × 3.2 mm
Particle, pore size 5 µm, 100 Å

19. First, ensure the **Method** and **Column** selected for your device are correct. You can edit them by clicking on their respective **names**.

20. Then, click **Quick Start** and the system will start running.



21. To see the chromatograph, click on the **Sample title** that appears at the top of the page.



22. As shown above, one can stop the Injection by clicking the **Stop** button.

23. If you are viewing the chromatograph, you can also stop the Injection by clicking **Stop**.



24. When doing manual injections, it is highly recommended that you wash the valve and loop to clear out any leftover sample and mobile phase.

25. See the Alltesta™ Mini Analyzer Set Up Guide for the recommended washing method.

Gradient

1. In some cases, to properly separate compounds, you will need 2 chemically similar mobile phases with variations only in concentration (for example, ACN/H₂O/H₃PO₄ at 10%/90%/.03% and at 50%/50%/.03%). In this case, you will need 2 pumps, with each one drawing from a different mobile phase. To set this up in HPLC.Cloud, click on the [Settings](#) button on your device's page as you did when you were initially adding component devices.

Example Alltesta

Stream LC T70I-A3LF fv1.242 • Online



2. For “Pump B”, select the 2nd pump you would like to use.

Settings

Pumps

Pump

Auto

Pump B (gradient)

HP00000003



Injection

☐ Combine injection and fraction collector

Autosampler

Auto

Tray

auto

Valve loop

Off

Detection

Detector

Auto

☐ LED saving mode

☐ Analog signal

Phase

Volume of A container, mL

250

Volume of B container, mL

250

Waste

Volume of waste container, mL

250


Save

- The system will now be in “Gradient Pump” mode and show both pumps (A and B).

Example Alltesta

Stream LC T70I-A3LF fv1.245 Online

Settings



Method EtOH/H₂O/H₃PO₄ - 10/90/.3%

Run time 6.5 min Injection 5.0 µL in 360 sec
A phase MeCN/H₂O - 50/50 Sample intake at 100 µL/min
B phase Flow rate 0.50 mL/min
Pressure min. 100 psi reach in 30 sec

- To set up the Gradient, click “Edit” on the Method section (above).

5. In the “Gradient” option now available, click “Add Row”.
6. Edit the % of A and B and the time you want the instruction to occur. You can add as many or as few (minimum 1) instructions as you prefer.

Method

Search and use preset

Method name

EtOH/H2O/H3PO4 - 10/90/.3%

General

Pump

Gradient

Autosampler

Thermostat

General

Run time, min

6

Comment

Pump

Mobile phase

EtOH/H2O - 10/90%

Buffer

H3PO4 - 0.1%

Flow rate, mL/min

0.50

Pressure min, psi

0

Pressure reach time, s

Gradient

Time, min	A, %	B, %
+ Add row		

Save pressure

ON

7. If you change the time for a row, it will automatically move to be in chronological order (two instructions cannot have the same start time, so they will appear in red until one is changed).

→ Pump

Mobile phase

EtOH/H₂O - 10/90%

Buffer

H₃PO₄ - 0.1%

Flow rate, mL/min

0.50

Pressure min, psi Pressure reach time, s

0

Gradient

Time, min	A, %	B, %	
0.0	0	100	X
2.0	66	34	X
4.0	34	66	X
6.0	100	0	X

+ Add row

8. After editing the % value of A or B, the site will automatically enter the reciprocal concentration.

→ Pump

Mobile phase

EtOH/H₂O - 10/90%

Buffer

H₃PO₄ - 0.1%

Flow rate, mL/min

0.50

Pressure min, psi Pressure reach time, s


0

Gradient

Time, min	A, %	B, %	
0.0	0	100	X
2.0	34	66	X
4.0	66	34	X
6.0	100	0	X


+ Add row

9. Click the “Save” button once you're done editing the method.
10. You can now start your injection, either by going via the “Samples” tab or via “Quick Start.”
11. The system will then control the flow rates of each pump in order to follow the instructions set in the Method.
12. Open the chromatograph for the current sample and click on the “**Data**” bar (Pressure, UV, etc), then the “**G**” to show the Gradient profile.



Import

1. To import externally collected chromatographic data, click on the 3 dots opposite the Sample's name.
2. This will bring up a menu with four different options; click "**Import.**"



3. This will open a menu to select the method and column used, as well as the ability to select a file to import.

Create injection from CSV

Method

EtOH/H₂O/H₃PO₄ - 10/90/.3%

Column

Primesep SB_Example

CSV file

Select File

Upload

4. Once the file is selected, click "**Upload**".
5. The imported file should now appear as a new injection for that Sample.

Fraction Collection

1. The Alltesta also supports Fraction Collection, which allows you to inject, detect, and collect your separated and purified samples.
2. The Fraction Collection function requires an additional Free-Standing Valve and some additional capillaries. The guide for proper set-up can be found in the [Alltesta™ Analyzer Fraction Collector Set Up](#). Once the valve and capillaries have been set-up correctly, continue with the following steps.
3. Click on the Settings button in the top-right corner of your device's home page.
4. In the “Injection” section, select the “Combine injection and fraction collector” box.

Settings

Pumps

Pump
Auto

Pump B (gradient)
Off

Detection

Detector
Auto

☒ LED saving mode

☐ Analog signal

Injection

Autosampler
Auto

Tray
48_6x8

☐ Combine injection and fraction collector

Valve loop
Off

Save

5. Ensure that the “Tray” is set to the proper size setting.
 - a. 48_6x8 corresponds to a 6x8 plate with 48 positions.
 - b. 96_8x12 corresponds to an 8/12 plate with 96 wells.
6. From the “Fraction Collector Valve” drop-down menu, select “Auto.”
 - a. You can also select the valve itself by selecting the number string that corresponds to your valve.

Settings

Pumps

Pump
Auto

Pump B (gradient)
Off

Detection

Detector
Auto

☒ LED saving mode
☐ Analog signal

Injection

Autosampler	Tray
Auto	48_6x8

☒ Combine injection and fraction collector


Fraction collector valve
Off

Save


7. Click on your method to edit it.
8. There should now be a “Fraction Collection” group of settings.
9. Set how long you want your initial delay to be, or choose from one of the pre-filled options.
10. Set how long you want each fraction to collect for.
11. Select which vial you want to start your collection in.
 - a. Note: If the amount of vials needed for your collection exceeds the available amount of vials, the software will deliver an error. For example, if you start your collection in F1 and have a 10 minute run time with 1 minute collections (and no delays), you will need 10 vials. However, by starting at F1, there are only 8 remaining vials (F1-F8).

ACCOUNT

1. Click the **Account** tab at the top of the screen. Here, you will be able to see your preet methods, columns, and account information.



2. When you click on this tab you will be automatically directed to the **Method Presets** page. You will be able to see your currently available methods, as well as create new ones.



3. In order to create a new method, click on the **New** button and fill out all the necessary information (below).

Create Method Preset

Method preset name

Method

- General
- Pump
- ⊗ Gradient
- ⊗ Autosampler
- ⊗ Thermostat

General

Run time, min

1

Comment

→ Pump

Mobile phase

Buffer

Flow rate, mL/min

1

Pressure min, psi

Pressure reach time, s

Save pressure

OFF

⊕ Injection

Sample amount, μ l

Injection time, s

Injection delay, s

1

5

Create

- Once you have filled out the information, click **Create** and your method will be saved.

- To view and add new columns, click on the **Columns** tab.

Columns [View archived](#) [+ New](#)

Name	Serial	Type	Lot	Pressure	Length × ID	Particle	Pore
Primesep D-46.100.0510	PDOTLZAW	Primesep D	C01-056A	1000 psi	100 mm × 4.6 mm	5 µm	100 Å
Primesep 200.150.32.0310		Primesep 200		1000 psi	150 mm × 3.2 mm	0.3 µm	100 Å
Obelisc R 21.150.0510		Obelisc R		1000 psi	150 mm × 4.6 mm	0.5 µm	100 Å
Primesep 100.150.46.0510		Primesep 100		1000 psi	150 mm × 4.6 mm	0.5 µm	100 Å

- In order to create a new column, click on the **New** button and fill out all the necessary information (below).

Create Column

Column name
Column

Type Length, mm

Serial Internal diameter, mm

Lot Particle size, µm


Max pressure, psi Pore size, Å

[Create](#)


- Once you have filled out the information, click **Create** and your column will be saved.
- Instead of filling out the column info manually, you can click the **Camera** icon in the serial box to open a QR Code Reader.

9. Hold up the QR Code on the column packaging to your camera and the system will fill in all the information automatically.

10. If you wish to logout, click on the **Info** tab and select **Logout**.



11. To enable a split, dual window view, click the box for **Allow split view**.



That's everything you need to get an understanding of how HPLC.Cloud works. It's extremely simple to navigate as well as teaching newcomers their own way around it.

Again, HPLC.Cloud is accessible via the cloud which means you can use any device at any time to control your lab as long as you have a connection to the internet.


www.sielc.com

www.hplc.cloud

www.newcrom.com

TROUBLESHOOTING

- Below are some tips and tricks to troubleshoot some common problems.
- Pump stuck but in working state.** The pump can sometimes get stuck during the Refill stage, as seen in the image below (Pump B is stuck).




- To fix this issue, simply unplug the pump's power, count to 5, then plug it back in again.
- The pump should refill properly once it restarts.

1. **Pump cannot reach or maintain minimum pressure.** Sometimes the system cannot reach or maintain the minimum pressure, and the system will abort the injection.
2. This is often caused by a leak or air in the pump.
3. First, ensure that all connections are tightened properly and that there is no liquid dripping from any connections.
4. If there is liquid, simply loosen the fitting, reset the capillary in the fitting, and screw it back in.
5. If there is no liquid, then there may be air in the pump.
6. As mentioned in the [Alltesta Analyzer Set Up](#), unscrew the PEEK plug currently in the purge port.




7. Screw in a fitting with a teflon capillary (translucent) and connect the capillary to a waste flask.
8. Ensure that the fitting for the pump's intake capillary is tightened completely and that the end of the capillary that is in the mobile phase flask is completely submerged and surrounded by liquid. If it's not, you may need a longer capillary or more mobile phase in the flask.

9. Manually turn the pump on by clicking “**Start**”.



10. If there is air in the pump, you should see bubbles in the tubing.
11. Let the pump continue to run until there are no more bubbles exiting the pump (You may need to let the pump run a complete cycle - 5.80 mL).
12. The pump should now be drawing from liquid only.
13. Unscrew the fitting and capillary and replace it with the original PEEK plug.

1. **Baseline drift in initial part of injection.** Sometimes there may be an initial baseline drift at the beginning of the scan.



2. To alleviate this issue, increase the “Injection Delay” time in the Method.

⊕ Injection

Sample amount, μ l	Injection time, s	Injection delay, s
10.0	10	
Needle depth, mm	Wash cycles	Intake rate, μ l/min
Shaking mode before injection	Shaking duration, s	

1. **Capillaries pop out.** The capillaries can pop out of their fitting when the pressure inside is too high (more than approximately 2000 psi).
2. First, reinsert the capillary and the fitting into the correct position.
3. Then, lower the flow rate of the method.

→ Pump

Mobile phase

EtOH/H₂O - 10/90%

Buffer

H₃PO₄ - 0.1%

Flow rate, mL/min

0.50

Pressure min, psi


0

Pressure reach time, s



- *A slower mobile phase will result in a lower internal pressure.*
- *A higher mobile phase will result in a higher internal pressure.*

1. **Inconsistent peak heights.** Sometimes you may get inconsistent peak heights for a particular sample or standard.
2. First, check that your sample is properly purified (for standards) or mixed (for solutions).
3. If that doesn't resolve the issue, (and you're using an Autosampler) then your injection time may be too short.
4. Simply increase the injection time in the Method.

 Injection

Sample amount, μl	Injection time, s	Injection delay, s
10.0	10	
Needle depth, mm	Wash cycles	Intake rate, $\mu\text{l}/\text{min}$
Shaking mode before injection	Shaking duration, s	

5. This will keep the needle in the sample for longer, giving the syringe pump enough time to collect the designated sample amount.
6. Alternatively, you can increase the intake rate, which will increase the speed at which the syringe pump collects the sample.

INDEX

Devices

- Edit
 - Click to add or delete a device from your Devices list.
 - **Delete**: Deletes the device.
 - **Add connector**: Opens the feature to add a new device.
 - XXXX-XXXX: Input device's Serial number.
 - **: Input device's key.
 - **Scan**: Opens the camera to scan a QR code and automatically add your device.
 - **Cancel**: Cancels and closes the **Add Connector** feature.
 - **Finish Editing**: Save settings.

Devices

 Edit

 **Example Alltesta**

Stream LC T70I-A3LF

 Online

+ Add connector

XXXX-XXXX

**

Add device

Scan

Cancel

- Settings
 - Allows user to select component devices:
 - Pump
 - Gradient Pump (if necessary)
 - Detector
 - Autosampler
 - Tray size
 - Device options include:
 - **Off**: Disconnects component device.
 - **Auto**: Automatically selects an available device.
 - **HP(P/D/A)00000000**: The device's name. Will have one of P/D/A depending on if device is Pump, Detector, or Autosampler.

Settings

Pumps

Pump

Auto

Pump B (gradient)

Auto

Detection

Detector

Auto

☒ LED saving mode

☐ Analog signal

Injection

Autosampler	Tray
Auto	48_6x8

Valve loop

Auto

Save

- Pump
 - **Start:** Manually starts the pump.
 - **Stop:** Manually stops the pump.
- Autosampler
 - **Quick Start:** Creates a new sample and starts an injection for the selected vial.
 - **Wash Needle:** Runs the “Wash Needle” sequence.
 - **Shake:** Shakes the needle (to remove any remaining solution).
- Detector
 - **Turn on LEDs:** Turns on all 4 LEDs.
 - **Calibrate:** Calibrates and “zeroes” the detector to a standard baseline signal.

Example Alltesta

Stream LC T70I-A3LF fv1.106 • Online



→ Pump

203938554752500300320048

✓ Ready

0.00 mL/min

Flow rate

0 psi

Pressure

27.75 °C

Temperature

■ 5.80 mL

▶ Start

■ Stop

⊕ Autosampler

HA00000000

✓ Ready

▶ Quick start

☰ Wash needle

≡ Shake

✱ Detector

HD00000000

✓ Ready

Red 630nm

Green 520nm

Blue 460nm

UV 275nm

☐ Turn on leds

⊙ Calibrate


- Valve Loop
 - For manual injections (ie not using an Autosampler).
 - **Quick Start**: Creates a new sample and starts an injection for the selected vial.
 - **Sample Loading mode**: The default valve mode. Also used for washing the loop.
 - **Sample Injection mode**: Valve mode for injecting the sample into the column.

Example Alltesta


Stream LC T70I-A3LF fv1.106

Online

Settings



- Method
 - **Edit/[Method Name]**: Hover over the column name and click. Opens the Method window to allow editing of the current Method. To save, click “**Save**.” To “undo” any immediate changes, click “**Revert**.”
 - **Search Bar**: Allows you to search and select preset Methods.



Method EtOH/H2O/H3PO4 - 10/90/.3%

Run time 5 min Injection 2.5 µl in 5 sec
 Mobile phase EtOH/H2O - 10/90%
 Buffer H3PO4 - 0.1%
 Flow rate 1.00 mL/min
 Pressure (save) min. 200 psi reach in 30 sec

Column Primesep SB_Example

Serial SBQLPLGL Length x ID 100 x 3.2 mm
 Type Primesep SB Particle, pore size 5 µm, 100 Å

Recent Runs

 Search

Name	Date
★ Caffeine ~ 19	08/31 10:40 AM
★ Caffeine ~ 18	08/31 10:40 AM
★ Caffeine ~ 17	08/31 10:40 AM
★ Caffeine ~ 16	08/31 10:40 AM
★ Caffeine ~ 15	08/31 10:40 AM
★ Caffeine ~ 13	08/31 10:40 AM
★ Caffeine ~ 12	08/31 10:40 AM
★ Caffeine ~ 10	08/31 10:40 AM
★ Caffeine ~ 9	08/31 10:40 AM
★ Caffeine ~ 8	08/31 10:40 AM

- **Run Time:** Run time for the scan (in minutes).
 - **Mobile phase:** Enter the mobile phase (and percentages of each) you will be using.
 - **Buffer:** Enter the buffer (and percentage) you will be using.
 - **Pressure min:** Set a minimum pressure for the system to reach before starting the injection (*optional*).
 - **Pressure reach time:** Set the amount of time the system has to reach the minimum pressure (*optional*).
 - **Save Pressure:** Turn on to save pressure data with the chromatograph (*recommended*).
 - **Sample Amount:** Sets how much of the sample (in μL) will be injected.
 - **Injection Time:** Sets how long (in s) the collected sample will have to be injected into the pump line.
 - **Injection Delay:** Sets a delay (in s) before injection begins.
 - **Needle Depth:** Sets (in mm) how deep the autosampler will plunge the needle into the sample vial.
 - **Wash Cycles:** Sets how many wash cycles the needle will run after each injection (includes shaking).
 - **Intake Rate:** Sets the speed (in $\mu\text{L}/\text{min}$) at which the syringe will collect the sample.
 - **Shaking mode before injection:** Determines the shaking pattern used.
 - **Shaking duration:** Sets (in s) how long the shaking will last.
- Below you'll find the suggested intake rate based on the sample amount.

Sample Amount (μL)	Suggested Intake Rate
1	10
5	50
10	100
25	250
50	500
100	1000

- Column
 - **[Column Name]:** Hover over the column name to show a drop-down menu, which allows you to select columns from your Cabinet.
 - **New:** Hover over the column name to show a “New” button. Allows you to create a new column and add it directly to your device. Click “Create” when you have input all relevant information.
 - **Edit:** Hover over the column name to show an “Edit” button Allows you to edit the current column information.

▶ Start
 Stop

UV 275nm

☐ Turn on leds

Calibrate

Method EtOH/H2O/H3PO4 - 10/90/.3%

Run time 5 min
 Mobile phase EtOH/H2O - 10/90%
 Buffer H3PO4 - 0.1%
 Flow rate 0.70 mL/min
 Pressure (save)..... min. 200 psi reach in 30 sec

Injection 2.5 µl in 5 sec

Column Primesep SB_Example

Serial SBQLPLGL
 Type Primesep SB

Length × ID 100 × 3.2 mm
 Particle, pore size 5 µm, 100 Å

- Recent Runs
 - Shows the most recent injections on your device.

Recent Runs

 Search

Name	Date
 Sample 6622  1	08/30 5:44 PM
 Caffeine  19	08/30 5:13 PM
 Coffee  2	08/30 5:04 PM
 Sample 6620  3	08/30 4:25 PM
 Sample 6620  2	08/30 4:19 PM
 Sample 6620  1	08/30 4:06 PM
 Sample 6619  1	08/30 3:58 PM
 Sample 6618  1	08/30 3:54 PM

- Search:** Click to search through previous runs on your device. Search by:
 - Type (Calibration, Sequence, Sample)
 - Keyword
 - Device (if using more than one device)
 - Time (Date Range)
 - Click the sample of interest to open its chromatograph.

Search

[View Archive](#)











Any type

 Keywords

Example Alltesta


Any time

☐ Search in description

Name	Sequence	Time
Monday, Aug 30		
 Sample 6622		5:44 PM
 Coffee		4:26 PM
 Sample 6620		4:06 PM
 Sample 6619		3:58 PM
 Sample 6618		3:54 PM
 Sample 6617		3:54 PM
 Energy Drink	 Caffeine Test	3:43 PM
 Tea	 Caffeine Test	3:43 PM

Samples

- Shows all recent Samples created in chronological order (newest to oldest).
- Search:** Click to search through all samples created. Search by:
 - Type (Calibration, Sequence, Sample)
 - Keyword
 - Search in Description:** enable to search through Sample descriptions as well.
 - Device (if using more than one device)
 - Time (Date Range)
- New:** Click to open a new Sample window. Click “Create” after inputting all relevant information.
- Show more:** Shows more Samples.



The screenshot shows the HPLC Cloud interface. At the top is a navigation bar with icons for HPLC Cloud, Devices, Samples, Calibrations, Sequences, and Cabinet. Below this is a large heading 'Samples' with a search bar and a '+ New' button. The main content area displays a table of samples, grouped by date. The table has columns for Name, Device, and Time. The samples listed are:

Name	Device	Time
Monday, Aug 30		
Sample 6622	Example Alltesta	5:44 PM
Coffee	Example Alltesta	4:26 PM
Friday, Aug 27		
Sample 6564	Research Alex's Alltesta	6:03 PM
Sample 6562	Research Alex's Alltesta	5:57 PM
Sample 6561	Research Alex's Alltesta	5:51 PM
Sample 6559	Research Alex's Alltesta	5:44 PM


Below the table is a 'Show more' link.

- Sample Chromatograph
 - Click on a tested Sample to open the Chromatograph for that test.
 - To edit the Sample name, click on the name, type in a new name, then click outside the text box to save it.
 - Injections**
 - All Injections for this sample are listed in numerical order.
 - Click on each number to browse between different Injections.


- Use the arrow keys on your keyboard to cycle between Injections.




- **New Injection:** Click to start a new Injection of this Sample. Ensure the correct Vial position, Device, Method, and Column are selected, then click “Start.”



- **3-dot Menus:** There are two 3-dot Menus, one that governs the Sample and one that governs the selected Injection.
 - **Sample Menu:**
 - *Import:* Allows you to import your own chromatograph data for a given sample.
 - *To Calibration:* Creates a Calibration based on the Injections done for that sample.
 - *Edit:*
 - *Archive:* Archives the entire Sample (including all Injections).




- **Injection Menu:**
 - *Edit:* Allows you to edit/create a description for the selected Injection.
 - *Archive:* Archives the selected Injection.




• Graphing Tools

- *Data:* Click to show available options
 - **Pressure:** Toggles on/off Pressure data.
 - **Red:** Toggles on/off Red LED data.
 - **Green:** Toggles on/off Green LED data.
 - **Blue:** Toggles on/off Blue LED data.
 - **UV:** Toggles on/off UV LED data.
 - **View All:** Turns on all LED data.




- *Wrench*: Click to enable graphing tools.




■ Peaks Editing


- *Add smart peak*: Click to add a peak and the software will match it properly with the baseline.




- *Add peak:* Click to add a custom peak.



- *Erase peak*: Click to erase a peak.



- *Find Peaks*: Click to let the software automatically determine and label peaks. You can delete erroneous peaks by clicking the “X” as shown below.




■ Baseline Editing

- **Add Baseline:** Click to open the Baseline tool. Every click on the graph will be used as a setpoint to create a baseline.
 - After desired setpoints are established, click “Finish.”
 - To undo a setpoint, click “Undo.”
 - To restart the drawing, click “Start Over.”




- *Erase Baseline*: Click to remove the drawn Baseline.




- *Find Baseline*: Click to let the software automatically determine the baseline.

■ Tools


- *Download:* Download the raw data.




- **Compare:** Click to add selected Injection to the Compare window. Click again to open the Compare window.




Compare



- *Favorite*: Click to add to Favorites.




- **Zoom In:** Click and drag over the area on the chromatograph you want to zoom in on.
- **Zoom Out:** Click on the graph (after zooming in) to zoom out to the original dimensions.



Calibrations

- Shows all recent Calibrations created in age-descending order (newest to oldest).
- **Search:** Click to search through all Calibrations created. Search by:
 - Type (Calibration, Sequence, Sample)
 - Keyword
 - **Search in Description:** enable to search through Calibration descriptions as well.
 - Device (if using more than one device)
 - Time (Date Range)
- **New:** Click to open a new Calibration window. Click “Create” after inputting all relevant information.



Calibrations

[+ New](#)

Name	Device	Time
Monday, Aug 30		
★ Caffeine	! New Example Alltesta	3:31 PM


Compound name

★


Description, optional

≡


Device



Method



Column



Use Ctrl+Enter to create



New Injection: Ensure the correct Vial position, Device, Method, and Column are selected, then click “Start.”

★ Caffeine

Calibration

9/1/2021

! New

Description

1 mg/mL

Device

Example Alltesta

Method


EtOH/H₂O/H₃PO₄ - 10/90/.3%

Column

Primesep SB_Example

Injections

+ New Injection



No injections for this calibration. Click **New Injection** to create an injection.

- The sample amount and unit can be edited after starting the injection.

★ Caffeine

Calibration

9/1/2021

! New

Description

1 mg/mL

Device

Example Alltesta

Method

EtOH/H₂O/H₃PO₄ - 10/90/.3%

Column

Primesep SB_Example

Injections

+ New Injection

1

Setup new injection

Description

Vial

Amount

Device and configuration

Device

Example Alltesta

• Online

Method


EtOH/H₂O/H₃PO₄ - 10/90/.3%

Column

Primesep SB_Example




Start

- **Amount:** Sets the amount of known sample in the sample Injection. Unit options are:
 - %
 - mg/mL
 - mol (moles)
 - mmol (millimoles)
- **Peak:** Displays the ratio of mAU (milli-Absorption Units) to time (min) at the peak's location.
- **Select:** When hovered over, sets that peak as the Calibration peak.




Sequences

- Shows all recent Sequences created in age-descending order (newest to oldest).
- **Search:** Click to search through all Sequences created. Search by:
 - Type (Calibration, Sequence, Sample)
 - Keyword
 - **Search in Description:** enable to search through Sequence descriptions as well.
 - Device (if using more than one device)
 - Time (Date Range)
- **New:** Click to open a new Sequence window. Click “Create” after inputting all relevant information.


<div>  Devices Samples Calibrations Sequences Cabinet </div>			
<h1>Sequences</h1>		<input type="text" value="Search"/>	+ New
Name		Device	Time
Monday, Aug 30			
 Caffeine Test	! New	Example Alltesta	3:42 PM
 Caffeine Test	! New	Example Alltesta	3:40 PM

Name




Sequence name


Description



Expected compounds, separate with semicolon




Device




Example Alltesta

Method



EtOH/H₂O/H₃PO₄ - 10/90/.3%


Use Ctrl+Enter to create



Create

Cabinet


- Shows Method Presets, Columns, and Account Info.
- Method Presets
 - **Search:** Search by Name only.
 - **New:** Click to open a new Method Preset window. Click “Create” after inputting all relevant information.
 - **View Archived:** Allows you to view deleted/archived Method Presets.











The screenshot displays the HPLC Cloud interface. At the top, a navigation bar includes icons for HPLC Cloud, Devices, Samples, Calibrations, Sequences, and Cabinet. Below this, a secondary bar shows Method Presets (selected), Columns, and Info. The main section is titled 'Method Presets' and includes a 'View archived' link, a search bar, and a '+ New' button. A table lists three method presets with columns for Name, Run time, Injection, Flow rate, Temp., and Mobile phase.

Name	Run time	Injection	Flow rate	Temp.	Mobile phase
15 min	15 min	3 µl / 1 s	1 mL/min	22 °C	ACN/H2O/TFA-50/50/0.25
5 min	5 min	5 µl / 2 s	1 mL/min	20 °C	ACN/H2O/TFA-50/50/0.25
30 sec	0.5 min	3 µl / 1 s	1 mL/min	21 °C	ACN/H2O/TFA-50/50/0.5


- Columns
 - **Search:** Search by Name only.
 - **New:** Click to open a new Column window. Click “Create” after inputting all relevant information.
 - **View Archived:** Allows you to view deleted/archived Columns.





 Devices
  Samples
  Calibrations
  Sequences
  Cabinet

 Method Presets
  Columns
  Info

Columns


 View archived

 Search

 New

Name	Serial	Type	Lot	Pressure	Length × ID	Particle	Pore
Primesep D-46.100.0510	PDOTLZAW	Primesep D	C01-056A	1000 psi	100 mm × 4.6 mm	5 µm	100 Å
Primesep 200.150.32.0310		Primesep 200		1000 psi	150 mm × 3.2 mm	0.3 µm	100 Å
Obelisc R 21.150.0510		Obelisc R		1000 psi	150 mm × 4.6 mm	0.5 µm	100 Å
Primesep 100.150.46.0510		Primesep 100		1000 psi	150 mm × 4.6 mm	0.5 µm	100 Å

- Info
 - Shows Login/Email info.
 - **Allow Split View:** Click to enable dual window mode.



Account Info

