


Tetra-alkylammonium Salts on a Thermo Scientific™ Acclaim™ Surfactant Column with Suppressed Conductivity Detection



Column: Thermo Scientific™ Acclaim™
 Surfactant 5 μ m, 4.6 \times 150 mm
 Pump: Thermo Scientific™ Dionex™
 ICS-3000 DP
 Mobile Phases: (A) Water
 (B) 100 mM formic acid
 (C) 70:30 acetonitrile:water (v/v)
 Flow: 1.0 mL/min
 Gradient Times: -5 0 10 14
 %A: 80 80 30 30
 %B: 15 15 15 15
 %C: 5 5 55 55
 Temperature: Dionex ICS-3000 DC at 30 °C
 Injection: Dionex ICS-3000 AS autosampler; 15 μ L
 Detector: Suppressed conductivity
 Suppressor: Thermo Scientific™ Dionex™ CSRS™
 ULTRA II 4 mm, external water mode,
 3.0 mL/min, 44 mA
 Peaks: 200 μ g/mL each
 1. Tetrabutylammonium chloride
 2. Tetrapentylammonium bromide
 3. Tetrahexylammonium bromide
 4. Tetraheptylammonium bromide

24092

Tetra-alkyl ammonium compounds are especially difficult to measure because of their lack of a useful UV chromophore. Ordinary HPLC columns don't cleanly resolve these compounds, and usually the analyte peaks tail severely. The Acclaim Surfactant column provides good peak shapes and resolution for these salts. Suppressed conductivity can detect as little as 15 ng on column.