


Ulcer medication drugs (Cimetidine, also known as Tagamet) generally work against ulcers by actually turning off the pathway that produces stomach acids and other gastrointestinal fluids. This application was developed for general interest. If desired this mixture can also be separated at higher pH.

Highlights

- Agilent ZORBAX StableBond columns offer good peak shape for basic compounds at low pH.
- StableBond columns are extremely stable at low pH. This is also true for short chain, more-polar bonded phases like ZORBAX SB-CN.



Conditions

Column: ZORBAX SB-CN (4.6 x 150 mm; 3.5 μ m) (Agilent P/N: 863953-905) Mobile Phase: 90% 25 mM NaH_2PO_4, Monobasic (pH 3.0 with H_3PO_4): 10% Acetonitrile Inj. Vol.: 2 μ L Flow: 1 mL/min Temp.: 35°C Det.: UV (220 nm)



Robert Ricker is an application chemist based at Agilent Technologies, Wilmington, Delaware.

For more information on our products and services, visit our website at: www.agilent.com/chem

Copyright[©] 2002 Agilent Technologies, Inc. All Rights Reserved. Reproduction, adaptation or translation without prior written permission is prohibited, except as allowed under the copyright laws.

For Research Use only. Not for use in diagnostic procedures.

Information is subject to change without notice.

Printed in the USA April 25, 2002 5988-6419EN



Agilent Technologies