

Analysis of Polymer Extracts by Gas Chromatography-High Resolution Time-of-Flight Mass Spectrometry With Electron and Chemical Ionization

David E. Alonso, Joe Binkley, and Jeff Patrick | LECO Corporation, St. Joseph, MI USA

Introduction

Comprehensive profiles of polymer extracts can provide valuable information for the quality control of products or identification of adulterated materials. Extractables and leachables may result from incomplete polymerization, surface residues, secondary packaging of components, or product degradation (Figure 1).




Figure 1: Potential Extractables and Leachables

Electron ionization mass spectrometry (EI-MS) can help characterize extract components through library comparisons provided adequate databases are available. Often library hits are ambiguous due to mass spectral similarity of isomeric or homologous compounds. Matters are complicated further by the inability to determine molecular formulas because of the absence of parent ions in mass spectra of labile compounds. High quality spectral data and mass accuracy values (<1.0 ppm) facilitated confident characterization of residual monomers, oligomers, additives, and impurities in samples.

Characterization was facilitated through a combination of hard and soft ionization methods and high resolution time-of-flight mass spectrometry (HRT, Figure 2).




Figure 2: Pegasus® GC-HRT (top) and Operational Modes (bottom)

In this study, three samples (foam, rubber, and siloxane polymer) dispersed in thermoplastic polyurethane (TPU) polymer were analyzed using LECO's Pegasus GC-HRT. The combination of electron impact and chemical ionization high resolution time-of-flight mass spectrometry (EI and CI-HRT) was vital for unambiguous identification of polymer extractables.

Experimental

Sample Preparation

Polymers were extracted using common organic solvents (e.g., isopropanol, acetonitrile, etc.).

The procedure used for extraction of the Siloxane/TPU polymer is illustrated below (Figure 3).




Figure 3: Sample Preparation – Siloxane/TPU Polymer

Instrument Parameters

GC Parameters

Agilent 7890 and 7693 Auto Sampler
 Restek Rxi-5Sil MS (30 m, 0.25 mm ID, 0.25 mm df)
 Column Type:
 Injection: 0.5 mL, Split 3:1 (CI Splitless); Inlet Temp. 250°C
 Oven: 60°C (2 min) to 190°C at 10°C/min (1 min) to 250°C at 50°C/min (10 min)
 Carrier Gas: He, Constant Flow (1.00 mL/min)

MS Parameters

LECO Pegasus GC-HRT
 Spectrometer:
 Ion Sources: LECO EI, CI
 Source Temp.: 250°C (CI 225°C)
 Transfer Line Temp.: 300°C
 Spectral Acquisition: 6 spectra/second
 Mass Range (m/z): 50–650 (CI 180–1400)
 Calibration: PFTBA (Internal)
 CI Reagent Gas: 5% Ammonia in Methane
 Folded Flight Path: HR (R = 25,000 FWHM)

Analysis Workflow




Figure 4: Polymer Analysis Workflow; EI and CI-HRT Analysis

Results (Rubber Extract)

The Analytical Ion Chromatogram (AIC) in Figure 5 shows a representative set of compounds in rubber extract: Radical inhibitors (e.g., butylated hydroxytoluene (BHT), 2,2'-methylene-bis(4-ethyl-6-*t*-butyl)phenol), and oligomers such as 1-isopropenyl-2,2,4,4-tetramethylcyclohexane (C₁₃H₂₄) and 1,1,5,5-tetramethyl-2-(1-methylethenyl)-3-(2,2,4-trimethylpentyl) cyclohexane (C₂₁H₄₀). The mass accuracy values for these compounds ranged from -0.44 to 0.44 ppm. Peak true (Deconvoluted) mass spectra for BHT and the C₁₃H₂₄ oligomer are displayed in Figure 6.




Figure 5: AIC, Rubber Extract




Figure 6: Peak True Mass Spectra for BHT and a C₁₃H₂₄ Oligomer in Rubber Extract

Results (Foam Extract)

An AIC of major components in foam extract are shown in Figure 7. Two monomers 2,4- and 2,6-toluene diisocyanate (2,4- and 2,6-TDI), and the radical inhibitors 6-*tert*-butyl-2,4-dimethylphenol and butylated hydroxytoluene were easily identified through spectral similarity and formula searches of their EI-HRT data. Mass accuracy values for these compounds ranged from -0.63 to 0.11 ppm. Identification of the lactones and catalytic foam agent, 2,2'-dimorpholinodiethyl ether, required a combination of EI-HRT and CI-HRT data analysis, as well as scientific database and literature searches. Peak True CI-HRT data for these compounds are displayed in Figure 8.




Figure 7: AIC and XIC Expansion Showing Diisocyanates 2,4- and 2,6-TDI in Foam Extract




Figure 8: CI-HRT Data of a Dilactone (A), Dimorpholino Diethyl Ether (B), and Tetralactone (C) in Foam Extract

Results (Siloxane/TPU Extract)

EI-HRT analysis of the extract resulted in the AIC shown in Figure 9. The extract contained methylated cyclosiloxane homologs (D4-D10) with empirical formula C₂H₄O₂Si, BHT, diisocyanate monomers and some lactones. Cyclosiloxanes D4-D8 were easily identified using spectral similarity matches and accurate mass data for fragment ions [M-CH₃]⁺; however, D9, D10, and the series of compounds labeled DX (D11-D17) gave similar spectra as shown by two of the homologs later determined to be D14 and D15 (Figure 10).




Figure 9: EI-HRT, AIC of Siloxane/TPU Extract




Figure 10: EI-HRT, Peak True Mass Spectra of D14 (A) and D15 (B)

CI-HRT data was critical for identification of the entire series of methylated cyclosiloxanes. An XIC of the CI-HRT data for D11-D17 is shown in Figure 11. Peak True CI-HRT data for D14 and D15 confirms the formulas for these homologs (Figure 12). Table 1 lists the formulas and mass accuracy values (Ave. = 0.43 ppm) for the major components in the sample.




Figure 11: CI-HRT, XIC of D11-D17




Figure 12: CI-HRT, Peak True Mass Spectra of D14 (A) and D15 (B)

Table 1: CI-HRT, Compounds in Siloxane/TPU Polymer

| Compound | Formula | Species | Expected Ion m/z | Observed Ion m/z | Mass Delta (m/z) | Mass Accuracy (ppm) |
|--------------|----------------|----------|------------------|------------------|------------------|---------------------|
| D4 | C8H24O4Si4 | H Gain | 297.08244 | 297.08219 | -0.00026 | -0.86 |
| D5 | C10H30O5Si5 | H Gain | 371.10123 | 371.10109 | -0.00014 | -0.37 |
| D6 | C12H36O6Si6 | NH4 Gain | 462.14657 | 462.14607 | -0.00050 | -1.08 |
| D7 | C14H42O7Si7 | NH4 Gain | 536.16536 | 536.16526 | -0.00010 | -0.19 |
| BHT | C15H22O | H Gain | 219.17434 | 219.17438 | 0.00003 | 0.15 |
| D8 | C16H40O8Si8 | NH4 Gain | 610.18416 | 610.18423 | 0.00007 | 0.12 |
| D9 | C18H54O9Si9 | NH4 Gain | 684.20295 | 684.20293 | -0.00002 | -0.03 |
| D10 | C20H60O10Si10 | NH4 Gain | 758.22174 | 758.22170 | -0.00004 | -0.05 |
| D11 | C22H66O11Si11 | NH4 Gain | 832.24053 | 832.24025 | -0.00028 | -0.34 |
| MDI | C15H10N2O2 | H Gain | 251.08150 | 251.08131 | -0.00020 | -0.78 |
| D12 | C24H42O12Si12 | NH4 Gain | 906.25532 | 906.25510 | -0.00023 | -0.25 |
| D13 | C26H78O13Si13 | NH4 Gain | 980.27811 | 980.27941 | 0.00129 | 1.32 |
| D14 | C28H84O14Si14 | NH4 Gain | 1054.29690 | 1054.29750 | 0.00059 | 0.56 |
| D16 | C30H90O15Si15 | NH4 Gain | 1128.31570 | 1128.31558 | -0.00012 | -0.11 |
| D16 | C32H96O16Si16 | NH4 Gain | 1202.33449 | 1202.33514 | 0.00066 | 0.55 |
| D17 | C34H102O17Si17 | NH4 Gain | 1276.35328 | 1276.35264 | -0.00064 | -0.50 |
| Tetralactone | C20H30O8 | H Gain | 401.21699 | 401.21690 | -0.00009 | -0.24 |
| | | NH4 Gain | 418.24354 | 418.24322 | -0.00033 | -0.78 |

Summary

The GC-HRT's high resolution, quality spectral data, and excellent mass accuracy values were very useful for the confident identification of monomers, oligomers, and additives in polymers. A combination of EI and CI-HRT was shown to facilitate confident identification of extractables and leachables from different types of polymer systems.