


Fat-Soluble Vitamin Separation Using a Thermo Scientific™ Acclaim™ Phenyl-1 Column



Column: Thermo Scientific™ Acclaim™ Phenyl-1,
3 μ m
Dimensions: 3.0 \times 150 mm
LC System: Thermo Scientific™ Dionex™ UltiMate™
3000 RSLC
Mobile Phase: Methanol:water 90:10 (v:v)
Flow Rate: 0.50 mL/min
Temperature: 30 °C
Inj. Volume: 5 μ L
Detection: Diode array; UV at 220 nm
Peaks:

1. trans-Retinol acetate 100 μ g/mL
2. Vitamin D₃
3. Vitamin D₂
4. δ -Tocopherol
5. γ -Tocopherol
6. α -Tocopherol
7. Vitamin E acetate
8. Vitamin K₂
9. Vitamin K₁

27979

Fat-soluble vitamins are a diverse group of micronutrients that often come in multiple forms. They are used as ingredients in fortified foods and dietary supplements. Because of the complexity and variety of fat-soluble vitamins and the sample matrices they occur in, it is beneficial to have generic methods that can separate different fat-soluble vitamins within a single analysis. The Acclaim Phenyl-1 column is a specially designed reversed-phase column with high aromatic selectivity, high hydrophobic retention, and unique and complementary selectivity. In this example, the Acclaim Phenyl-1 column is used in a classic reversed-phase application. This stationary phase offers a different selectivity compared to other reversed-phase HPLC columns.