

Rapid Characterization of Alkaloids using Probe ESI Q-TOF LCMS-9050 in OAD-MS/MS

Hidenori Takahashi, Kaoru Nakagawa, Mami Okamoto, Yuta Miyazaki, Yohei Arao, Tetsuo Iida
Shimadzu Corporation, Kyoto, Japan

1. Introduction to OAD-MS/MS

- While low-energy CID-MS/MS is one of the most effective fragmentation techniques for structural analysis, it may not be ideally suited for the analysis of certain isomers.
- Several novel fragmentation techniques have been proposed to complement low-energy CID-MS/MS.

Table 1. Example of proposed novel fragmentation techniques

Electron-based fragmentation

EIEIO, ECD(Electron Capture Dissociation) by Zubarev et al. (1996)

Anion-based fragmentation

ETD(Electron Transfer Dissociation) by Syka et al. (2004)

Photon-based fragmentation

IRMPD (Infrared), UVPD (Ultraviolet), BRID (Blank body infrared)

- We have introduced **neutral radical-based** fragmentation techniques to structural analysis of biomolecules, peptides and lipids, since 2016.

Neutral radical-based *Takahashi et al., *Anal. Chem.* 2018, 90, 12, 7230.

Charge-neutral radical-induced dissociation is available in both positive and negative ion modes!

O• OAD (Oxygen Attachment Dissociation)-MS/MS
 $[M+H]^+ + O^\bullet \rightarrow [M+H+O]^+ \rightarrow \text{fragments}$

H• HAD (Hydrogen Abstraction Dissociation)-MS/MS
 $[M+H]^+ + H^\bullet \rightarrow M^{\bullet+} + H_2 \rightarrow \text{fragments}$




Fig. 1. Shimadzu LCMS-9050 (Q-TOF) with OAD unit.

2. Introduction to a direct ionization of PESI

- Probe Electro Spray Ionization (PESI) is one of the direct ionization techniques. Fig. 2 shows the scheme of the PESI system.




Fig. 2. Schematic diagram of the DPiMS™ QT system.

3. PESI-OAD Synergistic Lipid Analysis

Check WP019

- OAD-MS/MS reveals double-bond (C=C) positions not accessible with CID-MS/MS. Atomic oxygen selectively oxidizes and cleaves at C=C.





Fig. 3. Typical OAD spectrum of a model lipid of PC (18:1).

4. PESI-OAD Synergistic Alkaloids Analysis

- Atomic oxygen selectively oxidizes carbon atoms adjacent to nitrogen in nitrogen-containing heterocycles, leading to OAD-specific fragmentation distinct from CID, as shown in Fig. 4 (m/z 450.2127).

Scheme 1. Proposed OAD mechanism of nitrogen-containing heterocycles






Fig. 4. OAD and CID spectra comparison of reserpine.


Lycorine $[M+H]^+$ 288.123

OAD specific fragment ion m/z 154.086



Galantamine $[M+H]^+$ 288.159

OAD specific fragment ion m/z 72.081



Solanine $[M+H]^+$ 868.505

OAD specific fragment ion m/z 150.127




Fig. 5. Example of OAD specific fragmentation for several alkaloids.

5. Conclusions

- We have successfully developed the LCMS-9050 system integrating OAD and PESI.
- OAD provides unique structural information on nitrogen-containing heterocycles, which is distinct from CID, particularly observed in alkaloids.
- OAD specific fragment ions observed in alkaloids offer the potential for rapid characterization of alkaloids.