

A COMBINED WORKFLOW FOR IN-DEPTH CHARACTERIZATION OF CYSTEINE-CONJUGATED ANTIBODY DRUG CONJUGATES

Liuxi Chen¹, Robert Birdsall¹, Henry Shion¹, Ying Qing Yu¹, Frank W. Kotch², April Xu³, Thomas J. Porter⁴, ¹Asish Chakraborty, and Weibin Chen¹
¹Waters Corporation, Milford, MA; ²Pfizer Bioprocess Research & Development, Pearl River, NY; ³Pfizer Analytical Research & Development, Pearl River, NY; ⁴Pfizer Analytical Research & Development, Andover, MA

OVERVIEW





Fig. 1. Positional isoforms of cysteine-conjugated

INTRODUCTION

The structural complexity and intrinsic heterogeneity of antibody drug conjugates (ADCs) impose a prominent analytical challenge to current characterization methods. In this study, we present an analytical strategy of integrating several orthogonal methods for the in-depth characterization of critical quality attributes (CQAs) of ADCs including drug-to-antibody ratio (DAR), drug distribution, drug conjugation sites and conjugate site occupancy.

The following analytical methodologies were incorporated to facilitate a workflow for the in-depth characterization of the ADCs.

- HIC-UV and SEC-LC/MS at intact protein level were performed for the automated determination of DAR and drug loading distribution using an integrated informatics platform for streamlined data acquisition, processing and reporting.
- An on-line 2D LC/MS (HIC/RPLC) approach was applied to elucidate the positional isomers of the ADCs.
- Subunit analysis of ADCs (e.g. HC, LC, 1/2Fc, Fab, or Fd) were also carried-out using RPLC-UV or RPLC-MS, which enabled a rapid comparison of the ADC samples for batch assessment.

METHODS

- HIC/UV, SEC-LC/MS, RP-LC/MS

Instrumentation
 LC: Waters ACQUITY H-Class Bio
 MS: Waters Xevo G2-S QToF

HIC (LC/UV)
Columns
 Waters ProteinPak Hi Res HIC

Native SEC-LC/MS
Columns
 ACQUITY UPLC Protein BEH SEC Column, 200Å, 1.7 µm, 4.6 mm X 150 mm

RP-LC/MS
Columns
 ACQUITY UPLC Protein BEH C4 Column, 300Å, 1.7 µm, 2.1 mm X 50 mm

MS Conditions
 Capillary: 3kV; Sample Cone voltage: 150 v; Source Temp: 500°C; Desolvation Temp: 350 °C; Desolvation Gas Flow: 800L/h

- 2D LC/MS (HIC/RPLC)

Instrumentation
 LC: Waters ACQUITY H-Class Bio with 2D Technology
 MS: Xevo G2 QToF
Columns
 Waters Protein Pak Hi Res HIC
 ACQUITY UPLC Protein BEH C4, 300Å, 1.7 µm, 2.1 mm X 50mm

- Peptide mapping

Instrumentation
 LC: Waters ACQUITY H-Class Bio
 MS: Xevo G2-XS QToF
Columns
 Waters ACQUITY UPLC CSH C18 Column, 130Å 1.7 µm, 2.1 mm X 100 mm

MS Conditions
 Capillary: 3kV; Sample Cone voltage: 120 v; Source Temp: 120°C; Desolvation Temp: 250 °C; Desolvation Gas Flow: 600L/h

RESULTS

Intact Analysis




Fig. 2. Cysteine-conjugated ADC analysis using HIC. Drug distribution was determined for three different samples with increasing drug load.

Fig. 3. Deconvoluted intact mass spectra for cysteine-conjugated ADCs from native SEC-LC/MS after deglycosylation.

2D-LC/MS (HIC/RPLC) - Positional Isomers Determination




Fig. 4. Heart-cut fractions of A) DAR 4, B) DAR 6a, and C) DAR 6c were performed from individual HIC separations of cysteine-conjugated ADCs. A reversed phase gradient of each cut produced up to 3 peaks representing subunits of the positional isomers. Deconvolution of each peak resulted in unambiguous identification of the isoform for each fraction.

Subunit Analysis





Fig. 5. RPLC-UV chromatogram overlays of the reduced cysteine-conjugated ADCs

Possible Subunit isoforms after enzyme treatment



RP-LC/MS




Fig. 5. Total ion chromatogram of reverse phase separated sub units. The subunit structures for peaks 1-7 were shown above. Cys-ADCs samples were treated by FabULOUS and then reduced by DTT.

Peptide Mapping




Fig. 6. Tryptic peptide mapping MS^E chromatogram of cys-conjugated ADC (Moderate). Heavy chain T21 peptides with two conjugation sites are shown as an example. Unconjugated T21 (1), T21 with 1 conjugation site (2 and 3), and T21 with 2 conjugation sites (4) are indicated on the chromatogram.

Chain	Pep#	Peptide sequence	Modifier	Drug Occupancy ratio
Light	1:T2	VTTTCR	ADC_cys	7.3%
Light	1:T11	SGTASVCLLNINFPYR	ADC_cys	1.2%
Light	1:T18	VYACEVTHQGLSSPVTK	ADC_cys	2.1%
Light	1:T20	GEC	ADC_cys	100.0%
Heavy	2:T11	AEDTAVVYCAR	ADC_cys	1.9%
Heavy	2:T15	STSGTAAALGCLVK	ADC_cys	1.9%
Heavy	2:T20	SCDK	ADC_cys	100.0%
Heavy	2:T21	THTCPPCPAPEAAGAPSVLFPPKPK	ADC_cys, CAM	5.9%
Heavy	2:T21	THTCPPCPAPEAAGAPSVLFPPKPK	ADC_cys, CAM	4.8%
Heavy	2:T21	THTCPPCPAPEAAGAPSVLFPPKPK	ADC_cys x2	24.6%
Heavy	2:T23	TPEVTCVVVDVSHEDPEVK	ADC_cys	1.5%
Heavy	2:T37	INQVSLTCLVK	ADC_cys	3.5%
Heavy	2:T42	WQQGNVFCSSVMHLEALHNHYTQK	ADC_cys	1.4%

Table 2. List of cys-conjugated peptides observed in the moderate loading sample. Drug occupancy ratio = MS intensity of conjugated/ (MS intensity of unconjugated + conjugated peptides)

CONCLUSIONS

- DAR values and drug loading distributions for cysteine-conjugated ADCs are automatically acquired from HIC-LC analysis and from native SEC-LC/MS analysis, and the results show excellent agreement.
- 2D-LC/MS provides unambiguous identification of positional isomers in cysteine-conjugated ADCs.
- LC/MS^E identifies 13 conjugation sites with drug occupancy ratio calculated.

References

- Details on the HIC-UV, SEC-LC/MS and RP-LC/MS analysis: 61st ASMS conference, poster number TP236
- TP236 Details on 2D LC/MS analysis: 61st ASMS conference, poster number T2265