

Jun Xiang Lee<sup>1</sup>, Sue Ann Lee<sup>2</sup>, Timothy Yan Ann Lim<sup>3</sup>, Wan Tung Liw<sup>1</sup>, Jie Xing<sup>1</sup>, Zhaoqi Zhan<sup>1</sup>

<sup>1</sup> Application Development & Support Centre, Shimadzu (Asia Pacific) Pte Ltd, 79 Science Park Drive #02-01/08, Singapore 118264; <sup>2</sup>School of Physical and Mathematical Sciences, Nanyang Technological University, <sup>3</sup>Department of Chemistry, Faculty of Science, National University of Singapore

## Overview

- A novel method is developed for screening and quantifying six very less-polar FTOHs and FOSEs in textiles on LCMS-8050 with ESI in negative MRM mode
- LC/MS/MS quantitation method with calibration range, linearity, LOD/LOQ and repeatability was established with mixed standards
- Matrix effect and recovery were evaluated for a procedure of MeOH extraction of textile sample with sonication at RT, matrix effect observed at 31~96% and recovery at 63% ~ 143%.

## 1. Introduction

Perfluorinated compounds (PFCs) represent a large collective of chemicals, which have been used widely in textile, leather and other consumer products for their unique water and oil repellent properties. However, PFCs including PFAS (per- and polyfluoroalkyl substances) also become environmental and health concerns. Some components like perfluorooctanesulfonic acid (PFOS) and perfluorooctanoic acid (PFOA) are regulated as restricted substances in textiles by industrial standards and regulations<sup>[1]</sup>. LC/MS/MS and GC-MS have used in analyses of various PFCs or PFAS in various environmental samples and consumer products. A LC/MS/MS method was reported previously for simultaneously analysis of 24 PFCs in textiles<sup>[2]</sup>. Perfluorooctane sulfonamides (FOSEs) and fluorotelomer alcohols (FTOHs) belong to less-polar PFCs and may degrade to form PFOS and PFOA. Normally, they are analysed using GCMS with chemical ionisation (CI) due to their non-polar properties<sup>[3,4]</sup>. In this study, an alternative method using LC/MS/MS was developed for the analysis of four FTOHs and two FOSEs in textiles to respond the recent demands in testing analysis for textile and consumer products.

## 2. Experimental

Four FTOH standards, i.e., 2-perfluorobutyl ethanol (4:2 FTOH), 2-perfluorohexyl ethanol (6:2 FTOH), 2-perfluorooctyl ethanol (8:2 FTOH) and 2-perfluorodecyl ethanol (10:2 FTOH) were purchased from Apollo scientific, while two FOSE standards, 2-(N-methylperfluoro-1-octanesulfonamido)-ethanol (N-MeFOSE) and 2-(N-ethylperfluoro-1-octanesulfonamido)-ethanol (N-EtFOSE) were purchased from Wellington Laboratories. Stock standard solutions were prepared individually in methanol. The six compounds studied were mixed together and a calibration series was prepared using acetonitrile and water as diluent. A white textile sample was cut into smaller pieces and 1 gram of the sample was weighed into a 50mL polypropylene centrifuge tube. 20mL of pure methanol was subsequently added. The samples were sonicated for 30 minutes, followed by centrifugation at 10,000 rpm for 5 min. The supernatant was filtered using 0.22µm Nylon filter. 2 mL of supernatant was blown to dryness under a gentle stream of nitrogen gas. The obtained sample was re-constituted with 1 mL of diluent (A CN : water (50/50)) and analysed on a Shimadzu LCMS-8050, a triple quadrupole system with heated ESI, and a Shim-Pack GIST C18 column.

Table 1. Analytical conditions of four PFOH and two FOSE on LCMS-8050 with ESI

Column	Shim-Pack GIST C18 2µm, (2.1 mm I.D. x 100 mm L)	Interface	ESI (heated)
Flow	0.3 mL/min	MS Mode	Negative mode
Mobile phase	A : 5mM Ammonium Acetate in water B : Acetonitrile	CID gas	Argon, 230 kPa
Oven Temp	40 °C	Block Temperature	100 °C
Injection vol	20 µL	DL Temperature	100 °C
Elution mode	B% : 30% (0.0 min) → 95% (15.0 to 18.5 min) → 30% (18.6 to 22 min)	Interface Temp.	250 °C
		Nebulizing Gas Flow	Nitrogen, 3.0 L/min
		Drying Gas Flow	Nitrogen, 3.0 L/min
		Heating Gas Flow	Zero Air, 17.0 L/min

## 3. Results and Discussion

### 3.1 MRM transitions and optimization

The six very less-polar compounds (4:2 FTOH, 6:2 FTOH, 8:2 FTOH, 10:2 FTOH, N-MeFOSE and N-EtFOSE) are ionised effectively with ESI to form acetate adduct ion [M+CH<sub>3</sub>COO]<sup>-</sup> in negative mode (Figure 1). The ions tend to fragment to form acetate ion or remain intact due to the extreme stability of the structures. With an automated MRM optimisation program in LabSolutions, MRM transitions and respective CE were obtained for two transitions, [M+CH<sub>3</sub>COO]<sup>-</sup> > [CH<sub>3</sub>COO]<sup>-</sup> and [precursor]<sup>-</sup> > [precursor]. These two MRMs for each compound were used, the former for quantitation and the latter for confirmation to set up an LC/MS/MS method.




Figure 1: MS Spectra and precursors of individual FTOHs and FOSEs in negative ESI mode on LCMS-8050

### 3.2 Establishment of LC/MS/MS method for detection and quantification

Fig.2 shows the MRM chromatograms of the six compounds, which are well separated and sharp peaks using a C18 column and gradient elution program (Table 1). Calibration curves of the six compounds were constructed in different concentration ranges due to their differences in sensitivity. FOSEs exhibited higher sensitivities than FTOHs. Linear calibration curves were established in ranges of 1~50 ng/mL for FTOHs (except for 4:2 FTOH) and 0.02~1 ng/mL for FOSEs, respectively. Excellent linearity was obtained for all compounds (R<sup>2</sup> > 0.999) (Figure 3).




Figure 2: MRM chromatograms of four FTOHs (1 ng/mL each) and two FOSEs (0.02 ng/mL each) in a mixed standard solution on LCMS-8050.




Figure 3: MRM calibration curves of four FTOHs and two FOSEs in pure solvent on LCMS-8050

Table 2. MRM-based calibration curves and sensitivity of four FTOHs and two FOSEs on LCMS-8050

No	Compound Abbr.	Formula	CAS No	RT (min)	Quantifier MRM	Reference MRM	R <sup>2</sup>	Range (ng/mL)	LOQ (ng/mL)	LOD (ng/mL)
1	4:2 FTOH	C <sub>8</sub> H <sub>5</sub> F <sub>9</sub> O	2043-47-2	6.67	323.2>59.0	323.2>323.2	0.9999	5~50	3.7	1.2
2	6:2 FTOH	C <sub>8</sub> H <sub>5</sub> F <sub>13</sub> O	647-42-7	9.37	423.1>59.1	423.1>423.1	0.9999	2~50	1.7	0.6
3	8:2 FTOH	C <sub>10</sub> H <sub>5</sub> F <sub>17</sub> O	678-39-7	11.57	523.2>59.0	523.0>523.0	0.9997	1~50	0.8	0.3
4	10:2 FTOH	C <sub>12</sub> H <sub>5</sub> F <sub>21</sub> O	865-86-1	13.60	623.2>59.1	623.2>623.2	0.9994	1~50	0.5	0.2
5	NMeFOSE	C <sub>11</sub> H <sub>5</sub> F <sub>17</sub> NO <sub>2</sub> S	24448-09-7	12.13	616.1>59.1	616.1>616.1	0.9995	0.02~1	0.005	0.002
6	NEtFOSE	C <sub>12</sub> H <sub>5</sub> F <sub>17</sub> NO <sub>2</sub> S	1691-99-2	12.81	630.1>59.1	630.1>630.1	0.9999	0.02 ~ 1	0.002	0.001

The details of the quantitation method including sensitivity are summarized in Table 2. Refer to Fig.3, LOD and LOQ values were determined with the lowest level standard using LabSolutions, with S/N >= 3 and 10, respectively. The S/N of each compound was calculated according to ASTM method with an interval of 0.5 min for the whole time range. For the four FTOHs, quantitation limit (LOQ) and detection limit (LOD) are at 0.5~3.7 ng/mL and 0.2~1.2 ng/mL, respectively. Much higher sensitivities were achieved for the two FOSEs as compared to FTOHs, with LOQ and LOD at 0.002~0.005 ng/mL and 0.001~0.002 ng/mL, respectively.

Repeatability of the method was evaluated at two concentrations, (i) 2 ng/mL for FTOHs and 0.04 ng/mL for FOSEs and (ii) 5 ng/mL for FTOHs and 0.1 ng/mL for FOSEs. The results were tabulated into Table 3. Good repeatability was obtained for all the 6 compounds with peak area RSD < 10% at both concentrations.

### 3.3 Matrix effect and recovery of FTOHs and FOSEs in textiles

A simple extraction procedure of textile sample with ultrasonication at RT was employed. The matrix effect (M.E.) and recovery (R.E.) of the method of the six compounds were evaluated at three concentrations: (i)10 ng/mL for FTOHs and 0.2 ng/mL for FOSEs, (ii)20 ng/mL for FTOHs and 0.4 ng/mL for FOSEs and (iii)50 ng/mL for FTOHs and 1 ng/mL for FOSEs. The results are shown in Table 4.




Figure 4: Individual MRM chromatograms of four FTOHs at 1 ng/mL and two FOSEs at 0.02 ng/mL in mixed standards

Table 3. Repeatability for FTOHs and FOSEs (n=6) of standards

No	Compound	Conc. (ng/mL)	Area RSD (%)	Conc. (ng/mL)	Area RSD (%)
1	4:2 FTOH	2.0	ND	5.0	6.7
2	6:2 FTOH		6.1		1.7
3	8:2 FTOH		7.5		3.1
4	10:2 FTOH	0.04	4.4	0.1	9.6
5	NMeFOSE		3.8		2.1
6	NEtFOSE		4.2		1.1

Table 4. Matrix effect and recovery of FTOHs and FOSEs spiked in white textile with MeOH extraction (n=3)

No	Compound	Conc. (i)		Conc. (ii)		Conc. (iii)	
		M.E	R.E	M.E	R.E	M.E	R.E
1	4:2 FTOH	71.1	107.6	95.7	68.7	83.8	68.7
2	6:2 FTOH	59.0	97.0	75.9	62.8	70.8	67.2
3	8:2 FTOH	46.2	91.0	70.1	70.9	67.8	74.5
4	10:2 FTOH	30.5	104.1	45.1	69.0	56.2	71.7
5	N-MeFOSE	56.8	142.9	73.5	106.9	67.9	110.5
6	N-EtFOSE	59.4	135.3	75.0	98.2	70.6	107.7




Figure 5: MRM chromatograms of spiked sample (ii) (a) Blank textile extract (b) Pre-spiked textile extract of 20 ng/mL each FTOH, 0.4 ng/mL each FOSE (c) Post-spiked textile extract of 20 ng/mL each FTOH, 0.4 ng/mL each FOSE.

Across all three concentrations, it was observed that the textile extracts resulted in ion suppression for all the six PFCs. Most PFCs displayed matrix effect ranging from 46.2% to 95.7%.10:2 FTOH was more significantly affected by the interferences, resulting in matrix effect of 30.5%~56.2% across. As shown in Figure 5, significant interference peaks appeared for every compound. This indicates that further clean up step may be needed for the MeOH extract before analysis to quantitate trace levels of FTOHs and FOSEs by ESI-LC/MS/MS method.

Recovery of the sample preparation was also evaluated. Overall, good recovery was obtained for all targets across the three spiked concentrations. The four FTOHs achieved a recovery of 62.8 ~ 107.6%, while the two FOSEs achieved a recovery of 98.2 ~ 142.9%.

## 4. Conclusions

A MRM based LC/MS/MS method for quantitation of very less-polar PFCs, i.e., two FOSEs and four FTOHs in textiles was developed on LCMS-8050. Linear calibration curves were established with MRM transition of acetate adduct ion to acetate ion in ESI negative mode. High sensitivity were achieved with mixed standards, The LOQs and LODs for four FTOHs are in the range of 0.5~3.7 ng/mL and 0.2 ~ 1.2 ng/mL, respectively. While much lower LOQs and LODs for two FOSEs were obtained at 0.002~0.005 ng/mL and 0.001~0.002 ng/mL, respectively. Ion suppression due to interferences in textile MeOH extract are present, resulting in matrix effect of 30.5~95.7%. Good recovery for the sample preparation was obtained with 62.8~142.9%. However, further study is needed in sample clean-up method to reduce interferences and ion suppression to achieve high sensitivity of the method.

## References

- Renner, R. PerfluorinatedSources Outside and Inside. Environ. Sci. Technol., 2004, 38(5), 80A
- Jun Xiang Lee, Zhe Sun, Jie Xing, Jia Yi Tan and Zhaoqi Zhan, A Fast LC/MS/MS Method for High Sensitivity Determination of Twenty-Six Perfluorocompounds in Textiles, ASMS poster session, TP 585
- Lahey, C. M.; Ting, G. W. E.; Yeong, H. X. C.; Jackie, J.; Chong, C. M.; Loo, L. C. Fast and Accurate Analysis of Fluorotelomer Alcohols and Acrylates using Triple Quadrupole GC-MS/MS. ASMS poster session. 2016, ThP156.
- Yeong, H. X.; Olivia, S.; Lahey, C. M.; Ting, G. W. E.; Loo, L. C. Fast Screening and Quantitation of PerfluorinatedSources from Textiles using Chemical Ionization GC-MS. ASMS poster session. 2017, TP295.

**Disclaimer: The products and applications in this poster presentation are intended for Research Use Only (RUO). Not for use in diagnostic procedures.**