

TheReporter

Reprinted from Volume 15, No. 6, 1996

T296056

© 1999 Sigma-Aldrich Co.

For more information, or current prices, contact your nearest Supelco subsidiary listed below. To obtain further contact information, visit our website (www.sigma-aldrich.com), see the Supelco catalog, or contact Supelco, Bellefonte, PA 16823-0048 USA.

ARGENTINA • Sigma-Aldrich de Argentina, S.A. • Av. Pueyrredon 2446/50 • Piso 5-B • Buenos Aires 1119
AUSTRALIA • Sigma-Aldrich Pty. Ltd. • Unit #2, 14 Anella Avenue • Castle Hill NSW 2154
AUSTRIA • Sigma-Aldrich Handels GmbH • Hebbelplatz 7 • A-1110 Wien
BELGIUM • Sigma-Aldrich N.V./S.A. • K. Cardijnplein 8 • B-2880 Bornem
BRAZIL • Sigma-Aldrich Quimica Brasil Ltda. • Rua Sabara, 566-Conj. 53 • 01239-010 São Paulo, SP
CANADA • Sigma-Aldrich Canada, Ltd. • 2149 Winston Park Dr., Oakville, ON L6H 6J8
CZECH REPUBLIC • Sigma-Aldrich s.r.o. • Pobrezni 46 • 186 21 Praha 8
DENMARK • Sigma-Aldrich Denmark A/S • Vejlegaardsvej 65B • DK-2665 Vallensbaek Strand
FINLAND • Sigma-Aldrich Finland/YA-Kemia Oy • Teerisuonkuja 4 • FIN-00700 Helsinki
FRANCE • Sigma-Aldrich Chimie • Chromatographie Supelco • L'Isle d'Abeau Chesnes - B.P. 701 • 38297 Saint-Quentin Fallavier Cedex
GERMANY • Sigma-Aldrich Chemie GmbH • Geschäftsbereich Supelco • Grünwalder Weg 30 • D-82041 Deisenhofen
GREECE • Sigma-Aldrich (o.m.) Ltd. • 72 Argonafton Str. • 16346 Ilioupoli, Athens
HUNGARY • Sigma-Aldrich Kft. • Nagy Diófa u. 7., IV fl. • H-1067 Budapest
INDIA • Sigma-Aldrich Co. • Survey No. 31/1, Sitharamapalaya • Mahadevapura P.O. • Bangalore 560 048
IRELAND • Sigma-Aldrich Ireland Ltd. • Airton Road • Tallaght • Dublin 24
ISRAEL • Sigma Israel Chemicals Ltd. • Park Rabin • Rohovot 76100
ITALY • Sigma-Aldrich s.r.l. • Via Gallarate, 154 • 20151 Milano
JAPAN • Sigma-Aldrich Japan K.K. • Division Supelco • JL Nihonbashi Building • 1-10-15 Nihonbashi Horidome-cho, Chuo-ku • Tokyo 103
KOREA • Sigma-Aldrich Korea • Samhan Camus Annex, 10th Floor • 17-26 Yoido-dong, Yungdeungpo-ku • Seoul
MALAYSIA • Sigma-Aldrich (M) Sdn. Bhd. • 9-2, Jalan 2/128, Taman Gembira • Off Jalan Kuchai Lama • 58200 Kuala Lumpur • Selangor
MEXICO • Sigma-Aldrich Quimica S.A. de C.V. • Calle 6 North No. 107 • Parque Industrial Toluca 2000 • 50200 Toluca
NETHERLANDS • Sigma-Aldrich Chemie BV • Postbus 27 • 3330 AA Zwijndrecht
NORWAY • Sigma-Aldrich Norway • Sandakerveien 102 • N-0483 Oslo
POLAND • Sigma-Aldrich Sp. z o.o. • Szelagowska 30 • 61-626 Poznań
PORTUGAL • Sigma-Aldrich Quimica, S.A. • P.O. Box 131 • Sintra 2710
RUSSIA • Sigma-Aldrich Russia • TOO Techmedbiochem • Makarenko Str. 2/21 • Building 1, Flat 22 • Moscow 103062
SINGAPORE • Sigma-Aldrich Pte. Ltd. • 102E Pasir Panjang Road • #08-01 Citilink Warehouse • Singapore 118529
SOUTH AFRICA • Sigma-Aldrich (pty) Ltd. • CNR Kelly & Ackerman Streets • Southern Life Industrial Park Unit • Unit 16/17 • Jet Park 1459
SPAIN • Sigma-Aldrich Quimica, S.A. • Apt. Correos 161 • 28100 Alcobendas, Madrid
SWEDEN • Sigma-Aldrich Sweden AB • Solkraftsvägen 14C • 135 70 Stockholm
SWITZERLAND • Supelco Switzerland • Industriestrasse 25 • P.O. Box 260 • CH-9471 Buchs
UNITED KINGDOM • Sigma-Aldrich Company Ltd. • Supelco UK • Fancy Road, Poole • Dorset BH12 4QH
UNITED STATES • Supelco • Supelco Park • Bellefonte, PA 16823-0048 • Phone 800-247-6628 or 814-359-3441 • Fax 800-447-3044 or 814-359-3044 • email:supelco@sial.com

H

This article is archived from a past issue of The Supelco Reporter. Information in the article was appropriate at the time of publication, but product specifications, catalog numbers, and availability may have changed over time.

If you have questions about applying methodology described in this article to a current application, please contact our technical service chemists.

 **SUPELCO**

ISO 9001
REGISTERED

Supelco is a member of the Sigma-Aldrich family. Supelco products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for a particular use. Additional terms and conditions may apply. Please see the reverse side of the invoice or packing slip.

Comparison of SPB-PUFA and Omegawax Capillary Columns for FAMES Analysis by GC

L. Sidisky
Supelco, Bellefonte, PA USA

Omegawax columns, developed to meet the requirements of AOAC Official Method 991.39 for the analysis of fatty acids, are based on polyethylene glycol phase chemistry. A new capillary column, SPB-PUFA, is coated with a polyalkylene glycol stationary phase, which has a slightly lower polarity than the traditional PEG phase. In comparing these two columns, we found the SPB-PUFA column produces a different elution pattern and, in some analyses, gives unique separations.

The Association of Official Analytical Chemists (AOAC) Method 991.39, *Fatty Acids in Encapsulated Fish Oils and Fish Oil Methyl and Ethyl Esters*, recommends analysis using a column with a bonded polyglycol phase, based on Carbowax® 20M (1). The Omegawax™ 320 column was developed to meet these requirements. It applies a polyethylene glycol (PEG) phase chemistry and is the primary column used and recommended for the analysis of highly polyunsaturated FAME samples. The new SPB™-PUFA column is based on a polyalkylene glycol (PAG) stationary phase, providing a lower polarity than the Omegawax column.

This study compares the SPB-PUFA and the Omegawax columns in the analyses of two FAME standards — Supelco™ 37 Component FAME Mix, which mimics the fatty acid composition of a broad range of food samples, and PUFA No. 1 Mix, which is representative of marine fish oil samples. The 37 component mix was analyzed using temperature programmed runs (Figure A). The PUFA mix was analyzed at 200°C on the Omegawax column and at 210°C on the SPB-PUFA column (Figure B). Analyses were performed at other temperatures, but are not shown here.

Overall, the SPB-PUFA and the Omegawax 320 analyses show differences in the elution patterns of key FAMES, especially with the more highly unsaturated longer chain length FAMES. Comparing these analyses indicates that the less polar SPB-PUFA column provides a truer carbon chain length separation. All even carbon number FAMES elute according to carbon number and degree of unsaturation. On the Omegawax column, there is some chain overlap, as the C22:6n3 FAME elutes after the C24:0 FAME.



In order to quantify the difference in polarity of the two columns, we compared the equivalent chain length (ECL) values for a series of unsaturated FAMES (Table 1). Also listed are the absolute differences in the ECL values (Δ ECL), at 200°C and per °C. The Δ ECL/°C figures indicate the temperature effect on the selectivity of the columns. The values for the saturated FAMES are identical on each column, as, by definition, the ECL value for a saturated FAME is always equivalent to its carbon number.

Table 1 shows the Omegawax column is more polar than the SPB-PUFA column, as demonstrated by the larger ECL values for each

Figure A. 37 Component FAME Mix

Column: **SPB-PUFA**, 30m x 0.32mm ID, 0.20 μ m film
Cat. No.: **24323**
Column: **Omegawax 320**, 30m x 0.32mm ID, 0.25 μ m film
Cat. No.: **24152**
Oven: 50°C (hold 2 min), 4°C/min to 210°C
Carrier: 25cm/sec, set at isothermal temp.
Det.: (2 x 10⁻¹¹), 260°C
Sample Conc.: 50mg/mL
Inj.: 1 μ L of 37 Component FAME Mix, split 100:1, 250°C

Component	Weight (%)	Component	Weight (%)	Component	Weight (%)
1. C4:0	4	13. C16:1n7	2	25. C20:2n6	2
2. C6:0	4	14. C17:0	2	26. C21:0	2
3. C8:0	4	15. C17:1	2	27. C20:3n6	2
4. C10	4	16. C18:0	4	28. C20:3n3	2
5. C11:0	2	17. C18:1n9c	4	29. C20:4n6	2
6. C12:0	4	18. C18:1n9t	2	30. C20:5n3	2
7. C13:0	2	19. C18:2n6c	2	31. C22:0	4
8. C14:0	4	20. C18:2n6t	2	32. C22:1n9	2
9. C14:1n5	2	21. C18:3n6	2	33. C22:2n6	2
10. C15:0	2	22. C18:3n3	2	34. C23:0	2
11. C15:1	2	23. C20:0	4	35. C24:0	4
12. C16:0	6	24. C20:1n9	2	35. C22:6n3	2
				37. C24:1n9	2





796-0582, 0581

Reprinted from Volume 15, No. 6, 1996

Figure B. PUFA No. 1 Mix

Column: **SPB-PUFA**, 30m x 0.32mm ID, 0.20µm film
 Cat. No.: **24323**
 Column: **Omegawax 320**, 30m x 0.32mm ID, 0.25µm film
 Cat. No.: **24152**
 Oven: 210°C Isothermal (SPB-PUFA)
 200°C Isothermal (Omegawax 320)
 Carrier: 25cm/sec, set at isothermal temp.
 Det.: (2 x 10⁻¹¹), 260°C
 Sample Conc.: 50mg/mL
 Inj.: 1µL of PUFA No. 1 Mix, split 100:1, 250°C

1. C14:0	6. C18:1n7	11. C20:1n9	16. C20:5n3
2. C16:0	7. C18:2n6	12. C20:1n7*	17. C22:1n11
3. C16:1n7	8. C18:3n3*	13. C20:2n6	18. C22:1n9
4. C18:0*	9. C18:4n3	14. C20:4n6*	19. C22:5n3
5. C18:1n9	10. C20:1n11*	15. C20:4n3*	20. C22:6n3
			21. C24:1n9*



*Components not always present in this naturally-derived mix. Components not marked are guaranteed to be present in mix.

796-584, 0583

of the unsaturated FAMES. The ECL data also reveal that as the degree of unsaturation increases from monoenoic (C_x:1) to hexaenoic (C_x:6) for the same carbon number, the selectivity difference between the columns increases. The ΔECL /°C data show that the more unsaturated the FAME, the more temperature sensitive is its reaction with the column phase. Though not shown, the data generated by evaluating the columns at a series of isothermal temperatures show that the polarity of the columns increases as the temperature increases. This concurs with a previous study showing a temperature-based effect on column selectivity (2).

Table 1. ECL Values for SPB-PUFA and Omegawax 320 at 200°C

FAME	SPB-PUFA ECL	Omegawax ECL	ΔECL	SPB-PUFA ECL/°C	Omegawax ECL/°C
C14:1n5	14.33	14.43	0.10	0.001	0.001
C15:1	15.32	15.42	0.10	0.001	0.001
C16:1n7	16.21	16.30	0.09	0.002	0.002
C17:1	17.20	17.29	0.09	0.002	0.002
C18:1n9c	18.14	18.23	0.10	0.002	0.001
C18:1n9t	18.14	18.23	0.10	0.002	0.001
C18:2n6c	18.50	18.70	0.20	0.005	0.002
C18:2n6t	18.54	18.70	0.17	0.001	0.002
C18:3n6	18.72	19.00	0.28	0.004	0.003
C18:3n3	19.04	19.33	0.30	0.004	0.003
C20:1n9	20.11	20.20	0.09	0.002	0.002
C20:2n6	20.48	20.67	0.19	0.003	0.003
C20:3n6	20.66	20.94	0.28	0.004	0.005
C20:3n3	20.81	21.17	0.36	0.005	0.005
C20:4n6	21.02	21.32	0.30	0.004	0.004
C20:5n3	21.35	21.81	0.46	0.006	0.006
C22:1n9	22.10	22.19	0.09	0.002	0.002
C22:2n6	22.47	22.66	0.19	0.003	0.003
C22:5n3	23.33	23.82	0.49	0.006	0.006
C22:6n3	23.53	24.09	0.55	0.007	0.006
C24:1n9	24.09	24.18	0.09	0.002	0.002

The SPB-PUFA column's FAME elution pattern was found to be different from that of the Omegawax column; therefore, the SPB-PUFA column can be used as a confirmational tool. In some cases the SPB-PUFA column gives unique separations not possible on the Omegawax column, allowing it to be the primary column for some analyses.

Ordering Information:

Description	Cat. No.
SPB-PUFA Column 30m x 0.32mm ID, 0.20µm film	24323
Omegawax 320 Column 30m x 0.32mm ID, 0.25µm film	24152
Supelco 37 Component FAME Mix 10mg/mL each component in methylene chloride. Components listed in Figure A.	47885-U
PUFA No. 1 Mix 100mg. Components listed in Figure B.	47033

References

- Official Methods of Analysis of AOAC International, 16th Edition, 1995; Supplement, March 1996, AOAC International, Gaithersburg, MD, Chapter 41, p.21. *Fatty Acids in Encapsulated Fish Oils and Fish Oil Methyl and Ethyl Esters*,
- L.M. Sidisky and H.R. Ridley, HRC & CC 3 (1991) 191.
References not available from Supelco.

For more information on the 37 Component Mix, request publication 196907, *Comparison of 37 Component FAME Standard on Four Capillary GC Columns*.

Trademarks

Carbowax — Union Carbide Corp.
 Omegawax, SPB, Supelco — Sigma-Aldrich Co.

Fused silica columns manufactured under HP US Pat. No. 4,293,415.

