



Analysis of Flavors and Off-Flavors in Foods and Beverages Using SPME

Robert E. Shirey and Leonard M. Sidisky

Supelco, Supelco Park, Bellefonte, PA, 16823 USA

Introduction

SPME is a convenient, solventless extraction technique that can be used to extract analytes from both liquid and solid matrices. The use of SPME for the analysis of flavors and off-flavors in food and beverages is important.

In this presentation, sample types such as non-alcoholic and alcoholic beverages, candy, and fruits are analyzed for flavor composition. The detection of off-flavors from rancid oils and fats and methods for quantifying pyrazines in peanut butter, and caffeine in coffee are presented. The ability to detect trace (low ppt) levels of odors in water is also shown. Background information concerning the fiber types typically used for these analyses is given along with guidelines on how to select the appropriate fiber for a wide variety of applications.

Fig. 1 - Extraction Procedure for SPME





Fig. 2 - Desorption Procedure for SPME



Available SPME Fiber Coatings

Non-Polar Fibers

Polydimethylsiloxane (PDMS): 100 μ m, 30 μ m, 7 μ m

Polar Fibers

85 μ m Polyacrylate

65 μ m Carbowax[®]-divinylbenzene (CW-DVB)

50 μ m CW-Templated resin (CW-TPR) – HPLC only

Bi-Polar Fibers

65 μ m PDMS-DVB

65 μ m PDMS-DVB StableFlex[™]

75 μ m Carboxen[™]-PDMS

50/30 μ m DVB-Carboxen-PDMS StableFlex

Fibers for the Analysis of Flavors and Fragrance

Fiber Types	Types of Analytes	Concentration Range
75μm Carboxen-PDMS	most gases, volatiles	low ppt to high ppb
50/30μm DVB-Carboxen-PDMS	some gases, volatiles and semivolatiles	low ppt to high ppb
65μm PDMS-DVB	volatiles and semivolatiles	high ppt to low ppm
100μm PDMS	volatiles and semivolatiles	low ppb to high ppm
65μm Carbowax-DVB	volatile free acids, polar oxygenates	mid ppt to mid ppm

Comparison of SPME to other Extraction Techniques

Extraction Technique	Types of Analytes and Matrix	Conc.	Range	Ease of use
Static headspace analyzer	gases, volatiles & some semivolatiles, liquids & solids	wide	high	medium
Dynamic headspace P&T	some gases, volatiles liquids only	narrow	high	easy
SPE	semivolatiles liquids only	wide	low	hard
SPME	some gases, volatiles & semivolatiles liquids and solids	wide	low	easy

Fig. 3 - Volatiles in White Wine by GC/MS Using SPME

Sample: white wine + 25% NaCl
SPME Fiber: Carboxen™/PDMS
Extraction: 10 min headspace, 40°C
Desorption: 3 min at 290°C
Column: VOCOL™, 30m x 0.25mm ID, 1.5µm
Detector: GC/MS, Quadrapole, m/z = 31-240

1. Sulfur dioxide
2. Ethanol
3. Methyl formate
4. Acetic acid
5. Ethyl acetate
6. Isobutanol
7. Isopentanol
8. 2-Methyl-1-butanol
9. Ethyl butyrate
10. 2,3-Butanediol
11. Hexanol
12. Isoamyl acetate
13. Ethyl hexanoate
14. Hexyl acetate
15. Octanoic acid
16. Ethyl octanoate




Fig. 4A - Artificial Cherry Flavored Candy by SPME





Fig. 4B - Artificial and Natural Cherry Flavored Candy



Conditions for Analysis of Hard Candy by SPME

Sample: 0.5g candy in 5mL water in a 15mL vial
SPME Fiber: DVB-Carboxen™-PDMS StableFlex™
Extraction: headspace, 30 min at 40°C
Desorption: 270°C for 5 min

Column: Meridian MDN-5, 30m x 0.25mm x 0.25µm


Oven: 45°C (1.5 min) to 260°C at 4°C/min

Inj.: split or splitless with 0.75mm liner, 270°C

Det.: ion trap mass spectrometer, m/z = 33-400
at 0.6 sec/scan

Selected ions used for quantitation.

Fig. 5 - Volatile Aroma Compounds in Apple Fruit



Conditions for Analysis of Apple Aromas

Sample: “Mutsu” apple fruit, 300-400g in a 3 liter flask

SPME Fiber: 65µm PDMS/DVB

Sampling: headspace, 4 min, under a stream of N₂

Desorption: 250°C for 90 sec, then cryofocused at -100°C


GC Column: (5% phenyl) polydimethylsiloxane, 25m x 0.1mm ID, 0.34µm film

Oven: 40°C (1.5 min) to 250°C at 50°C/min

Carrier: helium, 1.5mL/min

Det.: mass spectrometer, m/z = 40-300

Fig. 6 - Milk Sample Off-Flavors by SPME-GC/MS



Conditions for Analysis of Milk Off-Flavors

Sample: 3g of 2% milk + 10 μ L IS (20 μ g/mL 4-methyl-2-pentanone) (9mL GC vial)

SPME Fiber: PDMS/CarboxenTM, 75 μ m film

Extraction: headspace, 15 min with constant stirring at 45°C

Desorption: 5 min, 250°C

Column: Supel-QTM PLOT, 30m x 0.32mm ID

Oven: 70°C (2 min) to 140°C at 6°C/min (2 min hold) then to 220°C at 6°C/min (5 min hold)


Carrier: helium, 35cm/sec

Inj.: splitless (closed 2 min)

Det.: GC/MS ion trap, m/z = 33-300

The exposure of unsaturated fatty acids to UV light can result in cleavage at the double bond. The resulting products are easily oxidized to form aldehydes such as hexanal and heptanal as shown in **Figure 6**. These components produce an off-flavor that is undesirable. The analysis of these by-products was traditionally done by purge and trap. Marsili of Dean Foods noted that SPME with the Carboxen-PDMS was not only as sensitive as purge and trap, but SPME provided a wider linear range compared to purge and trap. Also, SPME was suitable for detecting dimethylsulfide another off-flavor from oxidation of fats.

Fig. 7 - Analysis of Potato Chips



Identified Components in Rancid and New Potato Chips

1. Acetic acid
2. Pentanal
3. Butanoic acid
4. Propyl acetate
5. Methyl butyrate
6. Hexanal
7. Octane
8. Methyl hexanal
9. Hexanoic acid
10. Heptanone
11. Heptanal
12. Heptanoic acid
13. Octanal
14. Octanoic acid
15. Nonanone
16. Nonanal
17. Butyl hexanoate
18. Decanal
19. Undecanone
20. Pentyl hexanoate
21. Dodecanone
22. Methyl heptanol
23. Dodecanal

Conditions for Analysis of Chips

Extraction Conditions:

Fiber: DVB-Carboxen-PDMS StableFlex or 100 μ m PDMS

Sample: 3 grams of crushed potato chips in 15mL vial

Extraction: heated headspace, 65°C for 20 min

Desorption: 3 min at 250°C

GC/MS Conditions:

Column: SPB™-1 SULFUR, 30m x 0.32mm ID, 4.0 μ m film

Oven: 45°C (hold 1.5 min) to 250°C at 12°C/min (hold 10 min)

Carrier Gas: helium, 40cm/sec

Injection Port: splitless/split, closed for 2 min at 250°C

Detector: quadrupole mass spectrometer, m/z = 35-290 @ 0.6 sec/scan

Fig. 8 - Peppermint Oil in Chocolate Cookie Bar





Fig. 9 - Analysis of Peanut Butter Flavors by SPME

Sample: 5g peanut butter in 40mL vial
SPME Fiber: DVB-Carboxen™-PDMS StableFlex™
Extraction: headspace, 30 min at 65°C
Desorption: 5 min, 270°C
Column: SUPELCOWAX™ 10,
30m x 0.25mm x 0.25µm film
Oven: 40°C (5 min) to 230°C at 4°C/min
Inj.: splitless/split, closed 0.5 min,
270°C, with 0.75mm liner
Det.: ion trap mass spectrometer,
m/z = 30-350 at 0.6 sec/scan
Selected ions used for
quantitation.



Flavor Components in Peanut Butter

Some Volatile Components in Peanut Butter

1. Carbon disulfide
2. 3-Methylbutanal
3. Pentanal
4. Dimethyl disulfide
5. Hexanal
6. 4-Methyl-pentene-2-one
7. 1-Methyl pyrrole
8. Heptanal

Pyrazines in Peanut Butter

9. 2-Methyl pyrazine
10. 2,5-Dimethyl pyrazine
11. 2,3-Dimethyl pyrazine
12. 2-Ethyl pyrazine
13. 2,6-Dimethyl pyrazine
14. 2-Ethyl-6-methyl pyrazine

Pyrazines in Peanut Butter (contd.)

15. 2-Ethyl-5-methyl pyrazine
16. Trimethyl pyrazine
17. 2-Ethyl-3-methyl pyrazine
18. 2,6-Diethyl pyrazine
19. 2-Ethyl-3,5-dimethyl pyrazine
20. 2,3-Diethyl pyrazine
21. 2-Methyl-5-isopropyl pyrazine
22. 3-Ethyl-2,5-dimethyl pyrazine
23. 5-Methyl-2-propyl pyrazine
24. 2-Methyl-5-propyl pyrazine
25. 2-Ethenyl-6-methyl pyrazine
26. 3,5-Diethyl-2-methyl pyrazine
27. 2-Ethenyl-5-methyl pyrazine
28. 2-Methyl-6-cis propenyl pyrazine
29. 2-Allyl-5-methyl pyrazine

Quantitation of Pyrazines in Peanut Butter

$$\left(\frac{\text{Area counts}}{\text{spiked pb}} \right) - \left(\frac{\text{Area counts}}{\text{unspiked pb}} \right) \times \left(\frac{\text{g spiked pb}}{\text{g unspiked pb}} \right) = \left(\frac{\text{Area counts}}{\text{pyrazine spike}} \right)$$

Area counts (spiked pyrazine) = ng/g for each pyrazine

Pyrazines in ppb = $\frac{\text{ng/g} \times \text{area counts (unspiked p.b.)}}{\text{area counts (spiked pyrazine)}}$

Analytes	ppb
2-Methyl pyrazine	158
2,5-Dimethyl pyrazine	526
2,3-Dimethyl pyrazine	47
2,6-Dimethyl pyrazine	16

The roasting of peanut butter (PB) produces the formation of pyrazines from the Maillard reaction. The nutty flavor and aroma in PB are the result of the pyrazines. By heating the peanut butter to 65°C, the pyrazines are released from the fat and transferred into the headspace. The DVB-Carboxen-PDMS fiber was ideal for extracting the pyrazines along with some of the smaller flavor components as shown in **Figure 9**.

Quantitation of peanut butter can be accomplished by spike addition. An equal weight of a peanut butter sample was placed into 2 vials. One vial was spiked with a known weight of pyrazines. Both the unspiked and spiked vials of PB were extracted with the same fiber using identical conditions. The difference in area counts between the two samples provided the area counts of the spiked pyrazines. By determining the amount of the spiked pyrazines per gram of PB, the amount of each pyrazine could be determined in the unspiked PB. The results obtained are within published results for pyrazines in PB.

Fig. 10A - Analysis of Regular Coffee Grounds by SPME

Sample: 5g coffee grounds in 40mL vial
SPME Fiber: DVB/Carboxen™/PDMS StableFlex™
Extraction: headspace, 30 min at 65°C
Desorption: 270°C for 5 min

Column: SUPELCO WAX™ 10, 30m x 0.25mm x 0.25µm film

Oven: 40°C (5 min) to 230°C at 4°C/min

Inj.: splitless/split, closed 0.5 min, 270°C, with 0.75mm liner

Det.: ion trap mass spectrometer, m/z = 30-350 at 0.6 sec/scan

Selected ions used for quantitation.





Fig 10B - Analysis of Decaffeinated Coffee Grounds by SPME

Sample: 5g coffee grounds in 40mL vial
SPME Fiber: DVB/Carboxen™/PDMS StableFlex™
Extraction: headspace, 65°C, 30 min
Desorption: 5 min, 270°C

Column: SUPELCOWAX™ 10, 30m x 0.25mm x 0.25µm film
Oven: 40°C (5 min) to 230°C at 4°C/min
Inj.: splitless/split, closed 0.5 min, 270°C, with 0.75mm liner
Det.: ion trap mass spectrometer, m/z = 30-350 at 0.6 sec/scan
Selected ions used for quantitation.



Components in Coffee

1. 2-Methyl furan
2. 2-Butanone
3. 2-Pentanone
4. 3-Methyl butanal
5. 2,5-Dimethylfuran
6. 2-Acetyloxy-2-propanone
7. 2-Ethyl hexanol
8. Dimethyldisulfide
9. Phenol
10. Hexanal
11. 2-Methyl thiophene
12. n-Methyl pyrrole
13. 4-Methylphenol
14. 2-Ethyl pyrrole
15. Pyridine
16. Pyrazine
17. Methyl pyrazine
18. 4-Methyl thiazole
19. 3-Hydroxy butanone
20. Dimethyl phenol (isomer)
21. 1,2-Ethanediol, monoacetate
22. 2,5-Dimethylpyrazine
23. 2,3-Dimethylpyrazine
24. 2-Ethylpyrazine
25. 2,6-Dimethylpyrazine
26. 2-Ethyl-6-methylpyrazine
27. 2-Ethyl-5-methylpyrazine
28. Trimethylpyrazine
29. 2-Ethyl-3-methylpyrazine
30. 2,6-Diethylpyrazine
31. 2-Ethenylpyrazine
32. 2-Ethyl-3,5-dimethylpyrazine
33. Glycerol
34. 2,3-Diethylpyrazine
35. 2-Ethyl-3,6-dimethylpyrazine
36. 2-Furancarboxaldehyde
37. 2-Isopropenylpyrazine
38. 3,5-Diethyl-2-methylpyrazine
39. Furfural formate
40. 2-Furonyl ethanone
41. Methyl benzoylformate
42. Furanmethanol acetate
43. 5-Methyl-2-furancarboxaldehyde
44. Furanmethanol proprionate
45. Furfanyl furan
46. Pyridine methanol
47. 2-Methyl-5-propenylpyrazine
48. Furanmethanol
49. 3-Ethyl-4-methyl-2,5-furandione
50. Pyrazinecarboxamide
51. 2-Ethyl-3-hydroxy-4H pyran-4-one
52. 1-(2-Furanylmethyl)-pyrrole
53. 2-Methoxyphenol
54. 1-(1H-pyrrole-2-yl)-ethanone
55. 4-Ethyl-2-methoxy phenol
56. 3-Phenylpropenal or 2-Methylbenzofuran
57. 3,5-Dimethylbenzoic acid

Comparison of Caffeine Levels in Coffee and Extraction Type

Coffee and Extraction	Regular	Decaffeinated	%Decaffeinated
Grounds HS	1202079	207422	83%
Brewed Immersed	13623252	1567167	88%
Brewed HS	77431	8347	89%

1 hour extraction time with DVB-Carboxen-PDMS StableFlex Fiber
HS = headspace at 65°C

Fig. 11A - Odor Agents at 1ppt in Water by SPME-GC/MS

Sample: 30mL water containing MIB and geosmin at 1ppt and 25% NaCl in a 40mL vial, at 65°C

SPME Fiber: 2cm DVB/Carboxen™/PDMS

Extraction: heated headspace, 30 min, 65°C, with rapid stirring

Desorption: 3 min, 250°C, splitter closed

Column: Meridian MDN-5, 30m x 0.25mm x 0.25µm film

Oven: 60°C (1 min) to 250°C at 15°C/min

Det.: mass spectrometer, m/z = 75-180 at 0.6 sec/scan (quantitation ions 95 and 112)




Fig. 11B - 2,4,6-Trichloroanisole in White Wine by SPME

Sample: 10ppt of 2,4,6-TCA spiked in 12mL of white wine and 2.5g of NaCl

SPME Fiber: 100 μ m PDMS

Extraction: heated headspace, 30 min, 50°C, with rapid stirring

Desorption: 5 min, 250°C


Column: Meridian MDN-5, 30m x 0.25mm x 0.25 μ m film

Oven: 60°C (1 min) to 250°C at 15°C/min

Carrier.: helium, 30cm/sec, 60°C

Det.: Ion trap mass spectrometer, m/z = 160-220 at 0.6 sec/scan

Inj.: splitless/split, closed 2 min, 250°



G000174

Drinking water that comes from reservoirs may contain blue-green algae. This algae produces by-products that have a highly undesirable odors. These by-products, geosmin and methylisoborneol (MIB), produce a musky odor that is easily detected at 10 ppt by the human nose. In some cases, the threshold is less than 5 ppt. Even though these odors are not harmful, they can produce many customer complaints when detected. As a result, the water utilities monitor for MIB and geosmin at levels less than 5 ppt. **Figure 11A** shows the capability of heated headspace SPME and selected ion MS to detect these odor components at 1 ppt. A special 2 cm-SPME fiber is used to enhance sensitivity.

In wine, the bleaching of cork can cause the formation of 2,4,6-trichloroanisole. Like geosmin and MIB this by-product also has a low odor threshold around 10-20 ppt. Using headspace SPME this odor can be detected at 10 ppt in from wine as shown in **Figure 11B**.

CONCLUSIONS

- **SPME can be used to detect flavors in both solid and liquid foods.**
- **Both volatile and semivolatile compounds can be analyzed .**
- **SPME can easily detect a wide range of analytes with one fiber.**
- **Specificity can be obtained with different types of fibers.**
- **Quantification is possible with analyte addition.**
- **SPME can detect analytes at trace and high concentration levels in one sample.**