

# The use of Desorption ElectroSpray Ionisation with a novel heated transfer line for the analysis of lubricated surfaces

Eleanor Riches<sup>1</sup>, Emmanuelle Claude<sup>1</sup>, Caitlyn Da Costa<sup>1</sup>, Jeff Goshawk<sup>1</sup>, Gordon Jones<sup>1</sup>; 1. Waters Corporation, Stamford Ave., Altrincham Rd, Wilmslow, Cheshire, SK9 4AX, U.K.

## INTRODUCTION

Desorption ElectroSpray Ionisation (DESI) is a well-established direct ionisation technique for coupling with mass spectrometry. The ion plume is transferred from the sample region into the front of the mass spectrometer via a narrow, metal transfer line. In this work, we present the deployment of a novel heated transfer line (HTL) for the analysis of metal working fluid (MWF) deposited on machined metal components.

## METHODS

**Mass spectrometer:** Xevo™ G2-XS QToF with a DESI™ XS ion source (Waters, Wilmslow, UK)

**Acquisition mode:** Sensitivity

**Capillary voltage:** 0.7 kV (+ve/-ve)

**Cone voltage:** 40 V (+ve/-ve)

**Source temperature:** 100 °C

**Gas pressure:** 7 psi (~48 kPa)

**Acquisition range:** 50 – 1200 m/z

**Solvent:** 95:5 MeOH:H<sub>2</sub>O + 500 pg/μL Leucine Enkephalin for lock mass, flowing at 2 μL/min.

**Acquisition rate:** 800 μm/s

**Step size:** 0.2 mm (between acquisition lines)

**Acquisition approach:** profiling

**Samples:** 2 x machined metal aperture plates

**Protocol:** both aperture plates underwent a standard post-machining, cleaning process. One clean aperture plate was soaked in MWF (Blasocut – Blaser Swisslube, Switzerland) for 5 min., then allowed to drip-dry for 5 min. Both plates, clean (blank) and soaked, were mounted side-by-side on microscope slides and positioned on the DESI XS stage (Figure 1). The transfer line was heated step-wise via a heater controlled by an applied voltage (Table 1). Data were acquired using MassLynx™ v4.2 and processed using MassLynx v4.2 and UNIFI™ v1.9 (both Waters, UK).

## RESULTS & DISCUSSION




Figure 1. Aperture plates mounted on glass microscope slides and positioned on the DESI XS stage – clean, blank plate on the left and MWF soaked plate on the right.

Table 1. Applied voltages, resulting heater current, and displayed temperature being applied to the HTL – the same values were used for both positive ion and negative ion acquisitions

Voltage (V)	Current (A)	Displayed Temp. (°C)
0	0.000	26
2	0.682	43
4	1.361	98
6	2.023	191
8	2.676	305
10	3.310	426
12	4.039	563




Figure 2. TIC for five replicate lines across the blank and soaked plates at each set voltage from 0 V (ambient temperature) to 12 V (maximum temperature) with the positive ion y-axes (% Intensity) linked, and the negative ion y-axes (% Intensity) linked.

## CONCLUSIONS

- A higher temperature applied to the HTL showed clear benefit for analysing the deposited MWF, in both positive ion mode and negative ion mode.
- In positive ion mode, a polymeric component of the MWF was observed at higher temperatures that was not observed at ambient temperature. This polymer was not observed at all on the blank plate.
- The  $\Delta m/z$  of the polymer was 44, suggesting it is poly(ethylene glycol) (PEG).
- In negative ion mode, a similar increase in regularly spaced clusters of ions was observed. These clusters of ions were not observed at all on the blank plate.
- The  $\Delta m/z$  between the ion clusters was 14, suggesting these are different series of hydrocarbons. This is indicative of the base oil used in the MWF.
- The key components of the MWF would not have been observed using a transfer line at ambient temperature – full characterisation of the MWF was only possible using the novel HTL.




Figure 3. [A] Soaked plate, positive ion mirrored spectra with linked axes, 0 V compared with 2 V. [B] Soaked plate, positive ion mirrored spectra with linked axes, 0 V compared with 12 V. [C] Soaked plate, negative ion mirrored spectra with linked axes, 0 V compared with 2 V. [D] Soaked plate, negative ion mirrored spectra with linked axes, 0 V compared with 12 V.