

Advanced Native and Top-Down MSⁿ Workflows on the timsOmni Platform Enable Deep Structural Analysis of Membrane Protein Complexes

Jeremy Norris¹; Abraham O Oluwole^{2,3}; Mariangela Kosmopoulou⁴; Dodge Baluya¹; Ioanna Barla⁴; Athanasios Smyrnakis⁴; Dimitris Papanastasiou⁴; Weston Struwe^{3, 5}; Carol V Robinson^{2,3}

¹ Bruker Scientific LLC, San Jose, CA, US; ²Department of Chemistry, University of Oxford, Oxford, United Kingdom; ³Kavli Institute for Nanoscience Discovery, University of Oxford, Oxford, United Kingdom; ⁴Fasmatech Science and Technology, Athens, Greece; ⁵Department of Biochemistry, University of Oxford, Oxford, United Kingdom

Introduction

- Membrane proteins are essential for cellular function and major therapeutic targets, yet their structural and functional characterization remains challenging due to diverse oligomeric states, proteoforms, post-translational modifications, and interactions with lipids, cofactors, and partner proteins.
- Native mass spectrometry (MS), combined with top-down (TD) analysis, enables detailed characterization of membrane protein assemblies and subtle modifications such as lipidation.
- Here, we present advanced MSⁿ workflows on the timsOmniTM platform for comprehensive analysis of membrane proteins, their modifications, and functional complexes.

Methods

- Sample introduction:** Coated, open type nanospray emitters (biased at 1.3kV), 1 mm i.d. glass inlet at 8 mbar
- Desolvation/isCID:** Ion acceleration at 120 eV to 200 eV, 0.5 mbar
- Mass selection:** Quadrupole mass filter for charge states of $m/z < 4500$ Th; Q2 segment of the OmnitrapTM platform used to isolate larger ions, which can be accumulated to enhance signal intensity
- Protein dissociation and fragmentation:** Fig. 1.
- Data analysis:** OmniScapeTM for MSⁿ data.




Fig. 1 Native Top-Down MSⁿ Workflow Overview.

Results

For the heteropentameric β -barrel assembly machinery BamABCDE complex, we achieved systematic dissociation of the intact complex into individual subunits using in-source collision-induced dissociation (isCID) followed by MSⁿ analysis. This approach enabled deep sequencing and detailed characterization of N-terminal lipidation on each lipoprotein subunit (Fig. 2).





Fig. 2 The timsOmni platform enables the disassembly of BamABCDE intact membrane-protein complex, using psMS2 CID (isCID), for targeted subunit characterization. Two subunits were interrogated on the MS³ level: BamE (workflow in blue arrows) and BamD (workflow in yellow arrows). Ion enrichment of low-abundance lipidated fragment ions significantly improves signal-to-noise, enabling MS⁴ fragmentation. Comprehensive localization and structural elucidation of N-terminal lipidation was achieved through MS⁴ rCID, revealing the lipid-anchor architecture directly from protein-derived fragment ions.



Learn more about
Structural
Proteomics



Low-abundance b-type fragments carrying the lipid modification, generated by MS³ resonance CID (rCID) in Q2, were selectively isolated and enriched, substantially improving signal-to-noise and enabling an additional stage of fragmentation.

The analytical strength of this ion-enrichment strategy is demonstrated by MS⁴ rCID of the b_{12}^{2+} product ion from the BamE subunit, which allowed direct structural elucidation of the N-terminal lipidation. MS⁴ rCID produced an extensive series of neutral-loss fragments, including palmitic acid (C16:0), oleic acid (C18:1), and the corresponding 1-oleoyl-2-palmitoyl-glycerol species. Additional neutral loss of an S-linked diacylglycerol unit (oleoyl plus palmitoyl) indicated attachment of the glycerol moiety to the cysteine side chain. Loss of a hexadecyl group further confirmed N-terminal linkage of a second palmitic acid to the cysteine residue, enabling full structural characterization of the lipid anchor.

The MSⁿ workflow was subsequently extended to the remaining lipoprotein subunits. In particular, MS³ CID of BamD yielded a characteristic series of b-type ions exhibiting CH_2 and C_2H_4 losses and the same lipid-associated mass shift observed for BamE. The most intense species, $[\text{b}-\text{CH}_2]_{10}^+$, was further interrogated by MSⁿ and enabled the confident localization of the lipid on the N-terminal Cysteine. Additionally, lipid-related losses from both the precursor and the fragments are in agreement with the corresponding findings in the MSⁿ workflow for BamE.

Conclusion

- Native top-down MSⁿ workflows on the timsOmni platform, exploiting ion enrichment, enable localization and detailed structural elucidation of lipidation on a multimeric protein complex.

**Proteomics: Intact Proteins and
Top-Down Analysis**