

Introduction to MetID Solution Software for Structural Elucidation of Expected and Unexpected Metabolites

Technical Report vol.11



1. Introduction

Identification of metabolites is an important step in the process of drug development. Mass spectrometry (LC/MS) is routinely used for metabolite identification and structural analysis. Typically, the MS¹ measurement of the control sample and actual sample are compared, and the results of metabolite ion masses (m/z) detected in the actual sample are compared to the m/z of expected metabolites. Metabolite identification generally involves sample analysis of urine and blood matrices which are often complex and include many substances. Accordingly, extremely large numbers of peaks are generated. If the m/z detected in the actual sample is very close to the m/z of the expected metabolite, it becomes difficult to distinguish one from the other. Even time-of-flight or Fourier transform mass spectrometers, though capable of accurate mass analysis, can still produce false positives or fail to differentiate among some structurally similar metabolite candidates and the parent compound.

This report introduces the software “MetID Solution”, developed to automatically extract and structurally analyze metabolite candidates through comparison of not only MS¹ but MS/MS and MSⁿ spectral results, such as can be generated by an LCMS-IT-TOF ion trap, time-of-flight mass spectrometer.

2. Metabolite Search from MS¹ Analysis Results

Fig. 1 shows a pattern diagram depicting the metabolite search method (chromatogram comparison) employed using the results of MS¹ analysis.




Fig. 1: Metabolite Candidate Extraction by Comparison of MS¹ Chromatograms

As shown in Fig. 1, comparison of chromatograms involves a control sample and an actual sample. The Extracted Ion Chromatograms (XIC) as generated using a high resolution mass spectrometer are compared.

Although chromatogram comparison is a useful means of searching for metabolites, these results obtained from observation of MS¹ alone cannot reliably determine whether or not a substance is actually a metabolite. For this, additional mass spectrometry data is required.

3. Including Structural Information for Comparison

One method for determining whether or not a substance is a metabolite is to compare the MS/MS spectrum of each substance in the actual sample with that of the parent drug. Basically, if a parent compound is metabolized, there is a strong possibility that metabolites will maintain partial structural integrity with respect to the original compound. Thus, commonalities (fragment ions or neutral loss) would be expected in the results of MS/MS measurement of the parent compound and its metabolites.

With the LCMS-IT-TOF, which is capable of high-speed MS/MS analysis, a large number of MS/MS spectra can be obtained from a single analysis of an actual sample.

MetID Solution software automatically finds the precursor ions in the MS/MS spectra having commonalities with the parent compound, and generates a table of common fragment ions and neutral losses. This process is depicted in the pattern diagram below.




Fig. 2: Metabolite Candidate Extraction by Statistical Processing of MS/MS Spectra

4. MetID Solution Features

MetID Solution software is a tool that utilizes a dual approach to improve the efficiency and reliability of metabolite identification.

	Chromatogram Comparison	MS/MS Spectral Comparison
Data used	MS ¹	MS/MS (up to MS ⁿ)
Method	Differences are extracted by comparing control sample and actual sample chromatograms (XIC).	MS/MS spectra having commonalities with referenced parent compound MS/MS spectrum are extracted.
Features	Conceptually easy. Possible with any mass spectrometer.	Fast, highly reliable.
Advantages	Possible with any mass spectrometer. Applicable for non-targeted search (no automated reference).	Control sample is not required. Uses PLS statistics to help find metabolites quickly. High robustness with respect to analytical conditions.
Shortcoming	Depends greatly on data repeatability. It takes time to accurately compare the XIC for all mass regions. Requires control samples	

5-1. Measurement Example

A metabolite sample of verapamil ($C_{27}H_{38}N_2O_4$; MW 454.283), a compound in which the principle metabolites are known, was measured by LCMS-IT-TOF, and data processing was conducting using MetID Solution.




Fig. 3: Verapamil Structural Formula — A to E indicate cleavage sites.




Fig. 4: Verapamil MS/MS Spectrum — * indicates neutral loss.

5-2. Determination of Metabolites by Chromatogram Comparison




Fig. 5: Chromatogram Comparison Results Window

5-3. Determination of Metabolites by MS/MS Spectral Comparison

Points of commonality between the parent compound MS/MS spectrum and the precursor ion in the analyte are indicated with positive values in the horizontal axis direction.

In this example, four types of metabolites were detected. In addition, fragment ions in common and neutral loss information are automatically summarized in the table.




Fig. 6: MS/MS Spectral Comparison Results Window (partial)

5-4. Metabolite Structure Prediction by MS/MS Spectral Comparison

Looking at P#141 (m/z 441.2737, RT 14.100) as the most relevant ion indicated in the horizontal axis, the mass difference (-14.0151 Da) with respect to the parent compound suggests that it corresponds to the neutral loss of CH₂ (14.0157), or in other words, that demethylation occurred. In addition, signals associated with the

B neutral loss and the fragment ions corresponding to D and E shown in Fig. 3 were also found in the MS/MS spectrum of P#141. This result suggests that the structures of the B neutral loss part and those of D and E are consistent with the overall structure of the original compound and P#141.




Fig. 7: Predicted Structure for P#141

One of the principle metabolites of verapamil is norverapamil (C₂₆H₃₆N₂O₄; MW 440.267). This structure does not contradict the information obtained from the MS/MS spectrum of P#141.




Fig. 8: Structural Formula of Norverapamil

Reference: S. Yamaguchi et al., J. Mass Spectrom. Soc. Jpn., 55, 83(2007) et al

Founded in 1875, Shimadzu Corporation, a leader in the development of advanced technologies, has a distinguished history of innovation built on the foundation of contributing to society through science and technology. We maintain a global network of sales, service, technical support and applications centers on six continents, and have established long-term relationships with a host of highly trained distributors located in over 100 countries. For information about Shimadzu, and to contact your local office, please visit our Web site at www.shimadzu.com



SHIMADZU CORPORATION. International Marketing Division
3. Kanda-Nishikicho 1-chome, Chiyoda-ku, Tokyo 101-8448, Japan
Phone: 81(3)3219-5641 Fax: 81(3)3219-5710
URL <http://www.shimadzu.com>

The contents of this brochure are subject to change without notice.