

Enables the Brevis™ GC-2050 to be Controlled by Thermo Scientific™ Dionex™ Chromeleon™ 7

# Shimadzu GC Driver for Chromeleon 7




Brevis™ GC-2050



Nexis™ GC-2030

## Enables Control of the Brevis GC-2050, a System that Offers Uncompromised Analytical Performance in a Small Footprint


The compact Brevis GC-2050, with a system width of 350 mm, including the autoinjector (AOC™-30i), allows for even more efficient use of lab space. Compared to the same Nexis GC-2030 configuration, the system width has been reduced by approximately 35 %. Despite its space-saving design, the GC oven design does not require a dedicated column, and general capillary columns can be used. Dual-line analysis is also supported.



## Facilitates Easy Daily Maintenance with Support for the Easy sTop Function

The GC-2050 is equipped with a function (Easy sTop) that automatically lowers the GC inlet temperature and simplifies liner (consumable) replacement by simply pressing the MAINT. button on the upper right of the GC front panel. It can also be operated from the Chromeleon window, helping to ensure a convenient workflow in the lab.



Note: Operating the Easy sTop function from Chromeleon is possible with both the GC-2030 and GC-2050.



## Supports He Gas Saving and Hydrogen Carrier Usage

The GC-2030/GC-2050 can be equipped with an optional built-in hydrogen sensor to detect potential leaks early and maintain a safe standby mode. If hydrogen leakage increases, the main power is turned off to prevent accidents.

When the GC-2030/GC-2050 is equipped with the optional gas selector, the carrier gas can be switched using the analysis method or via a Chromeleon window operation. Switching to an alternative carrier gas after sequence analysis is complete minimizes helium consumption.



The hydrogen sensor monitors the inside of the GC oven.

## Product Lineup

Description	Version	Versions with Functionality Verified
Shimadzu GC Driver for Chromeleon 7	2.40	Chromeleon 7.2 SR1 - 7.2 SR5 / 7.2.8 / 7.2.9 / 7.2.10 / 7.3 / 7.3.1 / 7.3.2

Note: This product supports Nexis GC-2030, Brevis GC-2050 and HS-20/HS-10 system control and HS-20 standalone control. GC-2010 (Plus/Pro) and GC-2014 (c) can be controlled using the Thermo Fisher Scientific Inc. driver provided with Chromeleon 7.

## Controllable Hardware

The Shimadzu GC Driver for Chromeleon 7 supports control of the following units.

**GC Unit** Nexis GC-2030, Brevis GC-2050

**Options** AOC-30i, AOC-20i (Plus) autoinjector, AOC-20s (U) autosampler, HS-20 (NX)/HS-10 headspace sampler, dual injection system

Unit	Device name
Sample Injector	GC-2030 : SPL-2030, WBI-2030, OCI-2030 (NX), PTV-2030, SINJ-2030
	GC-2050 : SPL-U(1.0)
Detector	GC-2030 : FID-2030, TCD-2030, ECD-2010 Exceed, FPD-2030, FTD-2030, BID-2030, SCD-2030, PTC-2030
	GC-2050 : FID-U(1.0), FPD-U(1.0), ECD-2010 Exceed U
Advanced Flow Technology	GC-2030 : Backflush, detector splitting, detector switching, heart-cut system
Additional Temperature Controller	GC-2030 : Auxiliary temperature control unit
Additional Flow Controller	GC-2030 : APC (3 auxiliary channels), APC (1 auxiliary channel)
	GC-2050 : APC (3 auxiliary channels), APC (1 auxiliary channel)
Options	GC-2030 : Gas selector, Low-temperature control solenoid valve set CRG-2030, External equipment control relay PRG-2010 Plus, PRG Box
	GC-2050 : Gas selector, Low-temperature control solenoid valve set CRG-2030

The following unit can be controlled using the Thermo Fisher Scientific driver included with Chromeleon 7.

By using the Shimadzu GC Driver for Chromeleon in conjunction with standalone control of the HS-20 headspace, the system consisting of both the GC and HS-20 (NX) can be controlled.

**GC Unit** GC-2010 (Plus) and GC-2014

**Options** HS-20 (NX) headspace sampler, AOC-20i (Plus) autoinjector, AOC-20s (U) autosampler

Unit	Device name
Sample Injection Port	GC-2010 (Plus) : SPL-2010 (Plus), OCI / PTV-2010 (Plus)
	GC-2014 : SPL-2014, DINJ-2014
Detector	GC-2010 (Plus) : FID-2010 (Plus), TCD-2010 (Plus), ECD-2010 (Plus), FPD-2010 (Plus), FTD-2010 (Plus)
	GC-2014 : DFID-2014, SFID-2014, TCD-2014, ECD-2014
Additional Flow Controller	GC-2010 (Plus) : APC (3 auxiliary channels)
Low-Temperature Oven Controller	GC-2010 (Plus) : CRG-2010 low-temperature control solenoid valve unit

Note: Up to six Shimadzu GC units can be controlled from a single instrument server.

It is not recommended to connect other vendor instruments to the same instrument server. Please prepare a dedicated server for controlling the Shimadzu GC units.

Brevis, AOC and Nexis are trademarks of Shimadzu Corporation or its affiliated companies in Japan and/or other countries.

Thermo Scientific, Dionex, Chromeleon and Thermo Fisher Scientific are trademarks of Thermo Fisher Scientific Inc. and its affiliated entities.



Shimadzu Corporation

[www.shimadzu.com/an/](http://www.shimadzu.com/an/)

**For Research Use Only. Not for use in diagnostic procedures.**

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

Company names, products/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation, its subsidiaries or its affiliates, whether or not they are used with trademark symbol "TM" or "®".

Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The contents of this publication are provided to you "as is" without warranty of any kind, and are subject to change without notice. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication.