

Experience New Benchmarks with OpenLab™

Shimadzu LC Driver Ver. 2.20, Compatible with OpenLab CDS




■ Enabling Control of Nexera™ lite inert (LC-40i/SIL-20A/SIL-20AC)

Nexera lite inert (LC-40i, SIL-20A, SIL-20AC) systems in the Shimadzu UHPLC Nexera series are now controlled from Agilent OpenLab CDS and OpenLab EZChrom. The cutting-edge functions of the Nexera series can be experienced first-hand with OpenLab. These include the bioinert system, intended to limit the adsorption of sample metals by configuring sample flow lines with nonmetals. The Shimadzu LC Driver was developed to comply with the Agilent RC.net open standard. Shimadzu LC instrument parameters can be configured in the OpenLab instrument method editing window and then saved as an instrument method.



Status Window



Method Editing Window

■ Product Type

	Description	Remarks
1	Shimadzu LC Driver, OpenLab CDS Compatible, 1 License	To install the driver in an existing OpenLab system, please provide the version of applicable software and other relevant information in advance.
2	Shimadzu LC Driver, OpenLab EZChrom Compatible, 1 License	

■ ■ List of Compatible Instruments

Nexera Series

Unit	Description
Solvent delivery unit	LC-40i
Autosampler	SIL-20A / SIL-20AC
Chemical reaction tank	CRB-40
High-pressure flow line switching valve	FCV-BOX / FCV-S

■ ■ List of Compatible Instruments

Series	Major Unit
Nexera series	SCL-40 / LC-40D XR / SIL-40C XR / CTO-40C / SPD-40 / SPD-M40
Nexera / Prominence™	CBM-20A / LC-30AD / SIL-30AC / CTO-20AC / SPD-20A / SPD-M30A
i-Series LC-2050/LC-2060	LC-2050 / LC-2060
i-Series	LC-2030 / LC-2040 / LC-2030 Plus / LC-2040 Plus

* Up to 4 Shimadzu LC systems and up to 2 PDA detectors can be controlled from a single computer or acquisition server (such as AIC). However, if 1 or 2 PDA detectors are connected, the maximum number of LC systems is reduced to 3 or 2 respectively. Additionally, if the SPD-M30A is used for data acquisition at high sampling speeds of 100 Hz or faster, only a single PDA detector can be connected.

* A Shimadzu LC/GC system and another brand of LC system cannot be connected to the same computer or AIC at the same time. Provide a separate computer or AIC dedicated for the Shimadzu LC/GC system.

Nexera and Prominence are trademarks of Shimadzu Corporation or its affiliated companies in Japan and/or other countries.
OpenLab is a trademark or a registered trademark of Agilent Technologies, Inc. in the United States and/or other countries.

First Edition: February, 2023



Shimadzu Corporation
www.shimadzu.com/an/

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. Company names, products/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation, its subsidiaries or its affiliates, whether or not they are used with trademark symbol "TM" or "®". Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

© Shimadzu Corporation, 2023
3655-02307-PDFIK