

The CLAM-2000 Challenge: Checking the Accuracy and Reproducibility of Automatic Sample Preparation of Steroids and Glucocorticoids in Plasma

MSACL 2017 #B01

Jana Rykl (1), Alessandro Visco (1), Michael Hayoz (2),
Yolanda Aebi (2), Michael Deuber (2), Cedric Bovet (2),
Jean-Christophe Prost (2)
(1) Shimadzu Switzerland,
(2) Inselspital Bern

The CLAM-2000 Challenge: Checking the Accuracy and Reproducibility of Automatic Sample Preparation of Steroids and Glucocorticoids in Plasma

Introduction

Manual sample preparation of clinical samples is a time consuming work. The CLAM-2000 robot is fully automated sample preparation module for LCMS. It can be used with commercial kits and homebrew methods. So far it intended use is research only.

We selected two steroid hormones (Cortisol, Testosterone) and two glucocorticoid drugs (Dexamethasone,

Prednisolone) for closer investigation of the efficiency of automatic sample preparation.

Cortisol and testosterone are hormones routinely tested in new born screening. Dexamethasone and Prednisolone treat many inflammatory and autoimmune conditions , both substances are banned under WADA anti-doping rules.




Fig 1 CLAM-2000 configuration

Materials and Methods

Manual Sample Prep:

(1) To 100 μ L plasma sample (with spiked compounds) add 200 μ L precipitation mix (containing internal standard), (2) Mix for 60 sec on vortex mix, (3) Centrifugation 5 min at 12000 g, (4) Transfer

supernatant to autosampler vial and inject 20 μ L into LCMS-8045 triple quadrupole coupled to an Nexera-X2 UHPLC. Analysis of triplicates (same sample preparation) in MRM spectrum mode using internal standards for

Automatic Sample Prep using the CLAM-2000:

(1) Conditioning the filtration vial with 20 μ L of methanol, (2) 50 μ L plasma sample (with spiked compounds) added to filtration vial, (3) 100 μ L precipitation mix (containing internal standard) added to filtration vial, (4) Shaking for 60 sec, (5) Filtration for 90

sec, (6) Transfer flowthrough from collection vial to HPLC autosampler and inject 20 μ L into the LCMS-8045 triple quadrupole to an Nexera-X2 UHPLC. Analysis of triplicates (separate sample preparation) in MRM spectrum mode using internal standards for confirmation.

The CLAM-2000 Challenge: Checking the Accuracy and Reproducibility of Automatic Sample Preparation of Steroids and Glucocorticoids in Plasma

Results and discussion




Fig. 2 Chromatograms: left manual prep, right: automatic sample prep. The filtration step during the automatic CLAM sample preparation seems to give sharper peaks for plasma samples.

Surprisingly we found sharper peak shapes for the chromatograms of all compounds using the automatic sample preparation (Figure 2). The check of the reproducibility using control samples showed more variation for automatic sample preparation, this is due to the fact that the CLAM-2000 robot prepares real independent technical replicates.

Table 1: Reproducibility check using five injections of a control sample

	Cortisol manual Conc. [ng/mL]	Cortisol automatic Conc. [ng/mL]
Average	602.28	614.75
%RSD	4.35	14.83
Maximum	637.98	688.82
Minimum	559.60	459.17
Std. Dev.	26.19	91.15

	Dexamethasone manual Conc. [ng/mL]	Dexamethasone automatic Conc. [ng/mL]
Average	416.02	471.76
%RSD	13.07	25.71
Maximum	518.25	659.82
Minimum	370.08	342.86
Std. Dev.	54.36	121.29

	Testosterone manual Conc. [ng/mL]	Testosterone automatic Conc. [ng/mL]
Average	5.18	5.25
%RSD	15.44	11.20
Maximum	6.39	5.98
Minimum	4.30	4.66
Std. Dev.	0.80	0.59

	Prednisolone manual Conc. [ng/mL]	Prednisolone automatic Conc. [ng/mL]
Average	392.00	471.77
%RSD	12.61	18.27
Maximum	485.20	580.61
Minimum	349.55	384.32
Std. Dev.	49.44	86.19

The CLAM-2000 Challenge: Checking the Accuracy and Reproducibility of Automatic Sample Preparation of Steroids and Glucocorticoids in Plasma

The statistical Passing-Bablok test did not show a significant difference in the results between manual and automatic sample preparation as shown in Figure 3. We found the same LOQ and the same linearity for all four compounds. Additionally automated sample preparation showed less variation of area of the internal standards.




Fig 3 Passing-Bablok test
 Red dotted line: 95% Confidence Interval
 Black line: Linear regression of results CLAM (y-axis) versus manual sample prep (x-axis)

The CLAM-2000 Challenge: Checking the Accuracy and Reproducibility of Automatic Sample Preparation of Steroids and Glucocorticoids in Plasma

Summary

Automatic sample preparation using the CLAM-2000 robot for LCMS gives similar results compared to manual preparation, but reduces the risk of human error and also the consumption of reagents. The filtration step in the automatic sample preparation seem to improve chromatographic peak shape.

First Edition: November, 2017



Shimadzu Corporation
www.shimadzu.com/an/

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. Company names, products/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation, its subsidiaries or its affiliates, whether or not they are used with trademark symbol "TM" or "®".

Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.