

# Eliminating the Fear Factor

## Flame Ionization Detector

Agilent Technologies



# The Flagship GCs from Agilent



8860




8890



Intuvo


# Agilent 8860/8890 Series GC Features – Detectors



- AUX 1 Detector (TOP) can be TCD or FPD+
- AUX 2 Detector (SIDE) can be TCD, FID, or ECD

**YES!** This means you can have up to 4 detectors on your 8890 GC or 3 detectors on your 8860 GC!

# Agilent Intuvo GC Features – Detectors




# Common GC Detectors

- Flame Ionization Detector (FID)\*
- Thermal Conductivity Detector (TCD)\*
- Nitrogen-Phosphorous Detector (NPD)\*
- Electron Capture Detector (ECD)
- Flame Photometric Detector (FPD)
- Chemiluminescence Detectors (XCD)
- Mass Selective Detector (MSD)


Today's topic is  
the FID!

# GC Detectors


Let's focus on the FID



Linear Dynamic Range  
Response increases *proportionally*  
with increased quantity.



# Flame Ionization Mechanism




- The FID is a **DESTRUCTIVE**, mass sensing detector
- Cations generated in the flame are counted and produce the detector signal
- Analytes that have the greatest number of low oxidation state carbons produce the **LARGEST** signal

## SELECTIVITY

- The FID will respond to **ORGANIC** compounds (C-H Bonds)
- There is little or NO response for the following:
  - Inert Gases
  - $H_2O$
  - $CO_2$
  - $CO$
  - $N_2$
  - $O_2$
  - $CS_2$
  - Heavily halogenated compounds


# Agilent FID Theory of Operation

- H<sub>2</sub> / Air = Flame
- Sample is burned producing positively charged ions
- Positive Ions are attracted to the negative voltage on the collector
- Collector current is converted to FID output in the Electrometer





# FID Detector Cutaway



# FID Detector Bases


- 6890/7890 GC –  
Capillary optimized FID base  
OR  
Packed FID base with  
capillary adapter included
- Intuvo GC –  
Capillary optimized FID base only  
-no packed column capabilities on  
Intuvo
- 8860/8890 GC –  
Capillary optimized FID base only  
A packed column adapter is  
included

Capillary optimized fitting



Capillary columns  
only


Adaptable fitting



Adapters are available  
for capillary columns or  
1/8" Packed Columns

1/4" Packed columns  
can be installed  
without an adapter

# Agilent 6890/7890/8860/8890/Intuvo FID Jets



Adaptable FID Jets – 6890/7890		
Jet Type	Part#	Jet Tip ID
Capillary	19244-80560	0.29 mm 0.011 in.
Packed	18710-20119	0.47 mm 0.018 in.
Packed Wide Bore	18789-80070	0.79 mm 0.030 in.
High Temp	G1531-80620	0.47 mm 0.018 in.

Capillary-Optimized FID Jets – 6890/7890		
Jet Type	Part#	Jet Tip ID
Capillary	G1531-80560	0.29 mm 0.011 in.
High Temp	G1531-80620	0.47 mm 0.018 in.



Intuvo/8860/8890 FID Jet		
Jet Type	Part#	Jet Tip ID
Jet	G4591-20320	0.011 in.




**Agilent's NEW  
FID Jet Design!**


Backwards  
Compatible back  
to 6890 GC

# FID Setup - Column Installation - FID

## Column Installation — From Troubleshooting and Maintenance Manuals



# FID – Optimizing Gas Flows for Sensitivity



## Gas Flow Guidelines:

Hydrogen 30-35 ml/minute

Carrier + Makeup 30-35 ml/minute


Air – 400 ml/minute

Nitrogen is the recommended makeup gas

H<sub>2</sub>-to-Air ratio should be between 8-12%

Total Inert Gas should have ~1:1 ratio to H<sub>2</sub>

## Method Editor Menu in Software



AKA: Column + Fuel = Constant

Select if using **Hydrogen Carrier**

The detector sees Constant Hydrogen Flow if carrier flow changes.

AKA: Column + Makeup = Constant

Best for Column – “**Constant Pressure**” Mode

Makeup is ramped during the run to compensate for dropping carrier flow – Detector sees constant Column + Makeup flow

AKA: Constant Makeup & Fuel Flow

Best for “**Constant Flow**” mode on the Column (He or N2)

Makeup Gas and Hydrogen Control are independent of the Carrier and remain constant throughout the analysis

- It is best to operate the FID at 300°C or hotter – But always 20°C hotter than the highest GC Oven Temp

### Gas Flow Guidelines:


Hydrogen 30-35 ml/minute

Carrier + Makeup 30-35 ml/minute (≈1:1 ratio to H2)

Air – 400 ml/minute

- Optimizing Flows 1:1 ratio of Carrier + Makeup to H2 flow
  - Nitrogen is the recommended makeup gas

- Make sure you check the base of the FID for any chunks of graphite when changing the column or jet to prevent those chunks from getting into the jet and causing a partial blockage
- Check your Lit Offset. Default is 2pA and sometimes with very clean gases the baseline will drop below 2pA causing your GC to think the flame is not it.
- Use the right **jet** for your application
  - High oven temps, bleeding columns, or clogging jets often – try the High Temp FID Jet



Wider opening to help prevent clogging

# FID Troubleshooting

## Flame Won't Light - Jet Diagnostic Test – 7890 Series GC

BACK DETECTOR (FID)		
H2 flow	75.00	75.00
Air flow	0.131	off<
Makeup (N2)	0.135	off

Pass


BACK DETECTOR (FID)		
H2 flow	75.00	75.00
Air flow	0.131	off<
Makeup (N2)	7.315	off

Fail

- Set H2 Flow to 75 ml/minute
- Turn the Makeup Gas “OFF”
- Monitor Makeup Gas “Actual Flow”
  - Low flow reading/minimal change – Jet is clear!
  - High Flow reading/big change - indicates Jet Plugging...



# FID Jet Restriction Test – How does it Work?



# Jet Diagnostic – Built into 8860/8890/Intuvo GC

Now the GC has the ability to run the Jet Restriction test from the front panel and generate a pass or fail.

The screenshot shows the Agilent GC Diagnostics interface. The top navigation bar includes 'Method', 'Sequences', 'DA Express', 'Diagnostics' (selected), 'Maintenance', 'Logs', 'Settings', and 'Help'. Below this, there are sub-sections: 'Warnings And Errors', 'Diagnostic Tests' (selected), 'System Health Report', and 'Detector Evaluation Reports'. The main content area is divided into three sections: 'Aux Detector 1 (2)', 'Back Inlet (5)', and 'Front Detector (2)'. The 'Front Detector (2)' section contains a list of tests: 'Jet Restriction Test' (highlighted with a red box), and 'Leakage Current Test'. At the bottom, there is a 'STATUS: READY' indicator and a table with columns for 'Sequence', 'Method', 'Sample Name', and 'Est. Remaining' (showing 8.17).

## Front Detector : Jet Restriction Test


Test Passed

Ok will update the tests info in the System Health Report

A quick, easy, un-intrusive test that can be run from the GC or the browser interface.

# FID Troubleshooting

## Baseline Spiking




### Possible Causes:

- Contamination of the jet, collector, and/or the base by particulates, stationary phase, graphite
  - Clean and Sonicate or Replace the Jet
  - Thoroughly clean the collector and chimney assembly
  - Thoroughly clean the base of the FID – especially if ferrule chunks are visible
- Electronic Noise – much less likely
  - Vibration
  - Electrical – Grounding problems, Power line disturbances, Power Line conditioners
  - Strong RF Interference

# FID Troubleshooting

## Baseline Drift and Wander




### Possible Causes:

- Contaminated/Saturated Gas Supply Traps
- Buildup of high MW sample contamination in the GC Column
- Column Bleed
- Contamination in the Base of the Detector
- Gas Leaks
- Column Fitting Leaks

New!  
Smart  
Gas Clean  
Sensor



# Baseline Evaluation – Intelligent GCs



The screenshot shows the software interface for configuring a detector. The top navigation bar includes 'Method', 'Sequences', 'DA Express', 'Diagnostics', 'Maintenance', 'Logs', 'Settings', and 'Help'. The left sidebar lists various components, with 'Front Detector - FID' selected under 'Detectors'. The main area is titled 'Blank Evaluation Setpoints - FID' and contains the following settings:

- Enable Blank Evaluation  
Requires detector signal(s) to be saved on Signals tab.  
**Load Default Setpoints**
- Initial Baseline**
  - Minimum Value (pA)\*: 8
  - Maximum Value (pA)\*: 11
  - Noise (pA)\*: 0.5
- Final Baseline**
  - Minimum Value (pA)\*: 8
  - Maximum Value (pA)\*: 9
  - Noise (pA)\*: 1
- Baseline evaluation**
  - Total peak area (pA\*sec)\*: 100
  - Maximum peak height (pA)\*: 3
  - Start time (min)\*: 0.5
  - End time (min)\*: 7

At the bottom, a status bar shows 'STATUS: READY' and a table with the following data:

Sequence	Method	Sample Name	Est. Remaining
			8.17

- Blank Evaluation allows you to set up expected signal levels at the start and end of the run and confirms they are as expected.
- When you run a blank, your smart GC monitors the blank and returns either a log entry if the blank is not blank, or stops the sequence.

# FID Troubleshooting

## Measuring Flows



### Gather your flowmeter and adapter tube first!

- Set oven temp to Ambient (35C)
- Turn off column flow
- Shut off all detector gases
- Turn off FID flame
- Cool detector
- Connect adapter to detector by inserting into FID as far as possible
- Turn on one detector flow at a time allowing the gas to equilibrate in the flowmeter tubing before taking a measurement



A separate adapter (19301-60660) is supplied for the FID. Insert the adapter into the detector exhaust vent as far as possible. You will feel resistance as the adapter O-ring is forced into the detector exhaust vent. Twist and push the adapter during insertion to ensure a good seal.





## Cleaning Kit: G1531-67000



Includes Small Parts:

- Teflon Collector Insulators
- Silicone Rubber Base Gasket
- Ignitor
- Cleaning brushes for Collector
- Cleaning Wires for Jet

## Rebuilding Kit: G1531-67001

Entire Collector Assembly for  
Extremely corrosive applications

**Note:** The FID Jet must always be ordered separately




Now let's step through a complete FID  
Maintenance/Rebuild Procedure together!



# FID Maintenance Procedures

## Ignitor Replacement




Disconnect the Ignitor



Remove with 5/16" wrench


Install the new Ignitor with  
Copper Washer –  
Tighten with 5/16" wrench





# FID Maintenance Procedures

## Collector Removal



### Loosen 3-T20 Screws

A common mistake people make is NOT removing the three screws **FIRST**. This makes getting to the jet quite difficult and you can damage the interconnect assembly and spring. Make sure you remove the screws first.



Pull the Collector Assembly straight up

# FID Maintenance Procedures

## Collector Disassembly and Cleaning



**Unscrew the knurled nut**



**Remove the castle**



**Remove the collector and  
Upper/lower PTFE Insulators**




**Remove/replace the  
silicone rubber gasket**

# FID Maintenance Procedures

## Jet Removal, Collector Removed




FID Jet



Using a 1/4" Nut Driver,  
Loosen the Jet

Remove the Jet

Be careful not to damage the  
Interconnect Spring




# FID Maintenance Procedures

## Check for Graphite Chunks in the Detector Base

**Indicator** ....Was the Column hard to remove or to re-install?



Shine a flashlight up  
From the oven to check  
For impacted graphite



Clean out the Detector Base – Top down  
A syringe plunger works well



Look for graphite chunks  
falling into the oven




# FID Maintenance Procedures

## Jet Cleaning or Replacement



**Jet sealing surface – if this is Compromised, replace the jet**

**Agilent recommends replacing the Jet**



**If the jet is plugged you can clean with a .010" cleaning wire**

**Jet Cleaning Wire (0.010") 9301-0985**

# FID Maintenance Procedures



## Jet Installation



FID Jet



Interconnect  
Spring



Using a 1/4" Nut Driver,  
Tighten the FID Jet 1/6<sup>th</sup> of  
A turn past Hand Tight  
(For a New Jet only)


# FID Maintenance Procedures

## Re-install the Collector



Position the cleaned or replacement collector directly over the detector base and push straight down.

The interconnect spring will “pop” into place!



Tighten 3-T20 Screws evenly and firmly

Re-connect Ignitor electrical cable

# Intelligent GC Features – Self Guided Maintenance

The screenshot displays the Agilent CrossLab software interface. The top navigation bar includes tabs for Method, Sequences, DA Express, Diagnostics, Maintenance (highlighted), Logs, Settings, and Help. A modal dialog box titled "Front Detector : Replace FID Jet" is open, showing a "Cancel" button in the top right. The dialog content includes the text "The instrument will now verify and finalize the maintenance task." followed by a bulleted list: "Reset the Early Maintenance Feedback counter", "Restore active method", and "Reinstall top oven cover". At the bottom of the dialog, it indicates "Step 12 of 12" and "Establish a flow of carrier gas", with "Back" and "Next" buttons. The background interface shows a "Perform Maintenance" button and a list of parts with checkmarks: "Collector ass...", "Ignitor igniti...", and "Jet injections...". The bottom status bar shows "STATUS: MAINTENANCE MODE" and a table with columns: Sequence, Method, Sample Name, and Est. Remaining (8.17).

These on-board procedures take you through common maintenance operations, guiding you step-by-step, cooling the instrument zones that are necessary, and even updating logs and EMF counters when finished!

These are accessed by way of the “perform maintenance” link in the maintenance tab.





For more information regarding Agilent's Gas Chromatography, please visit the website at:  
<https://www.agilent.com/en/products/gas-chromatography>

Videos for Simple Maintenance and Theory can be found on the Agilent YouTube page:  
<https://www.youtube.com/user/agilent/>

Checkout Agilent Communities where you can get answers and share insights:  
<https://community.agilent.com/welcome>

# A Sample of Customer Education Courses – GC

## e-Learning



Agilent University Course Code	Course Name
GC-MULTI-2100e	Advanced Operation of the Multi Mode Inlet (MMI)
GC-0GEN-1012s	GC Inlets Theory and Operation
GC-0GEN-1013s	GC Detectors Theory and Operation
SI-7696-1100s	Agilent 7696A Sample Prep WorkBench Operation
GC-MULTI-1240zs	Making Productivity Happen: an Agilent GC eLearning Series
SI-7693-1100s	Agilent 7693A ALS Basic and Advanced Operation and Maintenance
SI-7697-2100fs	Running Start for Agilent 7697A Headspace Sampler

## Classroom



Agilent University Course Code	Course Name
GC-0GEN-2000V2	Practical Gas Chromatography – 4 Day (R1915A)
GC-7890-2100cV4	Agilent 7890A/B GC and OpenLAB ChemStation Operation – 5 Day (R1778A)
GC-7890-2200cV3	Agilent 7890A/B GC Maintenance & Troubleshooting – 4 Day (R1914A)
GC-9000-2101cV3	Agilent 7890/9000 GC with OpenLAB 2.3 Essential and Advanced Operation – 4 Days
GC-9000-2103c	Agilent Intuvo 9000 GC with OpenLAB CDS ChemStation Operation – 4 Days
GC-8890-2100c	Agilent 8890 GC Operation with OpenLAB CDS ChemStation Edition – 4 Days

Enroll at Agilent University: <http://www.agilent.com/crosslab/university/>