

# High Throughput Quantitative Analysis of Multi-mycotoxin in Beer-based Drinks using UHPLC-MS/MS

**ASMS 2012** WP27-580

Masayoshi Tamura<sup>1</sup>; Keiko Matsumoto<sup>2</sup>; Jun  
Watanabe<sup>2</sup>; Junko Iida<sup>2</sup>; Yasushi Nagatomi<sup>1</sup>; Naoki  
Mochizuki<sup>1</sup>

<sup>1</sup>Asahi Group Holdings, Ibaraki, JAPAN;

<sup>2</sup>SHIMADZU CORPORATION, Kyoto, JAPAN

## High Throughput Quantitative Analysis of Multi-mycotoxin in Beer-based Drinks using UHPLC-MS/MS

### Introduction

Mycotoxins often exist as contaminants in grains. To ensure consumer food safety, manufactures of food and beverages have to strictly manage risks from such contaminants. To maintain the high-quality of food standards it is therefore essential to rapidly determine the concentrations of hazardous mycotoxins in foods or beverages.

UHPLC-MS/MS offers the best combination of selectivity, sensitivity, and speed for detection of these compounds in

complex matrices. In this study, a high throughput method for the quantification of 14 mycotoxins in beers was developed. Highest sensitivity of analysis is crucial to food safety, additionally, autosampler and system carry over need to be monitored to ensure these factors do not become a problem. In these experiments elimination of carry over was investigated through novel rinse condition cycles of the UHPLC autosampler.




Fig. 1 Structure of mycotoxins

### Methods and Materials

14 mycotoxins (patulin(PAT), nivalenol(NIV), deoxynivalenol (DON), aflatoxin(AF) B1, B2, G1, G2, T-2 toxin(T-2), HT-2 toxin(HT-2), zearalenone(ZON), fumonisin(FM) B1, B2, B3 and ochratoxin A(OTA)) were determined by LC-MS/MS using a UFLC HPLC system coupled to a LCMS-8030 triple quadrupole mass spectrometer.

The MRM method of 14 mycotoxins was optimized on

each compound-dependent parameter and MRM transition (Q1/Q3). As a result, all compounds were detected with high sensitivity by ESI. AFB1, B2, G1, G2, T-2, HT-2, FMB1, B2, B3 and OTA were detected in positive mode. While PAT, NIV, DON, ZON were detected in negative mode. Ultra Fast Polarity Switching of 15 msec enabled simultaneous determination of the compounds in both modes.



Fig. 2 LCMS-8030 triple quadrupole mass spectrometer

#### High Speed Mass Spectrometer

- Polarity Switching  
15 msec
- Scanning Speed  
Max. 15000 u/sec

# High Throughput Quantitative Analysis of Multi-mycotoxin in Beer-based Drinks using UHPLC-MS/MS

## Results

### Method development for 14 mycotoxins

For UHPLC separation, various LC mobile phase conditions were examined. Tailing of fumonisins peaks were observed when only ammonium acetate was added in mobile phase. It was found that pH of a mobile phase effected peak shape of fumonisins. In order to reduce tailing of fumonisins, acetic acid was added in mobile phase B and

the gradient program was controlled to maintain high concentration of acetic acid when fumonisins were eluted. By controlling the concentration of acetic acid and ammonium acetate with gradient program, 14 mycotoxins were separated and detected excellently in 11 minutes (Fig. 3).

### Analytical Conditions for LC-MS/MS

#### HPLC: UFLC system

Column: TriartC18 100 mm×2.0 mm, 1.9 μm  
 Mobile phase A: 10 mM Ammonium acetate - Water,  
 B: 2% Acetic acid - Methanol  
 Flow rate: 0.4 mL/min  
 Gradient program: B conc.2%(0 min) - 55%(3 min) - 85%(7.0-8.0 min) - 2%(8.01-11 min)  
 Column temperature: Column temperature: 40 C

#### MS: LCMS-8030 triple quadrupole mass spectrometer

Ionization: ESI, Positive/Negative MRM mode Ion spray voltage: -3.5 kV




Fig. 3 14 mycotoxins analysis by LC-MS/MS (PAT, NIV, DON, HT-2, T-2, OTA, ZON, FM B1/B2/B3 50ppb, AF B1/B2/G1/G2 10 ppb)

#### MRM Transition

Mycotoxin	MRM transition
AF G1 (+)	329.05 > 243.05
AF G2 (+)	331.00 > 245.00
AF B1 (+)	313.00 > 241.05
AF B2 (+)	315.00 > 259.00
HT-2 (+)	442.00 > 263.05 ([M+NH4] <sup>+</sup> )
T-2 (+)	483.95 > 305.00 ([M+NH4] <sup>+</sup> )
OTA (+)	404.10 > 238.90
ZON (-)	317.15 > 273.00
NIV (-)	371.10 > 281.25 ([M+CH3COO] <sup>-</sup> )
DON (-)	355.10 > 295.15 ([M+CH3COO] <sup>-</sup> )
PAT (-)	153.10 > 109.20
FM B1 (+)	722.45 > 334.30
FM B2 (+)	706.45 > 336.25
FM B3 (+)	706.45 > 336.25

Each mycotoxin standard was analyzed at six concentration levels. Good linearity was observed in the calibration curves, and excellent sensitivity was achieved.

## High Throughput Quantitative Analysis of Multi-mycotoxin in Beer-based Drinks using UHPLC-MS/MS

Table 1 Linearity 14 mycotoxins

Mycotoxin	Range	Coefficient(r2)	Mycotoxin	Range	Coefficient(r2)
AF G1	0.4-20 ppb	0.999	ZON	2-100ppb	0.999
AF G2	0.4-20 ppb	0.999	NIV	2-100ppb	0.999
AF B1	0.4-20 ppb	0.999	DON	2-100ppb	0.997
AF B2	0.4-20 ppb	0.999	PAT	10-100ppb	0.999
HT-2	2-100 ppb	0.998	FM B1	2-100ppb	0.995
T-2	2-100 ppb	0.999	FM B2	2-100ppb	0.994
OTA	2-100 ppb	0.999	FM B3	2-100ppb	0.997

### Rinse condition for eliminating carry over

Carry over of fumonisins was initially observed using the general rinse condition, because fumonisins formed complexes with trace metal ions in the sample's flow path.

Probably, several carboxyl groups of fumonisins coordinated with metal ion (Fig. 4).




Fig. 4 Possible coordination interaction with metal ion

For eliminating carry over, rinse solvent and rinse method were examined. The performance of Nexera autosampler SIL-30AC, which can wash both inner and outer needle surfaces with 4 different solvents, was used.

It was thought that carboxyl groups of fumonisins may preferentially pair with hydrogen ions in the presence of low pH, therefore formic acid was added to rinse solvent. When investigating rinse methods, it was discovered the inner and outer rinse of needle reduced carry over more than the outer rinse of needle. Finally the modified rinse

solvent consisted of: 1% formic acid aq./methanol / acetonitrile / isopropanol (1/1/1/1).

To test the modified rinse cycle method, one injection of the 100ppb fumonisins standard solution was followed by one blank injection to check for carryover. Figure 5 shows chromatograms of the standards of FMB2 and B3, and the following blank injection. Low carry over was observed in the blank injection. It resulted from washing fumonisins adsorbed inside needle with the needle's inner and outer rinse method and the effective rinse solvent.

## High Throughput Quantitative Analysis of Multi-mycotoxin in Beer-based Drinks using UHPLC-MS/MS




Fig. 5 Carry over evaluation of fumonisins




Fig. 6 HPLC path of SIL-30AC

## Quantitative Analysis of 14 mycotoxins in beer-based drinks

Mycotoxins were extracted from samples and were purified with a solid phase extraction (SPE) cartridge. 20 commercial beers were analyzed by using this method. The calibration curves were assessed using beer samples spiked with mycotoxins. PAT, AFB1, B2, G1, G2, NIV, T-2

and ZON were not detected in any of the beer samples. Some of the tested samples were found to be contaminated with DON, HT-2, OTA, FMB1, B2, and B3 at concentrations of less than their respective LOQs (each 5 ppb).

Table 2 Mycotoxins detected in analyzed samples

Producing country	Concentration of mycotoxin/ppb (detected rate)					
	DON	HT-2	FMB1	FMB2	FMB3	OTA
Mexico (1sample)			< 5 (1/1)	< 5 (1/1)		
USA (1sample)	< 5 (1/1)					
China (1sample)	< 5 (1/1)					
Philippine (1sample)			< 5 (1/1)			
Australia (1sample)						
Japan (5sample)			< 5 (3/5)	< 5 (2/5)		

Producing country	Concentration of mycotoxin/ppb (detected rate)					
	DON	HT-2	FMB1	FMB2	FMB3	OTA
Hol I and (2samples)		< 5 (1/2)				< 5 (1/2)
Ireland (2sample)						
England (1sample)	< 5 (1/1)		< 5 (1/1)			< 5 (1/1)
Germany (1sample)						
Czech (1sample)						
Belguin (2samples)	6.7 (1/2)	< 5 (1/2)				
	< 5 (1/2)		< 5 (1/2)	< 5 (1/2)		

< 5 (less than 5ppb)

## Conclusions

- High throughput LC-MS/MS method for 14 mycotoxins was developed, and could be applied to the quantification of these compounds in beers.
- Carry over of fumonisins was eliminated by using both the needle's inner and outer rinse method with effective rinse solvent.
- Results from these experiments indicate that the health risk to consumers posed by intake mycotoxins in commercial beers is relatively low.

## Acknowledgements

Reagents were provided from Wako Pure Chemical Ind., (Osaka, Japan) with substantial cooperation.



Shimadzu Corporation  
[www.shimadzu.com/an/](http://www.shimadzu.com/an/)

For Research Use Only. Not for use in diagnostic procedures.  
The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.