

MS_n Analyses for Tryptophan-Conjugated ADC Mimic by Miniature MALDI Digital Ion Trap Mass Spectrometer (MALDI-DIT-MS)

ASMS 2019 ThP409

Hideharu Shichi¹, Shuichi Nakaya¹, Katsuya Maruyama²,
Kosuke Hosoi¹, Takashi Nishikaze¹, Koichi Kojima¹,
Kei Kodera¹, Sadanori Sekiya¹, Shinichi Iwamoto¹,
Kounosuke Oisaki², Motomu Kanai², Koichi Tanaka¹
1 Shimadzu Corporation. Kyoto, Japan
2 Graduate School of Pharmaceutical Sciences,
The University of Tokyo, Bunkyo, Japan

MSn Analyses for Tryptophan-Conjugated ADC Mimic by Miniature MALDI Digital Ion Trap Mass Spectrometer (MALDI-DIT-MS)

Overview

Analysis of modification sites of chemically modified antibody Using MALDImini-1 Compact MALDI Digital Ion Trap Mass Spectrometer

Introduction

Antibody drug conjugates (ADCs), in which drugs are conjugated to antibodies, combine the high selectivity of antibodies with the effects of small molecule drugs to make it more effective than conventional low molecular weight drugs. It appeared in the 2000s as a drug expected to have cancer effects.

When artificially binding another compound to a protein,

such as ADC, the binding degree of that compound and its binding site become one of the critical quality properties.. In this experiment, we created a pseudo-ADC in which a low molecular weight compound was artificially bound to a research standard antibody, and modified group binding site analysis was attempted using the compact MALDI-DIT mass spectrometry "MALDImini™ -1." (Fig.1)



Figure 1 MALDImini-1 MALDI DIT MS

MSn Analyses for Tryptophan-Conjugated ADC Mimic by Miniature MALDI Digital Ion Trap Mass Spectrometer (MALDI-DIT-MS)

Methods

Standard antibody (NISTmAb, Humanized IgGκ monoclonal antibody, RM8671) in which Me-fluorescein-ABNO has been modified to tryptophan residue according to Seki et al.[1] (Fig.2), and untreated standard antibody (10 μg each) respectively internally tryptic digests and desalted with a Ziptip ® μC18 chip were loaded onto a MALDI target

plate. Furthermore, after accumulating and drying the matrix solution (0.5 μL), MSn analysis was performed using the compact MALDI-DIT mass spectrometry "MALDImini-1". The matrix used was DHB (2,5-dihydroxybenzoic acid).

[1] Seki, Y. et al.: J. Am. Chem. Soc. 2016, **138**, pp 10798-10801.




Figure 2 methyl fluorescein-ABNO conjugation to an antibody

Results

The modified and unmodified antibodies were measured using MALDImini-1, and their mass spectra were compared. Although almost all ions were common to both, some of the ions were detected only with the modified antibody (m/z 2416.9, 2430.7, 2452.7, 2560.9) (Fig. 3). Therefore, we decided to further analyze m/z 2416.9 and m/z 2560.9 that were considered to be from different molecules among the ions obtained only with the modified antibody. First, MS/MS analysis was performed on m/z 2416.9. As a

result, a fragment ion (m/z 1677.7) was detected (Fig.4), which appeared to be a peptide backbone after the modification group was eliminated. And this ion showed the same m/z value as the ion detected in the unmodified antibody. The MS³ spectrum of the modified antibody and the MS / MS spectrum of the unmodified antibody showed the same fragment pattern. Furthermore, when these data were searched by Mascot MS/MS, it was confirmed that it is a peptide sequence (²⁷⁸FNWYVDGVEVHNAK²⁹¹) derived from the heavy chain of the antibody.

MSn Analyses for Tryptophan-Conjugated ADC Mimic by Miniature MALDI Digital Ion Trap Mass Spectrometer (MALDI-DIT-MS)





Figure 3 MS spectra of trypsinized antibody (blue : modified, red : unmodified)

Similarly, MS/MS measurement was performed for the m/z 2560.9 ion. As compared the MS³ measurement result for the peptide backbone (m/z 1808.0) by MS/MS measurement of the modified antibody, with the MS/MS measurement results of the ions having the same m/z value by MS measurement of the unmodified antibody, both showed similar fragment patterns, and further, it was

found that the peptide sequence (³⁰⁵VVSVLTVLHQDWLNGK³²⁰) derived from the heavy chain of the antibody was obtained by MS/MS ions search (Fig.5). This suggests that the chemical modification used here is present in the tryptophan side chain in the above two peptide sequences.

MSn Analyses for Tryptophan-Conjugated ADC Mimic by Miniature MALDI Digital Ion Trap Mass Spectrometer (MALDI-DIT-MS)



- [IGHG1_HUMAN](#) Mass: 36083 Score: 58 Matches: 1(1) Sequences: 1(1)
 Ig gamma-1 chain C region OS=Homo sapiens GN=IGHG1 PE=1 SV=1
 Check to include this hit in error tolerant search or archive report

Query	Observed	Mr(expt)	Mr(calc)	Delta	Miss	Score	Expect	Rank	Unique	Peptide	
<input checked="" type="checkbox"/>	1	1677.7000	1676.6927	1676.7947	-0.1020	0	58	7.5e-005	1	U	K.FNWYVDGVEVHNAK.T

Figure 4 MSⁿ spectra of m/z 2416.9

MSn Analyses for Tryptophan-Conjugated ADC Mimic by Miniature MALDI Digital Ion Trap Mass Spectrometer (MALDI-DIT-MS)




Figure 5 MSⁿ spectra of m/z 2560.9

Amino acid sequence of the Heavy chain

QVTLRESGPALVKPTQLTLTCTFSGFSLSTAGMSVGVIRQPPGKALEWLADIWDDKKHYNPGLKDRLTISKDTSKNQVVLK
VTNMDPADTATYYCARDMIFNFYFDVWVGQTTVTVSSASTKGPSVFPLAPSSKSTSGGTAALGCLVKDYFPEPVTVSWNSGA
LTSGVHTFPAVLQSSGLYSLSSVTVPSSSLGTQTYICNVNHKPSNTKVDKRVPEKSCDKTHTCPPAPPELLGGPSVFLFPPK
KDTLMISRTPEVTCVVDVSHEDPEVKFNWYVDGVEVHNAKTKPREEQYNSTYRVVSVLTVLHQDWLNGKEYKCKVSNKAL
PAPIEKTISKAKGQPREPQVYTLPPSREEMTKNQVSLTCLVKGFYPSDIAVEWESNGQPENNYKTPPVLDSDGSFFLYSKLTV
DKSRWQQGNVFCSSVMHEALHNHYTQKLSLSLSPGK

Figure 6 Amino acid sequence of the heavy chain and position of the conjugation sides

MSⁿ Analyses for Tryptophan-Conjugated ADC Mimic by Miniature MALDI Digital Ion Trap Mass Spectrometer (MALDI-DIT-MS)

Conclusions

We analyzed an Tryptophan-Conjugated ADC mimic using the miniature MALDI-DIT-MS focusing on position and number of the conjugation sites. MSⁿ function of the miniature MALDI-DIT-MS could be a useful tool for quality assessment of next generation ADC with tryptophan conjugated.

Disclaimer: The products and applications in this presentation are intended for Research Use Only (RUO). Not for use in diagnostic procedures. Not available in China.

First Edition: June, 2019