


## Thermo Freestyle *Easily* Sends MSMS Spectra to NIST Search



**Uvinul 3000**

**MW 214**

*(UV Absorber/Light Stabilizer)*

[Name:](#) Methanone, (2,4-dihydroxyphenyl)phenyl-

[Formula:](#) C<sub>13</sub>H<sub>10</sub>O<sub>3</sub>

[MW:](#) 214 [Exact Mass:](#) 214.062994 [CAS#:](#) 131-56-6 [NIST#:](#) 492528 [ID#:](#) 172010 [DB:](#) mainlib

## Summary of Results


- Pos and Neg MSMS Spectra easily sent to NIST Search in accurate mass mode for searching
- Ion polarity imported into NIST search via MSP file
- Shows precursor  $m/z$  which is very useful
- See typical fields shown in NIST reference spectra below

[Name:](#) 2,4-Dihydroxybenzophenone  
[Precursor type:](#) [M+H]<sup>+</sup>  
[Instrument type:](#) HCD  
[Collision energy:](#) NCE=35% 15eV  
[Precursor m/z:](#) 215.0703  
[Formula:](#) C<sub>13</sub>H<sub>10</sub>O<sub>3</sub>  
[MW:](#) 214 [Exact Mass:](#) 214.062994 [CAS#:](#) 131-56-6 [NIST#:](#) 1190859 [ID#:](#) 81863 [DB:](#) hr\_msms\_nist  
[Comment:](#) NIST Mass Spectrometry Data Center  
[Notes:](#) Consensus spectrum; Acetonitrile/Water/Formic acid; Vial\_ID=3369; m<sub>z</sub>\_diff=-0.0004; Nreps=25  
[Spectrum type:](#) MS2  
[Instrument:](#) Thermo Finnigan Elite Orbitrap  
[Sample inlet:](#) direct flow injection  
[Ionization:](#) ESI  
[Collision gas:](#) N<sub>2</sub>  
[Ion mode:](#) P  
[InChIKey:](#) ZXDDPOHVAMWLBH-UHFFFAOYSA-N [Non-stereo](#)  
[4 m/z Values and Intensities:](#)  
77.0383 1.5 C<sub>6</sub>H<sub>5</sub>=p-C<sub>7</sub>H<sub>6</sub>O<sub>3</sub>/-3.6ppm 18/25  
105.0384 0.14 C<sub>7</sub>H<sub>5</sub>O 0.14001/0.0 05/05


## Approach for Qualitative Analyses in Thermo Orbitrap with Processing in Freestyle

- I really like the Orbitrap approach for qualitative analysis of compounds
- Very easy to use and quickly obtain useful data
- Acquires in Pos/Neg Mode
- Diode Array data
- Automatically determines major ions to fragment
- Collects the MSMS data at several different energies and stores in trap
- Then collects spectrum

# Uvinul 3000 Works Best in Negative Ion Mode Because Phenol




# Instrument Detects Ions and Then Automatically Decides Ones for MSMS Analyses User Looks Through List and Picks the One of Interest to Plot Vs. Time




User Finds  $m/z$  213 in Negative Ion Mode and Selects from List  
**NOTE** Positive and Negative Ions in Same List!  
*LOTS* of Ions in List!

# After User Select, Ion Chromatogram for M-1 Ion at $m/z$ 213 Plotted






# Can Obtain Diode Array Spectrum Shown on Bottom






# Export to NIST Button MSP Format




# Search Results in NIST MS Search of NIST Reference Library

#	Src.	Name
1	A	cc_2023_192 #3102 RT: 10.03 AV: 1...
2	A	-ESI Product Ion (5.155 min) Frag=5...
3	L	Fentanyl Carbamate (EI, SWGDRU...
4	L	Fentayl Analog (EI, SWGDRUG libr...
5	L	Amphetamine Analog (EI, SWGDR...
6	L	Methamphetamine Analog (EI, SW...
7	L	Cocaine (EI, SWGDRUG library)
8	L	Cortisone (EI, SWGDRUG library)



#	Lib.	Match	R.Match	Prob. (%)	RI	PSS.Match	DBs	Name
1	hr	657	786	56.4	705	28	...	2,4-Dihydroxybenzopheno
2	hr	565	674	4.73	565	10	...	7-Aminodesacetoxycepha
3	hr	560	683	4.10	575	...	...	NP-016409 [M-H]->[M-H-
4	M	552	679	3.25	2...	629	...	propanedinitrile, 2-[2-[4-(di
5	M	549	637	2.99	1...	678	...	4-tert-Butylphenol mesylate
6	M	541	559	2.38	2...	727	...	3-Methoxy-N-(6-methyl-pyr
7	M	539	610	2.25	1...	638	...	Fumaric acid, 2-isopropylp
8	M	538	630	2.19	2...	719	...	Isonicotinic acid, 1-adama
9	hr	520	558	1.27	655	1 G	...	N-(4-Methoxyphenyl)-1,3-b
10	M	516	767	1.13	2...	529	...	1,2-Benzenediol, o-(4-meth
11	hr	508	755	0.90	514	16	...	Loganic acid [M-H]- IT-FT
12	hr	506	835	0.85	513	1 M	...	6-Bromo-1,3-benzothiazole
13	M	505	509	0.83	1...	614	...	Methanone, (4-methoxyph
14	hr	504	867	0.80	504	...	...	5-Bromo-6-fluoro-1H-indol
15	M	502	595	0.76	2...	713	...	Nicotinic acid, 1-adamanty
16	M	501	664	0.74	2...	588	...	4-Methoxybenzoic acid, 3-I
17	M	500	745	0.72	2...	546	...	2-(2-Adamantan-1-yl)-2-oxo

InLib = 401, Hit List



(TextFile) cc\_2023\_192 #3102 RT: 10.03 AV: 1 NL: 1.24E9

Plot/Text of Search Spectrum | Plot of Search Spectrum | Spec List

Name: cc\_2023\_192 #3102 RT: 10.03 AV: 1 NL: 1.24E9

Precursor m/z: 213.0554

MW: N/A ID#: 30 DB: Text File


Comment: T: FTMS - c ESI d Full ms2 213.0554@hcd45.00 [47.9733-239.8665]

32 m/z Values and Intensities:

49	7.00	51	1.00	52	1.00	65	8.00	67	2.00
83	1.00	86	1.00	91	403.00	92	30.00	93	1.00
107	4.00	108	1.00	119	5.00	121	8.00	130	1.00
135	350.00	136	33.00	137	3.00	141	9.00	142	2.00
143	10.00	144	1.00	153	52.00	154	9.00	167	2.00
169	143.00	170	22.00	171	11.00	213	999.00	214	177.00
215	13.00	228	4.00						


Synonyms:

no synonyms.



cc\_2023\_192 #3102 RT: 10.03 AV: 1 NL: 1.24E9 | Head to Tail MF=774 RMF=853 | 2,4-Dihydroxybenzophenone

Head to Tail | Side by Side | Difference | Subtraction



(hr\_msms\_nist) 2,4-Dihydroxybenzophenone [M-H]- HCD 110% P=213.1

Plot/Text of Hit | Plot of Hit

Name: 2,4-Dihydroxybenzophenone

Precursor type: [M-H]-

Instrument type: HCD

Collision energy: NCE=110% 23eV

Precursor m/z: 213.0557

Formula: C<sub>13</sub>H<sub>10</sub>O<sub>3</sub>

MW: 214 Exact Mass: 214.062994 CAS#: 131-56-6 NIST#: 1689742 ID#: 420983 DB: hr\_msms\_nist

Comment: NIST Mass Spectrometry Data Center

Ion mode: N

Instrument: Thermo Finnigan Elite Orbitrap

Ionization: ESI

Collision gas: N<sub>2</sub>

Sample inlet: direct flow injection

Spectrum type: MS2

Notes: Sprec=Consensus Nprec=25/25 Mz diff=-2.8nm Vial ID=3369 Data source=Met 2017 08~125

## Full MSP File Format Exported from Thermo Freestyle to NIST Search

**Note:** Includes Precursor m/z Which Is Very Useful in NIST Search Program

```
Name: cc_2023_192 #3102 RT: 10.03 AV: 1 NL: 1.24E9
PrecursorMZ: 213.0554
DB#: 161
Comments: T: FTMS - c ESI d Full ms2 213.0554@hcd45.00 [47.9733-239.8665]
Num Peaks: 32
49 7.00
51 1.00
52 1.00
65 8.00
67 2.00
83 1.00
86 1.00
91 403.00
92 30.00
93 1.00
107 4.00
108 1.00
119 5.00
121 8.00
130 1.00
135 350.00
136 33.00
137 3.00
141 9.00
142 2.00
143 10.00
144 1.00
153 52.00
154 9.00
167 2.00
169 143.00
170 22.00
171 11.00
213 999.00
214 177.00
215 13.00
228 4.00
```

# MS Interpreter Will Interpret Fragments and Show the Mass Error, e.g. $m/z$ 137, 0.1 ppm error!

2,4-Dihydroxybenzophenone [M+H]<sup>+</sup> - MS Interpreter

File Edit View Options Help

Maximum Rate: 90(chain) @ 137.0233 m/z

Formula Calculator

Mass + = 137.0233 40 ppm Loss = 78.047

Formula C13H10O3(H) Parent = 215.0703

1 Ion	O+E	RDB	Mass	Ppm	C	H	O
C7H5O3	Even	5.5	137.02332	0.1	7	5	3

m/z	exact m/z	formula	loss	type	rate	rel.rate	ion.rate	ppm	abund	O+E
137	137.02332	C7H5O3	C6H6	H-Displacement	90	100	78	0.1	909	Even

