

Application News

Gas Chromatography Mass Spectrometry

No.M266A

GC/MS Forensic Toxicological Database (Part 2)

Introduction

In recent years, abuse of stimulants and other illegal drugs, particularly hallucinogens, shows no signs of abating; crime and poisoning events due to medicines, agricultural chemicals, and psychotropic drugs are on the increase, presenting a growing and serious social problem. In crime laboratories and forensic research laboratories, these chemicals are analyzed using a gas chromatograph mass spectrometer (GC-MS).

The detection of these harmful drugs by GC-MS requires investigation to determine the optimum analytical conditions and data processing conditions specific to each of the drugs. To eliminate the need to conduct these cumbersome operations, Shimadzu has developed its own proprietary "GC/MS Forensic Toxicological Database," which contains critical information for more than 1000 drugs*, including

psychotropic drugs, narcotics, stimulants, and pesticides. With emphasis on harmful drugs of abuse, difficult-to-obtain standard samples of these drugs and their metabolites are registered in the library, making the library ideal for forensic toxicological analysis. In addition, the database contains a library of mass spectra, method files containing the names, quantitation and reference ions, standard mass spectra, and retention indices specific to each of these drugs of abuse.

Here we introduce quick and easy techniques for detecting and conducting semi-quantitation of harmful drugs using method files.

* Includes derivatized compounds.

Content of Forensic Toxicological Method Files

Fig. 1 shows a view of the information contained in the GC/MS Forensic Toxicological Database, including chemical names, quantitation and reference ions, standard mass spectra, and retention indices.

Furthermore, relative response factors are also included for specific toxic substances, which can be used to calculate rough estimates of their concentrations in real-world samples.




Fig. 1 Registered Information in Method File

Table 1 Analytical Conditions

Model	: GCMS-QP2010 Series		
Workstation	: GCMSsolution Ver. 2.5 or later		
Column	: DB-5 ms (30 m × 0.25 mm I.D. df = 0.25 μm) or Rxi®-5 Sil MS (30 m × 0.25 mm I.D. df = 0.25 μm)		
-GC-		-MS-	
Column Temp.	: 60 °C (2 min)-10 °C/min-320 °C (10 min)	Interface Temp.	: 280 °C
Carrier Gas	: He (Constant Linear Velocity Mode)	Ion Source Temp.	: 200 °C
Linear Velocity	: 45.6 cm/sec	Scan Interval	: 0.3 sec
Injection Method	: Splitless		

Peak Detection Using Mass Chromatograms and Predicted Retention Time

Fig. 2 shows the screen used to analyze the data obtained from analysis of a drug-spiked urine sample using this method. When conducting actual analysis, it is very difficult to confirm the target peak associated with the drug of abuse because it is buried in the total ion current chromatogram (TICC), as shown in the upper part of the figure.

Use of this method file allows easy detection of the drug of abuse even in such cases of difficult-to-pinpoint target drug chromatographic peaks. First, the retention time of the target substance is predicted using the GCMSSolution AART (Automatic Adjustment of Retention Time) feature, based on the retention

times of n-alkanes measured beforehand, in addition to the registered retention index for each of the toxic substances. Then, the target drug substance is detected from the mass chromatogram (shown at the lower left of Fig. 2) of the quantitation and reference ions in the vicinity of the predicted retention time. These operations can be performed using the automatic processing function of the GCMSSolution software. Moreover, the final result can be verified through comparison with the standard spectrum, which is registered in the method file; the standard spectrum is shown here in the middle tier.




Fig. 2 Screen Capture of Data Processing Using GCMSSolution

Semi Quantitation-Calculation

This database method file contains relative response factors which are used for calculating semi-quantitative values for specific toxic drugs. As shown

in Fig. 3, the drug concentration can be roughly calculated without generating a calibration curve.

ID#	Name	Ret.Time	m/z	Area	Height	Search	SI	S/N	Conc	Unit
302	Benzalkonium chloride art									
303	5-MeO-DALT									
304	5-MeO-DPT									
305	4-OH-DPT									
306	5-MeO-DIPT	21.355	114.00	2974	1706		32	7.72	0.00000	ng
307	Seftipiline									
308	Piroheptine									
309	Endosulfan sulfate									
310	Carbamazepine	21.692	193.00	1224959	647709		98	172722	9.99980	ng
311	5-OH-DIPT									
312	Maprotiline									
313	Oxazepam									
314	Acetazolamide									
315	Profenamine									
316	Mazindol									
317	Sertraline									
318	Levallorphan									
319	Melatonin									
320	RAZ-509 (Haloxazolam)									
321	T-caine	22.431	99.00	10596	1891		36	1.81	0.00000	ng
322	Cannabidiol									
323	Doxulepin									
324	Fludiazepam									
325	4-AcO-DIPT									
326	Propiverine									

Fig. 3 Calculated Result Using Semi-Quantitation

Summary

The method file in the GC/MS Forensic Toxicological Database contains optimized analytical conditions, so toxic drug substances can be analyzed immediately. In addition, this method file contains various types of

information for many toxic substances, permitting semi-quantitation values to be calculated from an automatically identified drug of abuse.



SHIMADZU CORPORATION. International Marketing Division

3. Kanda-Nishikicho 1-chome, Chiyoda-ku, Tokyo 101-8448, Japan Phone: 81(3)3219-5641 Fax: 81(3)3219-5710