

Moving beyond ddMS²: Improving annotation confidence in untargeted metabolomics using higher resolution MS and parallel ion trap experiments

Brandon Bills, Sunandini Yedla, Rahul Deshpande, Bashar Amer, Ralf Tautenhahn, Susan Bird, Vlad Zabrouskov. Thermo Fisher Scientific, San Jose, CA. 95134. USA

Abstract

Purpose: Enhance the annotation of elemental formulas by increasing Orbitrap resolution to resolve isotopic fine patterns. Improve annotation confidence through supplementary MSⁿ or breakdown curve data via parallel ion trap scans.

Methods: NIST1950 SRM plasma extract was analyzed using a Thermo Scientific™ Orbitrap IQ-X™ Tribrid™ mass spectrometer at multiple resolutions, with AcquireX sequences incorporating ddMS², ddMS² with ddMS³, and ddMS² with breakdown curves, to evaluate the impact on the confidence of compound annotations.

Results: Compounds analyzed at higher resolutions were found to have more detailed isotopic fine patterns which in turn allowed for differentiation of isobaric formulas. This was more pronounced at m/z values above 400 where the number of potential elemental formulas within 1ppm mass accuracy was higher. Annotation confidence was improved through the addition of data collected in the ion trap. Collecting MS³ data provided additional criteria to compare and match against Thermo Scientific™ mzCloud™ mass spectral library or through FISH scoring and allowed for the differentiation of certain compounds where the MS² results were similar but the MS³ spectra differed. Break down curve data collected within the ion trap provided additional data points that could be used to differentiate similar compounds that produced different fragment intensities at different energy levels.

Introduction

In metabolomics, samples contain both known and unknown compounds needing characterization. Mass spectrometry helps determine elemental formulas via precursor m/z and structural composition through fragmentation data. However, higher m/z values lead to more potential formulas, complicating accuracy. Isomers and isobars further challenge single-spectrum characterization. This work aims to utilize increased Orbitrap resolution for better elemental formula annotation using isotopic fine patterns and boost confidence by collecting additional MSⁿ or breakdown curve data through parallel ion trap scans.

Materials and methods

Sample Preparation

NIST SRM 1950 plasma from Millipore Sigma (Massachusetts) was extracted and reconstituted in a methanol:water solution.

Test Method(s)

1 µL of plasma extract was injected on a Thermo Scientific™ Hypersil GOLD™ HPLC column using a Thermo Scientific™ Vanquish™ Horizon UHPLC system. Solvent A was water with 0.1% formic acid and solvent B was methanol with 0.1% formic acid. The gradient ranged from 0%B to 98%B.

Analysis was carried out on an Orbitrap IQ-X Tribrid mass spectrometer. The sample was analyzed in MS1 mode at 60K, 120K, 240K, 500K, and 1 million resolution. Fragmentation data was collected using Thermo Scientific™ AcquireX intelligent data acquisition to iteratively analyze the sample for more complete compound coverage. Three sequences incorporating ddMS², ddMS² with ddMS³, and ddMS² with breakdown curves were collected on the same sample.

Data Analysis

Raw data was processed using Thermo Scientific™ Freestyle™ software to predict theoretical isotopic patterns and isobaric species and Thermo Scientific™ Compound Discoverer™ software along with Thermo Scientific™ Mass Frontier™ spectral interpretation software to evaluate plasma compounds for elemental composition prediction, annotation confidence, and potential fragment structures.

Figure 1. Vanquish Horizon UHPLC system and Orbitrap IQ-X Tribrid mass spectrometer



Results

Resolution and elemental formula prediction

Evaluating Increased resolution

Elemental composition is predicted based on accurate mass and resolution. As m/z increases, the number of isobaric elemental compositions increases, even within 1 ppm mass error. Table 1 shows the number of predicted a2 isotopic peaks for compounds at different resolutions and the number of theoretical isobars.

Table 1. Predicted a2 isotopic peaks and theoretical isobars at different resolutions and mass accuracies

Compound Information			Number of a2 isotopic fine peaks* at resolution					Theoretical isobars* within mass tolerance		
Name	m/z	Formula	60K	120K	240K	500K	1M	1 ppm	3 ppm	5 ppm
Cystamine	153.0515	C ₄ H ₁₀ N ₂	2	3	3	3	3	1	2	3
Phenylalanine	166.0863	C ₉ H ₉ NO ₂	1	3	4	5	5	1	2	3
Melatonin	233.1284	C ₁₂ H ₁₆ N ₂ O ₂	1	2	4	4	5	1	2	4
Biotin	245.0954	C ₁₀ H ₁₆ N ₂ O ₃	2	2	4	5	5	3	9	12
Palmitic Acid	255.2330	C ₁₆ H ₃₂ O ₂	1	1	3	3	4	1	1	1
Glucose-1-P	261.0370	C ₆ H ₁₂ O ₆ P	1	1	2	2	3	5	11	22
Maltose	341.1089	C ₁₂ H ₂₂ O ₁₁	1	1	2	3	4	6	21	32
3'-AMP	348.0704	C ₁₀ H ₁₄ N ₅ O ₇ P	1	2	3	5	6	9	31	54
Cholic Acid	407.2803	C ₂₆ H ₄₆ O ₅	1	1	1	3	4	4	10	15
Coenzyme A	768.1225	C ₂₁ H ₃₆ O ₇ P ₂ S	1	1	2	4	5	86	277	400*

*peaks over 1% relative abundance of the most intense isotope in that cluster resolved at half maximum

*Only considering the elements: H, C, N, O, P, S

Isotopic fine peaks, peaks from heteroatoms, can eliminate incorrect elemental composition predictions. For example, figure 2 shows the a2 isotopic peak for 496.3398 at resolutions from 60K to 1M. At 1M four distinct peaks show up.

Figure 2. The a2 isotopic pattern for a peak at 496.3398 measured at increasing resolutions.




Figure 3 shows a2 patterns for three predicted elemental compositions for 496.3398. C₂₄H₅₁O₇NP has a small (but still present) peak for ¹³C+¹⁵N, a more intense peak for ¹⁸O, and a small peak for ¹³C+²H, which is closest to the a2 pattern at 1M.

Figure 3. The a2 patterns for three predicted elemental compositions for 496.3398

