

## Reducing false annotations in lipidomics using Compound Discoverer 3.4 software

Rahul Ravi Deshpande, **Ralf Tautenhahn**, Bashar Amer and Susan S. Bird, Thermo Fisher Scientific, San Jose, CA, USA

### Abstract

**Purpose:** Confident annotation is essential for translating untargeted lipidomics data into meaningful biological information. Thermo Scientific™ Compound Discoverer™ 3.4 software with Thermo Scientific™ LipidSearch™ software was evaluated for its ability to reduce false annotations by utilizing some of its key features including the detection of adducts and in-source or unintentional fragmentation.

**Methods:** Lipid Standards and bovine liver total lipid extracts were used for evaluation. Separation was carried out either using analytical flow or nano flow UHPLC coupled. Thermo Scientific™ Orbitrap™ Astral™ mass spectrometer or Thermo Scientific™ Orbitrap™ Excedion™ Pro hybrid mass spectrometer were used for detection.

**Results:** Adducts and unintentional fragments combined form more than 50% of all the features detected irrespective of the chromatographic flow rate or the mass spectrometer used. The Orbitrap Excedion Pro MS with its new front end for ion transmission can significantly reduce the in-source or unintentional fragmentation.

### Introduction

False annotations and unintended or in-source MS1 fragmentation are critical challenges in lipidomics, affecting data accuracy and interpretation. False annotations arise when lipids are misidentified due to spectral similarities, database limitations, or incorrect adduct assignments. Unintentional MS1 fragmentation occurs when lipids break down within the ion source or the ion path of the mass spectrometer, generating fragment ions that may be mistaken for distinct lipid species. This complicates lipid identification and quantification, increasing false discovery rates. Here we show the utility of Thermo Scientific Compound Discoverer 3.4 software with LipidSearch software for grouping of adducts and unintentional MS1 fragments to avoid false annotations and for confident annotation of lipids.

### Materials and methods

#### Sample Preparation

Avanti Research™ SPLASH™ LIPIDMIX™ (mixture of internal standard) and liver total lipid extract were purchased from Avanti Lipids. 1:100 dilution of the liver extract was used as the stock concentration of matrix.

#### Instrument Method

Lipid standards and bovine liver lipid extracts were purchased from Avanti Research™. Separation was carried out on a Thermo Scientific™ Accucore™ C30 HPLC column connected to a Thermo Scientific™ Vanquish™ Horizon UHPLC system. Extracted Lipids were also analyzed using nanoLC on a Thermo Scientific™ EASY-Spray™ PepMap™ Neo UHPLC column connected to a Thermo Scientific™ Vanquish™ Neo UHPLC system. The run time for analytical flow was 30 Minutes and for the nanoflow was 15 minutes. Data was acquired on Orbitrap Astral MS and Orbitrap Excedion Pro MS. 20X amount was loaded on the analytical flow compared to the nanoflow.

### Compound Discoverer 3.4 Software

- Includes LipidSearch software node for simplified lipid identification
- QC for small molecule experiments using internal standards
- Automatically detects and groups MS1 fragments
- Connects to Thermo Scientific™ mzCloud™ mass spectral library 2.0 for improved compound identification

### LipidSearch software in Compound Discoverer software

- LipidSearch software provides automated identification of lipids with a comprehensive database, containing 1.5 million lipid ions and their predicted fragment ions
- LipidSearch software node can be added to any untargeted data processing workflow
- New workflow template: *Lipidomics\Untargeted Lipidomics using LipidSearch software*




Figure 1. LipidSearch software available as a workflow node in Compound Discoverer software with a new template for untargeted lipidomics.

### LipidSearch: Lipid Identification

- Database contains >1.7 million lipid ions and predicted fragment ions
- Fragmentation patterns are utilized to grade identification quality and provide all possible annotations
- Lipid adduct ions and MS<sup>2</sup> and MS<sup>3</sup> fingerprints are also included




Figure 2. Predicted fragmentation of a typical lipid included in the LipidSearch software database

Category	Definition	Condition	Definition
C0	Headgroup specific ion	H1	C0 and C01 (or more than 2 C0) (or C01 for all substituents) are assigned
C01	Class specific product ion	H2	One of C0 or C01 is assigned
C1	Substituent product ion	S1	C1 and C01 from all substituents are assigned
C2	Other non-specific ions	S2	More than one of C1 or C01 from substituents are assigned
		Other	C2 – headgroup and substituents not defined

Grade	Rule	Note:
A	H1 and S1	Grade A assigned based on product ions for one class specific, C01, one headgroup specific, C0, and all substituent ions, C1
B	(H2 and S1) or (H1 and S2)	Grade B assigned based on the lack of H1 or S1 conditions
C	H1 or S2	Grade C assigned based on the lack of H or S conditions
D	Other	Grade D assigned based on lack of headgroup/substituent ions

Figure 3. Rule based gradation of lipid annotation in LipidSearch software. The grade of a lipid is determined by comparing its fragmentation spectra with the predicted fragmentation. Higher grade implies more confident annotation.

Calc. m/z	Ion	Category	PC 16:0_20:4	PC 20:1_16:3	PS 16:0_23:3	PS 23:0_16:3
826.55036	C <sub>42</sub> H <sub>84</sub> NO <sub>8</sub> P	M+HCOO <sup>-</sup>	☑	☑	☑	☑
766.55923	M-CH <sub>3</sub>	C0	☑	☑	☑	☑
480.30957	UPC14(9)-CH <sub>3</sub>	C01	☑	☑	☑	☑
303.23295	20:4-H <sup>+</sup>	C1	☑	☑	☑	☑
259.24312	20:4-H-CO <sub>2</sub> <sup>-</sup>	C1	☑	☑	☑	☑
255.23295	16:0-H <sup>+</sup>	C1	☑	☑	☑	☑
224.05934	C <sub>18</sub> H <sub>34</sub> NO <sub>8</sub> P	C2	☑	☑	☑	☑
265.19617	16:3-H-CO <sub>2</sub> <sup>-</sup>	C1	☑	☑	☑	☑
168.04312	(P-CH <sub>2</sub> )-CH <sub>3</sub> -H	C0	☑	☑	☑	☑
	Rank		1	2	3	4
	ID score		1.88	0.86	0.31	0.01
	Grade		H1 and S1 → A	H2 → C	D	D

Figure 4. The gradation of isomeric lipid species based on the fragmentation spectra for confident annotation.




Figure 5. Annotation of fragments in the MS2 spectra by LipidSearch software to help improve annotation confidence. The above spectra is for TG(16:1\_18:1\_18:1)

### Adduct Ions

Adduct ions are formed by the interaction of a precursor ion with one or more atoms or molecules to form an ion containing all the constituent atoms for the precursor ion as well as the additional atoms from the associated atoms or molecules. Common lipid adducts include +Na, +K, -H<sub>2</sub>O etc. depending on the lipid class.

### In-source or Unintentional Fragmentation

In-source or unintentional fragmentation refers to the breakdown of ions into smaller fragments within the instrument's ion source or ion path. This phenomenon occurs before the ions are analyzed in the mass analyzer leading to complex MS1 spectra.

