

Poster Reprint

ASMS 2023
Poster number TP 732

Improvement Upon a Multi-Residue Method for Nitrosamine Analysis in Losartan Drug Product Using an Enhanced LC/MS/MS System

Bertram Lee¹, Linfeng Wu¹, Shan-An Chan² and Winnie Huang³

¹Agilent Technologies, Inc., Santa Clara, CA, USA


²Agilent Technologies, Inc., Taiwan, Republic of China

³Agilent Technologies, Inc., Guangdong, China

Since 2018, nitrosamine impurities have been recognized as byproducts produced in trace amounts during manufacture of pharmaceutical drugs. They are members of the “cohort of concern” classifying them as potentially genotoxic impurities and probable carcinogens if they exceed the maximum allowable intake of 96 ng/kg for a treatment with long-term intake. A recent announcement for the recall of some angiotensin II receptor blocker (ARB) drug products made nitrosamine impurities a focus for regulatory agencies.¹ Therefore, it is important to determine their levels in final drug products with a high level of sensitivity and confidence.

In this study, a comprehensive analysis of eight nitrosamine compounds was carried out on the new 6495 LC/TQ (G6495D) coupled with the 1290 Infinity II Bio LC system and APCI source. It is demonstrated that the following compounds can be determined at low levels:

- N-Nitrosodimethylamine (NDMA)
- N-Nitrosomorpholine (NMOR)
- N-Nitrosomethylethylamine (NMEA)
- N-Nitrosopyrrolidine (NPYR)
- N-Nitrosodiethylamine (NDEA)
- N-Nitrosopiperidine (NPIP)
- N-Nitrosodi-n-propylamine (NDPA)
- N-Nitrosodi-n-butylamine (NDBA)



Sample Preparation

Nitrosamine standards were spiked into solvent blank (10:90 methanol:water) at nine different concentration levels ranging from 0.0125 to 10 ng/mL. The preparation of active pharmaceutical ingredient (API) matrix was carried out following these steps: 100 mg losartan potassium drug product was dissolved in 2 mL of solvent (50:50 methanol:water), followed by sonication for 30 minutes. Then, samples were centrifuged at 12,000 rpm for 10 minutes. Supernatants were collected and then diluted with water at a ratio of 1:5. Nitrosamine standards were spiked into the prepared API matrix at concentrations ranging from 0.05 to 1 ng/mL.

Equipment

Sample separation was performed using the Agilent 1290 Infinity II Bio LC system consisting of the following modules:

- 1290 Infinity II Bio high-speed pump (G7132A)
- 1290 Infinity II Bio multisampler with thermostat (G7137A)
- 1290 Infinity II multicolumn thermostat (G7116B)

The LC system was coupled to the Agilent 6495 triple quadrupole LC/MS (G6495D) equipped with the Agilent APCI source (G1947B). Agilent MassHunter Workstation software 12.0 was used for data acquisition.

The Agilent 6495 triple quadrupole LC/MS system with 1290 Infinity II Bio LC



Methods

Liquid chromatography/mass spectrometry (LC/MS) conditions and parameters are provided in Tables 1 and 2. The multiple reaction monitoring (MRM) settings for the compounds are listed in Table 3.

| Parameter | Value | |
|---------------------|---|----|
| Column | Agilent Poroshell 120 EC-C18, 3.0 X 150mm, 2.7 μ m (p/n 693975-302) | |
| Sampler Temperature | 4 °C | |
| Mobile Phase A | ddH ₂ O+0.1% Formic acid | |
| Mobile Phase B | MeOH+0.1% Formic acid | |
| Flow Rate | 0.5 mL/min | |
| Injection Volume | 20 μ L | |
| Column Temperature | 40 °C | |
| Gradient Program | Time (min) | %B |
| | 0.0 | 5 |
| | 3.5 | 5 |
| | 7.0 | 45 |
| | 9.0 | 60 |
| | 11.0 | 60 |
| | 15.0 | 65 |
| | 16.0 | 90 |
| | 16.1 | 5 |
| | 20.0 | 5 |

Table 1. 1290 Infinity II Bio LC Method.

| Parameter | Value |
|----------------------------|---|
| Ion Source | Agilent APCI source |
| Polarity | Positive |
| Gas Temperature | 300 °C |
| Drying Gas Flow | 7 L/min |
| Nebulizer | 25 psi |
| APCI Vaporizer Temperature | 350 °C |
| Capillary Voltage | 1000 V |
| Corona Current | 4.0 μ A |
| Scan Type | MRM with time segments |
| Detector gain factor (+) | 10 |
| Ion Funnel Mode | Fragile |
| LC diverter to waste | 0 – 2 min; 12 – 13 min; 14 – 18 min (the remaining time diverter to MS) |

Table 2. 6495 LC/TQ Parameters.

| Compound name | Precursor m/z | Product m/z | Dwell (ms) | CAV (V) | CE (V) | Polarity |
|---------------|---------------|-------------|------------|---------|--------|----------|
| NDMA | 75 | 43 | 50 | 3 | 16 | + |
| NDMA | 75 | 58 | 50 | 3 | 10 | + |
| NMOR | 117 | 45 | 50 | 4 | 21 | + |
| NMOR | 117 | 87 | 50 | 4 | 11 | + |
| NMEA | 89 | 43 | 50 | 3 | 12 | + |
| NMEA | 89 | 61 | 50 | 3 | 10 | + |
| NPYR | 101 | 41 | 50 | 3 | 24 | + |
| NPYR | 101 | 55 | 50 | 3 | 19 | + |
| NDEA | 103 | 47 | 50 | 4 | 20 | + |
| NDEA | 103 | 75 | 50 | 4 | 12 | + |
| NPIP | 115 | 41 | 50 | 3 | 24 | + |
| NPIP | 115 | 69 | 50 | 3 | 12 | + |
| NDPA | 131 | 43 | 50 | 4 | 10 | + |
| NDPA | 131 | 89 | 50 | 4 | 16 | + |
| NDBA | 159 | 41 | 50 | 3 | 20 | + |
| NDBA | 159 | 57 | 50 | 3 | 12 | + |

Table 3. Detailed Multiple Reaction Monitoring Settings.

Results and Discussion

Calibration Curve Analysis

To evaluate the quantification performance of nitrosamines, the calibration curves of the eight nitrosamine compounds were analyzed with concentrations ranging from 0.0012 to 10 ng/mL in solvent.




Figure 1. Total MRM Chromatogram of all Compounds at 0.6 ppb.

| Name | Concentration | | | | | | | | LLOQ (ng/mL) |
|------|------------------|--------|---------|--------------|------------------|--------|---------|--------------|--------------|
| | 0.01 ng/mL (n=3) | | | | 0.04 ng/mL (n=3) | | | | |
| | Response | S/N | RSD (%) | accuracy (%) | Response | S/N | RSD (%) | accuracy (%) | |
| NDMA | 13035 | 48.8 | 1.30 | 103.0 | 46851 | 130.7 | 0.001 | 92.6 | 0.01 |
| NMOR | 3286 | 102.1 | 0.01 | 89.4 | 12779 | 357.7 | 0.009 | 86.9 | 0.01 |
| NMEA | 5121 | 136.7 | 0.02 | 93.8 | 20411 | 436.3 | 0.019 | 89.2 | 0.01 |
| NPYR | 6792 | 66.9 | 0.01 | 103.6 | 31451 | 317.3 | 0.030 | 91.4 | 0.01 |
| NDEA | 1752 | 149.51 | 0.05 | 81.3 | 6921 | 374.2 | 0.012 | 85.5 | 0.02 |
| NPPI | 14630 | 162.5 | 0.02 | 109.5 | 64325 | 858.2 | 0.008 | 93.4 | 0.01 |
| NDPA | 8745 | 305.6 | 0.02 | 95.8 | 34844 | 1316.8 | 0.002 | 90.7 | 0.01 |
| NDBA | 28666 | 586.8 | 0.02 | 108.4 | 46273 | 769.8 | 0.005 | 102.5 | 0.04 |

Table 4. Average Response, Signal to Noise (S/N), %RSD and Accuracy at 0.01 and 0.04 ng/mL using the Agilent 6495 LC/TQ.




Figure 2. Calibration Curves for the Eight Nitrosamine Compounds on the Agilent 6495 LC/TQ.

Quantification in Losartan Potassium Drug Matrix

To examine the analytical performance of nitrosamine impurities in pharmaceutical drug products, the nitrosamine standards were spiked into the losartan potassium drug extract at concentrations of 0.0015 to 1 ng/mL. The results show that all eight nitrosamine compounds could be quantified with high confidence at 0.03 ng/mL in drug matrix. The analytical sensitivities for all eight analytes were well below the acceptable nitrosamine content that is stated by regulator agencies.^{2,3}

<https://www.agilent.com/en/promotions/asms>

This information is subject to change without notice.

DE60229055

© Agilent Technologies, Inc. 2023
Published in USA, May 31, 2023

NDMA is commonly a challenging compound for nitrosamine impurity analysis. Figure 3 shows the representative chromatograms of NDMA at matrix blank, and at 0.03 and 0.12 ng/mL in 10 mg/mL drug with the 6495 LC/TQ.




Figure 3. Representation Chromatograms of NDMA in Losartan Potassium Drug Matrix at Low Concentrations

Conclusions

This poster demonstrates that the Agilent 6495 triple quadrupole LC/MS system shows increased analytical sensitivity of nitrosamine impurities at the low concentration levels specified by regulatory requirements. Ability to detect lower than 0.01 ng/mL for most nitrosamine impurities. The 6495 LC/TQ system along with this method can be used to quantify these impurities in different ARB drug products, with some changes in chromatographic conditions based on the elution pattern of the drug product.

References

- 1 U.S. Food and Drug Administration. FDA Updates and Press Announcements on Angiotensin II Receptor Blocker (ARB) Recalls (Valsartan, Losartan, and Irbesartan). Last revised March 2023.
- 2 Center for Drug Evaluation and Research, U.S. Food and Drug Administration, Control of Nitrosamine Impurities in Human Drugs: Guidance for Industry. Last revised February 2021. Control of Nitrosamine Impurities in Human Drugs | FDA.
- 3 Committee for Medicinal Products for Human Use, European Medicines Agency. Procedure under Article 5(3) of Regulation EC (No) 726/2004. Procedure number: EMEA/H/A-5(3)/1490.