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 CustomerSupport

NotReady
Pump OFF

Destination (IP Address)
192 . 168 . 200 . 93 Disconnect Connect

Maintenance > Hydration and Regeneration of the Suppressor

- ✓ Preparation for Hydration and Regeneration of the Suppressor
- ✓ Inspection of the Suppressor
- ✓ Hydration Operation
- ✓ Pre-Check Before Regeneration Operation
- 5 Regeneration Operation**

Regeneration Operation

This explains the procedure for regenerating the suppressor. Use sodium hydroxide solution for the cation suppressor ICDS-Ci, and sulfuric acid solution for the anion suppressor ICDS-Ai.

Chapter

- Remove the suppressor and turn off the settings.
[00:00-00:50](#)
- Replace with cleaning solvent and perform flow path cleaning.
[00:50-02:20](#)
- Replace with ultrapure water and perform flow path cleaning.
[02:20-03:00](#)

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


Nebulizing Gas Flow	3 L/min
Heating Gas Flow	15 L/min
Interface Temperature	200 °C
Desolvation Temperature	355 °C
DL Temperature	200 °C
Heat Block Temperature	300 °C
Drying Gas Flow	5 L/min
Probe Position	+1mm
Interface Voltage	-0.5kV
Focus Voltage	-1kV

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0.05 – 50 µg/L in vial

0.2 – 200 ng/L in water sample after SPE

Thank you!

Welcome to the Shimadzu family.
Get ready to experience New Benchmarks.

Company names, product/service names and logos used in this presentation are trademarks and trade names of Shimadzu Corporation or its affiliates, whether or not they are used with trademark symbol "TM" or "®".