

Poster Reprint

ASMS 2022
Poster number TP070

Quantitative screening of 68 Semi-volatile Organic Compounds residues in water with Automated Liquid-Liquid Extraction and large volume injection by GC-MS

Ge Meng¹; Lakshmi Krishnan²; Anastasia Andrianova²;
Zhicong Wang¹; Aaron Boice²

¹ Agilent Technologies, Shanghai , China

² Agilent Technologies, Santa Clara, CA

Introduction

The China Ministry of Ecological Environment (MEE) has released the updated standards of pollutants analysis in drinking water¹, regulating analytical methods for multiple types of semi-volatile organic compounds (SVOCs), including pesticides, nitrosamines, polychlorinated biphenyls (PCBs), polycyclic aromatic hydrocarbons (PAHs) and phthalates. Sample preparation is a significant aspect of the analysis. Most of the standard methods applies conventional preparation workflow such as liquid-liquid extraction (LLE) and Solid phase extraction (SPE), needing a lot of manpower and solvent consumption. To help operators obtain enhanced efficiency in water SVOCs analysis, a fully automation solution including preparation (automated LLE) and determination (GC-MS) was developed for 68 targets quantitative screening.



Experimental



Figure 1. CTC PAL RTC with Agilent 8890GC-5977B MSD.

CTC PAL RTC with vortex module was brought into this solution for automated LLE. Agilent 8890 GC with Multi-mode Injection (MMI) and 5977B MSD were utilized as analytic platform. 15 mL water sample was transferred into 20 mL headspace vial with 1.5 g NaCl. The following preparation would be performed automatically including: 1 mL hexane addition, vortex, phase separation, supernatant transfer, internal standard addition and homogeneous mixing. After preparation, 10 μ L sample was injected into MMI under room temperature for concentration. With specific script, determination of sample (N) and preparation of sample (N+1) can be perfectly overlapped. With this solution, operators can obtain simultaneous screening results of 68 SVOCs (31 analytes indicators in Chinese drinking water quality standard) within 1 hour.

Experimental

GC	8890		
Mode	MMI Splitless		
liner	Ultra Inert, splitless, single taper, glass wool, PN: 5190-2293		
Heater	60 °C (10 μ L injection) / 280 °C (1 μ L injection)		
MMI program	10 μ L injection		
	rate	temp	Hold time
		60 °C	2 min
	720 °C/min	280 °C	20 min
	720 °C/min	60 °C	7 min
	1 μ L injection		
	rate	temp	Hold time
		280 °C	
Flow	1.5 mL/min, constant flow		
Purge flow to split vent	10 μ L injection: 60 mL/min @ 2.5 min		
	1 μ L injection: 60 mL/min @ 0.8 min		
Oven program	rate	temp	Hold time
		40 °C	4 min
	10 °C/min	300 °C	6 min
MSD Transfer Line	280 °C		
Column	HP-5MS UI, 30 m x 250 μ m x 0.25 μ m, PN: 19191S-433UI		
MS	5977B		
MS Source	230 °C		
MS Quad	150 °C		
Solvent Delay	5 min		
Acquisition type	Select Ion Monitor (SIM)		

Figure 2. GC and MS system settings for 68 SVOCs automated analysis.




Figure 3. Left: automated vortex at 2000 r/m; Right: transfer 200 μ L supernatant automatically.

Results and Discussion

For each analyte, calibration curve was built with at least 6 concentration levels (0.1-200 µg/L). Linearity coefficients were established greater than 0.99 for all curves. Accuracy and reproducibility of the workflow were evaluated with spiked drinking water matrix (pure water, tap water and reservoir water) at low and high concentration levels (10-100 µg/L). More than 95% analytes had average recoveries (n=6) within 70-130%; relative standard deviation (RSD%) of all the targets was lower than 20% (n=6), showing excellent repeatability of method. The lower quantitation limits (LOQ) were determinate with 7 replicates of spiked samples (1 µg/L) and calculated as lower than 0.1 µg/L. For all the 31 indicator compounds of drinking water quality, method LOQ was lower than minimum requirement limitation (MRL) in government regulation. This workflow solution reduces ~95% of manual operation and organic solvent consumption, in the meantime, provides comparable data performance as conventional methods.




Figure 4. 68 target SVOC analytes with 4 surrogates and 6 I.S..

Conclusion

A full automated workflow for screening of 68 SVOCs (31 analytes indicators in Chinese drinking water quality standard) was done within 1 hour. It resulted in ~95% reduced manual operation and organic solvent consumption, while still maintaining efficient and good analytical performance.

References

¹Standards for Drinking Water Quality (GB 5749-2022).

<https://explore.agilent.com/asms>

This information is subject to change without notice.

DE79278477

© Agilent Technologies, Inc. 2022
Published in USA, May 20, 2022