

# Scale Up to the Drug Discovery Fast Lane with Superficially Porous Particle Preparative LC Columns

Lakshmi Subbarao  
Applications Scientist  
Small Molecule Preparative LC Columns  
November 30, 2021



## In today's webinar, you will learn:

- Characteristics of columns that drive throughput
- What is a superficially porous particle
- How columns packed with superficially porous particles can be beneficial to the purification workflow



## This webinar is intended for:

- Purification scientists
- Users of reversed phase preparative LC in small molecule drug discovery labs
- Chemists and lab managers trying to drive high resolution and throughput

## What is Prep LC?

Workflow and understanding what is important in discovery environments

# What is Prep LC and why do we do it?



## Goal of preparative LC

- Collect one or more compounds from a sample

## Why do we do it?

- To obtain a “pure” product
- For compound characterization

## The focus

- To get required purity (purity)
- To get required amount (yield)
- In the required timeframe (throughput)

## High Throughput



### Objectives

- Purifying many different samples
- Small amounts of material for preliminary bioactivity testing or further characterization
- Target – main component

### Workflow

- Not much time can be spent on method development – analytical screen to determine best column, then utilize focused gradient to improve separation (if necessary)
- Requires small amounts of many different samples, typically 1-2 injections per sample

### Applications

- Combinatorial chemistry libraries
- Reaction cleanup
- Open access purification for medicinal chemists

### Objectives

- Purifying one sample
- Significant amounts of material for more in-depth testing and characterization
- Target – main component or an impurity

### Workflow

- Considerable time is spent optimizing method to maximize yield, purity and minimize solvent consumption
- Many injections of the same sample

### Applications


- Purification of target compounds (API), Natural products, etc..
- Impurity isolation for Identification, structure elucidation and activity testing

## Bulk Purification




# Objectives of prep LC: purity, yield and throughput


- Purity
- Yield
- **Throughput**



Longer run times to have enough separation between peaks



- **Throughput**
- Yield
- Purity



Shorter run times, but collections contain impurities

You can optimize 2 of the parameters by sacrificing the third  
The most important parameter depends on the application


# Objectives based on purification environment

## High Throughput



Objectives
<ul style="list-style-type: none"><li>• Purifying many different samples</li><li>• Small amounts of material for preliminary bioactivity testing or further characterization</li></ul>

Priority
<ul style="list-style-type: none"><li>• Purity – Yes</li><li>• Throughput – Yes</li><li>• Yield – No</li></ul>




## Bulk Purification




Objectives
<ul style="list-style-type: none"><li>• Purifying one sample</li><li>• Significant amounts of material for more in-depth testing and characterization</li><li>• Optimizing method to maximize yield, purity and minimize solvent consumption</li></ul>

Priority
<ul style="list-style-type: none"><li>• Purity – Yes</li><li>• Yield – Yes</li><li>• Throughput - No</li></ul>



# Steps in developing a prep method





# Steps 1 and 2: high throughput purification

1. Check  
sample  
formulation

2. Screen  
column  
phases

3. Optimize  
separation

4. Determine  
max sample  
load

5. Determine  
prep column  
size

6. Scale  
conditions to  
prep column



7. Purify the  
compound(s)

## 1. Check sample formulation

- Samples are crude mixes and can have complicated matrices

## 2. Screen column phases

- Screening combinations of columns/mobile phases using a generic gradient
- Column chemistries are ones that have historically high success rate
- The goal is to identify a column/mobile phase combination that will have the highest chance of separation



# Step 3: Optimize separation (high throughput purification)

1. Check sample formulation

2. Screen column phases

3. Optimize separation



4. Determine max sample load

5. Determine prep column size

6. Scale conditions to prep column

7. Purify the compound(s)

- Running a focused gradient can tremendously improve the critical pair's separation with minimal method development
- Focused gradient – reducing the slope of the original gradient specific to the peak(s) of interest




## Generic gradient (Adjusted):

5% B: 0 to 2.0 min  
85% B: 10.0 min  
85% B: 10 to 11.5 min  
Slope: 10% B/min


## Focused gradient (Adjusted):

5% B: 0 to 2.0 min  
25% B: 8.5 min  
85% B: 8.6 to 11.6 min  
Slope: 3.1% B/min



Want more info? Visit <http://explore.agilent.com/preparative-lc-primer>

# Steps 4 and 5: high throughput purification



## 4. Determine max sample load


- Generally not done for high throughput work since only small amount of sample is required

## 5. Determine prep column size

- Instead, a column is chosen based on anticipated load. The most popular column IDs for discovery work are 21.2 mm and 30 mm ID.

Column ID	Flow Rate	Difficult Separation ( $1.2 < \alpha < 1.5$ )
21.2 mm	20 mL/min	20 – 70 mg
30 mm	40 mL/min	40 – 140 mg

# Steps 6 and 7: high throughput purification



## 6. Scale conditions to prep column

- Premade focused gradients
- A generic gradient of 10-50% B would have 4 focused gradients:  
10-20%, 20.1-30%, 30.1-40%, 40.1-50%
- User can select gradient based on peak's elution percentage

### Manual calculation

$$f_p = f_a \left( \frac{d_p}{d_a} \right)^2, V_p = V_a \left( \frac{d_p}{d_a} \right)^2$$


$f_a$  = Analytical flow rate

$d_p$  = Prep internal diameter

$d_a$  = Analytical internal diameter

$V_a$  = Analytical injection volume

### Automated calculation



## 7. Purify the compound(s)

- Generic trigger settings to accommodate most collections
- Always a clean waste bottle in case a compound is missed
- MS detection is typically used to minimize number of collections made

## High Throughput



### What's important?

- Collecting enough material (10-100 mg)
- High purity
- Minimal method development


### Ways to improve throughput:

- Screening (step 2) – analytical columns
- Size (step 5) - preparative columns


# Superficially Porous Particles

What they are and how they work


# Understanding performance by understanding band broadening



Why do some columns produce peaks like this?



...and other columns produce peaks like that?



**How do I get the best performance from my column?**

# Van Deemter curves

- Van Deemter curves provide a simple way to compare columns
- Uses plate height (measure of band broadening), which is inversely proportional to plate number (N)

## Experimental

- Sample
  - Void volume marker (e.g., uracil)
  - Neutral analyte (naphthalene, toluene, etc.)
- Isocratic method that produces a symmetrical peak
- Inject sample across a wide range of flow rates
- Plot plate height (H) against linear velocity (u)

## Equations

$$N = 5.54 \left( \frac{t_r}{W_{1/2}} \right)^2$$

L = length of the column  
 $W_{1/2}$  = peak width at half height  
 $t_r$  = retention time of analyte

$$H = L/N$$

$$u = L/t_0$$

L = length of the column  
 $t_0$  = retention time of void volume marker




# Van Deemter equation – significance of A, B and C terms

$$H = A + \frac{B}{u} + Cu$$

Where:


- **A term:** eddy diffusion
  - Flow path of analyte through stationary phase particles
  - Particle size, size distribution and packing quality
- **B term:** longitudinal diffusion
  - Diffusion in the mobile phase
  - Only significant at very low flow rates
- **C term:** mass transfer resistance
  - Analyte traveling in/out of particle
  - **Significant at mid to high flow rates**

## Van Deemter plot for 5 $\mu\text{m}$ TPP column




# Superficially porous particles

- Superficially porous particle (SPP) – Solid core and porous outer layer
- Solid core shortens diffusion path, resulting in a narrower chromatographic peak over totally porous particle (TPP) column




## van Deemter plots for SPP and TPP columns (similar particle size)



- SPP columns have smaller C terms
- Can run them at higher flow rates while maintaining column efficiency

# Agilent InfinityLab Poroshell 120 preparative LC columns



## Made with 4 $\mu\text{m}$ Poroshell particles

- Symmetrical particles with smooth surfaces ( $\downarrow$  A term)
- Narrow particle size distribution ( $\downarrow$  A term)
- Short analyte diffusion path in/out of the porous layer ( $\downarrow$  C term)




# Generic purification methods on SPP and TPP columns

## Separation of withaferin A in ashwagandha

- Discovery labs regularly purify large batches of samples.
- They require small amounts (10 to 100 mg) of high purity fractions for downstream workup and characterization
- Purification of bioactive components in natural products represents a similar challenge to that of drug candidates in crude mixtures – both have complicated matrices.
- Withaferin A (WFA) is the most bioactive withanolide in ashwagandha<sup>1</sup>.

Withaferin A<sup>2</sup>



General Run Conditions	
Instrument	Agilent 1290 Infinity II autoscale preparative LC system
Sample	Ashwagandha extract in 2:1 ethanol:water, 100 mg/mL
Mobile Phase	A: Water + 0.1% formic acid
	B: Acetonitrile + 0.1% formic acid
Injection Volume	1 mL filtered extract

	Optimized	Elevated
Agilent InfinityLab Poroshell 120 SB-C18 21.2 x 50, 4 µm		
Flow rate:	25 mL/min	37.5 mL/min
Gradient:	5 – 95% B in 15 min	5 – 95% B in 10 min
TPP C18, 19 x 150, 5 µm		
Flow rate:	17 mL/min	25.5 mL/min
Gradient:	5 to 95% B in 18 min	5 to 95% B in 12 min


Link to application note: [5994-3518EN](#)

1. M. H. Mirjalili et al. *Acta Chromatographica* 25(2013)4, 745-754

2. [https://www.sigmaaldrich.com/content/dam/sigma-aldrich/structure3/198/mfcd10687098.eps/\\_jcr\\_content/renditions/mfcd10687098-medium.png](https://www.sigmaaldrich.com/content/dam/sigma-aldrich/structure3/198/mfcd10687098.eps/_jcr_content/renditions/mfcd10687098-medium.png)

# WFA separation at optimal and 1.5x optimal flow rate


Optimal  
flow rate



- Separates WFA from impurities
- Best overall resolution

- Separates WFA from leading impurity
- Trailing impurity elutes in WFA tail

1.5x  
optimal  
flow rate





- Still WFA from impurities
- Best overall resolution

- Impurities have almost completely coeluted

# Comparing best methods

## 1.5x optimal SPP flow rate vs optimal TPP flow rate

- SPP column maintains separation at high flow rate using a generic gradient
- TPP column loses resolution between optimal and high flow rates
- The SPP method is 45% faster than the TPP method




# Integrating SPP columns into the purification workflow

## Overcoming challenges



= How?



## Instrument compatibility

- Running conventional SPP sizes may exceed system pressure limits
- System volume causes enough band broadening to negate SPP performance benefits


## Sample loading limitations

- SPP phases lack the surface area necessary for adequate sample loading




# Addressing challenges – pressure and particle size

Proper particle size is key



- 4 μm SPPs provide speed and efficiency advantage over 5 μm TPPs
- SPP benefits without the high backpressure generated from a smaller particle size

## Van Deemter plots of 4 μm SPP and 5 μm TPP columns



# Addressing challenges – system volume


Agilent 1290 Infinity II autoscale preparative LC system

5-65% ACN gradient

- Efficiency benefits of 4  $\mu\text{m}$  SPP preparative columns can be seen on systems with standard configuration/plumbing
- **Note:** regardless of particle type, proper consideration should be given to tubing size, based on the diameter of the column

## Paraben Mix (in DMSO):


- Caffeine
- Methyl Paraben
- Ethyl Paraben
- Propyl Paraben
- Benzyl Paraben



Flow rate: Agilent – 25 mL/min; TPP – 16 mL/min; Mobile phase: A – water + 0.1% formic acid; B – acetonitrile + 0.1% formic acid; gradient: 5 – 65% B in 18 min; injection volume: Agilent – 80  $\mu\text{L}$ ; TPP – 64  $\mu\text{L}$ ; Detection: UV 254 nm.

# Addressing challenges – sample loading

- Stationary phases with higher surface areas will generally have more retention (and more separation between peaks) than columns with lower surface areas
- Higher surface areas (~400 m<sup>2</sup>/g) are used for additional loading



Phase	Surface Area (m <sup>2</sup> /g)	Pore Size (Å)
Poroshell SB-C18 (4 μm SPP)	130	120
Vendor X C18 (5 μm TPP)	185	135
Vendor Y C18 (5 μm SPP)	200	100
Vendor Y C18	<b>400</b>	100
Pursuit XRs C 18	<b>440</b>	100



Higher surface area phases are ideal for loading




- Lower surface area phases offer fast separations
- Smaller particle sizes offer higher efficiency to resolve complicated matrices

# Loading on InfinityLab Poroshell 120 SB-C18


**Sample:** Sulfanilamide (A) + Sulfamethoxazole (B) in 50/50 acetonitrile/water

- Sample 1: 30 mg/mL A + 30 mg/mL B
- Sample 2: 60 mg/mL A + 60 mg/mL B
- Sample 3: 120 mg/mL A + 120 mg/mL B

<b>Preparative instrument</b>	Agilent 1290 Infinity II preparative LC system
<b>Preparative column</b>	InfinityLab Poroshell 120 SB-C18, 21.2 x 150 mm, 4 $\mu$ m
<b>Flow rate</b>	25 mL/min
<b>Mobile phase</b>	55/45 acetonitrile+ 0.1 formic acid/water + 0.1% formic acid
<b>Injection volume</b>	900 $\mu$ L
<b>Wavelength</b>	238 nm

  
**High Throughput**

Discovery labs typically purify between 10 – 100 mg at a time




**Peaks still retain shape with over 200 mg on column**

# Mass loading comparison at fast flow rates

## Analyte - toluene


- Comparing peak widths at 10% peak height
- Adjusted flow rate and mobile phase composition on TPP C18 to match  $t_0$  and  $k'$  of toluene on Poroshell column
- Injection volume – 400  $\mu\text{L}$ 
  - Determined by volume loading study (not shown) on Poroshell 120 column
- Toluene sample concentrations in 50/50 ACN/DMSO (mg/mL)
  - 25, 50, 100, 150, 200, 250, 300, 350



Mass on column – 4 mg to 140 mg

# Mass loading results and ROI

## Peak width at 10% height vs mass on column



## For 1,000 sample campaign:

	TPP Column	InfinityLab Poroshell 120
Number of samples	1,000	1,000
Number of injections per sample	3	1
Number of injections per campaign	3,000	1,000
Total time (analysis time x number of injections)	175 hours	58 hours
Total lab overhead (lab overhead x total time)	\$21,875	\$7,292
Total scientist cost (scientist cost x total time)	\$26,250	\$8,750
Total consumables cost (consumables cost x no of injections)	\$62,621	\$23,361
Total cost per campaign	\$110,746	\$39,402
Purification cost per sample	\$111/sample	\$39/sample
<b>Cost savings</b>	<b>\$71,344 or \$71/sample</b>	
<b>% savings</b>	<b>64%</b>	

**InfinityLab Poroshell column delivers a 64% savings and 728% ROI over TPP column**

		TPP Column	InfinityLab Poroshell 120
Return on Investment	Total cost per campaign + net cost of change	\$110,746	\$48,022
	Agilent column saves		\$62,724
	<b>ROI% (savings/net cost to change)</b>		<b>728%</b>

# Integrating SPP columns in purification workflow

Driving throughput

## High Throughput



### What's important?

- Collecting enough material (10-100 mg)
- High purity
- Minimal method development

### Ways to improve throughput:

- Screening (step 2) – analytical columns
- Size (step 5) - preparative columns



# Step 2: Improving screening by leveraging $L/d_p$

1. Check sample formulation

**2. Screen column phases**

3. Optimize separation

4. Determine max sample load

5. Determine prep column size

6. Scale conditions to prep column


7. Purify the compound(s)

- Preferred columns: 2.1 x 50 mm, sub-2  $\mu\text{m}$  or smallest available to further improve throughput
- Small columns and fast gradients (2-3 min) are combined to minimize method development time
- Methods can be aggressive - The goal is to identify a column/mobile phase combination that will have the highest chance of separation, not a final separation
- Columns within 10  $L/d_p$  units have similar resolving power




High Throughput

Screen with sub-2  $\mu\text{m}$  columns and then scale with  $L/d_p$  (Length/particle diameter)



# Step 5: high throughput purification



## 5. Determine prep column size

- Instead, a column is chosen based on anticipated load. The most popular column IDs for discovery work are 21.2 mm and 30 mm ID. **If you have sufficient loading on a smaller column, do you really need a large column?**

Column ID	Flow Rate	Difficult Separation ( $1.2 < \alpha < 1.5$ )
21.2 mm	20 mL/min	20 – 70 mg
30 mm	40 mL/min	40 – 140 mg

### Using a smaller column:

- **Reduces column costs (smaller columns are less expensive)**
- **Reduces solvent consumption and waste disposal**
- **Reduces fraction dry down time**

# Summary

High throughput environments can benefit from analytical and preparative columns that:



- can resolve samples in complex matrices
- can be run aggressively while still maintaining performance
- have the loading capability to purify required amounts of sample

With little to no change to instrumentation

Columns packed with superficially porous particles:


- Have higher efficiency than their TPP counterparts
- Maintain efficiency at high flow rates
- Meet the loading needs of discovery and high throughput environments

Using conventional analytical and preparative hardware



# Resources for support

- **Product webpage:** [www.agilent.com/chem/prepcolumns](http://www.agilent.com/chem/prepcolumns)
- **Product brochure:** InfinityLab Poroshell120 Preparative LC Columns [5994-3601EN](#)
- **Primer:** Principles and Practical Aspects of Preparative Liquid Chromatography [5994-1016EN](#)
  - Great overview that covers fundamentals and best practices
- **Application note:** Developing Fast Purification Methods [5994-3518EN](#)
- **Part numbers:**
  - [670050-702](#): InfinityLab Poroshell 120 HPH-C18, 21.2 x 50 mm, 4  $\mu$ m
  - [670050-902](#): InfinityLab Poroshell 120 SB-C18, 21.2 x 50 mm, 4  $\mu$ m
  - [670150-702](#): InfinityLab Poroshell 120 HPH-C18, 21.2 x 150 mm, 4  $\mu$ m
  - [670150-902](#): InfinityLab Poroshell 120 SB-C18, 21.2 x 150 mm, 4  $\mu$ m



# Agilent InfinityLab

